

DRG[®] Vitamin Biotin (BIO-4884)

Revised 6 Aug. 2010 rm (Vers. 3.1)

USA: **RUO**

Please use only the valid version of the package insert provided with the kit.

INTENDED USE

The *Vitamin Biotin* Kit is a microtiter plate test kit based on a microbiological assay which measures the total biotin content in serum.

The test kit contains all required reagents, e.g. standard, medium and microtiter plate coated with a specific microorganism, sufficient for 96 determinations including standard curves. An ELISA reader is required for evaluation of the biotin content.

For *in vitro* use only.

INTRODUCTION

In the late 1950s the discovery of covalently bound biotin with a coenzyme function was reported. Biotin is a growth factor present in minute amounts in every living cell. The fact that humans have a requirement for biotin has been most clearly shown in 2 situations that result in biotin deficiency: 1) prolonged consumption of raw egg whites and 2) parenteral nutrition without biotin supplementation in patients with short-gut syndrome.

Biotin plays an indispensable role in numerous naturally occurring carboxylation reactions. The daily biotin requirements for adults lie between 100 and 200 mg. Renewed interest in the role of biotin in human nutrition and therapy stems from the accumulating reports of various biotin-responsive syndromes.

Indications for a determination of biotin

- Defects in enzymes (e.g. genetic deficiency of biotinidase)
- Short gut syndrome
- Malnutrition

PRINCIPLE OF THE TEST

Serum samples are diluted and added into the microtiter plate wells coated with *Lactobacillus plantarum* which metabolizes biotin. The presence of biotin both in standards [STD] and samples gives a biotin-dependent growth response. After incubation at 37°C for 48 h, the growth of *Lactobacillus plantarum* is measured turbidimetrically at 610 - 630 nm (alternative at 540 - 550 nm) in an ELISA-reader. A dose response curve of absorbance unit (optical density, OD at 610 nm) vs. concentration is generated using the values obtained from standard. Biotin present in the samples is determined directly from this curve.

DRG® Vitamin Biotin (BIO-4884)

Revised 6 Aug. 2010 rm (Vers. 3.1)

USA: **RUO**

MATERIAL SUPPLIED

Label	Kit Components	Quantity
PLATE	One <i>Lactobacillus plantarum</i> -precoated microtiter plate, ready to use	12 x 8 wells
DIL	Water 30 mL	4 x
ASYMED	Biotin assay medium	4 x
STD	Biotin Standard	4 x
CTRL	Biotin Control	4 x
FOL	Cover plastic foil	4 x
FRA	Replacement holder for 96-well plates	1 x

MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber, 37 °C
- Water bath (90°C - 100°C)
- ELISA-Reader 610 - 630 nm (540 - 550 nm)
- Micropipette 20 - 200 µL
- Micropipette 100 -1000 µL
- Micropipette tips to deliver 20 - 200 µL and 100 -1000 µL, sterile
- Pipettes of 5 and 10 mL
- 1.5 - 2 mL reaction vials, sterile
- 0.2 µm sterile polyethersulfon filter with a sterile tip
- 15 mL centrifugal tubes, sterile (e.g. Falcon tubes)

PREPARATION AND STORAGE OF REAGENTS

- Store test kit / reagents at 2-8°C.
- Prepare reagents freshly and use immediately after preparation. Discard remaining unused reagents and waste in accordance with country, federal, state, and local regulations.
- Put unused reagents (standard, medium) in the test kit and store at 2-8°C.
- Take as many microtiter strips as needed from kit. Store unused strips in the original package bag at 2-8° C to prevent contamination or moisture exposure.
- No warranty can be given after the expiry date (see label of test package).
- To run assay more than once, ensure that reagents are stored at conditions stated on the label. Prepare only the appropriate amount necessary for each assay. The kit can be used up to 3 times within the expiry date stated on the label.

DRG[®] Vitamin Biotin (BIO-4884)

Revised 6 Aug. 2010 rm (Vers. 3.1)

USA: **RUO****PRECAUTIONS**

- As the test is based on a microbiological method, the general guidelines for sterile work must be observed as far as possible, (work in a sterile bench, PCR-Hood, use of sterile instruments or equipment).
- Water quality is extremely important. Only the water delivered with the test kit should be used for medium dilution, standard reconstitution as well as for sample preparation.
- It is essential to run a standard curve for each separate assay.
- It is recommended to run a duplicate standard curve as well as a sample analysis.
- If a higher dilution results in a higher measured value, inhibitors like antibiotics might be present.
- Reagents should not be used beyond the expiration date shown on kit label.
- By finishing the test, the used microtiter plates [PLATE] should be autoclaved.
- Signs for reagent damage: The highest standard should have an absorption higher than 0.6 Extinction units ($A_{630nm} > 0,6$)

SAMPLE PREPARATION**Notes**

- Donor serum is used for analysis.
- Original samples should be kept light-protected at 2–8°C until measurement. The samples are stable for 1 day at 2-8°C in the dark. For longer storage samples should be frozen at -20°C.
- Hemolytic samples may give erroneous results and should not be used for analysis. Lipemic samples should be centrifuged at 13 000 x g before assaying.
- Samples with visible amounts of precipitates should be centrifuged (5 min at 10000 g) prior to measurement and the resulting supernatant should be used in the test.

Sample dilution

Serum samples and control [CTRL] should be diluted **1:20** (= dilution factor) with water [DIL] from the kit prior to analysis:

50 µL sample + 950 µL water [DIL]

ASSAY PROCEDURE**Procedural notes**

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test.
- The assay should always be performed according the enclosed manual.

Test preparations

Take as many microtiter strips as needed from kit. Put unused strips in the original package bag, and return the remaining parts of the test kit to the refrigerator. Bring all necessary reagents to room temperature.

Water [DIL] for medium [ASYMED], Standard [STD], Control [CTRL] and sample

Push the lid up, pull it back to the rim of the glass and then remove the entire seal by turning.

Assay medium

The medium must be freshly prepared before the test.

- Take the dry bag out of medium vial [ASYMED] by tweezers, shake off and discard.
- Add 10 mL of water to the assay medium, securely close the bottle and shake well.
The amount is sufficient for 6 strips.
- Heat the bottle with medium in a water-bath at 95 °C for 5 min, while shaking well at least twice. It is important to make sure that the medium bottle is firmly closed at all times.
- Quickly cool the medium bottle to under 30 °C.
- Filter the medium sterilely with a 0.2 µm filter in a 15 mL centrifuge test tube.

Standard [STD]

Before the test freshly prepare the standard curve solutions:

- Open the bottle of standard [STD], place the screw-top lid upside-down on the work bench.
- Add x mL (x = see QS test kit data sheet) water [DIL] from the test kit to the standard bottle [STD], close the bottle and shake (= standard concentrate).
- Add water [DIL] into 6 sterile reaction vials (capacity 1.5 – 2.0 mL) and then pipet the standard concentrate to the vials. Prepare a standard curve using the following scheme:

DRG[®] Vitamin Biotin (BIO-4884)

Revised 6 Aug. 2010 rm (Vers. 3.1)

USA: **RUO**

Biotin [µg / l]	Water [DIL] [µL]	+	Standard [STD] [µL]	=	Total volume [µL]
Blank: 0	500	+	0	=	500
Standard 1: 0.02	450	+	50	=	500
Standard 2: 0.06	350	+	150	=	500
Standard 3: 0.10	250	+	250	=	500
Standard 4: 0.14	150	+	350	=	500
Standard 5: 0.18	50	+	450	=	500

Control [CTRL]

- The control must be freshly prepared before the test.
- Open the bottle of control [CTRL], remove seal. Dispose of screw-top lid and seal.
- Add 0.2 mL water [DIL] from the test kit to the control bottle [CTRL], close the bottle and dissolve by vortexing the bottle (= control).
- Treat the control afterwards as the sample is treated.
- Pipette 150 µL of the diluted control [CTRL] into each well. We recommend to run a duplicate.
- For the concentration of the Control [CTRL] please see Control specification.

Test Initiation

- Take as many microtiter strips as needed from the kit and put them in the second microtiter strip holder [Holder]. Store unused strips in the original package bag at 2-8° C to prevent contamination or moisture exposure.
- A medium solution is sufficient for 6 strips.
- Put 150 µL biotin assay medium [ASYMED] in the cavities.
- Add 150 µL standard [STD], control [CTRL], respectively, sample in the cavities. Pre-rinse the pipette tip with standard and sample solution respectively.
- Carefully seal the cavities with plastic foil [Cover]. Important: the cavities must be made airtight by pressing down with the hand!
- Keep at **37 °C for 48 hrs** in an incubator.

Measurement

- Securely press the foil [Cover] down with the hand.
- Upturn the plate [PLATE] onto a tabletop and shake the germination well.
- Turn the plate over again and carefully remove the foil [Cover], beginning with the lower, left corner and pulling diagonally backwards at an angle of 180°.
- Remove air bubbles in the cavities using a pipette tip or a needle.
- Read turbidity in an ELISA-Reader at E 610 - 630 nm (alternatively at 540 - 550 nm)

DRG[®] Vitamin Biotin (BIO-4884)

Revised 6 Aug. 2010 rm (Vers. 3.1)

USA: **RUO****Please note**

- After 48 hrs incubation time, the microtiter platter may be stored for a maximum of 48 hrs in the refrigerator before measuring the turbidity.
- To prevent time-loss through public holidays or weekends, the microtiter plate may also be evaluated after 60 hrs incubation.

REFERENCES / LITERATURE

1. Burtis C A, Ashwood E R (Eds): Tietz Textbook of Clinical Chemistry, 3rd Edition, 1999
2. Coronel F et al. (1991) Treatment of hyperlipemia in diabetic patients on dialysis with a physiological substance. Am J Nephrol 11(1):32-6

GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- For *in vitro* use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the results of the test. DRG can, therefore, not be held reliable for any damage resulting from this.