3 α -Diol G-RIA-CT

KIP0151

Read entire protocol before use.

3 α -Diol G-RIA-CT

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 5 α -Androstane-3 α -17 β -Diol-Glucuronide (3 α -Diol G) in serum and plasma.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource 3 α -Diol G-RIA-CT Kit
- B. **Catalog number :** KIP0151 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

5 α -Androstane-3 α -17 β -Diol-Glucuronide (3 α -Diol G) is a C19 steroid. It is produced mainly as a metabolite of testosterone and dihydrotestosterone (DHT). It is largely produced in target peripheral tissues such as the skin, especially around hair follicles. The stimulation by large amounts of the 3 α -Diol G, leads to excessive hair formation – notably conspicuous where hair is not normally present in women.

In recent years the interest in the measurement of this steroid has increased among clinical investigators studying women suffering from idiopathic hirsutism.

Among the steroids known to be precursors for 3 α -Diol G, namely dehydroepiandrosterone (DHEA), dehydroepiandrosterone sulphate (DHEAS), DHT androstenedione and testosterone, only 3 α -Diol G has been shown to increase with hirsutism and decrease with treatment. This correlation has also been demonstrated in patients with polycystic ovaries (PCO).

3 α -Diol G determinations have therefore proved to be useful as an indicator in a variety of ways including monitoring the progress of treatment of idiopathic hirsutism and women with polycystic ovaries.

Furthermore diabetic patients, both men and women under cyclosporine A therapy, have shown 3 α -Diol G levels above normal, a side effect resulting in the appearance of hair in previously hairless areas.

IV. PRINCIPLES OF THE METHOD

A fixed amount of ^{125}I labelled 3 α -Diol G competes with the 3 α -Diol G to be measured present in the sample or in the calibrator for a fixed amount of antibody sites being immobilized to the wall of a polystyrene tube. Neither extraction nor chromatography are required because of the high specificity of the coated antibodies. After a 2 hours incubation at RT on a shaker, an aspiration step terminates the competition reaction. The tubes are then washed twice with 2 ml of Working Wash Solution and aspirated again. A calibration curve is plotted and the 3 α -Diol G concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution
 Tubes coated with anti 3 α -Diol G	2 x 48	grey	Ready for use
Ag ^{125}I TRACER: ^{125}I Iodine labelled 3 α -Diol G (HPLC grade) in buffer with bovine casein (1%) and azide (<0.1%)	1 vial 53 ml 185 kBq	Red	Ready for use
CAL 0 Zero calibrator in bovine serum, thymol (<0.1%) and gentamycine (< 0.1%)	1 vial lyophilized	yellow	Add 2 ml distilled water
CAL N Calibrators 3 α -Diol G - N = 1 to 5 (see exact values on vial labels) in bovine serum, thymol (<0.1%) and gentamycine (<0.1%)	5 vials lyophilized	yellow	Add 0.5 ml distilled water
WASH SOLN CONC Wash solution (Phosphate Buffer)	1 vial 50 ml	brown	Dilute 25 x with distilled water (use a magnetic stirrer).
CONTROL N Controls - N = 1 or 2 in human serum and thymol (<0.1%) .	2 vials lyophilized	silver	Add 0.5 ml distilled water

Note : Use the zero calibrator for sample dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 100 μl , 500 μl and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Tubes shaker (700 rpm)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators** : Reconstitute the zero calibrator with 2.0 ml distilled water and other calibrators with 0.5 ml distilled water.
- Controls** : Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 24 volumes of distilled water to 1 volume of Wash Solution (25x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 24 hrs., storage at -20°C is recommended.
- Avoid successive freezing and thawing.
- Heparinized plasma yields 15 % lower results than serum :

$$Y (\text{Hep. plasma}) = 0.86 \times (\text{serum}) + 0.52 \quad r = 0.96 \quad n = 13$$
- EDTA plasma yields 25 % lower results than serum :

$$Y (\text{EDTA plasma}) = 0.74 \times (\text{serum}) + 0.52 \quad r = 0.97 \quad n = 13$$

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a standard curve for each run, do not use data from previous runs.

B. Procedure

- Label coated tubes in duplicate for each calibrator, sample, control. For the determination of total counts, label 2 normal tubes
- Briefly vortex calibrators, samples and controls and dispense 100 μl of each into respective tubes.
- Dispense 0.5 ml of ^{125}I Iodine labelled 3 α -Diol G into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently.
- Incubate for 2 hours at room temperature with continuous shaking (700 rpm).
- Aspirate the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate. Avoid foaming during the addition of the Working Wash solution.
- Wash tubes again with 2 ml Working Wash solution (except total counts) and aspirate.
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations, reject obvious outliers.
- Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

- Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B₀(%)) values for each calibrator point as a function of the 3 α -Diol G concentration of each calibrator point.
- Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is to be used with results, "4 parameters" function curve fitting is recommended.
- By interpolation of the sample (B/B₀(%)) values, determine the 3 α -Diol G concentrations of the samples from the reference curve.
- For each assay, the percentage of total tracer bound in the absence of unlabelled 3 α -Diol G (B₀/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

3α-Diol G	cpm	B/Bo (%)
Total count	35967	
Calibrator		
0.0 ng/ml	13912	100.0
0.2 ng/ml	11443	82.2
1.2 ng/ml	8385	60.3
6.0 ng/ml	4611	33.1
25.0 ng/ml	1796	12.9
75.0 ng/ml	960	6.9

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit
Twenty zero calibrators were assayed along with a set of other calibrators.
The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 0.05 ng/ml.

B. Specificity

Compound	Cross-Reactivity (%)
5α-Androstane-3α-17β-diolG	100.00
5α-Androstane-3α-17β-diol	10.69
5α-Androstane-3α-17β-diol-3Glucuronide	5.86
5α-Dihydrotestosterone Glucuronide	1.75
Progesterone	0.03
Testosterone glucuronide	0
Testosterone	0
11β-hydroxytestosterone	0
5α-Dihydrotestosterone	0
5β-Dihydrotestosterone	0
Cortisol	0
Dehydroepiandrosterone	0
Estrone	0
Androstenedione	0
5-Androstene-3β-17β-diol	0
17β-Estradiol	0

Note : this table shows the cross-reactivity for the anti 3α-Diol G.

C. Precision

INTRA-ASSAY PRECISION				INTER-ASSAY PRECISION			
Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A	20	2.64 ± 0.15	5.7	A	10	2.76 ± 0.17	6.4
B	20	10.11 ± 0.50	4.9	B	10	10.32 ± 0.70	7.2

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST			
Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum	1/1	50.70	50.70
	1/2	25.35	24.58
	1/4	12.67	11.96
	1/8	6.33	6.49
	1/16	3.16	3.26
	1/32	1.58	1.48
	1/64	0.79	0.62

Samples were diluted with the zero calibrator.

RECOVERY TEST			
Sample	added 3α-Diol G (ng/ml)	Recovered 3α-Diol G (ng/ml)	Recovered (%)
Serum	50	48.4	96.8
	20	21.3	106.0
	10	9.5	95.0
	5	4.9	98.0

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 20 minutes after the calibrator has been added to coated tubes.

TIME DELAY		
Serum (ng/ml)	0'	20'
Serum 1	2.92	2.90
Serum 2	9.60	9.59

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- To the best of our knowledge, no international reference material exists for this parameter.

XV. REFERENCE INTERVALS

Population		Absolute range (ng/ml)	Median (ng/ml)
Female	Premenopausal	0.3 - 7.9	1.9
	Postmenopausal	0.1 - 5.9	1.5
	Hirsute	1.6 - 9.3	4.6
Male		1.0 - 23.6	6.4

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.
This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves. All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

XVII. BIBLIOGRAPHY

1. RITTMASER R.S. *et al.* (1988)
Androstenediol Glucuronide Isomers in Normal Men and Women and in Men infused with labelled Dihydrotestosterone.
JCE & M., 66, 1, 212-216.
2. RAO P.N. *et al.* (1987)
Isolation and Identification of Androstenediol Glucuronide from Human Plasma.
J. Steroid Biochem., 28, 5, 565-569.
3. HORTON R. *et al.* (1982)
3 α ,17 β -Androstenediol Glucuronide in plasma.
J. Clin. Invest., 69, 1203-1206.
4. JOURA E.A. *et al.* (1996)
Serum 3 α -androstenediol glucuronide is decreased in nonhirsute women with acne vulgaris.
Fertility and Sterility, 66, 6, 1033-1035.
5. GILAD S. *et al.* (1994)
Assesment of 5 α -reductase activity in hirsute women: comparison of serum androstenediol glucuronide with urinary androsterone and aetiocholanolone excretion.
Clinical Endocrinology, 40, 459-464.
6. THOMPSON D.L. *et al.* (1990)
Androsterone Glucuronide is a marker of adrenal hyperandrogenism in hirsute women.
Clinical Endocrinology, 32, 283-292.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS μl	CALIBRATORS μl	SAMPLE(S) CONTROLS μl
Calibrators (0 to 5)	-	100	-
Samples, controls	-	-	100
Tracer	500	500	500
Incubation	2 hours at RT with continuous shaking (700 rpm)		
Separation	-	aspirate	
Working Wash solution		2.0 ml	
Separation		aspirate	
Working Wash solution		2.0 ml	
Separation		aspirate	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP0151	P.I. Number : 1701001/en	Revision nr : 110608/1
-------------------------------------	-----------------------------	---------------------------

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XVII. LITERATUR

1.

RITTMAS­TER R.S. *et al.* (1988)
Androstanediol Glucuronide isomers in normal men and women and in men infused with labelled Dihydrotestosterone.
JCE & M., 66, 1, 212-216.

2.

RAO P.N. *et al.* (1987)
Isolation and identification of Androstanediol Glucuronide from human plasma.
J. Steroid Biochem., 28, 5, 565-569.

3.

HORTON R. *et al.* (1982)
3 α -17 β -Androstanediol Glucuronide in plasma.
J. Clin. Invest., 69, 1203-1206.

4.

JOURA E.A. *et al.* (1996)
Serum 3 α -androstanediol glucuronide is decreased in nonhirsute women with acne vulgaris.
Fertility and Sterility, 66, 6, 1033-1035.

5.

GILAD S. *et al.* (1994)
Assessment of 5 α -reductase activity in hirsute women: comparison of serum androstanediol glucuronide with urinary androsterone and aetiocholanolone excretion.
Clinical Endocrinology, 40, 459-464.

6.

THOMPSON D.L. *et al.* (1990)
Androsterone Glucuronide is a marker of adrenal hyperandrogenism in hirsute women.
Clinical Endocrinology, 32, 283-292.

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1.

RITTMASER R.S. *et al.* (1988)

Androstanediol Glucuronide Isomers in Normal Men and Women and in Men infused with labelled Dihydrotestosterone.

JCE & M., 66, 1, 212-216.

2.

RAO P.N. *et al.* (1987)

Isolation and Identification of Androstanediol Glucuronide from Human Plasma.

J. Steroid Biochem., 28, 5, 565-569.

3.

HORTON R. *et al.* (1982)

3 α ,17 β -Androstanediol Glucuronide in plasma.

J. Clin. Invest., 69, 1203-1206.

4.

JOURA E.A. *et al.* (1996)

Serum 3 α -androstanediol glucuronide is decreased in nonhirsute women with acne vulgaris.

Fertility and Sterility, 66, 6, 1033-1035.

5.

GILAD S. *et al.* (1994)

Assesment of 5 α -reductase activity in hirsute women: comparison of serum androstanediol glucuronide with urinary androsterone and aetiocholanolone excretion.

Clinical Endocrinology, 40, 459-464.

6.

THOMPSON D.L. *et al.* (1990)

Androsterone Glucuronide is a marker of adrenal hyperandrogenism in hirsute women.

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1. RITTMASER R.S. *et al.* (1988)
Androstanediol Glucuronide Isomers in Normal Men and Women and in Men infused with labelled Dihydrotestosterone.
JCE & M., 66, 1, 212-216.

2. RAO P.N. *et al.* (1987)
Isolation and Identification of Androstanediol Glucuronide from Human Plasma.
J. Steroid Biochem., 28, 5, 565-569.
J. Clin. Invest., 69, 1203-1206.

3. HORTON R. *et al.* (1982)
3α,17β-Androstanediol Glucuronide in plasma.

4. JOURA E.A. *et al.* (1996)
Serum 3α-androstanediol glucuronide is decreased in nonhirsute women with acne vulgaris.
Fertility and Sterility, 66, 6, 1033-1035.

5. GILAD S. *et al.* (1994)
Assesment of 5α-reductase activity in hirsute women: comparison of serum androstanediol glucuronide with urinary androsterone and aetiocholanolone excretion.
Clinical Endocrinology, 40, 459-464.

6. THOMPSON D.L. *et al.* (1990)
Androsterone Glucuronide is a marker of adrenal hyperandrogenism in hirsute women.
Clinical Endocrinology, 32, 283-292.

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	<u>Used symbols</u>
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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