CORTISOL-RIA-CT

KIPI28000

Read entire protocol before use.

CORTISOL-RIA-CT

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human Cortisol in serum, plasma or urine without extraction.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource CORTISOL-RIA-CT Kit
- B. Catalog number :** KIP128000 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological activity

Cortisol is the major glucocorticoid produced and secreted by the adrenal gland. In response to different stimuli (diurnal rhythm, stress, low blood sugar concentration), the cerebral cortex stimulates the hypothalamus to release the CRF (corticotrophin releasing factor). CRF causes the release from the pituitary gland of ACTH (adrenocorticotrophic hormone). Glucocorticoids are then synthesized in response to ACTH.

Since Cortisol levels depend upon the interaction of the hypothalamus, pituitary and adrenal glands, measurement of urine cortisol levels can also aid in the differential diagnosis of the diseases states of these glands.

IV. PRINCIPLES OF THE METHOD

A fixed amount of ¹²⁵I labelled steroid competes with the steroid to be measured present in the sample or in the calibrator for a fixed amount of antibody sites being immobilized to the wall of a polystyrene tube. Neither extraction nor chromatography are required because of the high specificity of the coated antibodies. After 45 minutes incubation at 37°C, an aspiration step terminates the competition reaction. The tubes are then washed with 3 ml of wash solution and aspirated again. A calibration curve is plotted and the cortisol concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with anti cortisol (monoclonal)	2 x 48	yellow	Ready for use			
<table border="1"><tr><td>Ag</td><td>¹²⁵I</td></tr></table> TRACER: ¹²⁵ Iodine labelled cortisol in phosphate-citrate buffer with proteins, ANS, and azide (<0.1%)	Ag	¹²⁵ I	1 vial 52 ml 116 kBq	red	Ready for use	
Ag	¹²⁵ I					
<table border="1"><tr><td>CAL</td><td>0</td></tr></table> Zero Calibrator in human serum and azide (0.1%)	CAL	0	1 vial 1 ml	yellow	Ready for use	
CAL	0					
<table border="1"><tr><td>CAL</td><td>N</td></tr></table> Calibrators - N = 1 to 5 (see exact values on vial labels) in human serum and azide (0.1%)	CAL	N	5 vials 0.5 ml	yellow	Ready for use	
CAL	N					
<table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table> Wash solution (TRIS-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1"><tr><td>CONTROL</td><td>N</td></tr></table> Controls – 1 and 2 in human plasma with azide (<0.1%)	CONTROL	N	2 vials lyophilised	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note : Use the zero calibrator for sera dilutions. If needed, urines can be diluted in an ordinary phosphate saline buffer (PBS).

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- 1. Distilled water
- 2. Pipettes for delivery of: 25 µl and 500 µl (the use of accurate pipettes with disposable plastic tips is recommended)
- 3. Vortex mixer
- 4. Magnetic stirrer
- 5. Water bath at 37°C
- 6. 5 ml automatic syringe (Cornwall type) for washing
- 7. Aspiration system (optional)
- 8. Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. **Controls:** Reconstitute the controls with 0.5 ml distilled water.
- B. **Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, controls are stable for 7 days at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.

- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Keep serum, plasma or urine at 2-8°C for 1-2 days. Keep frozen for longer periods. Highly lipemic or hemolyzed samples must be discarded. Presence of fibrin filaments in the plasma can interfere with the assay. Use clear samples only
- Collect urine samples during 24 hours without preservative. Record the total volume.
- Avoid subsequent freeze-thaw cycles.

X. PROCEDURE

- A. **Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

- B. **Procedure**
 1. Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes
 2. Briefly vortex calibrators, controls and samples and dispense 25 µl of each into the respective tubes.
 3. Dispense 500 µl of ¹²⁵Iodine labelled Cortisol into each tube, including the uncoated tubes for total counts.
 4. Shake the tube rack gently by hand to liberate any trapped air bubbles.
 5. Incubate for 45 minutes at 37° +/- 2°C, in a water bath.
 6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
 7. Wash tubes with 3 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
 8. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
 9. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B0 (\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

3. Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B0(%)) values for each calibrator point as a function of the cortisol concentration of each calibrator point. Reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B0 (%)) values, determine the cortisol concentrations of the samples from the calibration curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled cortisol (B0/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Cortisol	cpm	B/Bo (%)
Total count	63729	
Calibrator		
0 µg/L	35476	100.0
10 µg/L	31047	87.5
30 µg/L	25851	72.9
100 µg/L	17182	48.4
300 µg/L	9399	26.5
1000 µg/L	4050	11.4

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was **0.9 µg/L**.

B. Specificity

The percentage of cross-reaction estimated by comparison of the concentration yielding a 50% inhibition are respectively:

Compound	Cross-Reactivity (%)
Cortisol	100.00
17-OH-Progesterone	5.60
DOC	0.74
Corticosterone	0.51
11-DOC	0.27
Prednisone	0.25
Estradiol	0.17
Prednisolone	0.10
Testosterone	0.08
Progesterone	0.00

Note: this table shows the cross-reactivity for the anti cortisol

C. Precision

INTRA-ASSAY PRECISION				INTER-ASSAY PRECISION			
Serum	N	Mean (µg/L)	CV (%)	Serum	N	Mean (µg/L)	CV (%)
A	5	1.10	6.2	A	8	1.95	8.7
B	5	3.70	5.2	B	8	5.48	11.5
C	5	28.30	7.7	C	8	36.85	15.1

D. Accuracy

DILUTION TEST		
Dilution (serum)	Theoretical Concent. (µg/L)	Measured Concent. (µg/L)
1/1	-	630.3
1/2	315.2	298.5
1/4	157.6	159.3
1/8	78.8	81.3
1/16	39.4	39.1
1/32	19.7	21.7
1/64	9.9	11.1

Samples were diluted with zero calibrator.

RECOVERY TEST			
Sample (µg/L)	Added Cortisol (µg/L)	Recovered Cortisol (µg/L)	Recovered (%)
3.0	30.0	17.7	93.2
3.0	100.0	48.7	105.7
3.0	300.0	149.2	101.5
3.0	1000.0	509.2	98.5

Conversion factor :
From µg/L to nmol/L : x 2.758
From nmol/L to µg/L : x 0.363

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Cortisol levels vary diurnally and as a function of clinical suppression manoeuvres. Morning samples results are generally 2 to 3 times higher than these obtained with afternoon samples. Morning samples generally read less than 5 µg/L following metyrapone testing or dexamethasone suppression.

Identification	Range (*) (µg/L)
Serum samples	
Morning	50 – 250 µg/L
Afternoon	25 -125 µg/L
24 h urine samples	7 – 96 µg/L

(*) The range is based on 2.5 % and 97.5 % percentiles

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.
This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.
This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.
All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.
Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.
The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.
All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.
Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.
Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. Abraham, GE, Ed. (1981)
Radioassay Systems in Clinical Endocrinology.
Marcel Dekker Inc, New York
2. Rittler M and Geisler D (1983)
Radioimmunoassay of urinary cortisol: comparison of the procedures with and without previous cortisol extraction.
Meth. Find. Exptl. Clin. Pharmacol. 5(6), 403-406.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE (S) CONTROLS µl
Calibrators (0 to 5) Samples, Controls Tracer	- - 500	25 - 500	- 25 500
Incubation	45 minutes at 37 +/- 2°C, in a water bath		
/Separation Working Wash solution Separation	-	Aspirate (or decant) 3.0 ml Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP128000	P.I. Number : 1700435/en	Revision nr : 150324/3
---------------------------------------	-----------------------------	---------------------------

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2. Rittler M and Geisler D (1983)
Radioimmunoassay of urinary cortisol: comparison of the procedures with
and without previous cortisol extraction.
Meth. Find. Exptl. Clin. Pharmacol. 5(6), 403-406.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
PREC AGENT	Precipitating Agent
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

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