ZenTech

# ELIZEN Proinsulin

DEUTSCH (de)

ENGLISH (en)

FRANÇAIS (fr)



ESPAÑOL (es)

**ZenTech s.a.**

Liège Science Park  
Avenue du Pré-Aily 10,  
4031 ANGLEUR, Belgium  
Tel. : + 32-(0)4-361.42.32  
Fax : + 32-(0)4-367.00.63

[info@zentech.be](mailto:info@zentech.be)  
[www.zentech.be](http://www.zentech.be)

| ISO15223                                                                          | SYMBOLE FÜR MEDIZINPRODUKTE          | MEDICAL DEVICES SYMBOL         | SYMBLES APPLIQUES AUX DISPOSITIFS MEDICAUX | SÍMBOLOS PARA APARATOS DE USO EN MEDICINA              |
|-----------------------------------------------------------------------------------|--------------------------------------|--------------------------------|--------------------------------------------|--------------------------------------------------------|
|  | Limitierung der Lagertemperatur      | STORAGE TEMPERATURE LIMITATION | LIMITES DE TEMPERATURES                    | LIMITES DE TEMPERATURA DE ALMACENAMIENTO               |
|  | Chargencode                          | BATCH CODE                     | NUMERO DE LOT                              | CÓDIGO DE LOTE                                         |
|  | Verwendbar bis                       | USE BY                         | DATE D'EXPIRATION                          | CONSUMIR ANTES DE                                      |
|  | Das Handbuch zu Rate ziehen          | CONSULT OPERATING INSTRUCTIONS | LIRE LES INSTRUCTIONS                      | CONSULTAR LAS INSTRUCCIONES DE MANEJO O FUNCIONAMIENTO |
|  | Für in-vitro-diagnostische Anwendung | IN VITRO DIAGNOSTIC DEVICE     | IN VITRO DIAGNOSTIC                        | DISPOSITIVO DE DIAGNÓSTICO IN VITRO                    |
|  | Hergestellt von                      | MANUFACTURED BY                | FABRIQUE PAR                               | FABRICADO POR                                          |
|  | Katalog Nr.                          | CATALOGUE NUMBER               | REFERENCE                                  | NÚMERO DE CATÁLOGO                                     |

|                                                                                     | SYMBOLE (empfohlen von der EDMA )    | SYMBOLS (EDMA recommendations)          | SYMBLES (recommandations EDMA)       | SÍMBOLOS (RECOMENDACIONES DE LA EDMA)          |
|-------------------------------------------------------------------------------------|--------------------------------------|-----------------------------------------|--------------------------------------|------------------------------------------------|
|    | Anzahl der Bestimmungen (96)         | number of determinations (96)           | Nombre de déterminations (96)        | Número de determinaciones                      |
|    | Kalibratoren                         | Calibrators                             | Calibrateurs                         | Calibradores                                   |
|   | Kontollserum                         | Control serum                           | Sérum de contrôle                    | Suero de control                               |
|  | Mikrotiterplatte                     | Microtiterplate                         | Plaque de microtitration             | Placa de microvaloración                       |
|  | Blockierungspuffer                   | Blocking buffer                         | Tampon de blocage                    | Tampón de bloqueo                              |
|  | HRP- Konjugat                        | HRP Conjugate                           | Conjugué HRP                         | HRP conjugado                                  |
|  | Chromogenes Substrat                 | Chromogenic substrate                   | Substrat chromogène                  | Sustrato cromogénico                           |
|  | Waschlösung zum 10-maligen verdünnen | Washing solution to be diluted ten-fold | Solution de lavage à diluer dix fois | Solución de lavado que hay que diluir 10 veces |
|  | Stoplösung                           | Stop solution                           | Solution d'arrêt                     | Solución de parada                             |

# DEUTSCH

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# ENGLISH

## IMMUNOENZYMOMETRIC ASSAY FOR THE QUANTATIVE DETERMINATION OF PROINSULIN IN HUMAN SERUM AND PLASMA.

### E-BX-96 - 96 Assays

#### FOR IN VITRO DIAGNOSTIC USE ONLY

#### 1. CLINICAL APPLICATIONS

Proinsulin (PI), a polypeptide of 9390 MW (86 amino acids), is synthesized in the  $\beta$  cells of the Islets of Langerhans in the pancreas(1). The sequence of proinsulin is highly conserved in mammalian species and is homologous with IGF-1 and IGF-2 (2). This protein is processed to C-peptide and insulin forms (3). Both are secreted in equimolar amounts into the blood. In normal individuals, proinsulin is present in the circulation in very low concentrations (typical basal values 2-6 pmol/L). The level of proinsulin can be a reflection of  $\beta$  cell status and a consequence of dysfunction of PI processing and/or secretion (4). Proinsulin like material is increased in clinical conditions as insulinoma, familial hyperinsulinemia, non-insulin dependent diabetes mellitus (5, 6, 7, 8).

#### 2. PRINCIPLE OF THE ASSAY (I)

Proinsulin kit is a sensitive two-site sandwich enzyme linked immunosorbent assay. The microtiter plates are coated with a monoclonal antibody (S2) specific for epitope at the C-peptide/insuline A chain junction. S2 is able to bind intact PI, des (31,32)-PI and split (32,33)-PI but not insulin, C-peptide and the other "des" and "split" forms.

Prior to use a blocking buffer is added to the allocated wells. An aliquot of patient sample is then added to the wells and incubated. The wells are then washed to remove unbound antibody and other serum compounds. In a second incubation time, an enzyme labelled anti-proinsulin antibody (S53) is incubated in the wells. This antibody is specific for the epitopes at insulin  $\beta$  chain/C-peptide junction. S53 is able to bind to intact PI, des (64,65)-PI and split (65,66)-PI but not insulin, C-peptide and other "des" and "split" forms. The combination of the two Mabs has the ability to detect only the intact human proinsulin.

After washing away any unbound enzyme labelled anti-proinsulin antibody, the enzyme activity is measured by adding a chromogen substrate. The intensity of colour development is proportional to the concentration of proinsulin in the patient sample.

(I) Paule Housa et al. "First direct assay for intact human proinsulin". Clinical Chemistry 44:7, 15414-1519 (1998)

#### 3. REAGENTS PROVIDED WITH THE KIT

- Reagents are sufficient for 96 wells.
- Store the kit at 2-8°C.
- The expiry date of each reagent is shown on the vial label.

**1. Coated microtiterplate :** 96 breakable wells coated with a proinsulin specific antibody . Keep unused wells at 2-8°C, protected from moisture in the provided aluminium bag and carefully sealed.  
Store at 2 – 8°C

**2. Proinsulin 0 calibrator :** 1 vial of lyophilised calibrator. *Reconstitute with 3 ml of distilled water.* Blue coded. After reconstitution, keep the calibrator at –20°C (freeze-defreeze cycle: max 2 times).  
Preservatives : Merthiolate.  
Store lyophilized at 2 – 8°C

**3. Proinsulin 1 - 5 calibrators:** 5 vials of lyophilised calibrator. *Reconstitute with 1 ml of distilled water.* Blue coded. After reconstitution, keep the calibrator at –20°C (freeze-defreeze cycle: max 2 times). For the exact value, refer to the data sheet included.  
**The calibrators are calibrated against the WHO 1<sup>st</sup> International Reference reagent IRR 84/611.**  
Preservatives : Merthiolate.  
Store lyophilized at 2 – 8°C

**4. Control sera :** 2 vials of lyophilised sample. *Reconstitute with 1 ml of distilled water.* Blue coded. After reconstitution, keep the control serum at –20°C (freeze-defreeze cycle: max 2 times)  
For the exact value, refer to the data sheet included.  
Preservative : Merthiolate.  
Store lyophilised at 2 – 8°C.

**5. HRP conjugate:** 1 vial ( 11 ml) of anti-human proinsulin conjugated to horseradish peroxidase (HRPO). Preservatives : Neomycin and Merthiolate.  
Store at 2 – 8°C.

**6. Blocking buffer :** 1 vial (1.5 ml) of murin IgG in phosphate buffer. Ready for use.  
Preservatives : Merthiolate.  
Store at 2 – 8°C.

**7. Washing solution :** 1 vial (40 ml) of buffer with Tween 20. Preservatives: Sulfate Streptomycin and Amphotericin. Bring the vial content to 400 ml (final volume) with distilled water. The diluted washing solution is stable for 6 months at 2 - 8°C.  
Store at 2 – 8°C

**8. Chromogenic substrate:** 2 vials (13ml) of Tetramethylbenzidine in citrate-phosphate buffer and DMSO. Ready for use.  
Store at 2 – 8°C

**9. Stop solution :** 1 vial (15ml) of 0.5M H<sub>2</sub>SO<sub>4</sub>. Ready for use.  
Store at 2 – 8°C

#### 4. KIT REAGENTS

| Reagents              | Quantity   | Physical state    |
|-----------------------|------------|-------------------|
| Wells                 | 96         | Ready for use     |
| Calibrator 0          | 1 x 3 ml   | Lyophilised       |
| Calibrators 1 – 5     | 5 x 1 ml   | Lyophilised       |
| Control Sera          | 2 x 1 ml   | Lyophilised       |
| HRP conjugate         | 1 x 11 ml  | Ready for use     |
| Blocking buffer       | 1 x 1.5 ml | Ready for use     |
| Washing Solution      | 1 x 40 ml  | Concentrated 10 x |
| Chromogenic substrate | 2 x 13 ml  | Ready for use     |
| Stop solution         | 1 x 15 ml  | Ready for use     |

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Adjustable, automatic micropipettes with disposable tips.
- Graduated cylinder.
- Aspiration pump or automated well washing device.
- Microtiterplate spectrophotometer capable of measuring absorbances within a 0-3.0 A interval at 450 and 405 nm.
- Millimetric graph paper.
- Orbital shaker adjustable at 150 rpm.
- Distilled H<sub>2</sub>O.

#### 6. WARNINGS AND PRECAUTIONS

**In order to obtain reproducible results, the following rules must be observed:**

- Do not mix reagents of different lots.
- Do not use reagents beyond their expiry date.
- Use thoroughly clean glassware.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples; for this purpose, disposable tips should be used for each sample and reagent.

**In order to avoid personal and environmental contamination, the following precautions must be observed:**

- Use disposable gloves while handling potentially infectious material and performing the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and

reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be decontaminated and disposed of, in accordance with established safety procedures.

- Disposable ignitable material must be incinerated; disposable non-ignitable material must be sterilized in autoclave for at least 1 hour at 121°C. Liquid wastes must be added with sodium hypochlorite at a final concentration of 3%. Let the hypochlorite act for at least 30 minutes. Liquid wastes containing acid must be neutralized with appropriate amounts of base, before treating with sodium hypochlorite.
  - Chromogen and Blocking Reagent should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
  - Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.
- Some reagents contain sodium azide as preservative; to prevent build-up of explosive metal azides in lead and copper plumbing, reagents should be discarded by flushing the drain with large amounts of water.

## 7. SPECIMEN COLLECTION AND PREPARATION

Both serum and plasma may be used in the assay, although values with serum are slightly lower than plasma. Avoid using hemolyzed, lipemic or bacterially contaminated sera.

Thoroughly mix thawed specimens before assay and avoid repeated freeze/thawing cycles which may cause loss of antibody activity and give erroneous results.

Samples may be stored at 2°-8°C for up to 24 hours. For long term storage, samples should be stored frozen at -20°C or lower.

## 8. ASSAY SCHEME

| Wells                                             | Blank  | Calibrator (0-5) | Control Sera | Samples |
|---------------------------------------------------|--------|------------------|--------------|---------|
| <b>Reagents</b>                                   |        |                  |              |         |
| Blocking buffer                                   | ---    | 10 µl            | 10 µl        | 10 µl   |
| Calibrator (0 - 5)                                | ----   | 100 µl           | ----         | ----    |
| Control Sera                                      | ----   | ----             | 100 µl       | ----    |
| Samples                                           | ----   | ----             | ----         | 100 µl  |
| - Incubate: 30' at R.T., orbital shaker (150 rpm) |        |                  |              |         |
| - Aspirate and wash: 3 X 300 µl.                  |        |                  |              |         |
| Conjugate                                         | -----  | 100 µl           | 100 µl       | 100 µl  |
| - Incubate: 60' at R.T., orbital shaker (150 rpm) |        |                  |              |         |
| - Aspirate and wash: 3 x 300 µl                   |        |                  |              |         |
| Chrom Substrate                                   | 200 µl | 200 µl           | 200 µl       | 200 µl  |
| - Incubate: 15' at R.T., orbital shaker (150 rpm) |        |                  |              |         |
| Stop solution                                     | 100 µl | 100 µl           | 100 µl       | 100 µl  |
| - Read: 450-405 nm and 620 nm as reference.       |        |                  |              |         |

## 9. ASSAY PROCEDURE

- Allow reagents and samples to warm up at room temperature.
- Mix samples by inversion before use.
- The samples are not diluted.

1. Allocate the wells of Microtiterplate for calibrators, control sera and samples
2. Pipette 10 µl of Blocking buffer directly into the bottom of the wells.
3. Pipette 100 µl of each calibrators, control sera and sample into the corresponding wells.
4. Incubate for 30 minutes at room temperature on an orbital shaker (150rpm) (NB. The incubation time begins after the last sample addition)
5. Wash the wells 3 times with 300 µl of diluted Washing solution. Aspirate all liquid from the wells.
6. Add 100µl of HRP conjugate into the wells.
7. Incubate for 60 minutes at room temperature on an orbital shaker (150rpm).
8. Wash the wells 3 times with 300 µl of diluted Washing solution. Aspirate all liquid from the wells.

9. Pipette 200 µl of Chromogenic substrate solution into the wells
10. Incubate 15 minutes at room temperature on an orbital shaker (150rpm), avoid direct light exposure (NB. The incubation time begins after the first TMB addition)
11. Pipette 100 µl of Stop solution into the wells
12. Read the absorbance of the wells (450, 405 and 620nm (as reference)). Reading must be completed within 20 minutes from the end of the assay.

## 10. CALCULATION OF RESULTS

In order to obtain a better sensitivity, the present method employs spectrophotometric reading at two wavelengths (450 and 405 nm). For all O.D. overflow at 450nm, multiplied the O.D. 405nm by the correction factor calculated by the ratio between O.D. 450nm and O.D. 405nm. To establish this factor, used a sample reading at 405nm and 450nm.


Draw a calibration curve on millimetric graph paper, by plotting the calibrators concentration (x-axis) against the absorbance obtained for each standard (y-axis, logit). Corresponding Proinsulin concentrations in pmol/l are obtained by interpolating the absorbances of each sample on the calibration curve.

## 11. EXAMPLE OF CALCULATION

The values shown below must be considered as an example and must not be used in place of experimental data.

| Description     | Concentration (pmol/l) | O.D   |
|-----------------|------------------------|-------|
| Calibrator 0    | 0                      | 0,040 |
| Calibrator 1    | 3,6                    | 0,255 |
| Calibrator 2    | 11,1                   | 0,746 |
| Calibrator 3    | 20,4                   | 1,391 |
| Calibrator 4    | 47,5                   | 3,116 |
| Calibrator 5    | 99,5                   | 6,129 |
| Control serum 1 | 14,0                   | 0,961 |
| Control serum 2 | 26,6                   | 1,845 |
| P1              | 7,6                    | 0,581 |
| P2              | 22,8                   | 1,620 |
| P3              | 25,9                   | 1,992 |
| P4              | 76,7                   | 5,352 |

## 12. TYPICAL CALIBRATION CURVE

### 13. NORMAL VALUES

The normal values determined are only indicative since they may be affected by various agents. We recommend that each laboratory establish its own normal range.

| Serum                      | n  | Mean $\pm$ SD (pmol/l) |
|----------------------------|----|------------------------|
| Non obese fasting patients | 16 | 2.56 $\pm$ 1.28        |
| Fasting IDDM               | 9  | 4.67 $\pm$ 3.85        |
| Non-fasting IDDM           | 13 | 16.25 $\pm$ 17.75      |

The regression analysis for a plasma/serum comparison equation is:

$$Y (\text{plasma value}) = 1.20 \times X (\text{serum value}) + 1.17$$

Correlation Coefficient  $r = 0.98$  ( $n=48$ )

### 14. PERFORMANCES OF THE ASSAY

#### SENSITIVITY

##### Analytical sensitivity

The sensitivity was calculated based upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero Calibrator (mean value + 2 S.D.). This dose is 0.6 pmol/l.

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#### SPECIFICITY

Patient samples may contain human anti-mouse antibodies (HAMA) which are capable of giving falsely elevated or depressed results with assays that utilize mouse monoclonal antibodies. This assay has been designed to minimize interference from HAMA-containing specimens with the use of a blocking HAMA interferences (blocking buffer). Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.

#### Cross-reactivity

The following peptides were tested and no cross-reactivity has been observed :

|                           |                 |
|---------------------------|-----------------|
| Human Insulin             | < 10 000 pmol/l |
| Human C-Peptide           | 50 000 pmol/L   |
| Des (31,32)- Proinsulin   | < 200 pmol/L    |
| Split (32,33)- Proinsulin | 5000 pmol/L     |
| Des (64,65)- Proinsulin   | 200 pmol/L      |
| Split (65,66)- Proinsulin | 1000 pmol/L     |

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