ZenTech

RIAZENco TG Ab

Deutsch	(de)
English	(en)
Español	(es)
Ελληνικά	(el)
Français	(fr)
Italiano	(it)
Nederlands	(nl)
Polski	(pl)
Portuges	(pt)
БЪЛГАРСКИ	(bg)



ZenTech s.a.

Liège Science Park
Avenue du Pré-Aily 10,
4031 ANGLEUR, Belgium
Tel. : + 32-(0)4-361.42.32
Fax : + 32-(0)4-367.00.63

info@zentech.be

www.zentech.be

ISO15223	SYMBOLS FÜR MEDIZINPRODUKTE	MEDICAL DEVICES SYMBOL	SÍMBOLOS PARA APARATOS DE USO EN MEDICINA	ΣΥΜΒΟΛΟ ΙΑΤΡΙΚΩΝ ΣΥΣΚΕΥΩΝ ΚΑΤΑ	SYMBOLS APPLIQUÉS AUX DISPOSITIFS MÉDICAUX	SIMBOLI DISPOSITIVI MEDICI	SYMBOL CONFORM MEDISCHE HULPMIDDELEN	SYMBOL URZĄDZENIA MEDYCZNEGO	SÍMBOLOS DE DISPOSITIVO MÉDICO	СИМВОЛИ НА МЕДИЦИНСКАТА АПАРАТУРА
	LIMITIERUNG DER LAGERTEMPORATUR	STORAGE TEMPERATURE LIMITATION	LIMITES DE TEMPERATURA DE ALMACENAMIENTO	ΠΕΡΙΟΡΙΣΜΟΣ ΘΕΡΜΟΚΡΑΣΙΑΣ ΦΥΛΑΞΗΣ	LIMITES DE TEMPERATURES	Conservare a	TEMPERATUURLIMIET	OGRANICZENIA TEMPERATURY PRZECHOWYWANIA	LIMITE DA TEMPERATURA DE ARMAZENAGEM	ОГРАНИЧАВАНЕ ТЕМΠΕΡΑΤУРАТА НА СЪХРАНЕНИЕ
	CHARGENCODE	BATCH CODE	CÓDIGO DE LOTE	ΚΩΔΙΚΟΣ ΠΑΡΤΙΔΑΣ	NUMERO DE LOT	Lotto	LOT NUMMER	KOD SERII	LOTE	ΠΑΡΤΙΔΕΝ НОМЕР
	VERWENDBAR BIS	USE BY	CONSUMIR ANTES DE	ΗΜΕΡΟΜ. ΛΗΞΗΣ	DATE D'EXPIRATION	Data di scadenza	HOUDBAAR TOT	ZUŻYC PRZED	UTILIZADO POR	СРОК НА ГОДНОСТ
	DAS HANDBUCH ZU RATE ZIEHEN	CONSULT OPERATING INSTRUCTIONS	CONSULTAR LAS INSTRUCCIONES DE MANEJO O FUNCIONAMIENTO	ΑΝΑΤΡΕΨΤΕ ΣΤΙΣ ΟΔΗΓΙΕΣ ΛΕΙΤΟΥΡΓΙΑΣ	LIRE LES INSTRUCTIONS	Consulti l'istruzione per l'uso	RAADPLEEG DE GEBRUIKSAANWIJZING	ZAPOZNAĆ SIĘ Z INSTRUKCJĄ	CONSULTE O MANUAL DE OPERAÇÕES	ВИЖТЕ ИНСТРУКЦИИТЕ ЗА ЕКСПЛОАТАЦИЯ
	FÜR IN-VITRO-DIAGNOSTISCHE ANWENDUNG	IN VITRO DIAGNOSTIC DEVICE	DISPOSITIVO DE DIAGNÓSTICO IN VITRO	ΔΙΑΓΝΩΣΤΙΚΗ ΣΥΣΚΕΥΗ IN VITRO	DIAGNOSTIC IN VITRO	Dispositivo medico-diagnostico in vitro	IN VITRO DIAGNOSTISCH PRODUCT	URZĄDZENIE DO DIAGNOSTYKI IN VITRO	DISPOSITIVO DE DIAGNOSTICO IN VITRO	УРЕД ЗА IN VITRO ДΙΑΓΝΩΣΤΙΚΑ
	HERGESTELLT VON	MANUFACTURED BY	FABRICADO POR	ΚΑΤΑΣΚΕΥΑΖΕΤΑΙ ΑΠΟ ΤΗΝ	FABRIQUE PAR	Prodotto da	FABRIKANT	PRODUCENT	FABRICADO POR	ΠΡΟΙΖΒΕΔΕΝΟ ОТ
	KATALOG NR.	CATALOGUE NUMBER	NÚMERO DE CATÁLOGO	ΑΡΙΘΜΟΣ ΚΑΤΑΛΟΓΟΥ	REFERENCE	Numero di catalogo	CATALOGUS NUMMER	NUMER KATALOGOWY	NÚMERO DO CATALOGO	КАΤΑΛΟЖΕΝ НОМЕР

	SYMBOLS (EMPFOHLEN VON DER EDMA)	SYMBOLS (EDMA RECOMMENDATIONS)	SÍMBOLOS (RECOMENDACIONES DE LA EDMA)	ΣΥΜΒΟΛΑ (Συστάσεις EDMA)	SYMBOLS (recommandations EDMA)	SIMBOLI (raccomandazioni EDMA)	SYMBOLS (EDMA-AANBEVELINGEN)	SYMBOLS (ZALECANE PRZEEZ EDMA)	SÍMBOLOS (RECOMENDAÇÕES EDMA)	СИМВОЛИ (ΠΡΕΠΟΡΨΚИ НА EDMA)
	Anzahl der Bestimmungen	number of determinations	Número de determinaciones	Αριθμός προσδιορισμών	Nombre de déterminations (Σ)	Numero di determinazioni	Inhoud voldoende voor "n" tests	Liczba oznaczeń	Número de determinações	БРОЙ ИЗΜΕΡΒΑΝΙΑ
	beschichtete Röhrchen	Coated tubes	Tubos recubiertos	Επιστρωμένα σωληνάκια	Tubes coatés	Tubi cotati	Gecoate buisjes	Probówki z warstwą przeciwciał	Tubos adsorvidos	ΠΟΚΡΙΤΙ ΕΠΡΥΒΕΤΚИ
	Radioactiver Tracer	Radioactive tracer	Trazador radiactivo	Ραδιενεργός ιχνηθέτης	Traceur radioactif	Tracciante radioattivo	Radioactieve tracer	znacznik izotopowy	Marcador radioativo	РАΔΙΟΑΚΤΙΒΕΝ ΤΡΕЙСЪР
	Kalibrator	Calibrator	Calibrador	Βαθμονομητής	Calibrateurs	Calibratore	Kalibrator	Kalibrator	Calibrador	КАΛΙΒΡΑΤΟΡΙ
	Kontrolle	Control	Control	Ορός ελέγχου	Contrôle	Controllo	Controle	Kontrola	Controle	ΚΟΝΤΡΟΛΑ
	Waschlotion zum x-maligen Verdünnen	Wash solution to be diluted n-fold	Solución de lavado que debe diluirse n veces	Διάλυμα πλύσης που πρέπει να αραιωθεί n φορές	Solution de lavage à diluer n fois	Soluzione di lavaggio da diluire in n volte	In n-voud te verdunnen wasoplossing	Roztwór odmywający do n-krotnego rozcieńczenia	Solução de lavagem a ser diluída n X	Изμиващият разтвор трябва да бъде разреден n-пъти.
	RADIOAKTIVES MATERIAL	RADIOACTIVE MATERIAL	MATERIAL RADIOACTIVO	ΡΑΔΙΕΝΕΡΓΟ ΥΛΙΚΟ	matériel radioactif	MATERIALE RADIOATTIVO	RADIOACTIEF MATERIAAL	MATERIAL PROMIENIOTWORCZY	MATERIAL RADIOACTIVO	РАΔΙΟΑΚΤΙΒΕΝ ΜΑΤΕΡΙΑΛ

DEUTSCH

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ENGLISH

RADIOIMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF THE ANTI-THYROGLOBULIN AUTOANTIBODIES IN HUMAN SERUM R-CI-100 - 100 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

1. CLINICAL APPLICATIONS

Thyroid autoimmunity is more frequently registered in women. Antibody prevalence in women increases with age, rising from approximately 10% at the age of 18-24 up to 30% at the age of 55-65 for TgAb and from 15% at the age of 18-24 up to 24% at the age of 55-65 for TMAb.

2. PRINCIPLE OF THE ASSAY

The present method is based on a competitive radioimmunoassay (RIA). During the incubation, the monoclonal anti-Tg antibody on solid phase competes with the sample/calibrator anti-Tg autoantibodies for the specific sites of the Iodine¹²⁵ labelled Tg (tracer). After aspiration and washing, the radioactivity in the tubes is measured in a gamma counter. The degree of binding will be inversely proportional to the sample anti-Tg autoantibody concentration.

The calibrators of RIAZENco TG Ab are calibrated against the TgAb international reference preparation MRC65/93.

3. REAGENTS PROVIDED WITH THE KIT

- The reagents are sufficient for 100 determinations.
- Store the kit and reagents at 2-8°C.
- The expiration date of each reagent is shown on the label.

- 1 - Radioactive Tracer:** 1 vial containing 22 ml of iodinated thyroglobulin, red color coded. Radioactivity content : ± 137 kBq per vial. Preservative: NaN_3 (<0.1%).
- 2 - Calibrators: Zero Calibrator:** 1 vial containing 4.2 ml. Preservative: NaN_3 (<0.1%). Ready for use.
Calibrators 1 to 5: 5 vials containing 0.7 ml of anti-Tg (human serum), at the following concentrations: 20, 60, 200, 1000, 2000 IU/ml. Preservative: NaN_3 (<0.1%).
- 3 - Controls:** 2 vials containing 0.7 ml of anti-Tg (human serum). **For the exact value, refer to the value written on the Quality Control data sheet.** Preservative: NaN_3 (<0.1%).
- 4 - Coated Tubes:** 100 tubes coated with mouse monoclonal anti-Tg antibodies. Unused tubes must be stored at 2-8°C, protected from moisture.
- 5 - Washing Solution (50 x concentrated):** 1 vial of 20 ml Tris-HCl buffer with detergent and preservative NaN_3 (<0.1%). Bring to 1000 ml with distilled water. The diluted washing solution is stable for 2 months at 2-8°C.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Plastic test tubes
- Test tube racks.
- Adjustable, automatic micropipettes with disposable tips.
- Vortex mixer.
- Graduated cylinder.
- Aspiration pump or automated washing device.
- Scintillation gamma counter.
- Distilled water.
- Orbital shaker adjustable at 150 rpm.

5. WARNINGS AND PRECAUTIONS FOR USERS

For *in vitro* diagnostic use.

Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.

Radioactive Material - Not for Internal or External Use in Humans or Animals.

	erwartete Werte (nmol/l)	gemessene Werte (nmol/l)	Wiederfindung (%)
S1	-	50,4	-
S1 + KAL 0	25,2	32,0	127,0
S1 + KAL 1	35,2	39,1	111,1
S1 + KAL 2	55,2	59,6	108,0
S1 + KAL 3	125,2	111,8	89,3
S1 + KAL 4	525,2	392,4	74,7
S1 + KAL 5	1025,2	1251,1	122,0
S2	-	102,9	-
S2 + KAL 0	51,5	63,9	124,1
S2 + KAL 1	61,5	64,2	104,4
S2 + KAL 2	81,5	84,7	103,9
S2 + KAL 3	151,5	133,0	87,8
S2 + KAL 4	551,5	367,2	66,6
S2 + KAL 5	1051,5	1261,3	120,0

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This radioactive material may be received, acquired, possessed, and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation there from to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations of each country.
Physical characteristics of ¹²⁵I: see end of instructions

- Safety Precautions:** The following precautions should be observed in handling radioactive material:
 - o Store radioactive materials in a designated area.
 - o Do not eat, drink, smoke or apply cosmetics where radioactive materials are being handled.
 - o Do not pipette by mouth.
 - o Wear gloves when handling radioactive materials and wash hands thoroughly afterwards.
 - o Cover working area with disposable absorbent paper.
 - o Wipe up all spills immediately and thoroughly and dispose of the contaminated materials as radioactive waste.
 - o Dispose of the liquid radioactive waste into the sanitary sewage system if permitted by the local regulations.

- CHEMICAL HAZARD Sodium Azide (NaN₃) Warning:** Some of the reagents in this kit contain sodium azide as a preservative. For all such reagents, the concentration of sodium azide is <0.09% w/w. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Dispose of all non radioactive reagents by flushing with large amounts of water through the plumbing system.

Risk Phrases

R 21/22 Harmful in contact with skin and if swallowed.

Safety Phrases

S-26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S-28.1 After contact with skin, wash immediately with plenty of water.

S-46 If swallowed seek medical advice immediately and show the container or label.

- POTENTIAL BIOHAZARDOUS MATERIAL Warning:** This kit may contain some reagents made with human serum or plasma. The serum or plasma used has been tested by an FDA-approved method and found to be non-reactive for HIV-1/2 Antibodies, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 3rd Edition 1993.

6. SPECIMEN COLLECTION

It is recommended to use serum. Highly lipemic or hemolyzed samples must be discarded. Keep samples at 2-8°C for 1 day; for longer periods it is advisable to freeze samples in aliquots at -20°C. Repeated freezing and thawing of samples should be avoided.

7. ASSAY PROCEDURE

- Bring all reagents and samples to warm up at room temperature.
 - Before use, mix samples by gentle agitation before use.
 - For all calibrators, a duplicate measure is recommended
- 1- Prepare plain tubes for Total Counts, and coated tubes for Calibrators, Samples and Controls.
 - 2- Pipette **20 µl** of each Calibrator, Control Serums and Samples into the corresponding tubes.
 - 3- Add **200 µl** of Radioactive Tracer into all tubes.
 - 4- Mix the tubes on vortex and incubate for **90 minutes** at room temperature on an orbital shaker set at 150 rpm.
 - 5- Aspirate the contents of each tube, except the tubes for Total Counts.
 - 6- Add **2 ml** of diluted washing solution to every tube, except tubes for Total Counts and aspirate thoroughly or decant the contents of all tubes on absorbent paper.
 - 7- Count the radioactivity bound to the tubes for 1 minute in a gamma counter. We suggest to control the background of the instrument before counting the assay. In order to avoid variations in the sensitivity of the system, the background must be reduced to a minimum or adjusted properly.

NOTE :

If high anti-Tg values (> 2000IU/ml) are expected for any patient sample, the original sample should be diluted. The dilution should be done using zero calibrator. **Do not use any buffers for this purpose.**

ASSAY SCHEME

Tubes	Total Activity	Calibrators	Controls	Samples
Reagent				
Calibrators	----	20 µl	----	----
Control	----	----	20 µl	----
Serums				
Samples	----	----	----	20 µl
Tracer	200 µl	200 µl	200 µl	200 µl

- Mix and incubate: 90 min R.T. shaking (150 rpm)
- Aspirate and wash: 1 x 2 ml
- Count

8. CALCULATION OF RESULTS

$$\text{Binding capacity} = \frac{B_o}{T} \% = \frac{B_o \text{ Cpm}}{T \text{ Cpm}} \times 100$$

Percent binding for Calibrators, Controls and Samples =

$$\text{Binding \%} = \frac{B}{B_o} \% = \frac{\text{Calibrator or Sample Cpm}}{B_o \text{ Cpm}} \times 100$$

Draw the calibration curve on log/lin graph by plotting the B/B_o% obtained for each calibrator (y-axis) against the relative concentration (x-axis). Calculate the B/B_o% of each sample and read the anti-Tg concentration in IU/ml, by interpolating on the calibration curve.

EXAMPLE OF CALCULATION

The values reported below must be considered as an example and may not be used in place of experimental data.

Description	Average cpm	B/B _o (%)	TG Ab conc. (IU/ml)
Total Activity (T)	57953	-	-
Cal 0	28090	100	0
Cal 1	24340	86.6	20
Cal 2	16927	60.3	60
Cal 3	9706	34.5	200
Cal 4	4100	14.6	1000
Cal 5	2048	7.3	2000
CONTROL 1	10843	38.6	158.2
CONTROL 2	5757	20.5	604.7
P1	13998	49.8	91.5
P2	6634	23.6	464.8
P3	4270	15.2	942.9

INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

9. REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

Lower than 30 IU/ml	Negative for anti-Thyroglobulin auto-antibodies.
Between 30 - 70 IU/ml	Borderline
Higher than 70 IU/ml	Positive for anti-Thyroglobulin auto-antibodies.

10. PERFORMANCES OF THE ASSAY

SPECIFICITY

No cross-reactions have been observed with anti-TPO autoantibodies present in human serum.

SENSITIVITY

Analytical sensitivity

The sensitivity was calculated based upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero Calibration (mean value - 2 S.D.). This dose is 6 IU/ml.

Functional sensitivity

The functional assay sensitivity is the lowest value which is measured with a precision of maximum 20% inter-assay variance. For the anti-thyroglobulin autoantibodies, this value is lower than 15 IU/ml

PRECISION

Precision was evaluated upon intra- and inter-assay variability, at different analyte concentrations.

Intra-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. %	N
A	72.6	±	6.0	8.3	20
B	345.7	±	24.7	7.1	20
C	889.9	±	51.0	5.7	20

Inter-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. %	N
A	90.3	±	8.4	9.3	9
B	390.0	±	49.8	12.8	9
C	985.4	±	104.9	10.6	9

ACCURACY

Accuracy of the method has been checked by the recovery and parallelism tests:

Recovery Test

Samples, mixed with equal volumes of each calibrator, were tested.

	Expected (IU/ml)	Measured (IU/ml)	Recovery (%)
S1	-	50.4	-
S1 + CAL 0	25.2	32.0	127.0
S1 + CAL 1	35.2	39.1	111.1
S1 + CAL 2	55.2	59.6	108.0
S1 + CAL 3	125.2	111.8	89.3
S1 + CAL 4	525.2	392.4	74.7
S1 + CAL 5	1025.2	1251.1	122.0
S2	-	102.9	-
S2 + CAL 0	51.5	63.9	124.1
S2 + CAL 1	61.5	64.2	104.4
S2 + CAL 2	81.5	84.7	103.9
S2 + CAL 3	151.5	133.0	87.8
S2 + CAL 4	551.5	367.2	66.6
S2 + CAL 5	1051.5	1261.3	120.0

Parallelism Test

A serum with high analyte concentration was tested at different dilutions with the Zero Calibrator.

Dilution	Expected (IU/ml)	Measured (IU/ml)	Recovery (%)
S1 undiluted	-	1109	-
1/2	554.5	586.7	105.8
1/4	277.3	275.0	99.2
1/8	138.6	143.0	103.2
1/16	69.3	82.0	118.3

NB. Due to the heterogeneity of autoantibodies, for some patient samples a nonlinear dilution is possible

This package insert is downloadable at www.zentech.be

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1. Prentice L.M. et al. Geographical distribution of subclinical autoimmune thyroid disease in Britain : A study using highly sensitive direct assays for autoantibodies to thyroglobulin and thyroid peroxidase. Acta Endocrinologica , 1990, 123, 493-8.
2. Jaume J.C. et al. Thyroid peroxidase autoantibody epitopic « fingerprints » in juvenile Hashimoto's thyroiditis : evidence for conservation over time an in families. Clin. Exp. Immunol., 1996, 104, 115-123.
3. Czarnocka B. et al. Immunochemical properties of hTPO. Thyroperoxidase and Thyroid Autoimmunity, Ed. P. Carayon, J. Ruf. 1990, 207, 59-67.
4. Czarnocka B. et al. Majority of thyroid peroxidase autoantibodies in patients with autoimmune thyroid disease are directed to a single TPO domain. Autoimmunity, 1996, 23, 145-154.
5. Takamatsu J. et al. Correlation of antithyroglobulin and antithyroid-peroxidase antibody profiles with clinical and ultrasound characteristics of chronic thyroiditis. Thyroid, 1998, 8, 1101-6.
6. Smyth P. et al. Serum thyroid peroxidase autoantibodies, thyroid volume, and outcome in breast carcinoma. J. Endocrinol. Metab., 1998, 83, 2711-6.
7. Gauna A. et al. Immunological aspects of Graves' disease patients in different clinical stages. J. Endocrinol. Invest., 1989, 12, 671-7.
8. Stagnaro-Green A. et al. Detection of at-risk pregnancy by means of highly sensitive assays for thyroid autoantibodies. JAMA, 1990, 264,1422-5.

Physikalisches Daten of ¹²⁵I
Physical characteristics of ¹²⁵I
Caracteristicas fisicas ¹²⁵I
φυσικά χαρακτηριστικά του ¹²⁵I
Caractéristiques physique de ¹²⁵I
Caratteristica fisica¹²⁵I
Fysieke kenmerken van ¹²⁵I
Charakterystyka fizyczna I-¹²⁵
Características físicas ¹²⁵I
Физични характеристики на ¹²⁵I

t_{1/2} = 59.9 Tagen, days, dias, ημέρες, jours, giorni, dagen, dni, días,
дена
Wichtigste Emissionnen
Main emissions
Emissionnes principales
κύριες εκπομπές
Emissions principales
Emissione principale
Hoofduitzendingen
Głowne rodzaje promieniowania
Emissões principais
Основнилъчения

	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25