ZenTech

RIAZENco Aldosterone

DEUTSCH	(de)
ENGLISH	(en)
ESPAÑOL	(es)
ΕΛΛΗΝΙΚΑ	(el)
FRANÇAIS	(fr)
ITALIANO	(it)
NEDERLANDS	(nl)
POLSKI	(pl)
PORTUGES	(pt)
БЪЛГАРСКИ	(bg)


ZenTech s.a.

Liège Science Park,
Avenue du Pré-Ailly 10,
4031 ANGLEUR, Belgium
Tel. : + 32-(0)4-361.42.32
Fax : + 32-(0)4-367.00.63

info@zentech.be

www.zentech.be

ISO15223	SYMBOLS FÜR MEDIZINPRODUKTE	MEDICAL DEVICES SYMBOL	SÍMBOLOS PARA APARATOS DE USO EN MEDICINA	ΣΥΜΒΟΛΟ ΙΑΤΡΙΚΩΝ ΣΥΣΚΕΥΩΝ ΚΑΤΑ	SYMBOLS APPLIQUES AUX DISPOSITIFS MEDICAUX	SIMBOLI DISPOSITIVI MEDICI	SYMBOL CONFORM MEDISCHE HULPMIDDELEN	SYMBOL URZĄDZENIA MEDYCZNEGO	SÍMBOLOS DE DISPOSITIVO MÉDICO	СИМВОЛИ НА МЕДИЦИНСКАТА АПАРАТУРА
	LIMITIERUNG DER LAGERTEMPORATUR	STORAGE TEMPERATURE LIMITATION	LIMITES DE TEMPERATURA DE ALMACENAMIENTO	ΠΕΡΙΟΡΙΣΜΟΣ ΘΕΡΜΟΚΡΑΣΙΑΣ ΦΥΛΑΞΗΣ	LIMITES DE TEMPERATURES	Conservare a	TEMPERATUURLIMIET	OGRANICZENIA TEMPERATURY PRZECHOWYWANIA	LIMITE DA TEMPERATURA DE ARMAZENAGEM	ОГΡΑΝΙΧΑΒΑΝΕ ΤΕΜΠΕΡΑΤУРАТА НА СЪХΡΑΝЕНИЕ
LOT	CHARGENCODE	BATCH CODE	CÓDIGO DE LOTE	ΚΟΔΙΚΟΣ ΠΑΡΤΙΔΑΣ	NUMERO DE LOT	Lotto	LOT NUMMER	KOD SERII	LOTE	ΠΑΡΤΙΔΕΝ ΝΟΜΕΡ
	VERWENDBAR BIS	USE BY	CONSUMIR ANTES DE	ΗΜΕΡΟΜ. ΛΗΞΗΣ	DATE D'EXPIRATION	Data di scadenza	HOUBAAR TOT	ZUŻYC PRZED	UTILIZADO POR	СРОК НА ГОДНОСТ
	DAS HANDBUCH ZU RATE ZIEHEN	CONSULT OPERATING INSTRUCTIONS	CONSULTAR LAS INSTRUCCIONES DE MANEJO O FUNCIONAMIENTO	ΑΝΑΤΡΕΨΤΕ ΣΤΙΣ ΟΔΗΓΙΕΣ ΛΕΙΤΟΥΡΓΙΑΣ	LIRE LES INSTRUCTIONS	Consulti l'istruzione per l'uso	RAADPLEEG DE GEBRUIKSAANWIJZING	ZAPOZNAĆ SIĘ Z INSTRUKCJĄ	CONSULTE O MANUAL DE OPERAÇÕES	ΒΙЖΤΕ ΙΝΣΤΡУКЦИЈИТЕ ЗА ЕКСПΛΟΑΤΑЦИЈА
IVD	FÜR IN-VITRO-DIAGNOSTISCHE ANWENDUNG	IN VITRO DIAGNOSTIC DEVICE	DISPOSITIVO DE DIAGNÓSTICO IN VITRO	ΔΙΑΓΝΩΣΤΙΚΗ ΣΥΣΚΕΥΗ IN VITRO	IN VITRO DIAGNOSTIC	Dispositivo medico-diagnostico in vitro	IN VITRO DIAGNOSTISCH PRODUCT	URZĄDZENIE DO DIAGNOSTYKI IN VITRO	DISPOSITIVO DE DIAGNOSTICO IN VITRO	УРЕД ЗА IN VITRO ΔΙΑΓΝΩΣΤΙΚΑ
	HERGESTELLT VON	MANUFACTURED BY	FABRICADO POR	ΚΑΤΑΣΚΕΥΑΖΕΤΑΙ ΑΠΟ ΤΗΝ	FABRIQUE PAR	Prodotto da	FABRIKANT	PRODUCENT	FABRICADO POR	ΠΡΟΙΖΒΕΔΕΝΟ ΟΤ
REF	KATALOG NR.	CATALOGUE NUMBER	NÚMERO DE CATÁLOGO	ΑΡΙΘΜΟΣ ΚΑΤΑΛΟΓΟΥ	REFERENCE	Numero di catalogo	CATALOGUS NUMMER	NUMER KATALOGOWY	NÚMERO DO CATALOGO	ΚΑΤΑΛΟЖΕΝ ΝΟΜΕΡ

	SYMBOLS (EMPFOLHEN VON DER EDMA)	SYMBOLS (EDMA RECOMMENDATIONS)	SÍMBOLOS (RECOMENDACIONES DE LA EDMA)	ΣΥΜΒΟΛΑ (Συστάσεις EDMA)	SYMBOLS (recommandations EDMA)	SIMBOLI (raccomandazioni EDMA)	SYMBOLS (EDMA-AANBEVELINGEN)	SYMBOLS (ZALECANE PRZEZ EDMA)	SÍMBOLOS (RECOMENDAÇÕES EDMA)	СИМВОЛИ (ΠΡΕΠΟΡЪΚИ НА EDMA)
	Anzahl der Bestimmungen	Number of determinations	Número de determinaciones	Αριθμός προσδιορισμών	Nombre de déterminations (Σ)	Numero di determinazioni	Inhoud voldoende voor "n" tests	liczba oznaczeń	Número de determinações	БРОЙ ИЗΜΕΡΒΑΝΙΑ
CAL	Kalibratoren	Calibrators	Calibradores	Βαθμονομητές	Calibrateurs	Calibratore	Kalibratoren	Kalibrator	Calibrador	ΚΑΛΙΒΡΑΤΟΡΙ
CONTROL	Kontrolle	Control	Control	Ορός ελέγχου	Contrôle	Controllo	Controle	Kontrola	Controle	ΚΟΝΤΡΟΛΑ
SORB CT	Beschichtete Röhrchen	Coated tubes	Tubos recubiertos	Επιστρωμένα σωληνάκια	Tubes coatés	Tubi cotati	Gecoate buisjes	Probówki z warstwą przeciwiścial	Tubos adsorvidos	ΠΟΚΡΙΤΙ ΕΠΡΥΒΕΤΚΙ
ALDO ^{125I}	Radioactiver Tracer	Radioactive tracer	Trazador radiactivo	Ραδιενεργός ιχνηθέτης	Traceur radioactif	Tracciante radioattivo	Radioactieve tracer	znacznik radioaktywny	Marcador radioativo	ΡΑΔΙΟΑΚΤΙΒΕΝ ΤΡΕΙΣЪР
RCNS x ml H ₂ O	Wiederherstellen mit x ml destilliertem Wasser	Reconstitute with x ml distilled water	Reconstituir con x ml de agua destilada	Ανασύσταση με x ml απεσταγμένου νερού	Reconstitué avec x ml d'eau distillée	Ricostituite con x ml di acqua distillata	Reconstitueer met x ml gedestilleerd water	Rozpuścić w x ml wody destylowanej	Reconstitua com x ml água destilada	ΡΕΚΟΝΣΤΙΤΥΡΑΙΤΕ С Χ ΜL ΔΕΣΤΙΛΙΡΑΝΑ ΒΟΔΑ
	RADIOAKTIVES MATERIAL	RADIOACTIVE MATERIAL	MATERIAL RADIOACTIVO	ΡΑΔΙΕΝΕΡΓΟ ΥΛΙΚΟ	MATERIAL RADIOACTIF	MATERIALE RADIOATTIVO	RADIOACTIEF MATERIAAL	MATERIAL PROMIENIOTWORCZY	MATERIAL RADIOACTIVO	ΡΑΔΙΟΑΚΤΙΒΕΝ ΜΑΤΕΡΙΑΛ

DEUTSCH

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ENGLISH

RADIOIMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF ALDOSTERONE IN HUMAN SERUM AND URINE

R-CW-100 - 100 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

1. CLINICAL APPLICATIONS

Aldosterone is a steroid hormone with molecular weight of 360 produced by the zona glomerulosa of the adrenal cortex.

The aldosterone assay is extremely important in the diagnosis of primary hyperaldosteronism (Conn's disease), in the so-called essential hypertension, in arterial pressure changes induced by the administration of abnormal amounts of certain types of food (licorice) or drugs (potassium/sodium), in nephrosis syndromes etc.

Blood levels of aldosterone follow also a circadian rhythm and are strictly related to states of rest or exercise, it is necessary to study blood aldosterone levels by means of stimulation and/or suppression dynamic tests.

2. PRINCIPLE OF THE ASSAY

The present method is based on a competitive radioimmunoassay (RIA). During the incubation, the sample/calibrator aldosterone competes with the aldosterone labeled with Iodine 125 (tracer) for the specific sites of the antiserum coated on the tubes. Separation is based on the use of antibody coated tubes, where the anti-aldosterone antiserum is fixed on the tube walls. After aspiration, the radioactivity in the tubes is measured in a gamma counter. The degree of binding will be inversely proportional to the sample/calibrator hormone concentration.

3. REAGENTS PROVIDED WITH THE KIT

- The reagents are sufficient for 100 determinations.
- Store the kit and reagents at 2-8°C.
- The expiration date of each reagent is shown on the label.

1 Radioactive Tracer: 1 vial (50 ml) of ¹²⁵I-Aldosterone in phosphate buffer. Radioactivity contents: 85 KBq. Preservative: NaN₃ (<0.1%).

2 Calibrators: 8 vials of Aldosterone in human serum. Preservative: NaN₃ (<0.1%). Lyophilized. Reconstitute the Zero Calibrator with 3 ml of distilled water and the Calibrators 1-6 with 2 ml of distilled water. After reconstitution, aliquot and freeze unused calibrators and control at -20°C. For the exact value, refer to the Quality Control data sheet.

3 Control: 1 vial of Aldosterone in human serum. Preservative: NaN₃ (<0.1%). Lyophilized. Reconstitute with 2 ml of distilled water. After reconstitution, aliquot and freeze unused control at -20°C. For the exact value, refer to the Quality Control data sheet.

4 Coated Tubes: 100 tubes coated with anti-Aldosterone antibody. Unused tubes must be stored at 2-8°C, protected from moisture.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Plastic test tubes
- Test tube racks.
- Adjustable, automatic micropipettes with disposable tips.
- Vortex mixer.
- Graduated cylinder.
- Aspiration pump or automated washing device.
- Scintillation gamma counter.
- Distilled water.
- Orbital shaker adjustable at 150 rpm.
- Hydrochloric acid.

5. WARNINGS AND PRECAUTIONS FOR USERS

For *in vitro* diagnostic use.

Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.

Radioactive Material - Not for Internal or External Use in Humans or Animals.

This radioactive material may be received, acquired, possessed, and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation there from to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations of each country. Physical characteristics of ¹²⁵I: see end of instructions

- Safety Precautions:** The following precautions should be observed in handling radioactive material:
 - o Store radioactive materials in a designated area.
 - o Do not eat, drink, smoke or apply cosmetics where radioactive materials are being handled.
 - o Do not pipette by mouth.
 - o Wear gloves when handling radioactive materials and wash hands thoroughly afterwards.
 - o Cover working area with disposable absorbent paper.
 - o Wipe up all spills immediately and thoroughly and dispose of the contaminated materials as radioactive waste.
 - o Dispose of the liquid radioactive waste into the sanitary sewage system if permitted by the local regulations.
- CHEMICAL HAZARD Sodium Azide (NaN₃) Warning:** Some of the reagents in this kit contain sodium azide as a preservative. For all such reagents, the concentration of sodium azide is <0.1% w/w. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Dispose of all non radioactive reagents by flushing with large amounts of water through the plumbing system.

Risk Phrases
R 21/22 Harmful in contact with skin and if swallowed.
Safety Phrases
S-26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S-28.1 After contact with skin, wash immediately with plenty of water.
S-46 If swallowed seek medical advice immediately and show the container or label.
- POTENTIAL BIOHAZARDOUS MATERIAL Warning:** This kit may contain some reagents made with human serum or plasma. The serum or plasma used has been tested by an FDA-approved method and found to be non-reactive for HIV-1/2 Antibodies, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 3rd Edition 1993.

6. SPECIMEN COLLECTION

The assay can be performed in serum or urine samples. Since Aldosterone shows a circadian rhythm we suggest to note down the sample collecting time. Highly lipemic or hemolyzed samples must be discarded. Keep samples at 2-8°C for 1-2 days; for longer periods it is advisable to freeze samples at -20°C. Repeated free zing and thawing of samples should be avoided.

Urine Samples

Urine collection.

Collect 24-hour urine by adding boric acid (10 g/l) as preservative. Record the volume and store at 2-8°C for 1-2 days; freeze at -20°C for longer time periods.

Urine Hydrolysis:

- 1 - Prepare uncoated tubes, one for each urine sample.
- 2 - Pipette 100 µl of urine sample.
- 3 - Add 1 ml of 0.1 M hydrochloric acid into each tube.
- 4 - Incubate for 15-20 hours at 30°C.
- 5 - Use 50 µl of this solution for the test.

7. ASSAY PROCEDURE

- Bring all reagents and samples to room temperature prior to use.
- Before use, mix the samples by inversion.
- For all calibrators, a duplicated measure is recommended.

- 1 - Prepare plain tubes for Total Activity (T) and Non-specific Binding (NSB) as well as coated tubes for Zero Calibrator (Bo), Calibrators (1-6), Control and Samples.

- 2 - Pipette 200 µl of each Calibrator, Control and Sample into each tube.
- 2a - **Urine Samples.** Pipette 50 µl of Hydrolyzed sample into each tube and add 150 µl of the Zero Calibrator.
- 3 - Pipette 200 µl of Zero Calibrator into the Non-specific Binding (NSB) tubes.
- 4 - Add 500 µl of Radioactive Tracer into all tubes.
- 5 - Mix test-tube rack manually. Do not use vortex.
- 6 - Incubate: **Procedure A:** 18-24 hours at room temperature; or **Procedure B:** 3 hours at room temperature with an orbital shaker set at 150 rpm.
- 7 - Carefully aspirate the incubation mixture from all tubes, except those for total activity, with an aspiration pump or decant by drying the edges of the tubes with blot-paper.
- 8 - Count the radioactivity in the tubes for 1 minute by using a gamma counter. We suggest to check the background of the instrument before counting the assay. In order to avoid variations in the sensitivity of the system, the background should be reduced to a minimum or adjusted properly.

ASSAY SCHEME

Tubes Reagent	Total Activity	NSB	Cali-brators	Control	Samples	Urine samples
Calibrator 0	----	200 µl	----	----	----	150 µl
Calibrators	----	----	200 µl	----	----	----
Controls	----	----	----	200 µl	----	----
Samples	----	----	----	----	200 µl	50 µl
Tracer	500 µl	500 µl	500 µl	500 µl	500 µl	500 µl

- Incubate: 18 – 24 h R.T. or 3 h R.T. shaking (150 rpm)
 - Aspirate
 - Count

8. CALCULATION OF RESULTS

Draw the calibration curve on logit-log or semi-log paper, with the calibration concentrations on the x-axis and the respective B/Bo % on the y-axis.

Calculate the B/Bo% for each sample and read the concentration by interpolating on the calibration curve to obtain the Aldosterone concentration in the tested samples, expressed in pg/ml.

Urine Samples.

The concentration read on the calibration curve must be multiplied by factor 44 to obtain the Aldosterone concentration in pg/ml. In order to express Aldosterone concentration in µg/24h:

$$\mu\text{g}/24\text{h} = \frac{\text{Concentration (pg/ml)} \times 44 \times \text{Urine excreted (liters)}}{1000}$$

EXAMPLE OF CALCULATION

The values reported below must be considered as an example and may not be used in place of experimental data.

Description	Average cpm	B/BO (%)	Aldosterone (pg/ml)
Total activity (T)	41999	-	-
NSB	154	-	-
CAL 0	19213	100	0
CAL 1	16439	85.4	25
CAL 2	14667	76.1	50
CAL 3	12268	63.6	100
CAL 4	8184	42.1	300
CAL 5	5844	29.9	600
CAL 6	3669	18.4	1500
CONTROL	11263	58.3	133
P1	14905	77.4	45.5
P2	9869	51.0	191
P3	7934	40.8	323

9. REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

Serum Aldosterone:	
- At rest:	10 - 160 pg/ml
- In motion:	35 - 300 pg/ml
Urinary Aldosterone:	2.8 - 30 µg/24h.

10. PERFORMANCES OF THE ASSAY

SPECIFICITY

The present method has shown the following cross-reactions: 100% with Aldosterone, 17.2% with 3β,5α Tetrahydroaldosterone, 0.12% with 3α,5β Tetrahydroaldosterone, less than 0.017% with Prednisolone, less than 0.003% with Cortisol, Corticosterone, Cortisone, 11-Deoxycortisol, 11-Deoxycorticosterone, 17-Hydroxyprogesterone, Progesterone and less than 0.0017% with Prednisone, Testosterone and Androstenedione.

The percentage of cross-reactivity is calculated by using the Abraham formula: $X/Y \times 100$, X and Y being the weight of the substance to be assayed and the weight of the interfering substance, both able to reduce the binding capacity by 50%.

SENSITIVITY

Analytical sensitivity

The sensitivity was calculated based upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero Calibrator (mean value - 2 S.D.). This dose is 1.4 pg/ml.

Functional sensitivity

The functional assay sensitivity is the lowest value which is measured with a precision of maximum 20% inter-assay variance. For the RIAZENco Aldosterone, this value is 25 pg/ml.

PRECISION

Precision was evaluated upon intra- and inter-assay variability at different analyte concentrations.

Intra-assay

Serum	Mean	±	S.D.	C.V.	N
	(pg/ml)	(pg/ml)	(pg/ml)	(%)	
1	42.3	±	5.8	13.7	20
2	177	±	6.7	3.8	20
3	278	±	14.8	5.3	20

Inter-assay

Serum	Mean	±	S.D.	C.V.	N
	(pg/ml)	(pg/ml)	(pg/ml)	(%)	
1	38.2	±	7.1	18.6	9
2	169	±	12.6	7.5	9
3	284	±	17.7	6.2	9

ACCURACY

Accuracy of the method has been checked by the recovery and parallelism tests.

Recovery Test

Samples, mixed with equal volumes of each calibrator, were tested.

Added	Expected (pg/ml)	Measured (pg/ml)	Recovery (%)
S1	-	97.2	-
S1 + CAL 0	48.6	41.6	85.6
S1 + CAL 1	61.1	69.4	113.6
S1 + CAL 2	73.6	88.8	120.7
S1 + CAL 3	98.6	121	122.7
S1 + CAL 4	199	249	125.1
S1 + CAL 5	349	455	130.4
S1 + CAL 6	799	933	116.8
S2	-	142	-
S2 + CAL 0	71.0	68.4	96.3
S2 + CAL 1	83.5	83.4	99.9
S2 + CAL 2	96.0	96.9	100.9
S2 + CAL 3	121	137	113.2
S2 + CAL 4	221	244	110.4
S2 + CAL 5	371	431	116.2
S2 + CAL 6	821	952	116.0

Parallelism Test

Serums with high analyte concentration was tested at different dilutions with the Zero Calibrator.

Dilution	Expected (pg/ml)	Measured (pg/ml)	Recovery (%)
S1 undiluted	-	128	-
1/2	64.1	61.6	96.1
1/4	32.1	30.9	96.3
1/8	16.0	14.9	78.4
S2 undiluted	-	351	-
1/2	176	161	91.5
1/4	87.8	75	85.4
1/8	43.9	36.3	82.7
1/16	22.0	25	113.6

This package insert is downloadable at www.zentech.be.

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- 1 - Bizollon C.A. et al. Solid-phase radioimmunoassay of plasma Aldosterone. Steroids, 1974, 23,809-21.
- 2 - Malvano R. et al. Methodological simplifications in radioimmunoassay of urinary Aldosterone. Clin. Chim. Acta, 1976, 66,331-43.
- 3 - Lijnen P. et al. Direct radioimmunoassay of plasma Aldosterone in normal subjects. Clin. Chim. Acta, 1978, 84, 305-14.
- 4 - Klumpp F. et al. Radioimmunologische Bestimmung von Aldosteron im Plasma. Z. Klin. Chem. Klin. Biochem., 1974, 12, 128-33.
- 5 - Pratt J.J. et al. Special problems in the radioimmunoassay of plasma Aldosterone without prior extraction and purification. Clin. Chim. Acta, 1978, 84, 329-37.
- 6 - St. Cyr M.J. et al. Quantitation of plasma Aldosterone by radioimmunoassay. Clin. Chem.,1972, 18, 1395-1402.
- 7 - Sulon J. et al. Radioimmunoassay for Aldosterone and Deoxycorticosterone in human plasma. Comparison of various antisera and determination of normal values. Clin. Chim. Acta, 1976, 72, 275-84.
- 8 - Young S.C. et al. Hypertensive complications in patients with primary aldosteronism. Current therapeutic Research, 1991, 50, 317-25.

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