
KH001-02 Anti AGEs Monoclonal Antibody (Clone No. 6D12), Peroxidase conjugated

| | | Application | |
|--------------------|--|----------------|----------------|
| Primary Source | - | WB | 0.25-5.0 µg/mL |
| Type | Monoclonal | IHC | 2.0 µg/mL |
| Immunogen | AGEs-BSA | ICC | Not tested |
| Raised in | Mouse | ELISA | 0.1-0.5 µg/mL |
| Myeloma | P3U1 | FCM | Not tested |
| Clone number | 6D12 | Neutralization | Not tested |
| Isotype | IgG1 | IP | Not tested |
| Source | Ascites | | |
| Purification notes | ProteinG | | |
| Cross Reactivity | every animal species | | |
| Concentration | 0.1 mg/mL | | |
| Contents (Volume) | 20 µg (200 µL/vial) | | |
| Label | peroxidase conjugated | | |
| Buffer | PBS [containing 2 % Block Ace as a stabilizer, 0.1 % Proclin as a bacteriostat] | | |
| Storage | Store below -20 °C. Once thawed, store at 4 °C. Repeated freeze-thaw cycles should be avoided. | | |

Immunohistochemical staining of renal proximal tubules and glomeruli in patients with diabetic nephropathy, using anti-AGEs antibody 6D12
Yamada, K. et al.,
Clinical nephrology, Vol.42, 354-361, 1994

Immunohistochemical staining of the early stage of human atherosclerotic lesions of the aorta with anti-AGEs antibody 6D12.
Kume, S. et al.,
American Journal of Pathology, Vol.147, 654-667, 1995

Note

Reaction of protein amino groups with glucose leads, through the early products such as a Schiff base and Amadori rearrangement products, to the formation of advanced glycation end products (AGEs). Recent immunological studies using anti-AGEs antibody (6D12) demonstrated the presence of AGEs-modified proteins in several human tissues: (i) human lens (nondiabetic and noncataractous), (ii) renal proximal tubules in patients with diabetic nephropathy and chronic renal failure, (iii) diabetic retina, (iv) peripheral nerves of diabetic neuropathy, (v) atherosclerotic lesions of arterial walls, (vi) β 2-microglobulin forming amyloid fibrils in patients with hemodialysis-related amyloidosis, (vii) senile plaques of patients with Alzheimer's disease, (viii) the peritoneum of CAPD patients, (ix) skin elastin in actinic elastosis, and (x) ceroid/lipofuscin deposits. These results suggest a potential role of AGEs-modification in normal aging as well as age-enhanced disease processes. This antibody named as 6D12 has been used to demonstrate AGEs-modified proteins in these human tissues, indicating potential usefulness of this antibody for histochemical identification and biochemical quantification of AGEs-modified proteins.

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Reference

- | | | |
|------------------------|--|---|
| 1 Horiuchi S. et al.: | Immunochemical approach to characterize advanced glycation end products of the Maillard reaction; Evidence for the presence of a common structure. | J. Biol. Chem. 1991 Apr 25;266(12):7329-32. |
| 2 Araki N. et al.: | Immunochemical evidence for the presence of advanced glycation end products in human lens proteins and its positive correlation with aging. | J. Biol. Chem. 1992 May 25;267(15):10211-4. |
| 3 Ikeda K. et al.: | N ϵ -(carboxymethyl) lysine protein adduct is a major immunological epitope in proteins modified with advanced glycation end products of the Maillard reaction. | Biochemistry. 1996 Jun 18;35(24):8075-83. |
| 4 Hammes H-P et al.: | Modification of vitronectin by advanced glycation alters functional properties in vitro and in the diabetic retina. | Lab Invest. 1996 Sep;75(3):325-38. |
| 5 Hemebring M. et al.: | Elimination of damaged proteins during differentiation of embryonic stem cells. | Proc Natl Acad Sci U S A. 2006 May 16;103(20):7700-5. Epub 2006 May 3. |
| 6 Chang KC. et al.: | Aminoguanidine prevents arterial stiffening and cardiac hypertrophy in streptozotocin-induced diabetes in rats. | Br J Pharmacol. 2006 Apr;147(8):944-50. |

WARNING AND PRECAUTION

- Not for diagnostic use. The safety and efficacy of product in diagnostic or other clinical uses has not been established.
- Harmful by inhalation, in contact with skin and if swallowed. Do not breathe dust. Avoid contact with skin and eyes.
- If contact with skin and eyes, wash all affected areas with large volume of water. If inhaled remove to fresh air. In severe case obtain medical attention.
- Wash hand thoroughly after handling the product.
- Do not use this product if container is broken or some contaminants are detected.
- When preserving the product, Close the container, ensure it does not fall aside or down.
- Dispose of the container and expired reagents in accordance with federal, state and local government regulations.
- Do not use the container and accessories of the product for other purpose.

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