

For research use only

# Anti Human Macrophage Surface Antigen Monoclonal Antibody (Clone No. AM-3K)

This anti-macrophage monoclonal antibody, AM-3K, was produced by using human alveolar macrophages as immunogen. AM-3K reacts intensely with most of macrophages in lymphoreticular organs and in many other organs and tissues. AM-3K also reacts with the macrophages in many pathological conditions. However, this antibody does not react with dendritic cell population, such as epidermal Langerhans cells, interdigitating cells in the paracortex of lymph nodes, nor follicular dendritic cells. Lymphocytes, granulocytes and freshly isolated blood monocytes are also negative. Reaction products for AM-3K were found on the cytoplasmic membrane of macrophages by immunoelectron microscopy.

In both cryostat sections and formalin-fixed paraffin sections, this antibody recognizes the antigen presenting on the cell surface membrane of tissue macrophages, but not monocytes or dendritic cells.

The molecular weights of the antigen recognized by AM-3K are 120 and 70 kDa.

Package Size  $50 \mu g (200 \mu L / vial)$ 

Format Mouse monoclonal antibody 0.25mg/mL

Buffer Block Ace as a stabilizer, containing 0.1% Proclin as a bacteriostat

Storage Store below –20°C

Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided


Clone No. AM-3K Subclass IgG1

Purification method The splenic lymphocytes from BALB/c mouse, immunized with human alveolar

macrophages, were fused to myeloma NS-1 cells. The screening of the hybridoma cells was performed on cryostat sections of human lung. The cell line (AM-3K) with positive reaction was grown in ascitic fluid of BALB/c mouse, from which the antibody was purified by Protein G affinity

chromatography.


Working dilution for immunohistochemistry:  $10 \mu$  g/mL, on frozen sections and paraffin sections. Antigen retrieval (microwave 10min, 0.01M citrate buffer, pH2.0) recommended.



Human liver (paraffin section):

Kupffer cells are positively stained

Takeya M., Second Department of Pathology, Kumamoto University School of Medicine, Kumamoto, Japan



Human lymph node (paraffin section):.

Macrophages in marginal sinus are positively stained.

Takeya M.,

Second Department of Pathology, Kumamoto University School of Medicine, Kumamoto, Japan

# Anti Human Macrophage Surface Antigen Monoclonal Antibody (Clone No. AM-3K)

## [Specificity]

Organ	reaction	
	positive	negative
Thymus	Macrophages in cortex	
	Macrophages in medulla	
Spleen	Red pulp macrophages	IDCs in PALS
	White pulp marginal zone macrophages	
Lymph nodes	TB macrophages in follicles	IDCs in paracortical areas
	Sinus macrophages	
Lungs	Alveolar macrophages	
Liver	Kupffer cells	
Skin	Dermal macrophages	Langerhans cells
Brain	Microglial cells	
Others		Renal tubules
		Blood monocytes

PALS=periarteriolar lymphatic sheath; TB=tingible body; IDCs=interdigitating cells

### [Interspecies reactivity]

Positive: Human, Monkey, Horse, Bovine, Pig, Goat, Dog, Cat, Rabbit, Guinea pig.

### [References]

- Zeng L., Takeya M., and Takahashi K. (1996) AM-3K, A novel monoclonal antibody specific for tissue macrophages and its applications to pathological investigation. *Journal of Pathology 178*: 207-214
- Yamate J., Yoshida H., Tsukamoto Y., Ide M., Kuwamura M., Ohashi F., Miyamoto T., Kotani T., Sakuma S., Takeya M. (2000) Distribution of cells immunopositive for AM-3K, a novel monoclonal antibody recognizing human macrophages, in normal and diseased tissues of dogs, cats, horses, cattle, pigs and rabbits. *Vet Pathol* 37(2): 168-176
- Zeng L., Takeya M., Ling X., Nagasaki A., Takahashi K. (1996) Interspecies reactivities of anti-human macropharge monoclonal antibodies to various animal species. *J Histochem Cytochem* 44(8): 845-853
- Frangogiannis NG, Burns AR., Micheal LH., Entman ML. (1999) Histochemical and morphological characteristics of canine cardiac mast cells. *Histochem J* 31(4):221-229

#### Supplier



3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013, JAPAN Tel +81(3) 5510-2347 Fax +81(3) 5510-0133 URL: http://www.sceti.jp/export/ e-mail: exp-pet@sceti.co.jp