GanglioCombi-*Light*

with Enzyme Label Mix IgG/IgM

anti-Ganglioside Autoantibodies
ELISA
(GM1, GD1b, and GQ1b)

EK-GCL 96 wells
(24 patient profiles; IgG/IgM-Mix)

Revision date: 2012-11-16

ENGLISH

INTENDED USE

BÜHLMANN GanglioCombi®-Light ELISA is designed for the quantitative determination of the clinically most relevant IgG and IgM autoantibodies directed against GM1, GD1b, and GQ1b in human serum (IgG/IgM Mix) (1-8).

PRINCIPLE OF THE ASSAY

BÜHLMANN GanglioCombi®-Light ELISA is based on the enzyme-immunometric assay technique. Gangliosides GM1, GD1b and GQ1b have stripwise been coated onto the wells of the microtiter plate. Calibrator, controls and patient sera are incubated for two hours in the microtiter wells and anti-ganglioside autoantibodies (Ab) present in the sample bind to the immobilized gangliosides. After washing off unbound substances, a horseradish-peroxidase (HRP) labeled antibody mix against human IgG and IgM is added to the wells and incubated for another two hours. Following a second washing step in which unbound antibody-enzyme reagent is removed, a substrate solution containing tetramethylbenzidine (TMB) is added to the wells. A blue color develops in proportion to the amount of anti-ganglioside autoantibodies bound to the microplate in the initial step. Color development is stopped by adding an acidic stop solution (diluted sulfuric acid) which turns the blue solution into yellow. The intensity of the color is measured at 450 nm.

The measured absorbance is proportional to the titer of anti-ganglioside antibodies present in a given sample. The titers of anti-ganglioside autoantibodies are expressed as ratios of the calibrator and can be assigned to titer categories.

BÜHLMANN GanglioCombi®-Light ELISA (EK-GCL) is designed to quantify autoantibodies of IgG and IgM isotype, both in the same test. (IgG/IgM-Mix-conjugate; 1 profile/patient). For further differentiation of IgG or IgM isotypes (individual IgG- and IgM-conjugate), refer to EK-GCL-GM (2 profiles/patient).

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Microtiter Plate precoated with gangliosides	12 x 8 wells	B-GCL-MP	Ready to use
Plate Sealer	3 pieces		
Wash Buffer Concentrate (10X) with preservatives	1 bottle 100 ml	B-GCO-WB	Dilute with 900 ml of deionized water
Incubation Buffer with preservatives	1 bottle 100 ml	B-GCO-IB	Ready to use
Calibrator Lyophilized with preservatives	1 vial	B-GCO-CA	Add 1.5 ml of Incubation Buffer
Negative, Low and Medium Control Lyophilized with preservatives	3 vials	B-GCO-CONSET	Add 1.5 ml of Incubation Buffer
Enzyme Label Mix Anti-human IgG and IgM Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 ml	B-GCO-ELGM	Ready to use
TMB Substrate TMB in citrate buffer with H ₂ O ₂	1 vial 11 ml	B-TMB	Ready to use
Stop Solution 0.25 M sulfuric acid	1 vial 11 ml	B-STs	Ready to use Corrosive agent

Table 1

STORAGE AND SHELF LIFE OF REAGENTS

Sealed / Unopened Reagents	
All sealed/unopened kit components are stable at 2-8°C until the expiration date printed on the labels.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the aluminium pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store for up to 4 months at 2-8°C.
Diluted Wash Buffer	Store for up to 4 months at 2-8°C.
Calibrator	Store for up to 1 month at 2-8°C. Do not freeze!
Controls	
Incubation Buffer	Store at 2-8° until expiration date printed on the labels.
Enzyme Label	
TMB Substrate	Store at 18-28°C.
Stop Solution	

Table 2

PRECAUTIONS

SAFETY PRECAUTIONS

- Both, Calibrator (B-GCO-CA) and Controls (B-GCO-CONSET) of this kit contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- Substrate and Stop Solution:** Substrate and Stop Solution: The Substrate TMB (B-TMB) contains Tetramethylbenzidine, hydrogen peroxide (H₂O₂) and dimethylformamide. The Stop Solution (B-STs) contains sulfuric acid (0.25 M). Each of those reagents is irritant to eyes, skin and mucous membranes. Avoid contact with Eyes, skin and cloths. Wear suitable protective clothing, gloves and eye protection. After contact with eyes or skin, wash immediately with plenty of water.
- Unused solution should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

Kit components

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- Residues in the microtiter plate** wells result from the production process. They are removed in the washing step (Assay procedure step 3) and do not affect the results.
- Steps 3-9:** Use cold (2-8°C) reagents for all these steps and keep them cold while pipetting.
- Steps 3, 6, 9:** Make sure that the wells are completely empty after the last washing cycle.
- Step 9:** Adjust TMB Substrate to room temperature (18-28°C) before using it.
- Step 11:** Shake microtiter plates during the incubation with substrate. Depending on the plate shaker, we recommend 400-600 rpm. The solution should be moved in the wells but must not spill over.
- If an **automated washer is used**, "plate mode" should be chosen so that dispensing is performed sequentially on all strips before aspirating.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.

- Microwells cannot be re-used.

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 20 µl, 100 µl and 1 ml pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the reconstitution of the wash buffer.
- Squeeze bottle for wash buffer or automatic microtiter plate washer.
- Blotting paper.
- Microtiter plate rotator.
- Microtiter plate reader for measurement of absorbance at 450 nm.

SPECIMEN COLLECTION AND STORAGE

- The procedure requires <0.1 ml of blood and <50 µl of serum, respectively.
- Refer to page 4 to learn about the interference of hemolyzed, lipemic or icteric samples.
- Collect blood into plain tubes (no anti-coagulant), avoid hemolysis, leave to clot for one hour, centrifuge for 10 minutes at approximately 1500 x g at room temperature (18-28°C), collect the serum.
- We recommend freezing aliquots of patient samples if you need to store samples in order to avoid repeated freezing/thawing.
- Store serum samples at ≤ -20°C up to 4 months. For long-term storage we recommend -70°C (samples are stable for >1 year). Frozen samples should be thawed and vortexed thoroughly prior to use.

ASSAY PROCEDURE

1. **Dilute** all patient samples to be investigated 1:50 with Incubation Buffer. Use 30 µl of serum + 1470 µl of Incubation Buffer. Mix by vortexing and leave diluted samples and reconstituted calibrator and controls for 30 minutes at 2-8°C prior to pipetting (refer to step 4 a-c).
2. **Prepare a plate-frame** with the required number of strips to test the patient samples. Reseal the remaining strips in the foil pouch together with the desiccant packs **immediately after usage**. Store refrigerated.

Note: Use cold reagents in steps 3 to 9.

3. **Wash** the coated wells twice using at least 300 µl of **cold** Wash Buffer per well. Empty the wells and tap the plate firmly onto blotting paper to remove remaining liquid completely.

Note: Immediately proceed to the next step (4a-e).

- 4a. Calibrator: Pipet 100 µl of the Calibrator into the well A1 (refer to Figure 1).
- 4b. Controls: Pipet 100 µl of the Medium Control into well B1, Low Control into well A2 and Negative Control into the well B2 (refer to Figure 1).

Note: If more than three strips per run are used, Calibrator and Controls can be tested in duplicates (see Figure 1).

- 4c. Patient serum: Pipet 100 µl of diluted patient serum 1 into the wells C1 - E1 (refer to Figure 1).
- 4d. Patient serum: Pipet 100 µl of diluted patient serum 2 into the wells F1 - H1 (refer to Figure 1)
- 4e. Pipet 100 µl of diluted patient sera 3-12 into the subsequent wells (refer to Figure 1).

5. Cover the plate with a Plate Sealer and **incubate** for 2 hours ± 5 minutes at 2-8°C (do not shake the plate).
6. Remove the Plate Sealer. Empty the wells and **wash** three times using at least 300 µl of **cold** Wash Buffer (2-8°C) per well. Empty the wells and strike the plate firmly onto blotting paper in order to remove washing buffer completely.
7. Add 100 µl of **Enzyme Label Mix** to all wells.
8. Cover the plate with a Plate Sealer and **incubate** for 2 hours ± 5 minutes at 2-8°C (do not shake the plate).
9. Remove the Plate Sealer. Empty the wells and **wash** three times using at least 300 µl of **cold** Wash Buffer (2-8°C) per well. Empty the wells and strike the plate firmly onto blotting paper.

Note: Let TMB Substrate Solution adjust to room temperature (18-28°C).

10. Add 100 µl of **TMB Substrate Solution** to each well.
11. Cover the plate with a Plate Sealer, incubate the plate on a **plate shaker** at 400-600 rpm incubate for 30 minutes ± 2 minutes at 18-28°C protecting the plate from direct light.
12. Add 100 µl of **Stop Solution** to all wells. Proceed to step 13. within 30 minutes.
13. Read the **absorbance** at 450 nm in a microtiter plate reader.

STANDARDIZATION

The Calibrator included in this kit has been calibrated against internal reference material. It has been adjusted to **200 % ratio**.

RESULTS AND CALCULATION

Calculation of results:

1. Record the absorbance (OD) at 450 nm for each well (Calibrator, Controls and patient samples).
2. Average the duplicate Calibrator and Control values (if available).
3. Results are expressed as ratio of the absorbance of samples and the (averaged) absorbance of the Calibrator

$$\% \text{ Ratio} = \frac{\text{absorbance of samples and Controls}}{\text{absorbance of Calibrator}} \times 200$$

Programs to calculate results as % ratio are available on most microplate readers.

Note: The results presented in Table 11 are examples. Calibrator and Controls must be used in each individual assay.

QUALITY CONTROL

A good understanding of this instruction for use is necessary to obtain reliable results. These will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following the instruction for use.

Since there is no control serum for anti-ganglioside antibodies commercially available, we recommend using a positive, and negative serum pool for internal quality control.

All controls must be within established confidence ranges (% ratio). The confidence ranges of the Controls are lot-specific and printed on the QC data sheet delivered with this kit.

Performance characteristics should be within established limits. If these characteristics are not in conformity with established limits and repetition excludes handling failures, check the following issues: i) Have all reagents, used in step 3-9, been kept at 2-8°C? ii) accuracy of the pipets, thermometers, and timers, iii) settings of ELISA washer and reader, iv) expiration date of the reagents v) storage and incubation conditions vi) color of the TMB Substrate Solution (should be colorless) vii) purity of the water.

PERFORMANCE LIMITATIONS

Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

PERFORMANCE CHARACTERISTICS

Intra-Assay Precision (Within-Run): 3.3 %. The intra-assay precision was calculated from results of 12 values of two IgM and/or IgG samples in a single run. The results in Table 12 were calculated as % ratio as described in "Results and Calculation".

Inter-Assay Precision (Run-to-Run): 10.5 %. The inter-assay precision has been determined by measuring two different serum samples with antibodies to 3 gangliosides in 20 different runs. The results in Table 13 have been calculated as % ratio as described in "Results and Calculation".

Detection limit (LoB): 12 Incubation Buffer replicates were assayed in a single run. The detection limit expressed as the % ratio of the calibrator was calculated to be ≤ 10 %.

Detection limit (LoQ): The inter- and intra-assay precision of serum samples between 9 to 307 % ratio was determined. The mean values and %CV were calculated for each sample. Figure 2 shows a polynomial X/Y plot of intra- and inter-assay precision vs. % ratio. The resulting curve was below the upper limit of 15% and 20% CV for intra and inter-assay precision, respectively, manifesting reliable results.

Linearity: The linear range of the test system was assessed according to CLSI guideline EP06-A. The system is linear in the diagnostic relevant range between 40 and 200% ratio. Results above 200% ratio are assessed clinically correct and can be diluted into the linear range by an additional 1:5/1:10 dilution.

Specificity: Different human serum samples containing specific anti-ganglioside IgM and/or IgG autoantibodies were incubated over night with the corresponding soluble antigen in different concentrations and subsequently tested in the BÜHLMANN GanglioCombi®-Light ELISA according to the assay procedure. Specificity of the autoantibody binding was demonstrated by inhibition with the corresponding antigen at concentrations between 1 and 100 µg/ml (data not shown).

INTERFERING SUBSTANCES

No interference is detected with the following substances up to the following concentrations: Triglycerides (Intralipid®): 3000 mg/dL; conjugated bilirubin: 60mg/dL; unconjugated bilirubin: 40 mg/dL and hemoglobin: 400 mg/dL.

EXPECTED VALUES AND CUT-OFF

In co-operation with the institutions mentioned below, a cut-off of 50% has been established. Values < 30% have been clearly classified as negative. The established titer categories are based on n=100 blood donors¹ (adult men and women between 18 to 70 years of age) and n=277 pathological samples². Sera were assayed for each of the three gangliosides (GM1, GD1b and GQ1b), according to the assay procedure. Results of normal blood donors are shown in Table 14.

Guidelines for the Use of Cut-off and Titer categories/ratios (%):

Category	Ratio (%)
Negative	< 30%
Grey zone	30-50%
Cut-off	50%
Positive	> 50-100%
Strongly Positive	> 100 %

¹) collected at the blood donation centre, University Hospital of Basel

²) received from Friedrich-Baur-Institute, Ludwig-Maximilians-University of Munich; Department of Neurology, University of Basel; Department of Neurology, University of Lyon.

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Figure 1 Microtiter plate set-up
IgG/IgM-Mix conjugate

EK-GCL												
EK-GCL IgG/ IgM Mix												
	1	2	3	4	5	6	7	8	9	10	11	12
Calibrator	CAL	CTRL	CAL	CTRL								
CTRL	Med	Neg	Med	Neg								
GD1b												
GQ1b												
GM1												
GD1b												
GQ1b												
GM1												

24 patients / Kit (1 MP / Kit)
minimal load: 4 patients / run

Table 11 Example of Results

B-GCO IgG/IgM Mix-ELGM	Absorbance (OD ₄₅₀)	Ratio [%]	Result
Calibrator	1.415		
	1.445		
Calibrator Avg.	1.430	200.0	
Medium Control	0.498	69.6	
	0.482	67.6	
Med. Control Avg.	0.490	68.6	
Low Control	0.195	27.2	
	0.191	26.8	
Low Control Avg.	0.193	27.0	
Negative Control	0.090	12.6	
	0.100	14	
Neg. Control Avg.	0.095	13.4	
Sample 1 GM1 (well D)	0.544	76	Positive
Sample 1 GD1b (well G)	0.745	104	Strongly positive
Sample 1 GQ1b (well H)	0.090	13	Negative

Table 12 Intra-Assay Precision (Within-Run)

Ganglioside	Mean [% Ratio]	SD [% Ratio]	CV [%]
GM1	85	2.4	2.9
	157	5.1	3.2
GD1b	111	4.0	3.6
	310	7.9	2.6
GQ1b	287	9.2	3.2
	83	3.8	4.5
Mean			3.3

Table 13 Inter-Assay Precision (Run-to-Run)

Ganglioside	Mean [% Ratio]	SD [% Ratio]	CV [%]
GM1	75	7.1	9.5
	183	27.8	15.2
GD1b	95	9.6	10.1
	307	26.8	8.8
GQ1b	92	7.5	8.2
	263	26.6	10.1
Mean			10.5

Figure 2 Precision Profile


Table 14 Blood Donors

Category	anti-Ganglioside-autoantibodies (%)
	GM1, GD1b, GQ1b
Negative	89-98
Greyzone	1-7
Positive	1-7
Strongly positive	0-1

Legend table 14: Distribution frequency of individual anti-Ganglioside-autoantibodies (%) tested with EK-GCL (IgG/IgM Mix) and classified into titer categories (n=100); calculated over the three gangliosides:

Negative : 89-98 % of blood donors

Grey zone: 1-7% of blood donors

Positive: 1-7% of blood donors

Strongly positive: 0-1 % of blood donors

Table description: cf. "Results and Calculations", "Performance Characteristics", and "Expected Values and Cut-Off" (pages 3/4)

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
1. Willison HJ and Yuki N: *Peripheral neuropathies and anti-glycolipid antibodies*. Brain 125, 2591-2625 (2002).
2. Latov N: *Antibodies to glycoconjugates in neurological disease*. Clin Aspects Autoimm 4, 18-29 (1990).
3. Pestronk A: Invited Review: *Motor neuropathies, motor neuron disorders, and antiglycolipid antibodies*. Muscle Nerve 14, 927-936 (1991).
4. Steck AJ: *Use of anti-glycoconjugate antibody assays in neuropathy*. Workshop presentation at the Nat. Meeting Am. Acad. Neurology, Seattle, May 6-13, 1995.
5. Steck AJ and Kappos L: *Changing concepts in inflammatory and paraneoplastic neuropathies*. Curr Neurol 16, 191-212 (1996).
6. Humbel RL and Schmit P: *Anticorps antigangliosides et neuropathies peripheriques*. Rev Med Liege 51, 368-375 (1996).
7. Quarles RH and Weiss MD: *Autoantibodies associated with peripheral neuropathy*. Muscle Nerve 22, 800-822 (1999).
8. Nobile-Orazio E: et al: *How useful are anti-neural IgM antibodies in the diagnosis of chronic immune-mediated neuropathies?* JNS 266, 156-163 (2008).

BÜHLMANN GanglioCombi® Light ELISA

TIME TO RESULT: 4.5 HOURS

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
MP	Microtiterplate Mikrotiterplatte Plaque de microtitration Microplaca Microplaca
BUF WASH 10X	Wash Buffer Concentrate (10x) Waschpuffer-konzentrat (10x) Concentré de tampon de lavage Tampón de lavado concentrado (x10) Tampón de lavado concentrado (x10)
BUF INC	Incubation Buffer Inkubations-Puffer Tampón d'incubation Tampón de incubación Tampón de incubación

Symbol	Explanation
CONTROL -	Negative Control Negativkontrolle Contrôle négatif Control negativo Control negativo
CONTROL L	Low Control Kontrolle tief Contrôle faible Control bajo Control bajo
CONTROL M	Medium Control Kontrolle mittel Contrôle moyen Control medio Control medio
CAL	Calibrator Kalibrator Calibrateur Calibratore Calibrador
EL MIX	Enzyme Label Mix Enzymmarker-Mix Marqueur enzymatique Mix Marcato enzimatico Mix Mezcla de marcadores enzimáticos
SUBS TMB	TMB Substrate TMB Substrat Substrat TMB Substrato di TMB Substrato de TMB
SOLN STOP	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de parada

Printing Date
2013-01-30