GHB

Gamma-Hydroxybutyric Acid

Enzymatic Assay

KK-GHB


Revision date: 2010-05-21

ENGLISH

INTENDED USE

The BÜHLMANN gamma-hydroxybutyric acid (GHB) kit has been designed for the direct and quantitative determination of GHB in human urine and serum. The assay can be adapted to clinical chemistry analyzers according to specific instrument protocols.

PRINCIPLE OF THE ASSAY

The GHB is converted to its metabolite, (SSA, succinic semialdehyde), by the action of GHB-specific recombinant enzyme (Rec. Enzyme) and oxidized nicotinamide adenine dinucleotide (NAD^+). The increase in absorbance at 340 nm resulting from the reduction of NAD^+ into NADH is proportional to the amount of GHB in the sample.

REAGENTS SUPPLIED AND PREPARATION

Reagents ¹⁾	Quantity	Code	Reconstitution
Incubation Buffer	1 vial 12 ml	B-GHB-IB	Ready to use
Cofactor lyophilized NAD^+	1 vial	B-GHB-CF	Add 5.6 ml of deionized water ⁴⁾
Enzyme lyophilized recombinant Enzyme	2 vials	B-GHB-E	Add 4.2 ml of deionized water ⁴⁾ , do not vortex
Calibrators ²⁾ lyophilized GHB	2 vials	B-GHB-CASET	Add 2 ml of deionized water ⁴⁾
Controls Low / High ³⁾ lyophilized GHB in human urine	2 vials	B-GHB-CONSET	Add 2 ml of deionized water ⁴⁾

Table 1

- ¹⁾ All reagents contain sodium azide (<0.1%) as preservative.
²⁾ After reconstitution the GHB concentration of the Calibrators is 10 and 100 mg/L, respectively.
³⁾ The Controls contain lot-specific amounts of GHB. Refer to the QC data sheet for actual concentrations.
⁴⁾ For reconstitution see chapter Procedural Notes.

STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
The unopened kit components are stable at 2 - 8 °C, except Cofactor and Enzyme that must be stored at -20 °C upon receipt. Do not use reagents after expiration of the kit.	
Opened / Reconstituted Reagents	
Incubation Buffer	Store for up to 2 months at 2 - 8 °C
Cofactor	
Enzyme	Store for up to 15 days per vial (1 month in total) at 2 - 6 °C; refer to Procedural Notes
Calibrators	Store for up to 2 months at 2 - 8 °C
Controls	

Table 2

WARNINGS AND PRECAUTIONS

- B-GHB-CASET and B-GHB-CONSET contain GHB at <100 mg/L. At this concentration GHB is considered as non-toxic.
- All reagents contain sodium azide at <0.1%. Sodium azide may react with copper or lead and form explosive compounds in metal drain lines. Upon disposal, flush with a large volume of water to prevent azide build up.
- Do not ingest.
- Use only deionized H_2O for reconstitution. See chapter Procedural Notes.

- Calibrators and Controls contain components of human urine. Thus the risk of transmitting pathogens of known or unknown origin cannot be completely excluded.
- Urine samples should be handled as potentially infectious.
- All products and samples containing human material should be handled in accordance with good laboratory practice using appropriate precautions.

MATERIALS REQUIRED BUT NOT PROVIDED

- Accurate pipeting devices (1, 5 or 10 ml pipettes)
- Test tubes/rack
- Clinical chemistry analyzer with an optical filter at 340 nm
- Deionized water
- 0.9 % saline solution

SPECIMEN COLLECTION AND STORAGE

- Depending on the dead volume of the analyzer, ~110 μl sample volume must be transferred into the instrument specific sample tubes. Less than 10 μl of urine or serum sample is needed to run the test. For details refer to the instrument manual.
- Urine: Collect urine (pH 5-8) without additives, keep the samples refrigerated, and have them filtered or centrifuged. Collect the supernatant.
- Serum: Refer to page 4 to learn about the interference of hemolytic, lipemic or icteric samples. Collect blood into plain tubes (no anti-coagulant), avoid haemolysis, leave to clot for one hour, centrifuge for 10 minutes at approximately 1500 x g at room temperature (18-28 °C), and collect the supernatant.
- Sample storage: Urine and serum samples can be stored at 2-8 °C for at least two weeks. For longer storage keep the samples frozen at -20 °C. GHB stability in samples has intensively been investigated (8).

PROCEDURAL NOTES

Reconstitution of Reagents:

- Calibrator and Control vials are under vacuum. Thus, it is mandatory to remove the caps slowly and carefully in order to prevent loss of material.
- Cofactor (B-GHB-CF), Calibrators (B-GHB-CASET) and Controls (B-GHB-CONSET): Vortex the vials for 30 seconds after reconstitution or use a suspension mixer for 15 minutes.
- Enzyme (B-GHB-E): **DO NOT VORTEX!** Use a suspension mixer for 15 minutes at RT until the lyophilized Enzyme has completely dissolved. Alternatively, gently swivel the vial after reconstitution, leave it 15 minutes at RT, and gently move it upside-down afterwards. **Do not leave the vial at RT for more than 15 minutes.** Reconstitute the Enzyme refrigerated at 2-6 °C, if you need to extend the reconstitution time.
- Important:** The Enzyme must be kept at 2-6 °C. If the temperature range of the clinical chemistry analyzer in use does not fulfil this requirement, **reagents should be recapped and stored in the refrigerator below 6 °C.** On board, a temperature above 6 °C may cause a loss of enzyme activity.

Thus the OD obtained in the end point measurement might decrease. Controls included in the kit can be used to monitor the assay performance. If the controls do not fit to the range indicated in the QC data sheet a recalibration of the standard curve is mandatory.

ASSAY PROCEDURE

GHB enzymatic assay is designed for automated analysis on clinical chemistry analyzers. It has to be programmed as a three reagent test. For details refer to the analyzer specific protocols.

The following protocol has been established for KoneLab 30 (Thermo) and is used here as a general example for the test procedure. For applications on other instruments, we will provide assistance.

1. 100 µl of Incubation Buffer (R1)
+ 8 µl of urine or serum sample (S)
+ 7 µl of deionized water
2. +50 µl of Cofactor (R2)
→ Incubate for 1.5-2 min at 37 °C
3. +85 µl of Enzyme (R3)
4. T_{0 min.}: Measurement at 340 nm (Blank)
→ Incubation for 5-6 min at 37 °C
5. T_{5-6 min.}: Measurement at 340 nm

RESULTS

Calibration

Use endpoint mode at 340 nm. The absorbance obtained after reading 1 (M1) and 2 (M2) is calculated using the absorbance differences (M2-M1). The standard curve will be calculated from the two calibrators using a linear regression mode. Refer to Table 11 and Figure 1 for typical results and standard curve. Depending on the instrument it might be necessary to read additionally at a secondary wavelength e.g. (700 nm).

Important: Calibration of the assay should be carried out **every 10 days**, when a **fresh vial of Enzyme** is reconstituted or when the controls are out of the confidence range.

Samples and Controls

The absorbance at 340 nm is recorded for the Calibrators and the GHB concentration for each Control and sample is calculated as indicated above.

STANDARDIZATION

The GHB Calibrators have been calibrated against commercially available GHB quantified with HPLC.

QUALITY CONTROL

- The values of the Low and High Control provided with the kit must be within the lot specific range indicated on the corresponding QC data sheet. Otherwise, the standard curve has to be re-calibrated.
- It is good laboratory practice to record the following data for each assay: kit lot number, reconstitution dates of kit components, results for Calibrators and Controls, and internal urine or serum pool.
- The precision and expected values of standard curve and controls should be within established limits. The confidence limits for the controls are lot-specific and printed on the QC data sheet delivered with the kit.
- Reliable results will only be obtained by using precise laboratory techniques (current GLP guidelines) and by accurate attention of this instruction for use.

LIMITATIONS

Samples exceeding 250 mg/L should be diluted 1:10 with 0.9% saline solution (1 volume sample + 9 volumes saline solution). The results have to be multiplied by 10.

If the blank of urine samples exceeds 1.0 OD at 340 nm, results may be falsely positive.

Positive results should be confirmed by alternative methods.

Substances and/or factors other than those investigated in the specificity study may interfere with the test and cause false results.

CALCULATION

Automated procedure: For calculation of the results use the endpoint mode with two calibrators (10 and 100 mg/L). Refer to the instrument manual for further details.

Dynamic range: 5 - 250 mg/L

PERFORMANCE CHARACTERISTICS

Assay performance characteristics have been determined on KoneLab 30.

Intra-assay precision: Urine: 0.8 - 4.6 %; serum: 1.0 - 3.8 %. The intra-assay precision has been determined by repeated measurements of one urine and serum sample, each spiked with 3 different concentrations of GHB (10, 50, 100 mg/L) in a single run. The mean value is 2.1 % in urine and 2.0 % in serum (Table 12).

Between-days precision: Urine: 2.4 - 7.4 %; serum: 2.7 - 8.3 %. The between-days precision has been determined according to NCCLS EP5-A2 by repeated measurements of two controls and one urine and serum sample, each spiked with 3 different concentrations of GHB (10, 50, 100 mg/L) in 2 runs per day over a period of 20 work days. The mean value is 4.0 % in urine and 4.7 % in serum (Table 13).

Dilution linearity: Urine 98 - 102 %; serum 100 - 110 %. Two urine and serum samples with elevated GHB concentration have been diluted with 0.9% NaCl solution (Table 14). The mean value is 100 % in urine and 105 % in serum. The dilution linearity (observed vs. expected) has been calculated in a concentration range between approximately 100 and 6 mg/L. It is linear between 5 and 250 mg/L.

Spiking recovery: Urine: 102 - 115 %; serum: 103 - 112 %. Five different urine and serum samples have been spiked with increasing amounts of GHB and analyzed in 3 runs. The mean value is 105 % in urine and 109 % in serum (Table 15).

Analytical sensitivity: 1.5 mg/L. The analytical sensitivity is dependent on the precision of the clinical chemistry analyzer used. By repeated measurements with 0.9% NaCl (n= 60) the imprecision of the analyzer was determined by calculating the mean +3 SD value and converted it into GHB mg/L (Table 16).

Functional sensitivity: Urine: 2.8 mg/L; serum: 4.5 mg/L. The functional sensitivity of serum and urine was determined by repeated measurements of 5 urine and serum samples at concentrations between 0 and 250 mg/L. The sensitivity at 15 % CV was determined to be 2.8 and 4.5 mg/L, respectively (Figure 2).

Specificity: The substances listed in Table 17 have been analyzed between 10 and 1000 mg/L to determine the enzyme specificity.

INTERFERING SUBSTANCES

Therapeutic drugs and drugs of abuse: The **therapeutic drugs** and **drugs of abuse** listed in Table 18 and Table 19 were evaluated according to the approved guideline for interference testing in clinical chemistry, EP7-A2 (7), on the KoneLab 30. No interference has been observed up to the listed concentrations.

The interference of **Ethanol** is shown in Figure 3. Up to 3‰ the measured GHB concentration is below 10 mg/L. 1 g/L Ethanol raises the GHB value by 3 mg/L.

Serum Indices: No interference is detected with the following substances up to the listed concentrations: **Triglycerides** (Intralipid® 275 mg/dL; equivalent to 7.7 mmol/L triglycerides), **conjugated bilirubin** (360 µmol/L; 30 mg/dL), **unconjugated bilirubin** (513 µmol/L; 30 mg/dL) or **haemoglobin** (3.1 mmol/L; 500 mg/dL) on KoneLab 30.

EXPECTED VALUES

The following reference range for urine and serum (97.5th percentile) has been determined from individuals, who did not consume GHB:

Urine: 97.5 th percentile:	10.6 mg/L
Median:	3.3 mg/L
n = 75	
Serum: 97.5 th percentile	4.0 mg/L
Median:	0.68 mg/L
n = 50	

METHOD COMPARISON

34 urine samples have been compared to a published IC-GHB method (6) (Figure 4):

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Table 12: Intra-Assay Precision

Sample spiked with GHB	Mean [mg/L]	SD [mg/L]	CV [%]
Urine + 10 mg/L	16.4	0.75	4.6
Urine + 50 mg/L	57.9	0.49	0.8
Urine + 100 mg/L	107.8	0.85	0.8
Mean			2.1
Serum + 10 mg/L	10.8	0.42	3.8
Serum + 50 mg/L	54.1	0.69	1.3
Serum +100 mg/L	104.4	1.05	1.0
Mean			2.0

Table 13: Between-Day Precision

Sample spiked with GHB	Mean [mg/L]	SD [mg/L]	CV [%]
Urine + 10 mg/L	12.2	0.9	7.4
Urine + 50 mg/L	52.6	1.21	2.3
Urine + 100 mg/L	100.4	2.36	2.4
Mean			4.0
Serum + 10 mg/L	11.9	0.99	8.3
Serum + 50 mg/L	55.0	1.71	3.1
Serum +100 mg/L	106.9	2.89	2.7
Mean			4.7
Control Low	14.6	0.77	5.3
Control High	75.1	1.43	1.9

Table 14: Dilution Linearity

Example 1 of 2	Dilution	Observed [mg/L]	Expected [mg/L]	O/E [%]
Serum1	1:1	279.9	299.0	94
	1:3	103.3	99.7	104
	1:6	53.3	49.8	107
	1:12	26.5	24.9	106
	1:24	13.0	12.5	104
	1:48	6.0	6.2	96
	1:96	2.9	3.1	92
Serum 1 mean				100
Serum 2 mean				110
Mean serum				105
Urine 1	1:6	580.5	580.5	100
	1:12	297.3	290.3	102
	1:24	149.9	145.1	103
	1:48	76	72.6	105
	1:96	38	36.3	105
	1:192	18	18.1	99
Mean urine 1				102
Mean urine 2				98
Mean urine				100

Table 15: Spiking Recovery

Example 1 of 5	Native [mg/L]	spiked with [mg/L]	Observed [mg/L]	Expected [mg/L]	O/E [%]
Serum1	1.6	10	11.9	11.6	100
		25	28.1	29.6	105
		50	54	51.6	106
		100	102	101.6	102
Mean serum 1					103
Mean serum 2					110
Mean serum 3					110
Mean serum 4					109
Mean serum 5					112
Mean serum					109
Urine1	6.1	10	16.4	16.1	102
		25	32.2	31.1	104
		50	59.5	56.1	106
		100	104.7	106.1	99
Mean urine 1					102
Mean urine 2					102
Mean urine 3					115
Mean urine 4					106
Mean urine 5					102
Mean urine					105


Table 16: Analytical Sensitivity (AS)

Mean (ΔOD)	SD (ΔOD)	Mean+3SD (ΔOD)	AS [mg/L]
0.00034	0.00056	0.00202	1.5

(n=60)

Figure 2: **Functional Sensitivity**

Precision profile for urine and serum

Table 17: **Enzyme Specificity**

Component/substrate	Substrate specificity
γ -hydroxyvaleric acid, (GHV)	<0.1 %
γ -butyrolactone, (GBL)	4.0 %
1,4-butanediol, (1,4-BD)	<0.1 %
α -hydroxybutyric acid, (AHB)	<0.1 %
β -hydroxybutyric acid, (BHB)	<0.1 %
γ -hydroxyvaleric acid, (GHV)	<0.1 %
succinic acid	<0.1 %
γ -valerolactone, (GVL)	<0.1 %

InterferencesTable 18 **Interference of therapeutic drugs**

Therapeutic drugs	No interference up to
Amikacin	53.2 μ mol/L
Caffeine	84.5 μ mol/L
Carbamazepine	75.4 μ mol/L
Cyclosporine	480 nmol/L
Digoxin	4.43 nmol/L
Ethosuximide	970 μ mol/L
Gentamicin	18.6 μ mol/L
Lithium	1.77 mmol/L
Lithium (Vitros)	2.25 mmol/L
Methotrexate	7.78 μ mol/L
Paracetamol	1.44 mmol/L
Phenobarbitone	223 μ mol/L
Phenytoin	74.9 μ mol/L
Primidone	58.8 μ mol/L
Salicylic acid	2.86 mmol/L
Theophylline	177 μ mol/L
Tobramycin	19.5 μ mol/L
Valproic acid	998 μ mol/L
Vancomycin	20.3 μ mol/L

Table 19 **Interference of Drugs of abuse**

Reagent	No interference up to
Amphetamines (d-Methamphetamines)	1250 ng/mL
Barbiturates (Secobarbital)	375 ng/mL
Benzodiazepines (Oxazepam)	390 ng/mL
Cannabinoids (11-Nor- Δ -9-THC-9-COOH)	65 ng/mL
Cocaine (Benzoyllecgonine)	500 ng/mL
LSD	1 ng/mL
Methadone	375 ng/mL
Opiates (Free Morphine)	2500 ng/mL
Phencyclidine	31 ng/mL
Propoxyphene (Norpropoxyphene)	375 ng/mL
Nortryptiline	375 ng/mL



Figure 3 **Interference of Ethanol**Figure 4 **Correlation between KK-GHB and IC-GHB**

Table description: “Calculation” (pg. 3), “Performance Characteristics” (page.4), and “Interpretation of Results” (page. 4).

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1. Bravo DT et al.: *Reliable, Sensitive, Rapid and Quantitative Enzyme-Based Assay for Gamma-Hydroxybutyric Acid (GHB)*. J Forensic Sci, Mar. 2004, Vol. **49**, No.2. Paper ID JFS2003165
2. Mari F et al.: *What constitutes a normal ante-mortem urine GHB concentration*. J of Forensic and Legal Medicine **16** (2009) 148-151
3. Andresen H et al.: *Liquid Ecstasy - A Significant Drug Problem*. Deutsches Ärzteblatt Int. 2008; **105**(36): 599-603
4. Haller Ch. Et al.: *GHB urine concentrations After Single-dose Administration in Humans*. J Anal Toxicol. 2006; **30** (6): 360-364
5. Elliott S.: *The many faces of gamma-hydroxybutyrate (GHB)* Drug Monitor Vol. 4 Issue 2, 2003
6. M.Jordi et al.: *GHB Determination with Ion Chromatography* Therap. Drug Monit. 25 (2003) 486
7. Clinical Laboratory Standards Institute: *Interference Testing in Clinical Chemistry*; approved guideline (EP7-A2), 2005
8. H. Andresen, et al.: *Gamma-Hydroxybutyrate in Urine and Serum: Additional Data Supporting Current Cut-Off Recommendations*, Forensic Sci. Int. (2010), doi:10.1016/j.forsciint.2010.03.035

Symbol	Explanation	Symbol	Symbol
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad	BUF INC	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone di incubazione Tampón de incubación
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo	COF	Cofactor Kofaktor Cofacteur Cofattore Cofactor
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote	ENZ	Enzyme Enzym Enzyme Enzime Enzime
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>	CAL A	Calibrator A Kalibrator A Calibrateur A Calibratore A Calibrador A
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos	CAL B	Calibrator B Kalibrator B Calibrateur B Calibratore B Calibrador B
	Consult Instructions for Use Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso	CONTROL L	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura	CONTROL H	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
	Upper limit of temperature Temperaturobergrenze Limite de température maximale Limite de temperatura maximale Limite de temperatura maximale		

Printing Date
2010-06-17