Quantum Blue[®] Calprotectin

**Quantitative
Lateral Flow Assay**

LF-CAL20-RD 20 tests

Revision date: 2010-06-02

ENGLISH

INTENDED USE

The BÜHLMANN Quantum Blue® Calprotectin Assay is an immunoassay designed for the quantitative determination of Calprotectin in human stool samples (1-5) in combination with the BÜHLMANN Quantum Blue® Reader. For laboratory use only.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of Calprotectin antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) highly specific for Calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the extracted and diluted stool sample. The Calprotectin/anti-Calprotectin gold conjugate binds to the anti-Calprotectin antibody coated on the test membrane (test line; test band) and the remaining free anti-Calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (control line; control band). The signal intensities of the test line and the control line are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

The BÜHLMANN Quantum Blue® Calprotectin Assay is standardized with the BÜHLMANN Calprotectin ELISA (order code: EK-CAL).

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cartridge vacuum-sealed in a foil bag	20 pieces	B-LFCAL-TC	
Extraction Buffer	1 bottle 120 ml	B-CAL-EX	Ready to use
Extraction Device	1 bag	B-CAL-RD20	20 Tubes 20 Base Caps 4x5 Spatulas
RFID Chip Card	1 piece	B-LFCAL-RCC	

Table 1

REAGENTS SUPPLIED UPON REQUEST

Controls for BÜHLMANN Quantum Blue® LF-CAL20-RD can be ordered separately.

Controls Low* / High*	2 vials	B-CAL-CON	Ready to use
	0.5 ml		

* Lot specific concentration of Calprotectin

Table 2

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8°C until the expiration date printed on the labels.

WARNINGS AND PRECAUTIONS/REAGENTS

Human Source Material: None of the reagents provided with this kit is of human origin.

Patient specimens should be handled as potentially biohazardous, taking appropriate precautions.

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortexing device for stool extraction
- Precision pipettes with disposable tips: 10-200 µl and 1 ml
- Centrifuge for Microtubes
- Eppendorf tubes for dilution of the extracts (1:16)
- Timer (optional)
- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)

- Gloves and Lab coat
- Soft tissues or blotting paper

SPECIMEN COLLECTION AND STORAGE

Using the Extraction Device included in the test kit, less than 100 mg of native stool is required.

Collect stool samples into clean tubes and store refrigerated at 2-8°C for up to 6 days.

Freezing of the sample may result in slightly increased Calprotectin values due to release from neutrophils present in the sample. Therefore, it is of advantage for longer storage, to keep the extracts at -20°C. The extracts are stable for at least 4 months.

Important: The sample must be collected without any chemical or biological additions in the collection device.

ASSAY PROCEDURE

The Assay Procedure consists of two individual parts:

- Extraction of stool samples
- Lateral Flow Assay/Reading

The Assay Procedure is described and illustrated in details below on pages 15ff.

PROCEDURAL NOTES

- All reagents and test samples must be equilibrated at room temperature (18-28°C) before starting the assay. The Quantum Blue® Reader must be switched on and programmed for the Calprotectin assay before starting the assay (see Quantum Blue® Reader Manual).
- In order to receive reliable and quantitative results it is important to homogenize the stool sample entirely in the Extraction Device. Avoid contamination with stool at the top of the tube.
- In order to receive reliable and quantitative results it is important to centrifuge the homogenized extracts. The centrifugation can be done either before or after the dilution step.
- Diluted samples should be used within several hours and cannot be stored for a longer time period.
- Undiluted extracts can be stored at -20°C for at least 4 months. It is important to centrifuge the extracts before storage (5 min at 3.000 x g). After centrifugation the supernatant must be transferred into a fresh storage tube.

QUALITY CONTROL

- A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by following this instruction for use accurately.
- Patient samples that are not properly handled may cause inaccurate results.
- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices, ii) expiration dates of reagents and iii) storage and incubation conditions.

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see Figures 1A and 1B). It is used as functional test control only and can not be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 12 minutes of incubation time (Figure

1A), the concentration of Calprotectin present in the stool sample is below the detection limit. If a Test Line (T) is detectable after 12 minutes of incubation time (Figure 1B), the Calprotectin concentration present in the stool sample is calculated by the Quantum Blue® Reader.

- If only the Test Line (T) is detectable after 12 minutes of incubation time (Figure 1C), the test result is invalid and the Calprotectin assay has to be repeated using another test cartridge.
- If neither the Control Line (C) nor the Test Line (T) are detectable after 12 minutes of incubation time (Figure 1D), the test result is invalid and the Calprotectin assay has to be repeated using another test cartridge.
- As the Quantum Blue® Reader allows a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a lot-specific threshold after 12 minutes of incubation time, the test result is also invalid and the Calprotectin assay has to be repeated using another test cartridge.

STANDARDIZATION AND INTERPRETATION OF RESULTS

- The Lateral Flow Assay is calibrated with the BÜHLMANN Calprotectin ELISA (order code: EK-CAL).
- The BÜHLMANN Quantum Blue® Reader uses a lot-specific standard curve to calculate the Calprotectin concentration. This lot-specific standard curve is generated with the median values ($n \geq 5$ measurements each) from at least 10 Calibration points and at least five different stool samples with known Calprotectin concentrations. The assay range is between 30 and 300 µg/g.
- For quantitative measurements, unknown samples reading above 300 µg/g should be further diluted (e.g. 1:10) with Extraction Buffer and assayed again according to the procedure. The resulting dilution factor must be multiplied by the measured concentration to obtain the final results.
- Calprotectin values should be used as supplementary data available to the physician in establishing a diagnosis.

EXPECTED VALUES

Estimation of faecal Calprotectin is a reliable and simple way to distinguish organic from functional gastrointestinal diseases.

The following data were established with the BÜHLMANN Calprotectin ELISA (order code: EK-CAL). In a clinical study, 401 symptomatic consecutive patients scheduled for colonoscopy were investigated (Table 11). Endoscopy examination showed 273 patients with functional diseases, whereas 128 patients had various organic diseases (colitis, Crohn's, ulceration, diverticulitis, polyps, adenomas, cancer, or infectious diseases).

ROC curve analysis (AUC: 0.935) resulted in a clinical cut-off at 50 µg/g. Applying this cut-off for the differentiation between organic and functional diseases, a clinical sensitivity and specificity of 84.4% and 94.5%, respectively can be reached. The suggested cut-off level for adults (<50 µg/g) can be used for children aged from 4 to 17 years regardless of sex (6). While faecal Calprotectin levels from adults and children are comparable, levels of newborns can be significantly increased (5).

Therefore, samples **higher than 50 µg/g can be regarded as positive** and patients must be further investigated by colonoscopy or other tests for organic inflammation. It is recommended to re-test samples measured in the grey zone from 40 to 60 µg/g.

PERFORMANCE CHARACTERISTICS

Inter-Assay Precision (Run-to-Run): 15.6 %. The inter-assay precision of the Quantum Blue® Calprotectin assay was calculated from 7 extracted stool samples within the measurable range and divided into aliquots. The aliquots were tested 20 times each, according to the assay procedure. The values are presented in Table 12.

Analytical Sensitivity; Limit of Detection: <3 µg/g. Neat Extraction Buffer was tested 20 times according to the assay procedure. Mean and standard deviation were calculated for the absorbance values (in mV). The minimal detectable dose of Calprotectin was calculated to be clearly below 3 µg/g by adding two standard deviations to the mean absorbance and intersecting this value with the generic standard curve for the corresponding lot of the test cartridge.

Functional Sensitivity; Limit of Quantification: <30 µg/g. Eleven different stool samples with values between 13.6 and 403.6 µg/g Calprotectin were assayed 20 times each according to the assay procedure. The %CV and the mean values were calculated for each sample. The limit of quantification corresponds to the concentration of Calprotectin with a precision below 20%CV. The resulting precision profile (Figure 2) allows the precise measurement within the range from 30 to 300 µg/g.

Dilution Linearity/Parallelism: 110 %. Four stool samples with elevated Calprotectin values were extracted according to the assay procedure. The extracts were diluted with Extraction Buffer and each dilution was subsequently assayed on two different test cartridges according to the assay procedure and averaged for each test point. The expected values were calculated from the observed value found with the first dilution. The results are presented in Table 13.

The Quantum Blue® Calprotectin assay is linear between 30 and 300 µg/g (Figure 3). Values higher than 300 µg/g should be diluted additionally (e.g. 1:10) for quantitative determination.

Spiking Recovery: 91 %. Four stool samples with low Calprotectin levels were spiked with different amounts from diluted human buffy coats containing high concentrations of Calprotectin. Each sample was measured before and after spiking on two different test cartridges according to the assay procedure. The results are presented in Table 14.

Method Comparison: $R^2=0.93$. 26 patients samples within the measurable range were analyzed using the Quantum Blue® Calprotectin assay and compared with the values obtained by the BÜHLMANN Calprotectin ELISA (order code: EK-CAL). The linear regression analysis results in a coefficient of correlation of $R^2 = 0.93$, a slope of 1.02 and an intercept of 9.2 µg/g. The linear regression is presented in Figure 4.

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
ENGLISH

DEUTSCH

Tools provided with the Extraction Device (from left to the right):


- Tube consisting of Body, Stopper Cap and Spiral Coil
- Spatula
- Base Cap with Sample Chamber and Tag

It is recommended to observe standard safety precautions when following the procedure below including wearing gloves.

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
1. Label the empty tube with Patient ID.
2. Carefully press the sample into the hollow cavity (Sample Chamber) in the Base Cap and level off the surface.

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
4. Remove the Stopper Cap at the top and Pipette 4 ml of the Extraction Buffer into the Tube and carefully recap the tube.

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5. Homogenize the sample for at least 1 minute using a vibration mixer (e.g. Vortex mixer) at maximum speed until no large particles can be seen.

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6. Sample preparation for measurement. Dilute extract in a fresh microtube 1:16, e.g. add 20 µl of the extract and 300 µl of Extraction Buffer. Vortex.
7. Centrifuge the diluted extract for 5 min at 3'000 x g*) and continue to the lateral flow assay procedure.

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*) Alternatively the stool can be settled for 10 minutes. The supernatant has to be used for the lateral flow assay.

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
ENGLISH

Make the Reader ready:
(see Appendix I or consult the Reader Manual).

- Press the ENTER Button for 1 second to turn on the Reader.
- Wait for the Self test and the Calibration check which are done automatically.
- Select <SCAN> and press ENTER.
- Select <Test Type> <CAL_720> and press ENTER.
- Select <User ID> and press ENTER.
- Select <SAMPLE ID> and press ENTER.
- Select <Next> and press ENTER.
- Confirm <Lot> ID and <Product> ID by selecting <Yes> and press ENTER.
<Product> ID must be <LF-CAL>.

Change lot-specific test parameters using the RFID Chip Card:
(mandatory if <Lot> ID displayed by the Reader is not identical with <Lot> ID printed on the label of the test cartridge or on the lot-specific QC data sheet delivered with the test kit).

Change <Lot> ID and <Product> ID by selecting <No> and press ENTER.

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

Place the RFID chip card on the Tag position of the RFID Reader. The new lot-specific test parameters will be automatically transferred into the Firmware of the Reader. The reader is now ready to be used with the new Lot of the test kit.

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Make the Calprotectin Cartridge ready:

Open the protective foil of the test cartridge by either tearing or cutting the foil with a pairs of scissors. Label the test cartridge with Patient ID.

Open the tray of the reader and Place the test cartridge on it.

Attention:

- To avoid mechanical wear do not press on the clamp on the backside (hatched arrow) while fixing the test cartridge into the tray.

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
- Make sure that the sample loading port is located on the right side as indicated on the figure by the arrow.

1. Load 60 ul of diluted and centrifuged stool extract onto the sample loading port of the test cartridge with a precision pipette.

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2. Immediately after loading, start the Timer by pressing the ENTER Button of the Reader.

Attention: To avoid contamination of the optical device of the Reader close the tray only after about 30 seconds.

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3. After 12 minutes the test cartridge is automatically read and the test result calculated.^{*)}
4. Press ENTER Button to go back to the MAIN MENU.

Note: Results of the last 100 runs will be stored within the MEMORY MENU and can be recalled (see the Reader Manual).

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Alternatively, the test can be run in REMOTE MODE with additionally provided Software (LF Control). For instruction please refer to the manual of the software.

^{*)} Alternatively, the test cartridge can be measured (scanned) without the fixed test time of 12 minutes. For this, please select <Test Type> <CAL_0> and perform a "manual" timing of the 12-minutes incubation step.

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
SHORT DESCRIPTION OF THE QUANTUM BLUE® READER, THE RFID READER AND THE RFID CHIP CARD


Screen (Display)
Bildschirmanzeige
Ecran
Monitore
Monitor

Loading Insert for test cartridge (Tray)
Ladestation für die Testkassette
Chargement de l'insert pour la cartouche de test (plateau)
Sede di caricamento della cartuccia (vano portacartuccia)
Punto de carga del cartucho de prueba (bandeja)

ENTER Button (central)
Zentraler ENTER-Knopf (in der Mitte)
Bouton ENTER (au milieu)
Pulsante ENTER (centrale)
Botón ENTER (central)

Piloting Buttons (top, bottom, left, right)
Navigierknöpfe (oben, unten, rechts, links)
Boutons de commande (haut, bas, gauche, droite)
Pulsanti di comando (superiore, inferiore, sinistro, destro)
Botones de dirección (arriba, abajo, izquierda, derecha)

Tag Position of RFID Reader (place RFID chip card here)
Etikettenposition des RFID Readers (Position für die RFID Chipkarte)
Emplacement Tag du lecteur RFID (placer la carte à puce RFID à cet endroit)
Posizione tag dell'RFIDlettore (collocare qui la carta chip RFID)
Posición de lectura de etiquetas del dispositivo RFID (sitúe aquí la tarjeta con chip RFID)

RFID Chip card
RFID Chipkarte
Carte à puce RFID
Carta chip RFID
Tarjeta chip RFID

Connecting Plug to the I/O port(s) of the Quantum Blue® Reader
Verbindungskabel zum I/O Port des Quantum Blue® Readers
Raccorder la fiche au(x) port(s) I/O du lecteur Quantum Blue® Reader
Presca di connessione alla(e) porta(e) I/O del lettore Quantum Blue®
Conexión del enchufe al puerto(s) de E/S del lector Quantum Blue®

Figure 1A


Figure 1B


Figure 1C


Figure 1D


Table 11

Clinical Study

	Calprotectin (EK-CAL)	Lactoferrin
n	401	391
cut-off	50 µg/g	3 µg/ml
Sensitivity	84.4%	74.2%
Specificity	94.5%	91.0%
PPV	87.8%	79.3%
NPV	92.8%	88.4%
LR+	15.4	0.17
LR-	8.25	0.28

Table 12

Inter-Assay Precision

Sample	Mean [µg/g]	SD [µg/g]	CV [%]
S2	32.9	5.4	16.4
S3	32.2	6.3	19.5
S4	62.2	9.5	15.4
S5	89.9	13.4	14.9
S6	138.0	18.0	13.0
S7	185.0	31.0	16.5
S8	271.0	37.3	13.8
Mean			15.6

Figure 2:

Precision Profile


Table 13

Dilution Linearity/Parallelism

Sample	Dilution	Observed [µg/g]	Expected [µg/g]	O/E [%]
S 4	1:800	208	--	--
	1:1600	90	104	86
	1:3200	52	52	100
	1:6400	24	26	93
	Mean			93
S 5	1:800	149	--	--
	1:1600	75	75	101
	1:3200	53	37	141
	1:6400	24	19	129
	Mean			124
S 7	1:1600	224	--	--
	1:3200	126	112	112
	1:6400	60	56	107
	1:12800	30	28	107
	Mean			109
S 8	1:3200	254	--	--
	1:6400	150	127	118
	1:12800	71	63	112
	1:25600	28	32	88
	Mean			106
Mean				108

Figure 3:

Linearity Plot


Table 14

Spiking Recovery

Sample	µg/g	spiked with [µg/g]	Observed [µg/g]	Expected [µg/g]	O/E [%]
P1	52.5	25	62	78	79
		50	154	103	150
		100	143	152	94
		200	239	252	95
		265	272	318	86
	Mean				101
P2	81	25	107	106	101
		50	120	131	91
		100	162	181	89
		200	265	281	94
		265	357	346	103
	Mean				96
P3	56	25	77	81	94
		50	120	106	113
		100	110	156	71
		200	186	256	73
		265	272	321	85
	Mean				87
P4	160	25	153	185	82
		50	191	210	91
		100	190	260	73
		200	305	360	85
		265	337	425	79
	Mean				82
Mean					91

Figure 4:

Method Comparison


Table description: cf. "Standardization and Interpretation of Result", "Expected Values" and "Performance Characteristics" (pg. 3).

REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

1. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease*. Gastroenterol **123**, 450-460 (2002)
2. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality*. Lancet **356**, 1783-4 (2000)
3. Hessian PA and Fisher L: *The heterodimeric complex of MRP-8 (S100A8) and MRP-14 (S100A9). Antibody recognition, epitope definition and the implications for structure*. Eur J Biochem **268**, 353-63 (2001).
4. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease*. Gut **47**, 506-513 (2000).
5. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease*. Inflamm Bowel Dis **12(6)**, 524-34 (2006)
6. Fagerberg UL et al. : *Colorectal inflammation is well predicted by fecal calprotectin in children with gastrointestinal symptoms*. J Pediatr Gastroenterol Nutr **40**, 450-5 (2005)

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso

Symbol	Explanation
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Test Cartridge Test Kassette Cartouche de testi Cartuccia di rilevazione Cartucho de prueba
	Extraction Buffer Extraktions-Puffer Tampon d'extraction Tampone di estrazione Tampón de extracción
	Extraction Device Extraktionsbesteck Dispositif d'extraction Sistema di estrazione Dispositivos de extracción
	RFID Chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip

Printing Date
2010-06-11

Not for sale in the USA, due to current patent restrictions