Quantum Blue[®] Calprotectin High Range

**Quantitative
Lateral Flow Assay**

LF-CHR20 20 tests

Revision date: 2011-04-11

ENGLISH

INTENDED USE

BÜHLMANN Quantum Blue® Calprotectin High Range is an immunoassay designed for the quantitative determination of elevated Calprotectin concentrations in human stool samples (1-5) in combination with the BÜHLMANN Quantum Blue® Reader. For laboratory use only.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of Calprotectin antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) being highly specific for Calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the extracted and diluted stool sample. The Calprotectin/anti-Calprotectin gold conjugate binds to the anti-Calprotectin antibody coated on the test membrane (test line; test band) and the remaining free anti-Calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (control line; control band). The signal intensities of the test line and the control line are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cartridge vacuum-sealed in a foil bag	20 cartridges	B-LFCAL-TC	
Extraction Buffer	1 bottle 120 ml	B-CAL-EX	Ready to use
Chase Buffer	1 bottle 100 ml	B-CHR-CB	Dilution buffer; ready to use
RFID Chip Card	1	B-CHR-RCC	

Table 1

REAGENTS & MATERIAL SUPPLIED SUPPLEMENTARY

Fecal Extraction Devices described below are not delivered with the kit and either of them **has to be ordered** with the kit .

Smart-Prep	50 devices	B-CAL-RD
Schebo® Quick-Prep™	50 tubes	B-CAL-SOFI
	Filled with 1.3 ml Extraction Buffer	

Table 2

Controls **can be ordered** separately.

Controls Low* / High*	2 vials, 0.5 ml	B-CHR-CONSET	Ready to use
------------------------------	--------------------	--------------	--------------

Table 3

* Lot specific concentration of Calprotectin

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8°C until the expiration date printed on the labels.

WARNINGS AND PRECAUTIONS

Human Source Material: None of the reagents provided with this kit is of human origin.

Patient specimens should be handled as potentially biohazardous, taking appropriate precautions.

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex mixer for stool extraction
- Precision pipettes with disposable tips: 10-100 µl and 3 ml
- Centrifuge
- 5 ml polypropylene or polystyrene tubes for dilution of the extracts
- Timer (optional)
- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Gloves and Lab coat
- Soft tissues or blotting paper

SPECIMEN COLLECTION AND STORAGE

Using the Extraction Devices which can be supplied upon request, less than 1 g of native stool is required.

Collect stool samples into clean tubes and store refrigerated at 2-8°C for up to 6 days.

Freezing of the sample may result in slightly increased Calprotectin values due to release from neutrophils present in the sample. Therefore, it is of advantage for longer storage, to keep the extracts at -20°C. The extracts are stable for at least 4 months.

Important: The sample must be collected without any chemical or biological additives in the collection device.

ASSAY PROCEDURE

The Assay Procedure consists of three steps:

1. Extraction of stool samples:

The extraction is described in the instruction for use delivered with the respective extraction devices.

2. Sample processing:

- Dilute the extracts of the stool samples prior to using them for analysis 1:150 with Chase Buffer (eg. 20 µl sample and 3000 µl buffer).
- Vortex.
- Centrifuge the diluted extract for 5 min at 3'000 x g* and continue with the lateral flow assay procedure. Alternatively, the stool extract can be settled for 10 minutes.
- The supernatant has to be used for the lateral flow assay (refer to pages 18ff).

3. Lateral Flow Assay Procedure

The assay procedure is described and illustrated in details below on pages 18ff.

PROCEDURAL NOTES

- All reagents and test samples must be equilibrated at room temperature (18-28°C) before starting the assay. The Quantum Blue® Reader must be switched on and programmed for Calprotectin High Range (CHR_0 and CHR_900, respectively) before starting the assay (see Quantum Blue® Reader Manual).
- In order to receive reliable and quantitative results it is important to homogenize the stool sample entirely in the extraction device. Avoid contamination with stool at the top of the tube.
- Diluted samples (stored at 2-8 °C) should be used within 24 hours and cannot be stored for a longer time period.
- Undiluted extracts can be stored at -20 °C for at least 4 months. It is important to centrifuge the extracts before storage (5 min at 3'000 x g). After centrifugation the supernatant must be transferred into a fresh storage tube.

QUALITY CONTROL

- A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by following this instruction for use accurately.
- Patient samples that are not properly handled may cause inaccurate results.
- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices, ii) expiration dates of reagents and iii) storage and incubation conditions.

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see Figures 1A and 1B). It is used as functional test control only and can not be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 15 minutes of incubation time (Figure 1A), the concentration of Calprotectin present in the stool sample is below the detection limit. If a Test Line (T) is detectable after 15 minutes of incubation time (Figure 1B), the Calprotectin concentration present in the stool sample is calculated by the Quantum Blue® Reader.
- If only the Test Line (T) is detectable after 15 minutes of incubation time (Figure 1C), the test result is invalid and the Calprotectin High Range assay has to be repeated using another test cartridge.
- If neither the Control Line (C) nor the Test Line (T) are detectable after 15 minutes of incubation time (Figure 1D), the test result is invalid and the Calprotectin High Range assay has to be repeated using another test cartridge.
- As the Quantum Blue® Reader allows a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a threshold after 15 minutes of incubation time, the test result is also invalid and the Calprotectin High Range assay has to be repeated using another test cartridge.

STANDARDIZATION AND INTERPRETATION OF RESULTS

- The Lateral Flow Assay is calibrated with the BÜHLMANN Calprotectin ELISA (order code: EK-CAL).
- The BÜHLMANN Quantum Blue® Reader uses a lot-specific standard curve to calculate the Calprotectin concentration. This lot-specific standard curve is generated with the median values ($n \geq 10$ measurements each) from at least 10 Calibration points obtained from different stool samples with known Calprotectin concentrations. The assay range is between 100 and 1800 µg/g.
- For quantitative measurements, unknown samples reading above 1800 µg/g should be further diluted (e.g. 1:5) with Chase Buffer and assayed again according to the procedure. The resulting dilution factor must be multiplied by the measured concentration to obtain the final results.
- For quantitative measurements, unknown samples reading below 100 µg/g should be re-tested in the BÜHLMANN Quantum Blue® Calprotectin assay (order code: LF-CAL20).

LIMITATIONS

- The reagents supplied with this kit are optimised to measure human Calprotectin in extracted stool samples.
- Fecal Calprotectin values should be used as supplementary data available to the physician in establishing a diagnosis.
- A Calprotectin test result obtained with the BÜHLMANN Quantum Blue® High Range Calprotectin rapid assay should be confirmed by measuring other biomarkers (such as C-Reactive Protein) or by repeating the Calprotectin test (with ELISA, order code EK-CAL, or using a second test cartridge) and comparing to the CDAI score and the clinical status of the patient, before any clinical and/or therapeutic decision is made.
- During therapy follow-up, a Calprotectin test result obtained with the BÜHLMANN Quantum Blue® Calprotectin High Range rapid assay should be interpreted by comparing the result with previous ones and in the context of the patient's clinical picture.

PERFORMANCE CHARACTERISTICS

Limit of Blank (LoB): <60 µg/g calprotectin. The LoB has been established in two independent runs using two different lots of Test Cartridges with a total of 60 blank values each using Chase Buffer as sample. The LoB has been calculated as the average of the two runs in accordance with CLSI protocol EP17-A.

Limit of Detection (LoD): 95 µg/g calprotectin. The LoD has been established with two stool extracts with concentrations of 62 and 128 µg/g. The samples were measured with three different lots of Test Cartridges in 2 runs per day over 10-20 days. The SD values were determined and therefrom the LoD has been calculated in accordance with CLSI protocol EP17-A. The data presented are the average of three different lots of Test Cartridges.

Limit of Quantification (LoQ):

Lower LoQ: <100 µg/g calprotectin.

Upper LoQ: >1800 µg/g calprotectin.

The LoQ has been established with 11 stool extracts at concentrations between 24 and 2682 µg/g. Samples were measured with three different lots of Test Cartridges in 2 runs per day over 10-20 days. The Total Precision CV was determined, averaged for each sample and limit of 30% CV was applied.

Linearity: Serial dilution of four independent samples showed linearity between measured value and dilution factor within the indicated measuring range of 100-1800 µg/g ($R^2=0.985-0.996$) As an example, the medians of 4 replicates for each dilution of sample 1 are shown in Figure 2. The average difference between measured and expected concentrations was **3%** (recovery 103 %).

Recovery: 84-125%. Four stool extracts containing low calprotectin concentrations were spiked with different amounts from a highly concentrated calprotectin stock solution. Each sample was measured before and after spiking on three test cartridges according to the assay procedure. The three replicates were averaged and the mean values are presented in Table 16.

Repeatability: 21% CV. The repeatability of the Quantum Blue® Calprotectin High Range assay was calculated from

4 extracted stool samples containing between 128 and 917 µg/g calprotectin. Each sample was tested according to the assay procedure in one run of 20 replicates set up within 10 minutes. The mean values of 4 different lots of Test Cartridges are presented in.

Between-day Precision: <10% CV. The between-day precision of the Quantum Blue® Calprotectin High Range assay was calculated according to CLSI protocol EP5-A2 from six stool extracts between 128 and 1324 µg/g calprotectin tested as two replicates per day over 10-20 days. The mean values of three different lots of Test Cartridges are presented in Table 17.

Total Precision: 25% CV. The total precision of the Quantum Blue® Calprotectin High Range assay was calculated according to CLSI protocol EP5-A2 from six stool extracts between 128 and 1324 µg/g calprotectin tested in two runs per day (1 in the morning, 1 in the afternoon) over 10-20 days. The mean values of three different lots of Test Cartridges are presented in Table 17.

Inter-lot Precision: 21.1% CV. The inter-lot precision of the Quantum Blue® Calprotectin High Range assay was calculated from 4 extracted stool samples containing between 128 and 917 µg/g calprotectin. Each sample was tested according to the assay procedure in three runs of 20 replicates each set up within 10 minutes using three different lots of Test Cartridges (n=60 replicates each). The inter-lot data of the 4 samples varied from 17.0 to 22.9% CV.

Method Comparison: $R^2 = 0.86-0.93$. 54 stool extracts from patients with suspected IBD exhibiting calprotectin concentrations within the indicated linear range of the Quantum Blue® Calprotectin High Range assay were analyzed according to the assay procedure using two different lots of Test Cartridges (Lot A and Lot B) in two runs each and compared with the values obtained by the BÜHLMANN Calprotectin ELISA, extended range (order code: EK-CAL). The correlations are presented in Table 18 and one example is illustrated in Figure 3.

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Get the Reader ready:



(see Appendix I or consult the Reader Manual).

- Press the ENTER Button for 1 second to turn on the Reader.
- Wait for the Self test and the Calibration check which are done automatically.
- Select <SCAN> and press ENTER.
- Select <Test Type> <CHR_900>* and press ENTER.
- Select <User ID> and press ENTER.
- Select <SAMPLE ID> and press ENTER.
- Select <Next> and press ENTER.
- Confirm <Lot> ID and <Product> ID by selecting <Yes> and press ENTER.
<Product> ID must be <LF-CHR>.


Change lot-specific test parameters using the RFID Chip Card:

This is mandatory, if <Lot> ID displayed by the Reader is not identical with <Lot> ID printed on the label of the test cartridge or on the lot-specific QC data sheet delivered with the test kit.

Change <Lot> ID and <Product> ID by selecting <No> and press ENTER.

Place the RFID chip card on the Tag position of the RFID Reader. The new lot-specific test parameters will be automatically transferred into the Firmware of the Reader. The reader is now ready to be used with the new Lot of the test kit.

Make the Calprotectin Cartridge ready:

Open the protective foil of the test cartridge by either tearing or cutting the foil with a pairs of scissors. Label the test cartridge with Patient ID.

Open the tray of the reader and place the test cartridge on it.

Attention:

- In order to avoid mechanical wear, do not press on the clamp on the backside (hatched arrow) while fixing the test cartridge into the tray.
- Make sure that the sample loading port is

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
located on the right side as indicated on the figure by the arrow.

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1. Load 80 µl of diluted and centrifuged stool extract onto the sample loading port of the test cartridge with a precision pipette.

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2. Immediately after loading, start the Timer by pressing the ENTER Button of the Reader.

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Attention: To avoid contamination of the optical device of the Reader close the tray only after about 30 seconds.

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3. After 15 minutes the test cartridge is automatically read and the test result calculated.^{*)}
4. Press ENTER Button to go back to the MAIN MENU.

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Note: Results of the last 100 runs will be stored within the MEMORY MENU and can be recalled (refer to the Reader Manual).

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Alternatively, the test can be run in REMOTE MODE with additionally provided Software (LF Control). For instruction please refer to the manual of the software.

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^{*)} Alternatively, the test cartridge can be measured (scanned) without the fixed test time of 15 minutes. For this, please select <Test Type> <CHR_0> and perform a "manual" timing of the 15-minutes incubation step.

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
SHORT DESCRIPTION OF THE QUANTUM BLUE® READER, THE RFID READER AND THE RFID CHIP CARD


Screen (Display)
Bildschirmanzeige
Ecran
Monitore
Monitor

Loading Insert for test cartridge (Tray)
Ladestation für die Testkassette
Chargement de l'insert pour la cartouche de test (plateau)
Sede di caricamento della cartuccia (vano portacartuccia)
Punto de carga del cartucho de prueba (bandeja)

ENTER Button (central)
Zentraler ENTER-Knopf (in der Mitte)
Bouton ENTER (au milieu)
Pulsante ENTER (centrale)
Botón ENTER (central)

Piloting Buttons (top, bottom, left, right)
Navigierknöpfe (oben, unten, rechts, links)
Boutons de commande (haut, bas, gauche, droite)
Pulsanti di comando (superiore, inferiore, sinistro, destro)
Botones de dirección (arriba, abajo, izquierda, derecha)

Tag Position of RFID Reader (place RFID chip card here)
Etikettenposition des RFID Readers (Position für die RFID Chipkarte)
Emplacement Tag du lecteur RFID (placer la carte à puce RFID à cet endroit)
Posizione tag dell'RFIDlettore (collocare qui la carta chip RFID)
Posición de lectura de etiquetas del dispositivo RFID (sitúe aquí la tarjeta con chip RFID)

RFID Chip card
RFID Chipkarte
Carte à puce RFID
Carta chip RFID
Tarjeta chip RFID

Connecting Plug to the I/O port(s) of the Quantum Blue® Reader
Verbindungskabel zum I/O Port des Quantum Blue® Readers
Raccorder la fiche au(x) port(s) I/O du lecteur Quantum Blue® Reader
Presca di connessione alla(e) porta(e) I/O del lettore Quantum Blue®
Conexión del enchufe al puerto(s) de E/S del lector Quantum Blue®

TABLES AND FIGURES/ TABELLEN UND FIGUREN/ TABLES ET FIGURES/
TABELLE E FIGURE/ TABLAS E FIGURAS

Figure 1


Table 17

Precision

Sample	µg/g	Repeatability [CV]	Between Day [CV]	Total [CV]
S5	128	21.4%	5.4%	26.3%
S6	201	22.9%	6.5%	23.8%
S7	340	nd*	0.0%	26.0%
S8	622	16.5%	2.8%	24.8%
S9	917	24.1%	9.0%	22.4%
S10	1324	nd*	9.4%	26.9%

*nd = not determined

Figure 2

Linearity


Table 18

Correlations

Comparisons	n	R ²	Bias	Slope
Lot A, run 1 vs. EK-CAL	53	0.86	51	0.86
Lot A, run 2 vs. EK-CAL	52	0.87	76	0.86
Lot B, run 1 vs. EK-CAL	53	0.91	66	0.95
Lot B, run 2 vs. EK-CAL	51	0.86	94	0.97
Mean Lots A+B vs. EK-CAL**	54	0.91	70	0.91
Mean Lot A vs. mean Lot B	54	0.93	38	1.09

**Graph shown in Figure 3


Table 16

Recovery

Sample	µg/g	spiked with [µg/g]	Observed [µg/g]	Expected [µg/g]	O/E [%]
S1	122	100	257	222	116
		200	376	322	117
		500	845	622	136
		1000	1750	1122	156
		1500	1656	1622	102
	Mean				125
S2	161	100	320	261	123
		200	412	361	114
		500	931	661	141
		1000	1316	1161	113
	Mean				123
S3	163	100	313	263	119
		200	266	363	73
		500	629	663	95
		1000	1018	1163	88
		1500	1406	1663	85
	Mean				92
S4	184	100	270	284	95
		200	300	384	78
		500	671	684	98
		1000	841	1184	71
		1500	1329	1684	69
	Mean				84
Mean					105

Figure 3

Method Comparison

REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

1. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease*. Gastroenterol **123**, 450-460 (2002)
2. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality*. Lancet **356**, 1783-4 (2000)
3. Hessian PA and Fisher L: *The heterodimeric complex of MRP-8 (S100A8) and MRP-14 (S100A9). Antibody recognition, epitope definition and the implications for structure*. Eur J Biochem **268**, 353-63 (2001).
4. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease*. Gut **47**, 506-513 (2000).
5. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease*. Inflamm Bowel Dis **12(6)**, 524-34 (2006)

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso

Symbol	Explanation
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
TC	Test Cartridge Testkassette Cartouche de testi Cartuccia di rilevazione Cartucho de prueba
BUF EX	Extraction Buffer Extraktionspuffer Tampon d'extraction Tampone di estrazione Tampón de extracción
BUF CB	Chase Buffer Laufpuffer Tampon de chase Tampone di chase Tampón de chase
RCC	RFID Chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip

Printing Date
2011-05-23

Not for sale in the USA, due to current patent restrictions