ACE direct

Angiotensin Converting Enzyme REA

RK-ACD 100 tests

Revision date: 2007-09-27

ENGLISH

INTENDED USE

The BÜHLMANN ACE radioenzymatic assay (REA) is intended for the direct and quantitative determination of angiotensin converting enzyme (ACE) activity in serum (1-3).

PRINCIPLE OF THE ASSAY

Angiotensin converting enzyme (ACE) catalyzes the conversion of angiotensin I to angiotensin II. The enzyme also mediates the cleavage of a synthetic substrate, ³H-hippuryl-glycyl-glycine, into ³H-hippuric acid and the glycyl-glycine dipeptide (4). After acidification, the tritiated hippuric acid is separated from unreacted substrate by extraction with scintillation cocktail and measured in a beta counter.

This test has been standardized with the BÜHLMANN ACE Colorimetric kit (order code: KK-ACE) according to the method described in references 5 and 6.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Tracer ³ H-hippuryl-glycyl-glycine in HEPES buffer (pH=8.0)	1 vial 11 ml	B-ACD-TR	Ready to use
Calibrator¹⁾ lyophilized human serum	1 vial	B-ACD-CA	Reconstitute with 1 ml of deionized water and swirl gently to dissolve
Zero Calibrator lyophilized human serum free of any ACE activity	1 vial	B-ACD-0CA	Reconstitute with 1 ml of deionized water and swirl gently to dissolve
Control Normal /High¹⁾ lyophilized human serum	2 vials	B-ACD- CONSET	Reconstitute with 1 ml of deionized water and mix well before use
HCl Solution 1N hydrochlorid acid	1 vial 6 ml	B-ACD-HA	Ready to use
Scintillation Powder²⁾ Premixed	1 vial	B-ACD-SCP	Dissolve with 120 ml xylene and 40 ml of ethyl acetate, mix well

Table 1

¹⁾ Lot specific ACE activity in human serum matrix. Refer to the additional QC data sheet for exact concentrations.

²⁾ The following **non-xylene** scintillation fluid is an acceptable substitute for the scintillation powder: Ultima Gold F (PerkinElmer no. 6013179). Mix 120 ml of Ultima Gold F with 40 ml of ethyl acetate to obtain the scintillation reagent. Store in brown glass bottle for up to 6 months at 18-28°C.

STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Store at 2-8°C until expiration date printed on the label HCl Solution and Scintillation Powder can be stored at 18-28°C.	
Opened / Reconstituted Reagents	
Tracer	Stable at 2-8°C until expiration date marked on the label
Calibrator	Stable for 2 months at -20°C
Zero Calibrator	
Controls	
HCl Solution	Stable at 18-28°C until expiration date printed on the label
Scintillation Powder	Stable for 6 months at 18-28°C if stored in a brown glass bottle

Table 2

WARNINGS AND PRECAUTIONS

This kit contains radioactive material which does not exceed 188 kBq (5.0 µCi) of tritium (³H). The receipt, acquisition, possession, use and transfer are subject to the local regulations. Concerning the proper precautions for the handling and disposal of kit reagents, radioactive material, radioactive waste and patient specimens, we highly recommend to first consult the special local regulations of your country.

Reagents Containing Human Source Material: The Calibrator (B-ACD-CA) the Zero Calibrator (B-ACD-CA0)

and the Controls (B-ACD-CONSET) of this kit contain components of human origin. Each serum used in the preparation of the kit components was tested by a FDA-approved method and found negative for HBV surface antigen, HCV and HIV1/2 antibodies. Although those methods are highly accurate, there is no guarantee that this material cannot transmit Hepatitis or AIDS. *Therefore, all patient specimens and kit components should be handled as if capable of transmitting infections.* All products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions.

MATERIALS REQUIRED BUT NOT PROVIDED

- 10 µl, 50 µl, 100 µl and 1000 µl precision pipettes with disposable tips.
- Disposable 2 ml polypropylene tubes with tight screw-caps for the assay (preferably Sarstedt Minutubes no. 72.693).
- Cylinders and brown glass bottle for the preparation of the scintillation cocktail.
- Double distilled or deionized water.
- Xylene p.a. (e.g. Merck no. 8681) or, alternatively, Ultima Gold F scintillation fluid (PerkinElmer no. 6013179).
- Ethyl acetate p.a. (e.g. Merck no. 9623)
- 10 ml scintillation vials.
- Vortex mixer.
- Water bath set at 37°C.
- Beta counter.

SPECIMEN COLLECTION AND STORAGE

Serum: Since **EDTA inhibits** the ACE activity, only serum specimen must be used for the determination of ACE activity: Collect sufficient blood (at least 1 ml) by venipuncture into a plastic syringe or a venipuncture tube without any anticoagulant. Let stand the blood for 2-3 hours at 18-28°C, then centrifuge at 2-8°C and collect the serum. Freeze the serum specimen at -20°C or lower if not assayed within 5 days. ACE activity in sterile serum is stable for up to 30 days at 2-8°C and 6 months at -20°C or lower.

Lipemic sera do not interfere in the ACE REA assay.

Note: Many authors (e.g.7,8) have reported that ACE activity in patients under captopril medication increases as a function of the duration of storage. To prevent any loss of ACE inhibition, sera of patients under captopril medication should be assayed within 4 hours after blood collection. For storage up to 12 days, these specimen should be kept at -70°C (dry ice). This phenomenon does not occur in sera of patients treated with other ACE inhibitors such as enalapril (8).

Other specimen: Extraction procedures for urine and tissue specimen are described in references 9 to 11.

ASSAY PROCEDURE

1. Prepare 6 minitubes in duplicate for the blank (zero calibrator), the calibrator and the control serum. Prepare additional tubes in duplicate for each patient specimen.

Note: Do label neither tubes, nor caps. Any mark would reduce the counting efficiency!

2. Pipet 10 µl of the "zero" calibrator, calibrator, controls and patient samples, respectively, into the corresponding tubes.
3. Add 100 µl of tritiated substrate (Tracer) to each minitube. Vortex all tubes for at least 5 seconds.
4. Incubate all tubes for 60 minutes ± 5 minutes in a water bath at 37°C.
5. Add 50 µl of 1 N hydrochloric acid solution to each tube.

Vortex.

- Pipet 1.5 ml of scintillation cocktail into the tubes. Screw the caps strongly and vortex for 5 seconds.

Note: Let the minitubes stand for at least 1 hour before counting.

- Place each tube into a scintillation vial and count for 5 minutes in a beta counter.

RESULTS

Calculation: Record the cpm for all tubes (blank, calibrator, controls, samples) and calculate the mean cpm for each pair of tubes. Subtract the mean assay blank from the respective mean of each pair of tubes:

net counts = average cpm - average blank cpm

Calculate the corresponding ACE activity of the sample or control by dividing the net counts of the sample or control by the net calibrator counts and multiply by the ACE activity of the calibrator, A_s , (indicated on the calibrator vial):

$$\text{ACE activity} = \frac{\text{net counts for sample}}{\text{net counts for calibrator}} \times A_s$$

Definition: One unit of ACE activity is defined as the amount of enzyme required to release one μmole of hippuric acid per minute and per liter of serum at 37°C as determined by the colorimetric assay:

$$1 \text{ ACE unit} = \frac{1 \mu\text{mol hippuric acid}}{\text{min} \times \text{L}} = 1 \text{ U/l}$$

Standardization: The values obtained with this assay are identical to the values received with the standard colorimetric assay kit from BÜHLMANN LABORATORIES AG (order code: KK-ACE).

This standardization has the advantage that results obtained with either the colorimetric method (KK-ACE), the kinetic method (KK-ACK) or the radioenzymatic method (RK-ACD) are intercomparable.

QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the QC Data Sheet added to the kit. The use of at least one additional internal pool serum is recommended for internal quality control. There are also some commercial control sera available (e.g. Sigma A-6040 and A-7040, respectively). Due to methodological differences in analyzing the control sera, the values reported on the corresponding data sheets must not necessarily correspond to the data being received with the BÜHLMANN ACE REA assay.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents iii) storage and incubation conditions iv) purity of water.

PERFORMANCE CHARACTERISTICS

Intra-Assay Precision: 2.2%. The intra-assay precision was calculated from the results of 20 pairs of values from four different serum samples in a single run. The values are presented in Table 12.

Inter-Assay Precision: 4.3%. The inter-assay precision was calculated from four different serum samples and two controls tested in 20 different runs. The values are presented in Table 13.

Dilution Linearity: 115%. Three human serum samples have been diluted with Zero Calibrator (ACE free serum) and subsequently assayed according to the assay procedure. The values are presented in Table 14.

Analytical Sensitivity: 0.4 U/l. Twenty Zero Calibrators in duplicates have been determined according to the assay procedure. The mean and the standard deviation (SD) were calculated from the cpm values obtained. To determine the minimal detectable ACE activity the addition of 2 SD to the mean results in a cpm value which than has been calculated as described in "Results"

Functional Sensitivity: 2.0 U/l. Five different serum samples with mean values between 1.8 and 16.3 U/l were each tested 10 times in duplicates as an intra-assay. The coefficient of variation (CV) and the mean values were calculated for each sample. The functional least detectable dose (FLDD) was observed at a CV of 10%

Specificity: A serum sample of defined ACE activity has been inhibited with different dilution of angiotensin I (3.8×10^{-3} - 3.8×10^{-7} mol/l) and EDTA (1.3×10^{-3} - 1.3×10^{-7} mol/l) for two hours at $18-28^\circ\text{C}$ and subsequently measured according to the assay procedure. The results are presented in Table 15.

EXPECTED VALUES AND CUT-OFF

The serum ACE activity strongly depends on the genotype of the patients investigated (12). Therefore, reference values depend on the genetic pattern of the donors and differences can be explained by the frequency of the three genotypes within group of donors on investigation.

From a group of 80 adult (age 20-70 years) normal blood donors from Switzerland, a normal reference range (2.5-97.5th Percentile) has been determined to be (cf. Table 16):

26-60 U/l

In a study with 50 sarcoidosis patients, serum ACE levels were within a range of 45-135 U/l (2).

Serum ACE levels in individuals under 18 years of age are substantially higher and more variable than in adults. Those should be established separately (3).

METHOD COMPARISON

80 sera from apparently healthy donors (see above) were analyzed with the three different assays from BÜHLMANN AG. The results from ACE direct REA, ACE colorimetric (KK-ACE) and ACE kinetic (KK-ACK) were correlated:

ACE direct = $0.65 \times \text{ACE colorimetric} + 13.7 \text{ U/l}$;
 $r = 0.94$; $r^2 = 0.89$

ACE direct = $0.75 \times \text{ACE kinetic} + 6.9 \text{ U/l}$;
 $r = 0.96$; $r^2 = 0.93$

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APPENDIX I

REFERENCES/ LITERATURREFERENZEN/ RÉFÉRENCES/ RIFERIMENTI/ REFERENCIAS

1. Lieberman J: *Enzymes in sarcoidosis-angiotensin converting enzyme*. Clin Lab Med **9**, 745-755 (1989).
2. Studdy P R and Bird R: *Serum angiotensin converting enzyme in sarcoidosis – its value in present clinical practice*. Ann Clin Biochem **26**, 13-18 (1990).
3. Bénétteau-Burnat B et al: *Serum angiotensin converting enzyme in healthy and sarcoidotic children: comparison with the reference interval for adults*. Clin Chem **36**, 344-346 (1990).
4. Neels H M et al: *Improved micromethod for assay of serum angiotensin converting enzyme*. Clin Chem **28** (6), 1352-1355 (1982).
5. Lieberman J.: *Elevation of serum angiotensin converting enzyme (ACE) level in Sarcoidosis*. Am J Med **59**, 36-72 (1975)
6. Hurst P L and Lovell-Smith C.J.: *Optimized assay for serum angiotensin converting enzyme activity*. Clin Chem **27**, 2048-2052 (1981).
7. Kamoun P P et al: *Measurements of angiotensin converting enzyme in captopril-treated patients*. Clin Chim Acta **118**, 333-338 (1982).
8. Lieberman J and Zakria F: *Effect of captopril and enalapril medication on the serum ACE test for sarcoidosis*. Sarcoidosis **6**, 118-123 (1989).
9. Baggio B et al: *Increased urine angiotensin I converting enzyme activity in patients with upper urinary tract infection*. Clin Chim Acta **109**, 211-218 (1981).
10. Kato I et al: *Increased urinary excretion of angiotensin converting enzyme in patients with renal disease*. J Clin Chem Clin Biochem **20**, 473-476 (1982).
11. Yokoyama M et al: *Angiotensin converting enzyme in human prostate*. Clin Chim Acta **100**, 253-258 (1980).
12. Biller H, Zissel G, Müller-Quernheim J et al.: *Genotype-corrected reference values for serum angiotensin-converting enzyme*. Eur Respir J **28**, 1085-90 (2006).

Table 11 Example of Results

	cpm	ACE activity (U/l)	CV (%)
Blank	348		
Blank	365		
Blank Avg	357		3.4
Calibrator	4078	52.5	
Calibrator	4229	52.5	
Calibrator Avg	4153	52.5	2.6
Control NORMAL	2461	29.1	
Control NORMAL	2463	29.1	
Control NORMAL Avg	2462	29.1	0.1
Control HIGH	6776	88.8	
Control HIGH	6923	90.8	
Control HIGH Avg	6849	89.8	1.5
Sample 1	3729	46.6	
Sample 1	3573	44.5	
Sample 1 Avg	3651	45.6	3.0
Sample 2	1260	12.5	
Sample 2	1254	12.4	
Sample 2 Avg	1257	12.5	0.3

$$\text{Sample 1} = \frac{3651-357}{4153-357} \times 52.5 = 45.6 \text{ ACE Units}$$

Table 12 Intra-Assay Precision (Within-Run)

Sample	n	Mean (U/l)	SD (U/l)	CV (%)
1	20	52.8	1.6	3.0
2	20	37.4	0.7	1.9
3	20	29.1	0.7	2.5
4	20	23.2	0.4	1.5
Mean				2.2

Table 13 Inter-Assay Precision (Run-to-Run)

Sample	n	Mean (U/l)	SD (U/l)	CV (%)
Normal Control	20	25.9	1.1	4.2
High Control	20	58.8	2.4	4.0
5	20	52.6	2.1	4.0
6	20	38.9	2.2	5.6
7	20	28.9	1.2	4.3
8	20	23.5	0.9	3.7
Mean				4.3

Table 14 Dilution Linearity/Parallelism

Sample	Dilution	Observed (U/l)	Expected (U/l)	O/E (%)
9	1:1	54.2	—	—
	1:2	31.5	27.1	116
	1:4	15.8	13.6	117
	1:8	7.7	6.8	114
	1:16	4.2	3.4	124
	1:32	2.2	1.7	130
10	1:1	63.6	—	—
	1:2	35.5	31.8	112
	1:4	17.1	15.9	108
	1:8	8.6	8.0	108
	1:16	4.2	4.0	106
	1:32	2.4	2.0	121
11	1:1	56.8	—	—
	1:2	31.3	28.4	110
	1:4	15.8	14.2	111
	1:8	8.2	7.1	115
	1:16	4.1	3.6	115
	1:32	2.2	1.8	124
Mean				115

Table 15 Specificity

Inhibitor	Initial ACE activity [U/l]	Inhibitor concentration [mol/l]	ACE activity after inhibition [U/l]	Inhibition [%]
Angiotensin I	21.7	3.8×10^{-7}	21.5	1
		3.8×10^{-6}	21.5	1
		3.8×10^{-5}	14.7	32
		3.8×10^{-4}	3.9	82
		3.8×10^{-3}	1.2	94
EDTA		1.3×10^{-7}	21.4	1
		1.3×10^{-6}	21.5	1
		1.3×10^{-5}	21.4	1
		1.3×10^{-4}	11.3	48
		1.3×10^{-3}	2.8	87

Table 16 Normal values

	Adults
n	80
Age (years)	20-70
Mean (U/l)	41.1
SD (U/l)	10.0
Mean ± 2SD	21.1 – 61.0
Median (U/l)	40.1
IQR (U/l)	14.6
2.5-97.5th Percentile (U/l)	25.5 – 60.5

Table description: cf. “Results” and “Performance Characteristics” (page 3)

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RADIOENZYMATIC PROCEDURE							
Conical minitubes in duplicate	Standard, Control, Sample (μ l)	Tracer (μ l)		1N HCl (μ l)		Scintillation Cocktail (μ l)	
Blank	10	100		50		1500	Vortex Place each tube in a scintillation vial
Calibrator	10	100		50		1500	
Control NORMAL	10	100		50		1500	Let the tubes stand for at least 1 hour before counting
Control HIGH	10	100		50	Vortex all tubes	1500	
Patient Serum	10	100	Incubate all tubes for 60 $\pm$ 5 min at 37°C in a waterbath	50		1500	Count each vial for 5 minutes in a beta counter

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Radioactive Material Radioaktives Material matériel radioactif materiale radioattivo material radiactivo

Symbol	Explanation
TR	Tracer Tracer Traceur Elemento tracciante Trazador
CAL	Calibrator Kalibrator Calibrateur Calibratore Calibrador
CAL 0	Zero Calibrator Null Kalibrator Calibrateur Zero Calibratore Zero Calibrador Zero
CONTROL N	Control Normal Kontrolle Normal Contrôle Normal Controllo normale Control normal
CONTROL H	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
SOLN HCL	Hydrochloric Acid Salzsäure Acide chlorhydrique Acido cloridrico Ácido hidroclórico
SCP	Scintillation Powder Scintillationspulver Poudre de scintillation Polvere di scintillazione Polvo de centelleo

Printing Date
2007-10-11