VITAMIN B₆

Pyridoxal-5'-Phosphate ³H-REA

RK-VB6 100 tests

Revision date: 2008-10-29

INTENDED USE

The BÜHLMANN Vitamin B₆ radio-enzymatic assay (REA) is intended for the quantitative *in vitro* diagnostic determination of Pyridoxal 5'-Phosphate (PLP, vitamin B₆) in EDTA-Plasma (1-4).

PRINCIPLE OF THE ASSAY

³H-tyrosine is decarboxylated by the vitamin B₆-dependent enzyme tyrosine apodecarboxylase (Y-apoDC) from *Streptococcus faecalis* to ³H-tyramine. The activity of tyrosine apodecarboxylase is quantitatively dependent on the amount of pyridoxal 5'-phosphate present in the reaction mixture. The ³H-tyramine thus produced is selectively extracted in the scintillation cocktail (the excess of ³H-tyrosine remains in the aqueous phase) and can be measured by liquid scintillation counting.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Incubation Buffer Acetate buffer	1 vial 60 ml	B-VB6-IB	Ready to use
Calibrator Lyoph. pyridoxal 5'-phosphate (PLP) ¹⁾	3 vials	B-VB6-CA	Reconstitute with 1 ml of Incubation Buffer
Control Low/Normal ²⁾ lyophilized human serum	3 vials	B-VB6-CONSET	Reconstitute with 1 ml of Incubation Buffer
Enzyme Buffer	1 vial 6 ml	B-VB6-EB	Ready to use
Enzyme lyoph. Y-apoDC from <i>S. faecalis</i>	1 vial	B-VB6-E	Reconstitute with 5.5 ml of Enzyme Buffer
Tracer ³ H-tyrosine	1 vial 5.5 ml	B-VB6-TR	Ready to use
Borate Buffer	1 vial 21 ml	B-VB6-BB	Ready to use
Scintillation Powder ³⁾	1 vial 0.96 g	B-VB6-SCP	Dissolve in a mixture of 80 ml ethyl acetate and xylene each

Table 1

¹⁾ Lot specific amounts of PLP in rabbit serum. Refer to the additional QC data sheet. **Do not store reconstituted Calibrators**, but reconstitute a fresh Calibrator vial each time an assay is performed. After reconstitution, proceed immediately to the Calibrator dilution.

²⁾ Lot specific amount of PLP in human serum. Refer to the additional QC data sheets. **Do not store reconstituted Low and Normal Control**, but reconstitute a fresh control vial each time an assay is performed.

³⁾ The following **non-Xylene** scintillation fluid is an acceptable substitute for the Scintillation Powder: Ultima Gold F (PerkinElmer # 6013329). Mix 80 ml of Ultima Gold F with 80 ml of ethyl acetate to obtain the scintillation reagent. Store in a **brown** glass bottle at 18-28°C for up to 6 months.

STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Calibrators, Controls, Enzyme and Enzyme Buffer are to be stored at -20°C. Store the Incubation Buffer and the Tracer at 2-8°C. Borate Buffer and Scintillation Powder may be stored at 18-28°C. Do not use past kit expiration date printed on the labels.	
Opened / Reconstituted Reagents	
Incubation Buffer	Store at 2-8°C
Tracer	
Controls	Do not Store reconstituted
Calibrator	
Enzyme	Stable for 3 months at -20°C. Aliquot if repeated use is expected.
Enzyme Buffer	Store at -20°C
Borate Buffer	Store at 18-28°C
Scintillation Powder	Store at 18-28°C in a brown glass bottle.

Table 2

WARNINGS AND PRECAUTIONS

This kit contains **radioactive material** which does not exceed 518 kBq (14 µCi) of Tritium (³H) for the RK-VB6. The receipt, acquisition, possession, use and transfer are subject to the local regulations. Concerning the proper precautions for the handling and disposal of kit reagents, radioactive material, radioactive waste and patient specimens, we highly recommend to first consult the special local regulations of your country.

Reagents Containing Human Source Material: The Controls (B-VB6-CONSET) of this kit contain components of human origin. Each serum used in the preparation of the kit components was tested by a FDA-approved method and found negative for HBV surface antigen, HCV and HIV1/2 antibodies. Although those methods are highly accurate, there is no guarantee that this material cannot transmit Hepatitis or AIDS. *Therefore, all patient specimens and kit components should be handled as if capable of transmitting infections.* All products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions.

MATERIALS REQUIRED BUT NOT PROVIDED

- 25 µl, 50 µl, 500 µl and 1000 µl precision pipettes and multipipette with disposable tips.
- 5 ml volumetric pipette.
- Cylinders and brown glass bottle necessary for preparation of scintillation cocktail.
- Ethyl acetate p.a. (e.g. Merck no. 9623).
- Xylene p.a. (e.g. Merck no. 8681).
- Disposable polypropylene, polystyrene or glass tubes for preparation of calibrator and sample dilutions.
- Disposable 2 ml polypropylene tubes with tight screw-caps for the assay (preferably Sarstedt minitubes no. 72.693).
- Water bath set at 37°C.
- Scintillation vials.
- Beta-counter.

SPECIMEN COLLECTION AND STORAGE

The patient must not have a radioisotope scan prior to the collection of a blood sample.

A minimum volume of 0.5 ml of blood is recommended for duplicate determinations. Draw blood into an EDTA venipuncture tube. Centrifuge for 15 minutes at 1000 x g and 2-8°C immediately after collection. **Avoid exposure to light.** After centrifugation, collect the plasma in suitable tubes and store at -20°C or less if not assayed immediately. PLP in plasma will remain stable for up to 2 years if stored at -20°C or less and protected from light. Avoid repeated freeze-thaw cycles.

ASSAY PROCEDURE

Protect calibrators, controls and diluted plasma samples from exposure to direct light during the entire assay.

Preparation of Calibrators: In order to prepare an entire standard curve, 1:2 serial dilutions of the PLP Calibrator stock solution are required as follows:

1. Label five tubes A through E.
2. Pipet 500 µl of Incubation Buffer into tubes A – D.
3. Pipet 500 µl of PLP stock Calibrator into tubes E and D, vortex tube D.
4. Transfer 500 µl from tube D to tube C, vortex tube C.
5. Serially transfer 500 µl from tube C to tube B and then from tube B to tube A, vortex following each transfer.

Table 3 illustrates a complete standard curve prepared from 1:2 serial dilutions of a stock Calibrator Solution (example

with a lot-specific value of 192 nmol/l) and the resulting values of the individual working calibrators:

Tube	Blank	A	B	C	D	E
Final Dilution	Inc. Buffer	1:16	1:8	1:4	1:2	Stock Calibrator
Conc. [nmol/l]	0	12	24	48	96	192

Table 3

Note: The Calibrator solutions contain an effective concentration of PLP which is 20-fold lower than indicated in the table above for equalizing the 1:20 sample dilution (see below: Dilution of specimens). Therefore, sample and control values can be directly read from the standard curve without recalculation.

Dilution of Specimens: Patient samples must be diluted 20-fold just prior to performing the assay: Pipet 50 µl of EDTA plasma into a disposable Polypropylene, Polystyrene or Glass tube. Add 950 µl of Incubation Buffer and vortex.

Note: Controls must not be diluted.

Radioenzymatic Assay (REA)

1. Prepare eight suitable screw-cap microtubes (Sarstedt minitubes no. 72.693 or equal) in duplicate for blank, calibrators and controls. Prepare additional tubes in duplicate for patient samples.

Note: Do not label any of these tubes or caps. Any marks on the tubes or caps will quench the Scintillation counting efficiency.

- 2a. Pipet 25 µl of Incubation Buffer into the blank tubes.
- 2b. Pipet 25 µl of standard dilutions A to E into the corresponding tubes.
- 2c. Pipet 25 µl of Low and Normal Controls into the corresponding tubes.
- 2d. Pipet 25 µl of each diluted patient sample into the corresponding tubes.
3. Add 25 µl of Incubation Buffer to all tubes.
4. Add 50 µl of Enzyme solution to all tubes. Vortex.
5. Incubate all tubes for 30 ± 1 minutes in a water bath set at 37°C.
6. Add 50 µl of (³H)-tyrosine Tracer solution to all tubes. Cover the tubes and vortex.
7. Incubate all tubes for 60 ± 1 minutes in a water bath set at 37°C.
8. Pipet 200 µl of Borate Buffer into each tube.
9. Immediately add 1.5 ml of scintillation cocktail to all tubes. Screw caps on tightly.
10. Vortex each tube for at least 10 seconds. Incubate all tubes for 2 hours ± 5 min in a water bath set at 37°C to complete the extraction. Again, vortex each tube for at least 10 seconds.
11. Place each tube into an individual scintillation vial. Label the caps of the scintillation vials. Equilibrate scintillation vials for at least 60 minutes at 18-28°C in the dark.
12. Count each vial for 2 minutes in a Beta-counter.

RESULTS

Standard Curve: Record the counts per minute (cpm) for the blank and for the standard tubes. Calculate the mean cpm for each pair of tubes. Subtract the mean assay blank (blank tubes) from the respective mean of each pair of tubes:

$$cpm_{net} = cpm_{average} - cpm_{average\ blank}$$

Prepare a log-log graph paper and plot the net cpm (vertical axis) versus the pyridoxal 5'-phosphate (PLP) concentration in nmol/l (horizontal axis) for each of the calibrators. Draw the best fitting curve or calculate the standard curve using a four parameter algorithm.

Samples and Controls: Record the mean cpm for each pair of tubes. Subtract the mean assay blank (blank tubes) from the respective mean of each pair of tubes:

$$cpm_{net} = cpm_{average} - cpm_{average\ blank}$$

Locate the net cpm value of the Controls and samples on the vertical axis of the standard curve, draw a horizontal line intersecting the standard curve and read the pyridoxal 5'-phosphate (PLP) concentration in nmol/l from the horizontal axis. Since the 1 in 20 sample dilution is already considered by the standard labeling indicated, sample and control values can be read directly from the standard curve without recalculation. Samples which have been further diluted, because their concentrations were above the highest standard, should be multiplied by the additional dilution factor. Samples that fall below the sensitivity of the assay (<2.6 nmol/l) should be reported as "non-detectable".

Standardization: The BÜHLMANN PLP standard consists of a rabbit serum pool which is calibrated in at least 10 independent runs against a weighed-in preparation of PLP. The PLP material is checked for purity and proper concentration by UV scanning spectroscopy immediately before using it for calibration.

See Table 16 and Figure 1 for examples of results and standard curves. *These results and standard curves are for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the QC Data Sheet added to the kit. Because no international pyridoxal-5'-phosphate reference preparations are available the use of at least one additional pool serum is recommended for internal quality control.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents iii) storage and incubation conditions iv) purity of water.

LIMITATIONS AND TROUBLE SHOOTING

- Since PLP is light-sensitive (5), PLP Calibrator and Controls must be protected from light. It is also recommended to perform the vitamin B₆ REA protected from exposure to direct light.
- Do not label neither calibrator, control and sample tubes, respectively, nor corresponding caps. Any marks on the tubes or caps will reduce (quench) the scintillation counting efficiency. Labeling the caps of scintillation vials is allowed.
- Counting time should be sufficient to prevent statistical counting error: e.g., accumulation of 2000 cpm will yield 2.2% counting error, 10000 cpm will yield 1% counting error.
- PLP values should be used as supplementary data available to the physician in developing a diagnosis.

PERFORMANCE CHARACTERISTICS

Intra-Assay Precision (Within-Run): 4.6%. The intra-assay precision was calculated from the results of 15 pairs of values from each sample in a single run (n=3). The values are presented in Table 17 as nmol/l of PLP.

Inter-Assay Precision (Run-to-Run): 9.2%. The inter-assay precision was calculated from the results of 10 pairs of values (n=3) in 10 consecutive runs (over 15 days). The values are presented in Table 18 as nmol/l of PLP.

Specificity: Several compounds were substituted in the vitamin B6 assay and showed the cross-reactivities presented in Table 19.

Spiking Recovery: 99.8%. Two EDTA plasma samples containing a low (sample 1) and a normal (sample 2) endogenous PLP concentration, respectively, were spiked with equal volumes of synthetic PLP solutions and subsequently analyzed according to the assay procedure. The results are presented in Table 20.

Dilution Linearity/Parallelism: 101.4%. Two human EDTA plasma samples containing high concentrations of endogenous PLP were diluted serially with Incubation Buffer and subsequently analyzed according to the assay procedure. The results are presented in Table 21.

Analytical Sensitivity: 2.6 nmol/ml. Ten zero standard replicates were assayed in a single run. The minimum detectable concentration of PLP was calculated to be **2.6 nmol/l** by adding two standard deviations to the mean Blank cpm and intersecting this value with the standard curve obtained in the same run. Since the 1:20 dilution of plasma samples has already been taken into consideration for the calculation the effective sensitivity of the assay is 0.13 nmol/l.

Method Comparison: 100 EDTA plasma samples were split and analyzed by the BÜHLMANN Vitamin B₆ REA and by High-Performance Liquid Chromatography (HPLC). A linear regression analysis of the results is presented in Figure 2. The range of values was 8.0-146.5 nmol/l (mean = 39.8) for the REA and 0-117.5 nmol/l (mean = 29.1) for HPLC, respectively (6).

The REA values are some 30% higher than the HPLC values. This might be caused by the retention of PLP in the protein precipitate during the deproteinization step proceeding HPLC analysis.

EXPECTED VALUES

In a study using EDTA plasma samples from 123 apparently healthy donors supplied by the Swiss Red Cross Center in Basel, the **normal values** for pyridoxal 5'-phosphate (PLP) displayed in Table 22 were found.

These normal ranges should be used as guidelines only. It is recommended that each laboratory establishes its own expected range for its patient population.

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Figure 1: Example of a Standard Curve


Table 17: Intra-Assay Precision (Within-Run)

Sample	Mean Value [nmol/l]	SD [nmol/l]	CV [%]
Low	16.0	0.5	3.0
Normal	47.6	2.2	4.6
Elevated	69.7	4.4	6.3
Mean			4.6

Table 18: Inter-Assay Precision (Run-to-Run)

Sample	Mean Value [nmol/l]	SD [nmol/l]	CV [%]
Low	18.5	1.9	10.1
Normal	49.4	4.8	9.6
Elevated	92.9	7.3	7.9
Mean			9.2

Table 19: Specificity

Compound	Concentration	Cross reactivity
Pyridoxamine-phosphate	1 µmol/l 1 nmol/l	2.2 <0.5
Pyridoxamine	1 µmol/l 1 nmol/l	<0.5 <0.5
Pyridoxine-HCl	1 µmol/l 1 nmol/l	<0.5 <0.5
Pyridoxal-HCl	1 µmol/l 1 nmol/l	<0.5 <0.5

Table 20: Spiking Recovery

Sample	endogen PLP [nmol/l]	PLP added [nmol/l]	observed value (O) [nmol/l]	expected value (E) [nmol/l]	Recovery (O/E) [%]
1	16.0	10	23.3	26.0	90
	16.0	20	31.1	36.0	86
	16.0	40	46.5	56.0	83
	16.0	80	96.4	96.0	100
2	44.0	10	57.4	54.0	106
	44.0	20	65.2	64.0	102
	44.0	40	89.1	84.0	106
	44.0	80	155.0	124.0	125
Mean					99.8

Table 21: Dilution Linearity/Parallelism

Sample	Dilution	observed (O) [nmol/l]	expected (E) [nmol/l]	recovery O/E [%]
3	1:20	101.1	---	---
	1:40	47.9	50.6	95
	1:80	24.8	25.3	98
	1:160	12.3	12.6	98
	1:320	5.5	6.3	87
4	1:20	240.6	---	---
	1:40	143.3	120.3	119
	1:80	66.0	60.2	110
	1:160	29.7	30.1	99
	1:320	15.4	15.0	103
	1:640	7.8	7.5	104
Mean				101.4

Figure 2: Method Comparison


Table 22: **Expected Values**

Gender	Women	Men
n	54	69
Age (yrs, mean $\pm$ SD)	43.9 $\pm$ 13.9	48.4 $\pm$ 12.2
Mean Value [nmol/l]	60.8	57.8
Median Value [nmol/l]	38.5	51.3
Range [10-90 th Percentile]	21.6 - 104.1	29.8 - 102.4
Lowest Value [nmol/l]	11.5	20.9
Highest Value [nmol/l]	296.0	138.7

Table description: cf. "Results" (page 3), "Performance Characteristics" and "Expected Values" (page 4).

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RADIOENZYMATIC ASSAY PROCEDURE									
Minitubes in Duplicate	Standard, Control, Sample (μ l)	Incubation Buffer (μ l)	Enzyme Solution (μ l)		Tracer (μ l)		Borate Buffer (μ l)	Scinti-lation Cocktail (ml)	
Blank	--	50	50	Vortex and incubate at 37°C for 30 minutes ($\pm$ 1 min)	50	Cover tubes, vortex and incubate at 37°C For 60 minutes ($\pm$ 1 min)	200	1.5	Cap tubes tightly, vortex for 10 seconds and incubate at 37°C for 2 hours ($\pm$ 5 min) ↓ Vortex for 10 seconds ↓ Place tubes into scintillation vials, equilibrate for 1 hour in the dark and count tubes for 2 minutes in a beta counter
Standard A	25	25	50		50		200	1.5	
Standard B	25	25	50		50		200	1.5	
Standard C	25	25	50		50		200	1.5	
Standard D	25	25	50		50		200	1.5	
Standard E	25	25	50		50		200	1.5	
Low Control	25	25	50		50		200	1.5	
Normal Control	25	25	50		50		200	1.5	
Sample (1:20)	25	25	50		50		200	1.5	

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Upper limit of temperature Temperaturobergrenze Limite supérieure de température Limite superiore di temperatura Limite superior de temperatura
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Radioactive Material Radioaktives Material matériel radioactif materiale radioattivo material radiactivo

Symbol	Explanation
BUF INC	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone d'incubazione Tampón de incubación
CAL	Calibrator Kalibrator Calibrateur Calibratore Calibrador
CONTROL L	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
CONTROL N	Control Normal Kontrolle Normal Contrôle Normal Controllo normale Control normal
BUF ENZ	Enzyme Buffer Enzym-Puffer Tampon enzyme Tampone enzimatico Tampón enzimático
ENZ	Enzyme Enzym Enzyme Enzima Enzima
TR	Tracer Tracer Traceur Elemento tracciante Trazador
BUF BOR	Borate Buffer Borat Puffer Tampon borate Tampone borato Tampón borato
SCP	Scintillation Powder Szintilations-Pulver Poudre de scintillation Polvere di scintillazione Polvo de centelleo

Printing Date
2008-10-29