

# Human IGF-I ELISA

## (Insulin Like Growth Factor-I)

Cat. No.: RMEE20



### TECHNICAL FEATURES+APPLICATIONS


- High Specificity for IGF-I
- Correct measurement of IGF-I in **non-extracted** samples
- No physical separation of IGF-I from IGF-binding proteins required
- Elimination of interference by IGF-binding proteins through excess IGF-II
- 98,7 % recovery of recombinant IGF-I leads to correct absolute values
- Precise measurement of very low IGF-I levels: **sensitivity** of **0.09 ng/ml**
- Small sample volume requirement, thus ideal for young patients
- Inter- and Intra-assay variance: 6.8 and 6.7%

### INTRODUCTION

Insulin-like growth factors (IGF) I and II play a pivotal role in regulating the proliferation, differentiation and specific functions of many cell types (1-3). IGF-I is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 daltons (5). Its major regulators are growth hormone (GH) and nutrition (6), although its production in specific tissues is affected by a multitude of tropic hormones and other peptide growth factors. In contrast to many other peptide hormones, IGFs are avidly bound to specific binding proteins (IGFBP). The seven classes of IGFBPs which

are known at present (7,8,22) either bind IGF-I and IGF-II with similar affinities or show a preference for IGF-II (9,10).

A major problem of IGF-I measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-I/IGFBP-3 complexes during the assay incubation. Depending on the ratio IGF-I to IGFBP in the sample interference comes up (see example Figure 1):



**Figure 1.** Interference of IGFBP in IGF-I measurements. Known concentrations of IGF-I were assayed in the presence of 0.5 ng (left) or 5 ng (right) hIGFBP-3 by a conventional (□) and by the IGFBP-blocked assay (\*).

Therefore, various techniques were applied to physically separate IGF-I from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2,12,13). These techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries. The most widely used method is the acid-ethanol extraction (13,14) with a recovery of only 70-80 % of IGFBP-bound IGF-I as a result of co-precipitation. The absolute results of such an extraction are therefore false low (15). The extraction removes the IGFBPs only insufficiently and leads to reduction in sensitivity of the assay due to pre-dilution of the samples by the extraction procedure. Furthermore, the remaining IGFBP may still interfere in the assay. In addition, the acid-ethanol extraction is ineffective in specimens other than serum or plasma (e.g. cell culture media), in which determination of IGF-I is already difficult enough due to the fact that IGFBPs are frequently present at large excess.

To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation is not required before measurement.

## Clinical Significance

There are apart from GH, a number of variables that influence serum IGF-I. Decreased levels are found in states of malnutrition/malabsorption, hypothyroidism, liver disease, untreated diabetes mellitus, chronic inflammatory disease (1,6), malignant disease or polytrauma. High levels, on the other hand, are likely to be present in precocious puberty or obesity. Crucially important to the correct interpretation of IGF-I measurements is the relationship between age and IGF-I levels. It is certainly inadequate to use a common cut-off point to define "normal" levels for all age groups, particularly in children and adolescents.

Due to its GH-dependence, determination of serum IGF-I was shown to be a useful tool in diagnosis of growth disorders, especially with regard to GH deficiency (GHD) or acromegaly (6,16-19,23,24). The major advantage of IGF-I determination compared to GH determination is its stable circadian concentration; therefore a single measurement is sufficient. Hence IGF-I determination should be the first in a series of laboratory test. Clearly normal levels would then rule out disturbances of the GH-IGF-I-axis. Low levels, i.e. close to or below the age-related 5th percentile would indicate the necessity of further diagnostic efforts. Subnormal levels of IGF-I would be evidence for reduced GH secretion, if other causes of low serum IGF-I (e.g. malnutrition or impaired liver function) can be ruled out. For differentiation of healthy short children without GH deficiency and children with "classical" GH deficiency, the 0.1st percentile proved to be an appropriate cut-off point, especially after the age of eight. However, IGF-I levels of short children not suffering from GHD may nevertheless lie between the 0.1st and 5th percentile (19). In contrast, acromegaly is characterized by pathologically elevated IGF-I levels, which apparently reflect the severity of the disease better than GH-levels (17,18,20).

## Scientific Use

IGF-I is present in low concentrations in various body fluids and in conditioned cell culture media of many cell lines. However, the determination of IGF-I in these specimens is particularly difficult due to the presence of IGFBPs usually in excessive amounts. This explains why conventional assays, in which IGFBPs are not removed, result in incorrect IGF-I values, which reflect more the present amount of IGFBP rather than the exact concentration of IGF-I (Figure 1.) (15,21). The low IGF-I concentrations require often additional efforts after the extraction procedure to concentrate the extract for obtaining a satisfactory sensitivity. The IGFBP-blocked IGF-I ELISA avoids these problems and allows the simple, correct and sensitive IGF-I determination in numerous samples at the same time.

## INTENDED USE


This ELISA kit is suitable for the scientific and diagnostic measurement of IGF-I in human serum or plasma, other human body fluids or conditioned media of human cell lines. Due to the high cross-reactivity with IGF-I from other mammalian species, it can also be used as a **heterologous assay** for determination of IGF-I in **primates, cattle,**

pig, sheep, horse, donkey, goat, dog, cat, rabbit and guinea pig, however for rat, mouse and chicken derived samples the kit is not suited.

## PRINCIPLE

The ELISA for IGF-I RMEE20 is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. The IGF-I in the sample binds to the immobilized first antibody on the microtiter plate, the biotinylated

and Streptavidin-Peroxidase conjugated second specific anti-IGF-I-Antibody binds in turn to the immobilized IGF-I. In the closing substrate reaction the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-I-level of the samples.



In order to dissociate IGF-I from the IGFBPs, the samples must be diluted in an acidic buffer (Figure 2). The diluted samples are then pipetted into the wells, by this the pH-value will be neutralized. After neutralization of the samples, the excess IGF-II occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of resulting free IGF-I. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized.

Due to the extremely low cross-reactivity of the IGF-I antibody with IGF-II, the excess of IGF-II does not disturb the interaction with IGF-I.

The test runs like a conventional ELISA using a Streptavidin-Peroxidase-Enzyme Conjugate.

**Figure 2:** Principle of the IGFBP blocked IGF-I ELISA

## PERFORMANCE CHARACTERISTICS AND VALIDATION

The standards of the ELISA RMEE20 are **recombinant IGF-I** in concentrations of **1,5,15,30 and 50 ng/ml**, respectively the assay range covers –at recommended normal sample dilution- the **range from 21 to 1050 ng/ml**. By varying the sample dilution this can be adapted to the special individual requirements.

### Sensitivity

The **analytical sensitivity** of the ELISA RMEE20 yields **0.09 ng/ml** (2 SD of zero standard in 17fold determination).

The **Inter-** and **Intra-Assay variation** coefficients are less than 6.8% and 6.7 % respectively. Exemplary determinations are representend in the Table 1 and 2.

**Table 1: Inter-Assay-Variation**

	Mean value (ng/ml)	Standard deviation	VC %
<b>Sample 1</b>	174	11.79	6.79
<b>Sample 2</b>	494	11.11	2.25
<b>Sample 3</b>	142	8.68	6.11

**Table 2: Intra-Assay-Variation (n=18)**

	Mean value (ng/ml)	Standard deviation	VC %
<b>Sample 1</b>	144.8	9.63	6.65
<b>Sample 2</b>	140.79	7.15	5.08
<b>Sample 3</b>	138.02	7.86	5.69

### Linearity

Dilution:	Sample 1 (recalculated, ng/ml)	Dilution:	Sample 2 (recalculated, ng/ml)
1:10	137.2	1:10	439.1
1:20	133.5	1:20	500.2
1:40	133.6	1:40	499.2
1:80	134.6	1:80	490.5
1:160	134.4	1:160	494.5
1:320	135.7	1:320	526.4
		1:640	463.69
<b>AV / SD / VC%</b>	134.8 / 1.4 / 1.04	<b>AV / SD / VC%</b>	487.6/ 28.2 / 5.79

AV = Average Value , SD = Standard Deviation; VC = Coefficient of Variation

### Recovery

In different human sera the **recovery** was on average **98.67%** of the hypothetical expected amount.

The calibration of the IGF-I ELISA with regard to the WHO International Reference Standard preparation of IGF-I, **NIBSC Code 87/518** (25-27), yields a conversion factor of 1.075. The results must be multiplied by 1.075 to express IGF-I measurements according to the present international standard preparation (its effective content by respective unit definition).

The **recovery** of the International Standard 87/518 in ELISA yielded on average 93%, with the values between 88.6-99.3%.

### Validation

The IGF-I ELISA was developed in accordance with the IGF-I RIA. The clinical validation of the radioimmunoassay was achieved by determining the IGF-I levels in a large number of normal children and adults, normal short statured children without GH deficiency, girls with Ullrich-Turner syndrome, children with Silver-Russell syndrome, patients with GH deficiency, children with familial tall stature, Sotos syndrome, patients with acromegaly, and children with precocious puberty

## SPECIMEN COLLECTION, PREPARATION, AND STORAGE

Serum samples as well as Heparin and EDTA-Plasma samples are suited.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although IGF-I levels were found to be unaffected by few cycles in our experiments.

The high sensitivity of the assays allows IGF-I determinations in small sample volumes, which is limited by pipetting accuracy rather than the amount of IGF-I. For most of the determinations (Serum- or Plasma-Samples and no extreme values expected) the **dilutions of 1:10 to 1:50** in **Sample Buffer PP** should be suited. If necessary –depending on the expected IGF-I level- it is possible to use higher or lower dilution in **Sample Buffer PP**.

Generally, a dilution of **1:21 is very well suited** for serum or plasma samples. IGF-I concentration in other body fluids or cell culture supernatants however could differ strongly.

Suggestion for dilution protocol:

Please pipette 200 µl **Sample Buffer PP** in PE/PP-Tubes (application of a multi-stepper is recommended in larger series); subsequently add 10 µl serum or plasma samples (dilution 1:21). After mixing use 20 µl of this dilution per determination.

## REAGENTS PROVIDED

1)	<b>MTP</b>	<b>Microtiter plate</b> , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with anti-human IGF-I Antibody, packed in a laminate bag.
2)	<b>CAL</b>	<b>Standards A-E</b> , lyophilised, contain recombinant human IGF-I. Standard values are <b>between 1- 50 ng/ml</b> (1, 5, 15, 30 and 50 ng/ml) IGF-I and have to be reconstituted in <b>250 µl (each) in Sample Buffer PP</b> . After using store the reconstituted standards in the original flasks as soon as possible at –20°C. When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. 20 µl per well are used in the assay
3)	<b>BUF</b> X	<b>Dilution Buffer VP</b> , 25 ml, ready for use, please use for the dilution of Antibody Conjugate AK and Enzyme Conjugate EK.
4)	<b>BUF</b> X	<b>Sample Buffer PP</b> 25 ml, ready for use, please use for reconstitution of Standards and Control Serum KS and for dilution of samples and Control Serum KS.
5)	<b>Control</b>	<b>Control Serum KS</b> , 100 µl, lyophilised, contains human serum and has to be reconstituted in <b>100 µl Sample Buffer PP</b> . The reconstituted Control Serum KS must be stored in the original flask as soon as possible at –20°C after using. When using anew, please thaw it rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. The IGF-I target value concentration and the respective range is given on the vial label. The dilution of the Control Serum KS should be according to the dilution of the respective samples.
6)	<b>Ab</b>	<b>Antibody Conjugate AK</b> , 120 µl, 101fold Concentrate, contains the biotinylated anti-IGF-I antibody. Dilute before use <b>1:101 in Dilution Buffer VP</b> (e.g. 100 µl AK + 10 ml VP) and use 80 µl for each well in the assay.
7)	<b>CONJ</b>	<b>Enzyme Conjugate EK</b> , 140 µl, 101fold Concentrate, contains horseradish-peroxidase conjugate to streptavidin. Dilute before use <b>1:101 in Dilution Buffer VP</b> (e.g. 120 µl AK + 12 ml VP) and use 100 µl for each well in the assay.
8)	<b>WASHBUF</b> <b>20x</b>	<b>Washing Buffer (WP)</b> , 50 ml, 20 X concentrated solution.
9)	<b>SUBST</b>	<b>Substrate (S)</b> , 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H <sub>2</sub> O <sub>2</sub> Tetramethylbencidine.
10)	<b>H<sub>2</sub>SO<sub>4</sub></b>	<b>Stopping Solution (SL)</b> , 12 ml, ready for use, 0.2 M sulphuric acid, Caution acid!
11)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

## MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes (100 and 200 µl) Micropipettes and multichannel pipettes with disposable plastic tips

Distilled or Deionized water for dilution of the Washing Buffer (WP)

Vortex-mixer

Device to aspirate the standards and the samples from the wells (recommended because of the potential danger of infection by human samples)

Timer (120 min. range)

Reservoirs (disposable)

Plate washer and plate shaker (recommended)

Calibrated Micro plate reader ("ELISA-Reader") with filter for 450 and 620nm (or  $\geq 590$  nm)

Foil welding device for laminate bags (recommended)

## REAGENT PREPARATION

Room temperature incubation means: Incubation at 20 - 25°C.

The microtiter plate and all reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

The Standards **A – E** and **Control Serum KS** are reconstituted with the **Sample Buffer PP** provided in the Kit. It is recommended to keep the reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Use the **Dilution Buffer VP** for the dilution of **Antibody** and **Enzyme Conjugate concentrates AK and EK**.

The shelf life of the components after opening is not affected, if used appropriately. Store the unused seal stripes of the microtiter plate together with the desiccant at 2-8°C. Reconstituted Components (**Standards A – E** and **Control Serum KS**) should be stored at -20°C (or below). Freezing extends the expiry at least 2 months. When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

The 1:20 diluted Washing Buffer **WP** is only 4 weeks stable. Please dilute only according to daily requirements.

Before use, all kit components should be brought to room temperature. **Precipitates, possible in buffers, should be dissolved before use through mixing and warming.**

The **Substrate Solution S**, stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbencidine, is photosensitive – store and incubate in the dark.

When performing the assay, the Standards **A-E**, Control Serum **KS** and the samples should be pipetted as fast as possible (e.g., 15 minutes). To avoid distortions due to differences in incubation times the diluted Enzyme Conjugate **EK** as well as the succeeding **Substrate Solution S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution SL** should be added to the plate in the same order as the Substrate Solution **S**



## STORAGE CONDITIONS

The microtiter plate wells and all undiluted reagents are stable until the expiry date if stored in the dark at 2-8°C.

Store the unused seal strips and microtiter wells together with the desiccant at 2° to 8°C.

The Substrate Solution (S), stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbencidine, is photosensitive – store and incubate in the dark.

Store the unused seal stripes of the microtiter plate together with the desiccant at 2-8°C. Reconstituted components (**Standards A – E** and **Control Serum KS**) should be stored at -20°C (or below). Freezing extends the expiry at least 2 months. When using the Standards or Control Serum KS anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

## WARNINGS AND PRECAUTIONS

**For in-vitro diagnostic use only. For professional use only.**

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought **to room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature WILL affect the absorbance** readings of the assay. However, values for the patient samples will not be affected.

Reagents from Kits with different lot numbers should not be mixed. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of human and/or animal origin. Source of human serum for the Control Serum provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibody. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

**Stop Solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)**

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37      Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps. Some reagents contain  $\leq 0.095\%$  Kathon and/or ProClin 950 as preservatives. In case of contact with eyes or skin, flush immediately with water.

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.**

R20/21/R22	Harmful by inhalation, in contact with skin and if swallowed
R36/37/38	Irritating to eyes, respiratory system and skin
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1	After contact with skin, wash immediately with plenty of water
S36/37	Wear suitable protective clothing and gloves

**General first aid procedures:**

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

## **ASSAY PROCEDURE**

NOTES: All determinations (Standards, Control Serum and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

It is recommended to use the microtiterplate shaker, if this is not available, please use the alternative procedure (see below).

1) Add **80  $\mu$ l 1:101 diluted Antibody Conjugate AK** in all wells used

2) pipette in positions A1/2 **20  $\mu$ l Sample Buffer PP**

3) Pipette in positions B1/2 **20 µl of the Standard A (1 ng/ml)**

Pipette in positions C1/2 **20 µl of the Standard B (5 ng/ml)**,

Pipette in positions D1/2 **20 µl of the Standard C (15 ng/ml)**,

Pipette in positions E1/2 **20 µl of the Standard D (30 ng/ml)**,

Pipette in positions F1/2 **20 µl of the Standard E (50 ng/ml)**.

To control the correct accomplishment of the assay **20 µl** of the 1:21 (or in respective dilution ratio of the samples) in Sample Buffer diluted **Control Serum KS** can be pipetted in positions G1/2 .

Pipette **20 µl** each of the diluted samples (e.g. dilute 1:21 with Sample Buffer **PP**) In the rest of wells, according to your requirements.

4) Cover the wells with sealing tape and incubate the plate for **1 hour** at **room temperature** (shake at  $\geq 350$  rpm)

5) After incubation aspirate the contents of the wells and wash the wells 3 times **250 µl Washing Buffer WP** / well.

6) Following the last washing step pipette **100 µl** of the **1:101 diluted Enzyme Conjugate EK** in each well.

7) Cover the wells with sealing tape and incubate the plate for **30 Minutes** at **room temperature** (if possible shake  $\geq 350$  rpm).

8) After incubation wash the wells 3 times with Washing Buffer **WP** as described in step 5.

9) Pipette **100 µl** of the **Substrate Solution S** in each well.

10) Incubate the microtiter plate for **15 minutes in the dark** at **room temperature**.

11) Stop the reaction by adding **100 µl Stopping Solution SL** to all wells.

12) Measure the absorbance within **30 minutes at 450 nm (Reference filter  $\geq 590$  nm; e.g. 620 nm)**.

## ALTERNATIVE PROCEDURE

### Without Shaking

If there is no microtiter plate shaker available, the sample might not be mixed optimally within the well. In this case, it is recommended to dilute the samples **1:21 in Sample Buffer PP**. Then, further dilute these prediluted samples **outside of the wells** 1:5 in the 1:101 prediluted Antibody Conjugate **AK**. Subsequently add 100 µl of this external diluted sample per well. The following incubation without shaking should then take place at room temperature (20-25°C) as described in the assay procedure, however the incubation time has to be prolonged to 2 hours (instead of 1 hour). Please perform all other steps of the assay procedure according to the assay protocol given above.

## CALCULATION OF RESULTS

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.20, and the absorbance of standard E should be greater than 1.00.

Samples, which yield higher absorbance values than **Standard E**, are beyond the standard curve, for reliable determinations such samples should be retested at a higher dilution.

### Establishing the Standard Curve

The standards provided contain the following concentration of recombinant hIGF-I:

Standard	A	B	C	D	E
ng/ml	1	5	15	30	50

- 1) Calculate the **mean absorbance** (MA) value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance (MA) of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis on semi-log paper (lin-log).
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **Non-linear regression**, a higher-grade polynomial or four parametric logistic (4-PL) lin-log curve fit are suitable for the evaluation. A good fit is provided with cubic spline, 4 Parameter Logisitics or Logit-Log.
- 5) The concentration in ng/ml of the samples can be calculated by multiplication with the respective dilution factor.

## EXPECTED NORMAL VALUES

IGF-I levels are highly age-dependent in children, less so in adults until the age of about 60. The normal ranges in various age groups, which are log-normally distributed, are given in Table 5 by percentiles. Between 8 and 19 years of age, values are given for boys and girls separately, because the pubertal peak usually occurs approximately 2 years earlier in girls. A graphic presentation is shown in Figures 3, 4 and 5. A major problem for the interpretation of

IGF-I values arises from the fact that short stature is often due to developmental delay rather than any metabolic or endocrine disorder (constitutional delay of growth and adolescence). The sharp rise in IGF-I levels during puberty may therefore cause some uncertainty as to whether or not it would be appropriate to relate measured values to chronological age. It is recommended to take the pubertal stage into account (Table 4 and Figure 6) to get a more complete picture of this situation.

## LIMITATIONS OF PROCEDURE

IGF-I levels depend to a great degree on GH secretion. Diminished IGF-I values, however, do not prove GH deficiency, because a number of other factors can influence the plasma concentration of IGF-I and must therefore be taken into account in order to make a correct interpretation. IGF-I levels decrease during fasting (more than 1 day), as a result of malnutrition, malabsorption, cachexia, impaired hepatic function, or in hypothyroidism and untreated diabetes mellitus. They may also be low in chronic inflammatory disease and malignancies. IGF-I levels are high in states of accelerated sexual development. In clinical situations with hyperprolactinemia or in patients with craniopharyngioma, normal levels may be observed despite GH deficiency. In late pregnancy, IGF-I levels are moderately elevated.


## Appendix

**Table 4:** Normal range of serum IGF-I levels (ng/ml) at different pubertal stages according to Tanner. Because no significant difference between boys and girls is observed, both sexes are combined. Only children and adolescents between 7 and 17 years of age are included.


Pubertal Stage	Percentiles			
	0,1.	5.	50.	95.
1	61	105	186	330
2	105	156	298	568
3	113	196	352	631
4	171	268	431	693
5	165	263	431	706

**Table 5:** Serum levels of IGF-I in healthy subjects at various ages. Individuals between 8 and 19 years of age were classified according to gender, as the pubertal peak occurs almost 2 years earlier in girls than in boys


		Percentiles													
Age		0,1	1	5	10	20	30	40	50	60	70	80	90	95	99
0-2 y.		13	20	28	34	43	50	58	66	75	87	102	128	156	220
2-4 y.		20	29	40	48	59	68	77	87	98	111	129	159	189	260
4-6 y.		26	36	50	59	73	85	96	108	122	138	160	196	233	320
6-7 y.		34	46	62	72	87	99	111	124	138	155	176	212	248	332
7-8 y.		45	60	78	90	107	121	134	148	163	181	205	243	281	364
8-9 y.	boys	54	71	90	102	119	133	146	160	175	192	214	250	284	362
	girls	55	75	99	115	137	156	174	193	214	239	271	324	376	496
9-10 y.	boys	63	82	102	115	133	148	162	176	191	209	232	269	304	379
	girls	68	89	114	130	152	170	187	205	224	247	276	323	369	469
10-11 y.	boys	77	96	117	130	148	162	176	189	203	220	241	274	305	370
	girls	81	106	134	153	178	199	219	239	261	287	321	374	426	539
11-12 y.	boys	85	106	129	144	163	179	194	209	225	244	267	304	339	413
	girls	91	123	160	185	220	248	276	305	337	374	424	503	581	758
12-13 y.	boys	88	112	141	159	184	204	223	243	264	289	321	371	419	525
	girls	116	155	201	231	274	309	342	377	415	460	519	614	707	914
13-14 y.	boys	111	143	179	203	235	261	286	311	339	371	412	477	540	677
	girls	163	207	256	287	329	364	395	428	463	504	556	637	716	884
14-15 y.	boys	140	182	229	260	303	337	370	404	441	484	539	625	691	896
	girls	193	236	284	314	353	385	414	443	474	510	556	628	713	832
15-16 y.	boys	176	221	269	299	340	372	402	433	466	504	552	626	697	849
	girls	187	231	279	309	350	382	412	442	474	512	559	632	700	845
16-17 y.	boys	178	221	267	296	335	366	395	424	455	491	537	607	673	814
	girls	183	225	270	298	336	366	394	422	452	486	530	597	660	792
17-18 y.	boys	173	207	243	265	294	317	337	358	380	405	436	484	527	618
	girls	176	210	246	268	297	320	341	362	384	409	441	488	533	624
18-19 y.	boys	167	201	235	256	285	307	327	347	368	393	423	469	512	600
	girls	167	199	233	254	281	302	322	341	362	385	414	458	499	583
19-20 y.		158	189	220	240	265	285	304	322	341	363	391	433	471	550
20-30 y.		72	92	115	130	150	167	182	198	215	235	261	302	340	425
30-40 y.		68	87	109	123	142	158	173	188	204	223	248	287	324	404
40-50 y.		64	82	103	116	135	150	164	178	194	212	235	272	310	385
50-60 y.		60	77	97	110	127	142	155	169	184	201	224	260	292	369
60-70 y.		55	72	91	103	120	134	147	161	176	193	215	251	282	362
70-80 y.		25	35	47	55	67	78	88	98	110	124	142	173	207	276
>80 y.		21	30	40	47	58	67	76	85	95	108	125	153	184	245




**Fig. 3:** Age-dependent normal range of serum IGF-I levels in girls



**Fig. 4:** Age-dependent normal range of serum IGF-I levels in boys



**Fig. 5:** Age-dependant normal range of serum IGF-I levels in adults



## LITERATUR

1. Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem.*25:49-115
2. Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
3. Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors*. New York: Elsevier.

4. Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology*. 112:2215-2217.
5. Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem*. 253:2769-2276.
6. Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab*. 13:113-143.
7. Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh)*. 121:751-752.
8. Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab*. 74:1215-1216.
9. Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem*. 261:8754-8760.
10. Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J*. 8:2497-2502.
11. Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest*. 60:648-657.
12. Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med*. 109:355-363.
13. Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol*. 128:347-357.
14. Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab*. 51:781-788.
15. Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
16. Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr*. 109:428-433.
17. Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N Engl J Med*. 301:1138-1142
18. Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest*. 68:1321-1330.
19. Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence*. Mannheim: J + J Verlag; 102-117.
20. Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab*. 55:147-153.



21. Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh)*.118:374-380.
22. Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. *J Clin Endocrinol Metab* Vol 82, 4:1301-1303
23. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Horm Res* 54:60-68
24. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Horm Res* 55:155-124
25. Address NIBSC:Blance Lane, South Mimms, Potters Bar, Hertford EN 6 3 QG, Great Britain
26. Bristow AF, Gooding RP, Gaines-Das RE 1990 The international reference reagent for insulin-like growth factor-I. *Journal of Endocrinology*. 125:191-197.
27. Quarmby V, Quan C, Ling V, Compton P, Canova-Davis E 1998 How much insulin-like growth factor I(IGF-I) circulates? Impact of standardization on IGF-I assay accuracy. *J Clin Endocrinol Metab*. Vol.83, 4:1211-1216.
28. Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGFBPs. *Growth Regulation* 4 (Suppl. 1):11-19

---

## Distributed by

### **BioVendor - Laboratorní medicína a.s.**

CTPark Modrice

Evropska 873

664 42 Modrice

CZECH REPUBLIC

Phone: +420-549 124 185

Fax: +420-549 211 460

e-mail: [info@biovendor.com](mailto:info@biovendor.com)

<http://www.biovendor.com>

---

## European Union:

### **BioVendor GmbH**

Im Neuenheimer Feld 583

D-69120 Heidelberg, GERMANY

Phone: +49-6221-433-9100

Fax: +49-6221-433-9111

e mail: [roth@biovendor.com](mailto:roth@biovendor.com)

e mail: [abels@biovendor.com](mailto:abels@biovendor.com)

## USA, Canada and Mexico:

### **BioVendor, LLC**

1463 Sand Hill Road, Suite 227

Candler, NC 28715, USA

Phone: +1-828-670-7807, +1-800-404-7807

Fax: +1-828-670-7809

e mail: [info-us@biovendor.com](mailto:info-us@biovendor.com)

## China - Hong Kong Office:

### **BioVendor Laboratories Ltd.**

Room 4008, Hong Kong Plaza, No.188

Connaught Road West

Hong Kong, CHINA

Phone: +852 28030523

Fax: +852 28030525

e-mail: [lu@biovendor.com](mailto:lu@biovendor.com)

## China – Mainland Office:

### **BioVendor Laboratories Ltd.**

Room 2405, YiYa Tower, TianYu Garden, No.150

Lihe Zhong Road, Guang Zhou, CHINA

Phone: +86-20-87063029

Fax: +86-20-87063016

e-mail: [jhsu@biovendor.com](mailto:jhsu@biovendor.com)

---

## Manufactured by



Aspenhastr. 25

D-71770 Reutlingen

Germany

<http://www.mediagnost.de>

e-mail: [contact@mediagnost.de](mailto:contact@mediagnost.de)

---