S100 ELISA

Immunosorbent assay for quantitative measurement of
S100B protein in human serum

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1. INTENDED USE

This product is intended FOR *IN VITRO* DIAGNOSTIC USE ONLY.

Immunoassay for the in vitro quantitative determination of S100 in human serum. S100 is intended for use as an aid in the management of patients suffering from malignant melanoma.

2. PRINCIPLE OF THE TEST

The S100 ELISA is a two-site, one-step, enzyme linked immunosorbent assay. In the assay calibrators, controls and unknown samples react simultaneously with 2 solid phase capture antibodies and a detector antibody conjugated with horseradish peroxidase (HRP) during the incubation in the microtiter wells. After a washing step a TMB chromogen (Tetramethylbenzidine) is added and the reaction is allowed to proceed for 15 minutes. The enzyme reaction is stopped by adding a Stop Solution and the absorbance is measured at 450 nm.

3. KIT CONTENTS

The kit contains reagents for 96 determinations.

3.1 Anti-S100B - coated strips: ready to use reagent

1 x 96 microtiter wells coated with 2 monoclonal mouse anti-S100B antibodies. Contains BSA. Allow strips and holder to equilibrate to room temperature (20-25°C) in the foil pouch to protect wells from condensation. Once opened, reseal to assure proper closure in a desiccated environment.

NOTE: When the pouch containing strips is first opened, the desiccant indicator **MUST** be blue. If indication is not blue, do not use the strips.

3.2 HRP- conjugate: ready to use reagent

1 x 21 mL, blue-green liquid. Contains a monoclonal mouse anti-S100B antibody conjugated with HRP, BSA and 0.5% ProClin 300 as preservative.

3.3 Calibrators: lyophilized reagent

2 vials of each 6 calibrators. Reconstitute in 1.0 mL purified water. Calibrators consist of S100 bovine antigen. Calibrator S100 value is listed on the vial labels.

3.4 Controls: lyophilized reagent

2 vials of each 2 controls, Reconstitute in 1.0 mL purified water. Controls consist of S100 bovine antigen. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

3.5 Sample Diluent: ready to use reagent

1 x 10 mL, ready for use, colorless liquid. Contains 0.5% ProClin 300 as preservative.

3.6 Wash Buffer 10X: liquid solution

2 x 50 mL, a PBS-Tween concentrate. The 1X Wash Buffer is stable when stored at room temperature (20-25°C) for 5 days. **NOTE:** Do not use if buffer is turbid.

3.7 TMB Solution: ready to use reagent

1 x 16 mL, Buffered substrate and chromogen, colorless, 0.05% TMB (3,3', 5,5' Tetramethylbenzidine). Protect from light. **NOTE:** Do not use if blue precipitate is present.

3.8 Stop Solution: ready to use reagent

1 x 20 mL, colorless liquid. Contains 0.4N Sulfuric Acid. See Warnings and Precautions section.

4. SHELF LIFE AND STORAGE OF THE KIT

The kit is stable at 2-8°C until expiry date on the label. After opening, store all reagents at 2-8°C until expiration date on the label, except for calibrators and controls that can be frozen at -20°C, thawed and reused once. Do not use the kit after the expiry date.

5. WARNINGS AND PRECAUTIONS

This product is intended FOR *IN VITRO* DIAGNOSTIC USE ONLY.

Reagents Containing 3,3',5,5' Tetramethylbenzidine

Chromogen: This product contains 3,3',5,5' Tetramethylbenzidine (TMB) (0.05%) which has shown possible mutagenic effects in laboratory experiments.

6. TRACEABILITY OF STANDARDS AND CONTROLS

The calibrators and controls of the kit consist of S100B protein. This S100B protein originates from purified and well-characterized bovine brain material. Calibrators and Controls are calibrated against in-house reference S100B material.

7. MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes for 50, 100 and 150 μ L.
- Vortex mixer.
- Purified water.
- Plate Shaker: 800 rpm, amplitude 1.5 mm.
- Microplate washing device.
- Microplate ELISA Reader 450 nm reading capability.
- Calculation of results: Cubic Spline algorithm.

8. SPECIMEN COLLECTION, PREPARATION AND STORAGE

- Serum is recommended. Do not use EDTA- or sodium citrate plasma.
- Collect samples using standard procedures.
- During blood collection ensure that the sample is not contaminated with EDTA or sodium citrate. The blood should be processed and serum/plasma collected within 2–4 hr.
- Storage of serum samples at 2–8°C: 24 hours. For longer storage freeze to below –20°C. Samples are stable at –20°C for 6 months. Frozen sera must be gently but thoroughly mixed after thawing. Avoid repeated freeze-thaw cycles.
- Do not use samples, which are grossly lipemic or contaminated. Samples containing precipitates must be centrifuged before testing. Do not use heat-inactivated samples.
- High S100B samples (>5 μ g/L) can be diluted. S100 ELISA diluent should be used for dilution.

9. PREPARATION OF REAGENTS BEFORE USE

- Allow unopened reagents and samples to reach room temperature (20-25°C) before use.
- Calibrators: Reconstitute in 1.0 mL purified water. Let stand 20 minutes. Mix carefully.
- Controls: Reconstitute in 1.0 mL purified water. Let stand 20 minutes. Mix carefully.
- 10x Wash Buffer: Dilute 1:10 (1x Wash Buffer) with purified water.


10. STORAGE OF REAGENTS AFTER USE

- Reconstituted calibrators and controls are stable for 90 min at room temperature. If all of the calibrators and controls are not used immediately, they can be frozen at -20°C , thawed and reused once.
- All the other opened components should be stored at $2 - 8^{\circ}\text{C}$ in the original containers.
- The TMB substrate and stop solution, which are taken out of the original vials, **must be discarded**.
- Fresh 1x Wash Buffer should be prepared at least once every 5 days.

11. ASSAY PROCEDURE

- 11.1 Prepare all the reagents as described in the previous section. Mix the samples on a Vortex before pipetting. **Note:** The test must be started within 90 minutes of the reconstitution of the calibrators and controls.
- 11.2 Pipette 50 μ L of calibrators, controls and unknown samples into the wells.
- 11.3 Add 150 μ L conjugate to all wells.
- 11.4 Cover the plate and incubate for 2 hours on a plate shaker (800 rpm) at room temperature (RT).
- 11.5 Wash 3 times with 300 μ L of 1x Wash Buffer. **Note:** It may be necessary to decant sample and conjugate from plate manually prior to wash step if an automated plate washing system is used.
- 11.6 Add 100 μ L TMB substrate to all wells.
- 11.7 Cover the plate and incubate 15 ± 2 minutes on a plate shaker (800 rpm) at room temperature.
- 11.8 Stop the reaction by adding 100 μ L Stop Solution. Add the Stop Solution in the same order and speed, which was used for the TMB substrate.
- 11.9 Read the absorbance at 450 nm using a microplate reader within 15 minutes.

12. FLOW DIAGRAM OF ASSAY

13. QUALITY CONTROL

Each laboratory must include at least one high and one low control in every assay to monitor assay performance. The controls provided with the kit may be utilized for this purpose. The kit controls contain S100B protein at both high and low concentrations and have been assayed using the DiaSorin S100B kit to establish expected ranges. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. The controls should be treated as unknown samples and assayed in duplicate. Assays should be considered valid only if results meet the laboratories criteria for acceptability.

14. CALCULATION OF RESULTS

Interpretation of Results


The Cubic Spline algorithm should be used for calculation of results. When using a Microplate Reader with automatic data calculation refer to the users manual of the instrument.

This is an example only and should not be used in calculations.

Sample	ABS 450 nm	µg/L
CAL A	0.087	0.00
CAL B	0.154	0.12
CAL C	0.221	0.25
CAL D	0.339	0.49
CAL E	1.516	2.49
CAL F	2.692	4.46
CONTROL L	0.420	0.62
CONTROL H	1.113	1.65
Sample 1	0.109	0.40
Sample 2	0.389	0.59
Sample 3	1.223	2.03

Standard Curve

This is an example only and should not be used in calculations.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

Reference Range

Cut-off was determined to 0.15 µg/L (the 95%-ile of 100 blood donor samples). Each laboratory should, however, establish its own upper reference limit.

Accuracy

The accuracy of the results of this assay is directly related to the degree of care exercised during pipetting, vortexing, aspiration, and adherence to methodology and temperature requirements.

Measuring Range

The measuring range is up to 5 µg/L. Concentrations for high samples can be obtained by diluting with S100 ELISA Diluent and repeating the assay.

Precision

Serum samples at different concentration levels were evaluated running one assay per day over ten operating days. Intra-assay and Inter-Assay precision was estimated by analysis of variance (ANOVA). Within the range of concentrations from 0.18 to 4 µg/L, the within run imprecision is < 10% and the total imprecision is < 15%.

Sample µg/L	Intra-assay variation (CV%)	Inter-assay variation (CV%)
0.18	9.4	7.2
0.27	6.6	12.4
0.61	3.8	4.7
1.9	4.4	3.1
2.2	1.9	4.3
4.1	2.9	3.5

Sensitivity (LDL)

The detection limit is 0.03 µg/L (B0+ 3SD).

Dilution Linearity

Dilution linearity was assessed by diluting three clinical samples with S100B Sample Diluent and assaying. Dilution linearity demonstrates a functional assay range from 0.2 µg/L to 4.5 µg/L.

Sample	Dilution	Calculate Value (µg/L)	Recovery %
1	1:1	2.94	100
	1:2	3.03	103
	1:5	3.68	125
	1:10	3.29	112
2	1:1	3.08	100
	1:2	3.49	113
	1:5	3.80	123
	1:10	4.16	135
3	1:1	3.77	100
	1:2	3.76	100
	1:5	3.81	101
	1:10	3.74	99

Compounds and Cross-Reactivity

There were no significant differences observed in S100B concentration levels after the addition of heparin, heparin/protamine or propofol compared to the control at the concentrations tested.

Compounds Tested for Potential Interference

Drug Tested	Conc. Tested
Heparin	5 IU/mL
Heparin/Protamine	100 IU/mL/10 mg/mL
Propofol	80 µg/mL

Gao, F, Harris, DNF, Sapsed-Byrne, S. Publication Perfusion 1997; 12:163-165

Potential interfering endogenous substances were added at the following final concentrations: bilirubin, 0.125 µg/mL; triglycerides, 12.5 mg/mL; hemoglobin 1, 2.5 and 5 mg/mL. Each test samples was assayed in duplicate. Elevated hemoglobin may affect sample results.

Interference by Endogenous Substances

Potential Interferent	Conc. Tested	Mean Control Value	Mean Test Value	Difference	% Error
Bilirubin	125 µg/mL	0.083	0.090	0.007	8%
Triglycerides	12.5 mg/mL	0.084	0.083	0.001	2%
Hemoglobin	1 mg/mL	0.087	0.091	0.004	5%
	2.5 mg/mL	0.087	0.097	0.010	11%
	5 mg/mL	0.087	0.092	0.005	6%

16. LIMITATIONS OF THE PROCEDURE

High-Dose Hook

No high-dose hook effect is seen up to 3000 µg/L.

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1. Abraha, H.D., et al., Serum S-100 Protein: a Potentially useful Prognostic marker in cutaneous melanoma, **Br. J. Dermatol.**, 1997; 137: 381-385.
2. Abraha, H.D., et al., Serum S-100 Protein, Relationship to Clinical Outcome in Acute Stroke, **Ann. Clin. Biochem.**, 1997; 34: 366-370.
3. Blomquist, S., et al., The Appearance of S-100 Protein in Serum During and Immediately After Cardiopulmonary Bypass Surgery: A Potential Marker for Cerebral Injury, **J. Cardiothorac. Vasc. Anaesth.**, 1997; 11: 699-703.
4. Bonfrère, J.M.G., et al., The Luminescence Immunoassay S-100B: A Sensitive Test to Measure Circulating S-100B: Its Prognostic Value in Malignant Melanoma, **Br. J. Cancer.**, 1998; 77: 2210-2214.
5. Cho, K.M., et al., Immunohistochemical Study of Melanocytic Nevus and Malignant Melanoma With Monoclonal Antibodies Against S-100 Subunits, **Cancer**, 1990; 66: 765-771.
6. Djureen-Mårtensson, E., et al., Serum S-100B Protein as a Prognostic Marker in Malignant Cutaneous Melanoma, **J. Clin. Oncol.**, 2001; 19: 824-831.
7. Gaynor, R., et al. S-100 protein: a marker for human malignant melanomas, **Lancet** 1981; 1: 869-871.
8. Haglid, K., et al. An Immunological Study of Human Brain Tumours Concerning the Brain Specific Proteins S-100 and 14.3.2., **Acta Neuropathol.**, 1973; 24: 187-196.
9. Hauschild, A., et al., Predictive Value of Serum S-100B for Monitoring Patients with Metastatic Melanoma During Chemotherapy and/or Immunotherapy, **Br. J. Dermatol.**, 1999; 140: 1065-1071.
10. Hauschild, A., et al., Prognostic Significance of Serum S100B Detection Compared with Routine Blood Parameters in Advanced Metastatic Melanoma Patients, **Melanoma. Res.**, 1999; 9: 155-161.
11. Hauschild, A., et al., S-100B Protein Detection in Serum is a Significant Prognostic Factor in Metastatic Melanoma, **Oncology**, 1999; 56: 383-344.
12. Heizmann, C.W., et al., S100 Proteins: Structure, Functions and Pathology, **Front Biosci.**, 2002; 7: 1356-1368.
13. Henze, G., et al., Serum S-100 – A Marker for Disease Monitoring in Metastatic Melanoma, **Dermatol.**, 1997; 194: 208-212.
14. Hidaka, H., et al., Purification and Characterization of Adipose Tissue S-100B Protein, **J. Biol. Chem.**, 1985; 258: 2705-2709.
15. Ingebrigtsen, T., et al., Increased Serum Levels of Protein S-100 After Minor Head Injury: A Biochemical Serum Marker with Prognostic Value? **J. Neurol. Neurosurg. Psychiatry**, 1995; 59: 103-104.
16. Ingebrigtsen, T., et al., Serial S-100 Protein Serum Measurements Related to Early Magnetic Resonance Imaging after Head Injury, **J. Neurosurg.**, 1996; 85: 945-948.
17. Ingebrigtsen T et al. Management of Minor Head Injury: The Value of Early Computed Tomography and Serum Protein S-100 Measurements, **J. Clin. Neurosci.**, 1997; 4: 29-33.
18. Ishiguro, Y., et al., Determination of Three Enolase Iso- Enzymes and S-100 Protein in Various Tumors in Children, **Cancer Res.**, 1983; 43: 6080-6084.
19. Isobe, T., et al., S-100aa Protein is Present in Neurons of Central and Peripheral Nervous System, **J. Neurochem.**, 1984; 43: 1494-1496.
20. Jury, C.S., et al., Rising Levels of Serum S100 Protein Precede Other Evidence of Disease Progression in Patients with Malignant Melanoma, **Br. J. Dermatol.**, 2000; 143: 269-274.

21. Jönsson, H., et al., Significance of S-100 Release After Coronary Artery Bypass Grafting, **Ann. Thorac. Surg.**, 1998; 65: 1639-16442.
22. Kato, K. and Kimura, S., S-100 α (aa) Protein is Mainly Located in the Heart and Striated Muscles, **Biochim Biophys. Acta**, 1985; 842: 146–150.
23. Krähn, G., et al., S100b is a More Reliable Tumor Marker in Peripheral Blood for Patients with Newly Occurred Melanoma Metastases Compared with MIA, Albumin and Lactate-Dehydrogenase, **Anticancer Res.**, 2001; 21: 1311–1316.
24. Nakajima, et al., An Immunoperoxidase Study of S-100B Protein Distribution in Normal and Neoplastic Tissues. **Am .J. Pathol.**, 1984; 116: 497-503.
25. Rosén, H., et al., Increased Serum Levels of the S-100 Protein are Associated with Hypoxic Brain Damage after Cardiac Arrest, **Stroke**, 1998; 29: 473-477.
26. Schultz, E.S., et al., Clinical and Prognostic Relevance of Serum S-100B Protein in Malignant Melanoma, **Br. J. Dermatol.**, 1998; 138: 426-430.
27. Semba, R., et al., Purification of S-100 Protein from Rat Kidney, **Brain Res.**, 1987; 401: 9–13.
28. Stefansson, K., et al., S-100 Protein in Tissue Tumors Derived from Schwann Cells and Melanocytes, **Pathol.**, 1982; 106: 261–268.
29. Stefansson, K., et al., S-100 Protein in Human Chondro- Cytes., **Nature**, 1982; 295: 63–64.
30. Weiss, S., et al., Value of S-100 Protein in the Diagnosis of Soft Tissue Tumors with Particular Reference to Benign and Malignant Schwann Cell Tumors. **Lab. Invest.**, 1983; 49: 299–308.
31. Westaby, S., et al., Serum S-100 Protein: A Potential Marker for Cerebral Events During Cardiopulmonary Bypass. **Ann. Thorac. Surg.**, 1996 ; 61: 88–92.
32. Wunderlich, M.T., et al. The Early Neurobehavioural Out- Come neurobehavioural Outcome After Stroke is Related to the Release of Neurobiochemical Markers for Brain Damage. **Stroke**, 1999; 30: 1190-1195.
33. Zimmer, D., et al. The S-100 Protein Family: History, Function and expression, **Brain Res. Bull.**, 1995; 37: 417–429.

The data presented in this instruction booklet has been carefully compiled from our records and from the scientific literature and we believe them to be accurate and reliable. We do not, however, give any warranties or make any representations with respect hereto, nor is freedom from any patent to be inferred.

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
SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

	Português	Dansk	Svenska	Ελληνικά
	Conformidade Europeia	Europæisk øverensstømmeise	Europeisk øverensstømmeise	Ευρωπαϊκή Συμμόρφωση
	Prazo de validade	Udløbsdato	Utgångsdatum	Ημερομηνία Λήξης
	Fabricante	Producent	Tilverkare	Κατασκευαστής
	Consultar as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Συμβουλευτείτε τις Οδηγίες Χρήσης
IVD	Diagnóstico in vitro	In vitro diagnostik.	Diagnostik in vitro.	In vitro διαγνωστικό ιατροτεχνολογικό προϊόν.
	Limites de Temperatura	Temperaturgrænse	Temperatur- begrænsning.	Περιορισμοί θερμοκρασίας
SORB	Fase sólida (poços revestidos)	Fast fase (Coatede slanger)	Solid fas (Belagda strångar)	Στερεή φάση (/Επικαλυμμένες χάντρες)
DIL SPE	Diluyente da Amostra	Fortyndingsvæske, prøve	Provutspådningsmedel	Αραιωτικό δείγματος
Ab ¹²⁵ I	Anticorpo marcado com ¹²⁵ I	Tracer: antistof mærket med ¹²⁵ I	Spårelement: antikropp betecknad med ¹²⁵ I	Ιχνηθέτης: Αντίσωμα σημασμένο με ¹²⁵ I
CAL	Calibrador	Kalibrator	Kalibrator	Βαθμονομητής
CONTROL	Soro Controlo	Kontrolserum	Kontrollserum	Ορός μάρτυρα

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
+44(0) 1344 401 430
FAX: +44(0) 7884 050812

13534

34753
8/10

PRINTED IN U.S.A.