



## DRG<sup>®</sup> Human $\beta$ Amyloid 1 – 42 Colorimetric (COL-3995)



Revised 25 Aug. 2010 rm (Vers. 1.1)

### 1 NAME AND INTENDED USE

The DRG Human  $\beta$  Amyloid (hA $\beta$ 42) ELISA is for determination of hA $\beta$ 42 in samples (e.g., serum, tissue culture medium, etc.). The assay will recognize both natural and synthetic forms of hA $\beta$ 42. The anti-human A $\beta$ 42 antibody used in this kit is capable of selectively detecting A $\beta$ 42 and not A $\beta$ 40/A $\beta$ 43 (4-6).

**This kit has been configured for research use only and is not to be used in diagnostic procedures.**

### 2 PRINCIPLE OF THE TEST

The DRG Human A $\beta$ 42 kit is a solid phase sandwich Enzyme Linked-Immuno-Sorbent Assay (ELISA). A monoclonal antibody specific for the NH<sub>2</sub> terminus of hA $\beta$  has been coated onto the wells of the microtiter strips provided. Samples, including standards of known hA $\beta$  content, control specimens and unknowns are pipetted into these wells, followed by the addition of a rabbit antibody specific for the 1-42 sequence of hA $\beta$ . Bound rabbit antibody is detected by the use of a horseradish peroxidase-labeled anti-rabbit antibody.

After removal of excess anti-rabbit antibody, a substrate solution is added, which is acted upon by the bound enzyme to produce color. The intensity of this colored product is directly proportional to the concentration of hA $\beta$  (1-42) present in the original specimen.

**Revised 25 Aug. 2010 rm (Vers. 1.1)**

### **3 WARNINGS AND PRECAUTIONS FOR USERS**

#### **3.1 General Remarks**

- This kit contains materials with small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. Upon disposal, flush drains with a large volume of water to prevent azide accumulation. Avoid ingestion and contact with eyes, skin and mucous membranes. In case of contact, rinse affected area with plenty of water. Observe all federal, state and local regulations for disposal.
- All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and Health Administration when handling and disposing infectious agents.

#### **3.2 Remarks for the Test**

1. When not in use, kit components should be refrigerated. All reagents should be warmed to room temperature before use.
2. Microtiter plates should be allowed to come to room temperature before opening the foil bags. Once the desired number of strips has been removed, immediately reseal the bag and store at 2 - 8°C to maintain plate integrity.
3. Samples should be collected in pyrogen/endotoxin-free polypropylene tubes.
4. Samples should be frozen if not analyzed shortly after collection. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely, microfuge and mix well with diluent prior to analysis.
5. When analyzing samples, add a protease inhibitor cocktail with AEBSF (a serine protease inhibitor) and prepare the standard dilutions using the same diluent as used with the biological samples. Serine proteases can rapidly degrade A $\beta$ peptides, thus using AEBSF (water soluble and less toxic than PMSF) at a 1 mM final concentration is very helpful. Keep samples on ice until ready to apply to plate.
6. When possible, avoid use of badly hemolyzed or lipemic sera. If large amounts of particulate matter are present, centrifuge prior to analysis.
7. Standards can be dispensed into a reservoir and transferred to the plate using a multi-channel pipette. Be sure to aspirate all previous standard solution prior to adding the next standard (in ascending order) to the reservoir.
8. It is recommended that all standards and samples be run in duplicate. Prepare standards and samples in 5 mL polypropylene tubes. We recommend diluting samples 1:2 to 1:5 for A $\beta$ 42.
9. Sample matrix has a dramatic impact on A $\beta$ recovery. To ensure accurate quantitation, the standard curves must be generated in the same diluent as the samples.
10. Avoid use of sodium azide, thymol or thimerosal as sample preservatives. These agents inhibit measurement of A $\beta$ peptide.
11. The *Standard/Sample Diluent* is a sterile fluid and should be dispensed and used to avoid contamination. Use immediately after opening or store frozen.
12. Samples that are >1000 pg/mL should be diluted with Standard/Sample Diluent.
13. When pipetting reagents, maintain a consistent order of addition from well-to-well. This ensures equal incubation times for all wells.
14. Do not mix or interchange different reagent lots from various kit lots.
15. Do not use reagents after the kit expiration date.
16. Read absorbances within 2 hours of assay completion.
17. We recommend covering plate with adhesive plate covers during incubation steps.



**DRG<sup>®</sup> Human  $\beta$  Amyloid 1 – 42 Colorimetric (COL-3995)**



**Revised 25 Aug. 2010 rm (Vers. 1.1)**

18. All residual wash liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. **Never** insert absorbent paper directly into the wells.
19. Because Stabilized Chromogen is light sensitive, avoid prolonged exposure to light. Also avoid contact between Stabilized Chromogen and metal, or color may develop.

Revised 25 Aug. 2010 rm (Vers. 1.1)

#### 4 KIT COMPONENTS

Reagents	Quantity 96 Test kit	Preparation and Storage
<b>hA<math>\beta</math>42 Standard</b> , synthetic peptide Refer to vial label for quantity and reconstitution volume.	1 vials	Lyophilized 2 – 8°C
<b>Standard/Sample Diluent</b> ; Sterile filtered.	1 x 60 mL	2 – 8°C
<b>Antibody -Coated Wells</b> 96 wells per plate. Plate pre-coated with mAb to NH <sub>2</sub> terminus of A $\beta$	1 plate	2 – 8°C
<b>Detection Antibody Diluent</b> . Contains 15 mM sodium azide.	1 x 11 mL	2 – 8°C
<b>Detection Antibody</b> . Rabbit anti-hA $\beta$ 42. Contains 15 mM sodium azide.	1 x 50 $\mu$ L	2 – 8°C
<b>Secondary Antibody (100x)</b> . Anti-Rabbit IgG (horseradish peroxidase labeled). Contains 3.3 mM thymol.	1 x 125 $\mu$ L	Concentrate 2 – 8°C
<b>Secondary Antibody Diluent</b> : Contains 3.3 mM thymol.	1 x 25 mL	2 – 8°C
<b>Wash Buffer (25x)</b>	1 x 100 mL	Concentrate 2 – 8°C
<b>Stabilized Chromogen</b> , Tetramethylbenzidine (TMB);	1 x 25 mL	2 – 8°C
<b>Stop Solution</b>	1 x 25 mL	2 – 8°C
<b>Plate Covers</b> , adhesive strips	3	2 – 8°C

**NOTE:** Once opened, standard, standard diluent, detection and secondary antibodies should be stored at -20 to -80°C for improved shelf life.

#### 5 MATERIALS REQUIRED BUT NOT SUPPLIED

- Standard Reconstitution Buffer [55 mM Sodium Bicarbonate Buffer (NaHCO<sub>3</sub>, ultrapure grade), pH 9.0].
- Calibrated adjustable precision pipettes, preferably with disposable plastic tips (a manifold multi-channel pipette is desirable for large assays), beakers and graduated cylinders.
- Reagent reservoirs for dispensing standards, antibody solutions and substrate.
- Plate washer: automated or manual (squirt bottle, manifold dispenser, etc.).
- Deionized or distilled H<sub>2</sub>O.



Revised 25 Aug. 2010 rm (Vers. 1.1)

- Microtiter plate reader capable of measurement at or near 450 nm.
- Shaking platform (for low to moderate shaking) with plate and tube vortex adapter, or a mini-orbital shaker.
- Disposable 12 x 75 mm polypropylene tubes for diluting standards and samples.
- 4-(2-aminoethyl)-benzenesulfonyl fluoride (AEBSF) or protease inhibitor cocktail containing AEBSF.
- Absorbent paper towels.
- Graph paper: linear (Cartesian), log-log, or semi-log, as desired.

## Assay Procedure

### 5.1 Directions for Washing

**Incomplete washing will adversely affect the test outcome.** All washing must be performed with *Wash Buffer* provided.

Washing can be performed manually as follows: completely aspirate the liquid from all wells by gently lowering an aspiration tip (aspiration device) into the bottom of each well. Take care not to scratch the inside of the well.

After aspiration, fill the wells with at least 0.4 mL of diluted wash solution. Let soak for 15 to 30 seconds, then aspirate the liquid. Repeat as directed under **Assay Procedure**. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

Alternatively, the wash solution may be put into a squirt bottle. If a squirt bottle is used, flood the plate with wash buffer, completely filling all wells. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

If using an automated washer, we suggest 4 wash cycles for all wash steps except the last wash step where 5 wash cycles are recommended.

### 5.2 Preparation and Storage of Reagents

This assay has been calibrated against the mass determination of highly purified native A $\beta$  where mass was corrected for peptide content by amino acid analysis.

#### 5.2.1 Preparation of Standard Reconstitution Buffer

Dissolve 2.31 grams of sodium bicarbonate in 500 mL of distilled or deionized water.

Add 2 N sodium hydroxide until pH is 9.0. Filter solution through a 0.2  $\mu$ m filter unit.

#### 5.2.2 Reconstitution and Dilution of hA $\beta$ 42 Standard

*Note:* Polypropylene tubes may be used for standard dilutions.

1. Remove the *hA $\beta$ 42 Standard* vial from storage and let equilibrate to room temperature (RT). Reconstitute the *hA $\beta$ 42 Standard* to 1.0  $\mu$ g/mL with Standard Reconstitution Buffer (55 mM sodium bicarbonate, pH 9.0). Refer to the standard vial label for instructions. Swirl or mix gently and allow vial to sit for 5 minutes at room temperature then transfer to ice for 90 minutes with occasional mixing. Briefly vortex prior to preparing standards.
2. Generation of the standard curve using the A $\beta$  peptide standards provided in the kit must be performed using the same composition of buffers used for the diluted experimental samples. For example, if brain extracts are diluted 1:10 with *Standard/Sample Diluent*, then the buffer used to dilute standards should be 90% *Standard/Sample Diluent* and 10% brain extraction buffer (including AEBSF at a final concentration of 1 mM).

Revised 25 Aug. 2010 rm (Vers. 1.1)

3. Add 0.100 mL of the reconstituted standard to a tube containing 0.900 mL *Standard/Sample Diluent* or as otherwise dictated by the treatment of experimental samples. Label as 100,000 pg/mL hA $\beta$ 42. Mix.
4. Add 0.100 mL of the 100,000 pg/mL standard to a tube containing 0.900 mL *Standard/Sample Diluent*. Label as 10,000 pg/mL.
5. Add 0.200 mL of 10,000 pg/mL standard to a tube containing 1.8 mL of *Standard/Sample Diluent*. Label as 1000 pg/mL.
6. Add 1.0 mL of *Standard/Sample Diluent* to each of 6 tubes labeled 500, 250, 125, 62.5, 31.25, 15.63, 0 pg/mL hA $\beta$ 42.
7. Make serial dilutions of the standard as described in the dilution table on the next page. Mix thoroughly between steps.

**Dilution table:**

Standard:	Add:	Into:
1000 pg/mL	Prepare by adding 0.2 mL of 10,000 pg/mL standard to 1.8 mL of Diluent Buffer.	
500 pg/mL	1.0 mL of the 1000 pg/mL std.	1.0 mL of the Diluent Buffer
250 pg/mL	1.0 mL of the 500 pg/mL std.	1.0 mL of the Diluent Buffer
125 pg/mL	1.0 mL of the 250 pg/mL std.	1.0 mL of the Diluent Buffer
62.5 pg/mL	1.0 mL of the 125 pg/mL std.	1.0 mL of the Diluent Buffer
31.25 pg/mL	1.0 mL of the 62.5 pg/mL std.	1.0 mL of the Diluent Buffer
15.63 pg/mL	1.0 mL of the 31.25 pg/mL std.	1.0 mL of the Diluent Buffer
0 pg/mL	1.0 mL of the Diluent Buffer	An empty tube

Remaining reconstituted hA $\beta$ 42 standard may be stored at -80°C for up to 4 months.  
Return the Standard Diluent Buffer to the freezer.

### 5.2.3 Preparation of Samples

Prepare one or more dilutions of each sample. These dilutions should be made in *Standard/Sample Diluent*, although the exact dilution must be determined empirically (e.g., 1:2 and 1:5 represent a reasonable range). This dilution must be performed because certain components in samples can interfere with the detection of the A $\beta$  peptides and/or A $\beta$  levels are too high. AEBSF should be added to the diluted samples and the standards at a final concentration of 1 mM in order to prevent proteolysis of the A $\beta$  peptides.

### 5.2.4 Preparation of Detection Antibody Solution

Prepare Detection Antibody Solution by thoroughly mixing contents of the Detection Antibody vial with the entire contents of the Detection Antibody Diluent (11 mL).

Revised 25 Aug. 2010 rm (Vers. 1.1)

### 5.2.5 Preparation of Secondary Antibody Solution

**Please Note:** The Anti-rabbit IgG-HRP 100x concentrate is in 50% glycerol. This solution is viscous. To ensure accurate dilution, allow Anti-rabbit IgG-HRP concentrate to reach room temperature. Gently mix. Pipette Anti-rabbit IgG-HRP concentrate slowly. Remove excess concentrate solution from pipette tip by gently wiping with clean absorbent paper.

1. Dilute 10  $\mu$ L of this 100x concentrated solution with 1 mL of Secondary Antibody Diluent for each 8-well strip used in the assay. Label as Anti-rabbit IgG-HRP Working Solution.

*For Example:*

# of 8-Well Strips	Volume of Secondary Antibody Concentrate (HRP 100x conc.)	Volume of Diluent
2	20 $\mu$ L solution	2 mL
4	40 $\mu$ L solution	4 mL
6	60 $\mu$ L solution	6 mL
8	80 $\mu$ L solution	8 mL
10	100 $\mu$ L solution	10 mL
12	120 $\mu$ L solution	12 mL

2. Return the unused Secondary Antibody (HRP 100x concentrate) to the refrigerator.

### 5.2.6 Preparation of Wash Buffer

Allow the 25x concentrate to reach room temperature and mix to ensure that any precipitated salts have redissolved. Dilute 1 volume of the 25x wash buffer concentrate with 24 volumes of deionized water (e.g., 50 ml may be diluted up to 1.25 liters, 100 ml may be diluted up to 2.5 liters). Label as Working Wash Buffer.

Store both the concentrate and the Working Wash Buffer in the refrigerator.

The diluted buffer should be used within 14 days.

## 5.3 Assay Procedure

**Be sure to read the *Warnings and Precautions for Users* section before carrying out the assay.**

**Note:** The protocol described below has been developed to provide sensitive and reproducible detection of A $\beta$  peptides. Although the assay uses a standard double antibody ELISA approach combined with colorimetric detection, the inherent complexity in handling/detecting A $\beta$  peptides must be recognized.

1. Reconstitute hA $\beta$ 42 Standard.
2. Prepare samples and standards with appropriate diluents.
3. Wash plate 4x with Working Wash Buffer, then pat dry on a paper towel. See DIRECTIONS FOR WASHING.
4. Add 100  $\mu$ L/well of A $\beta$  peptide standards, controls, and dilutions of samples to be analyzed. Incubate for 2 hours while shaking on an orbital plate shaker at RT **or** overnight at 4°C without shaking.
5. Wash plate 4x with Working Wash Buffer and pat dry on a paper towel. See DIRECTIONS FOR WASHING.
6. Pour Detection Antibody Solution into reservoir; immediately add 100  $\mu$ L/well of Detection Antibody Solution to the plate. Incubate at RT for 2 hours while shaking.

Revised 25 Aug. 2010 rm (Vers. 1.1)

7. Wash plate 4x with Working Wash Buffer, then pat plate dry on a paper towel. See DIRECTIONS FOR WASHING.
8. Pour Secondary Antibody Solution into reservoir; immediately add 100  $\mu$ L/well of Secondary Antibody Solution to the plate. Incubate plate at RT for 2 hours with shaking.
9. Wash plate 5x with Working Wash Buffer, pat dry on a paper towel. See DIRECTIONS FOR WASHING.
10. Add 100  $\mu$ L of Stabilized Chromogen to each well. The liquid in the wells will begin to turn blue.
11. Incubate for **30 minutes at room temperature and in the dark. Please Note: Do not cover the plate with aluminum foil or metallized mylar.** The incubation time for chromogen substrate is often determined by the microtiter plate reader used. Many plate readers have the capacity to record a maximum optical density (O.D.) of 2.0. The O.D. values should be monitored and the substrate reaction stopped before the O.D. of the positive wells exceed the limits of the instrument. The O.D. values at 450 nm can only be read after the *Stop Solution* has been added to each well. If using a reader that records only to 2.0 O.D., stopping the assay after 20 to 25 minutes is suggested.
12. Add 100  $\mu$ L of Stop Solution to each well. Tap side of plate gently to mix. The solution in the wells should change from blue to yellow.
13. Read the absorbance of each well at 450 nm having blanked the plate reader against a chromogen blank composed of 100  $\mu$ L each of Stabilized Chromogen and Stop Solution.  
Read the plate within 2 hours after adding the Stop Solution.

## 6 CALCULATION OF THE RESULTS

1. Plot on graph paper the absorbance of the standards against the standard concentration. (Optimally, the background absorbance may be subtracted from all data points, including standards, unknowns and controls, prior to plotting.) Draw the best smooth curve through these points to construct the standard curve. If using curve fitting software, the four parameter algorithm provides the best curve fit.
2. Read the hA $\beta$ 1-42 concentrations for unknown samples and controls from the standard curve plotted in step 1.  
**Multiply value(s) obtained for serum/plasma sample(s) by the appropriate factor to correct for the sample dilution in section 7.2.3** (Samples producing signals higher than that of the highest standard (1000 pg/mL) should be further diluted in *Standard/Sample Diluent Buffer* and reanalyzed, multiplying the concentration found by the appropriate dilution factor.)

### 6.1 Typical Data

The following data were obtained for the various standards over the range of 0 to 1000 pg/mL hA $\beta$ 42.

Standard hA $\beta$ 42 (pg/mL)	Optical Density (450 nm)
0	0.221
15.63	0.241
31.25	0.257
62.5	0.316
125	0.385
250	0.601
500	1.059
1000	1.864





**DRG<sup>®</sup> Human  $\beta$  Amyloid 1 – 42 Colorimetric (COL-3995)**



**Revised 25 Aug. 2010 rm (Vers. 1.1)**

## **7 LIMITATIONS OF PROCEDURE**

Do not extrapolate the standard curve beyond the 1000 pg/mL standard point; the dose-response is non-linear in this region and accuracy is difficult to obtain. Dilute samples >1000 pg/mL with Standard/Sample Diluent Buffer, reanalyze these and multiply results by the appropriate dilution factor.

The influence of various drugs, aberrant sera (hemolyzed, hyperlipidemic, jaundiced, etc.) and the use of biological fluids in place of serum samples have not been thoroughly investigated. The rate of degradation of native hA $\beta$  in various matrices has not been investigated. The immunoassay literature contains frequent references to aberrant signals seen with some sera, attributed to heterophilic antibodies. Though such samples have not been seen to date, the possibility of this occurrence cannot be excluded.

**This kit is for research use only. Not for human therapeutic or diagnostic use.**



**DRG® Human  $\beta$  Amyloid 1 – 42 Colorimetric (COL-3995)**



**Revised 25 Aug. 2010 rm (Vers. 1.1)**

**8 REFERENCES**

1. Cotman, C.W. (1997) The  $\beta$ -amyloid peptide, peptide self-assembly, and the emergence of biological activities. A new principle in peptide function and the induction of neuropathology. *Ann. N. Y. Acad. Sci.* 814:1-16.
2. Estus, S., et al. (1997) Aggregated  $\beta$ -amyloid protein induces cortical neuronal apoptosis and concomitant "apoptotic" pattern of gene induction. *J. Neurosci.* 17:7736-7745.
3. Bancher, C., et al. (1997) Mechanisms of cell death in Alzheimer's disease. *J. Neural Transm. Suppl.* 50:141-152.
4. Savage, M.J., et al. (1998) Turnover of  $\beta$ -amyloid protein in mouse brain and acute reduction of its level by phorbol ester. *J. Neurosci.* 18:1743-1752.
5. Howland, D.S., et al. (1998) Modulation of secreted  $\beta$ -amyloid precursor protein and  $\beta$ -amyloid peptide in brain by cholesterol. *J. Biol. Chem.* 273:16576-16582.
6. Citron, M., et al. (1999)  $\beta$ -secretase cleavage of Alzheimer's amyloid precursor protein by the transmembrane aspartic Protease BACE. *Science* Vol. 286 Oct. 22:735-741.
7. Andreasen, N., et al. (1999) Cerebrospinal fluid  $\beta$ -amyloid (1-42) in Alzheimer disease: differences between early- and late-onset Alzheimer disease and stability during the course of disease. *Arch. Neurol.* 56 (6):673