



Revised 6 Jul 2005

IVD

NOTE:	This Instructions-For-Use (IFU) can be used for:
-------	--

Catalog No.: RIA1479 With Membrane (10 plate) AND RIA1480 Without Membrane (10 plate)

INTRODUCTION

INTENDED USE

The DRG Phenylalanine Microwell Enzyme Assay (PHEA-MW EA) is designed for the quantitative determination of phenylalanine in neonatal blood spots. It is for use as an aid in screening newborns for elevated phenylalanine levels.

CLINICAL PHYSIOLOGY

Phenylketonuria (PKU) is a congenital disease resulting in elevated blood levels of phenylalanine and excessive excretion of phenylpyruvic acid in urine. The disorder is characterized by irreversible mental retardation, skin abnormalities and behavioral disturbances (1). It is one of the most common hereditary amino acid disorders and affects approximately one in 15,000 newborns (2).

PKU is caused by the absence or deficiency of the hepatic enzyme, L-phenylalanine hydroxylase. This enzyme is necessary for the metabolic conversion of phenylalanine, an essential amino acid, to tyrosine. The deficiency results in an increase in both phenylalanine and its deaminated metabolite, phenylpyruvic acid, in plasma, urine and cerebrospinal fluid (3-7).

Excessive accumulation of phenylalanine and related metabolites in the body, especially cerebrospinal fluid, interferes with proper development of the brain. The injury to the cerebral tissue begins within the second and third week of life and progresses with time, becoming maximal at eight to nine months. Brain damage and the resulting progressive mental retardation may be minimized if the newborn is placed on a low phenylalanine diet soon after birth (1, 8-10).

Screening for elevated levels of phenylalanine in neonates is important for early detection and treatment of PKU. Phenylalanine is frequently determined in plasma or whole blood during the first few days after birth (11). During the mid 1960's, the United States began a program requiring mass screening of newborn infants using a dried blood spot specimen (12,13). This has since become widespread practice throughout the world and has resulted in early detection and dietary management of PKU in infants (14-18). It is now generally accepted that the mental retardation associated with PKU can be prevented by a low phenylalanine diet, provided that the disease is detected in early infancy.

Syndromes other than PKU have been recognized to demonstrate high blood phenylalanine levels. These syndromes are referred to as atypical PKU or benign persistent hyperphenylalaninemia (PHP). While classical PKU is caused by a deficiency in phenylalanine hydroxylase, the atypical or variant forms of PKU are related to those enzymes affecting the production of tetrahydrobiopterin. Tetrahydrobiopterin is the cofactor required for the activity of phenylalanine hydroxylase. All patients with persistent hyperphenylalaninemia should be investigated to rule out the tetrahydrobiopterin-deficient forms of phenylketonuria (19).





Revised 6 Jul 2005



The DRG PHEA-MW EA is designed to aid in newborn screening for elevated phenylalanine levels. Following enzymatic oxidation of the phenylalanine in the blood-spot sample, a simple color reaction occurs demonstrating rapid and quantitative determination of phenylalanine levels in neonates.

CLINICAL APPLICATIONS

Newborn Screening: The DRG PHEA-MW EA is designed as an aid in screening newborns for PKU. Presumptive positive results should be confirmed by analysis with another method.

FACTORS AFFECTING NORMAL VALUES

- 1. Newborn Age: In patients affected with PKU, phenylalanine concentrations begin to rise progressively in the first 6-12 hours after birth (20). To clearly identify elevated phenylalanaine levels, it is recommended that sampling should occur after 24 hours of age. The American Academy of Pediatrics recommends that infants initially screened before 24 hours should be re-screened for PKU (18). Other determinates of when to take patient samples depend on national requirements. Please observe local established guidelines.
- 2. **Sample Age:** Although blood spot samples are stable for up to six (6) months under proper storage conditions, prompt analysis of samples for the presence of phenylalanine should take place according to screening facility protocol. Once the sample is collected it should be thoroughly dried (21).
- 3. Newborn Diet: A diet containing protein, such as breast milk or formula, is required for phenylalanine accumulation and subsequent identification of PKU in newborns (22, 23).

PRINCIPLE OF TEST

The DRG PHEA-MW EA is an enzyme assay which follows an acid extraction of the phenylalanine contained in a 3/16" (5 mm) blood spot. In this assay, a single blood spot is extracted in a special semi-porous membrane bottom plate. Following extraction, all 96 extracts are simultaneously transferred to a conventional microwell plate, by means of a vacuum manifold. A neutralizing buffer is added to each acid extract, followed by addition of a combination enzyme-substrate reagent. This rapidly oxidizes the phenylalanine to phenylpyruvate, reducing NAD to NADH⁺ in the process. The NADH⁺ reduces a colorless dye to form a colored end-product with an absorbance maximum at 570 nm. The absorbance read is directly proportional to the amount of phenylalanine in the sample.

NOTE: All references to phenylalanine concentrations are in whole blood units.

The DRG PHEA-MW EA is simple to perform, rapid and easily automated. The basic assay is diagramed below:





Revised 6 Jul 2005



		AMOUNT	PER KIT
COMPONENT	LABEL COLOR	Single Kit	Bulk Kit
Membrane Transfer Plate	Tan	2 plates	10 plates
Collection-Reaction Plate	Tan	2 plates	10 plates
PHE/GAL Blood Spot Standards	Green	1 card	3 cards
Trl-Level Blood Spot Controls	Green	1 card	3 cards
TCA Extraction Solution	Yellow	22 mL	110 mL
Phenylalanine Neutralizing Solution	Yellow	11 rnL	55 mL
Reagent A	Red	11 mL	55 mL
Phenylalanine Reagent B	Green	11 mL	55 mL

REAGENTS PROVIDED AND LABEL COLOR CODE

REAGENT DESCRIPTION AND PREPARATION

NOTE: Some reagents contain sodium azide which has a tendency to build up in lead or copper plumbing forming potentially explosive metal azides. Always flush large quantities of water through the plumbing after the disposal of these reagents.

A. MEMBRANE TRANSFER PLATE

96-well porous microwell plates. STORAGE: May be stored at room temperature or at 2-8°C with kit. STABILITY: Refer to kit expiration.







Revised 6 Jul 2005

B. COLLECTION-REACTION PLATE

96-well uncoated microwell plates. STORAGE: May be stored at room temperature or at 2-8°C with kit. STABILITY: Refer to kit expiration.

C. PHE/GAL BLOOD SPOT CALIBRATORS

Five standards are provided in mg/dL of whole blood. These standards are prepared in a bovine whole blood matrix containing sodium azide as a preservative and are spotted on Schleicher and Schuell's filter paper #903. Refer to foil package for concentrations.

STORAGE: 2-8°C. For long term storage greater than 7 days, store at -20°C. Reseal bag after use. STABILITY: Refer to expiration date on foil package.

D. TRI-LEVEL BLOOD SPOT CONTROLS

Three controls are provided. These controls are prepared in a bovine whole blood matrix containing sodium azide as a preservative and are spotted on Schleicher and Schuell's filter paper #903. Refer to foil package for concentration ranges.

STORAGE: 2-8°C. For long term storage greater then 7 days, store at -20°C. Reseal bag after use. STABILITY: Refer to expiration date on foil package.

E. TCA EXTRACTION SOLUTION

Trichloroacetic acid (TCA) in deionized water. STORAGE: 2-8°C. STABILITY: Refer to expiradon date on kit vial.

F. PHENYLALANINE NEUTRALIZING SOLUTION

Sodium carbonate (Na2CO3) in deionized water. STORAGE: 2-8°C. STABILITY: Refer to expiration date on kit vial.

G. REAGENT A

Substrate buffer with sodium azide (NaN3) as a preservative. STORAGE: 2-8°C. For long term storage greater than 7 days, store at -20°C. STABILITY: Refer to expiration date on kit vial.

H. PHENYLALANINE REAGENT B

Nicotinimide adenine dinucleotide (NAD) and phenylalanine dehydrogenase (PDH) in buffer with bovine calf serum and sodium azide (NaN3) as a preservative. PDH is a peptide produced by microbial fermentation with *Sporosarcina sp.* The enzyme contains no material of human, animal or plant origin. STORAGE: 2-8°C. For long term storage greater then 7 days, store at -20°C. STABILITY: Refer to expiration date on kit vial.

LIMITATIONS, PRECAUTIONS AND GENERAL COMMENTS

A. The reagents supplied in this kit are for In-Vitro Diagnostic Use Only.

B. Strict adherence to the protocol is advised to obtain reliable results. Any modifications made to the reagents or







Revised 6 Jul 2005



assay procedure are the responsibility of the user.

- C. A standard curve must be established for each run. A run should consist of a maximum of two consecutive plates. For larger assays, the timing of pipetting should be staggered to ensure uniform plate processing.
- D. Extrapolated values outside the range of standards (see foil package for concentration) are approximate and should be reported as "less than the value of the lowest standard" or "greater than the value of the highest standard".
- E. This assay is designed to be used with samples which are exclusively collected on SCHLEICHER and SCHUELL (S&S) FILTER PAPER #903. Changes in filter paper lots may affect patient results. All laboratories should record lot numbers of S&S #903 paper to monitor possible changes.
- F. The DRG PHEA-MW EA is to be used as an initial screening test for the quantitation of phenylalanine levels in whole blood spots. Individual screening laboratories are responsible to further test those patients with increased phenylalanine levels using documented confirmatory tests for PKU.
- G. Several variables not directly related to serum levels of phenylalanine may affect the determined value. These include hematocrit, moisture content of the paper, rate of blood deposition and lot of filter paper utilized. For these reasons, uniform collection techniques coupled with careful examination of the blood spots before analysis is recommended (24).
- H. A blood test to screen for elevated concentrations of phenylalanine should be performed on all newborn infants after the onset of feeding. For further information, please see recommendations from the American Academy of Pediatrics (18).
- I. Infants receiving special treatment or prolonged hospital stay should be tested in a timely manner after the onset of feeding. Feeding status must be monitored to ensure adequate protein intake in these infants (25).

SPECIMEN COLLECTION AND HANDLING

Infant screening programs may differ from one another in the amount of sample required. A heel stick sample collected from a 5/8" (16 mm) diameter blood spot pre-printed on Schleicher and Schuell Filter Paper #903 is suggested. The following summary is described in detail in NCCLS publication LA4-A4.

- A. Be sure the required information on the specimen card has been completed, including the infant's name, time and date of birth, and time and date of collection. Also indicate pre- or post-term, the infant's weight and feeding status and whether or not the infant was a twin.
- B. Collect the blood from the heel of an infant usually 24 to 72 hours postpartum. Sampling times may vary from center to center. See Limitations Section (items H and I).
- C. Wash the heel with soap and water and wipe dry. Swab area with alcohol and allow to air dry. Disinfectants containing iodine should be avoided.
- D. Puncture infant's heel with sterile lancet and wipe away first drop of blood. Make sure the tip of lancet is no longer than 2.5 mm. Allow another drop of blood of adequate volume to form, and gently touch the specimen





VD

Revised 6 Jul 2005



- E. Place the filter paper card horizontally on a clean surface and allow to air dry for at least 6 hours at room temperature. Avoid direct sunlight.
- F. Place each specimen in its own paper envelope and transport to the laboratory within 24 hours of drying.
- G. The receiving laboratory should store the sample at 2-8°C in a moisture proof environment shielded from direct light.
- H. Although phenylalanine has been reported to be stable in blood spots in excess of six (6) months, not all components are as stable and neonatal blood spot storage must be considered on an analyte by analyte basis.

EQUIPMENT AND REAGENTS REQUIRED

In addition to the materials provided with this kit, the following supplies are required. All items are available from DRG.

- 3/16" (5 mm) or 1/8" (3 mm) diameter paper punch.
- A. Plate reader able to read absorbance at 550-570 nm.
- B. Multi-channel and single channel micropipets calibrated to 50 μL and 100 $\mu L.$
- C. Rotary horizontal shaker.
- D. Vacuum manifold. (Requires a standard hose bench-top vacuum outlet or equivalent external vacuum source.)

ASSAY PROCEDURE

A. Assay Preparation

- 1. Prior to beginning assay, label each collection-reaction plate to match its complementary membrane transfer plate. Additionally, each laboratory should determine the time necessary for complete transfer of extract to the collection plate using the vacuum manifold. The time required should be between 30 and 60 seconds.
- 2. For a two plate assay (192 wells), transfer the entire contents of REAGENT A vial to the PHENYLALANINE REAGENT B vial. Transfer mixture between vial A and B a few times to ensure adequate mixing. Combine REAGENT A and PHENYLALANINE REAGENT B immediately prior to use. **Do not store or re-use.**
- 3. If less than 22.0 mL of REAGENT A and PHENYLALANINE REAGENT B is required, mix equal volumes of vial A and B for desired volume. Promptly store the remaining un-mixed reagents at 2-8°C, or at -20°C for long term storage.

B. Assay Steps

1. For the five (5) STANDARDS, three (3) CONTROLS, and variable unknowns: Punch one 3/16" (5 mm) or two 1/8" (3 mm) blood spot(s) into the appropriate wells of the membrane transfer plate. Be sure to leave an empty







IVD

well for the reagent blank.

- 2. Add 100 µL of TCA EXTRACTION SOLUTION to each well.
- 3. Rotate mix for 60 minutes at room temperature (R.T.).
- 4. Simultaneously transfer acid extracts to the empty collection-reaction plate using the vacuum manifold.
- 5. Add 50 μ L of PHENYLALANINE NEUTRALIZING SOLUTION to each well and shake gently by hand for 10 seconds.
- 6. Add 100 µL of pre-combined REAGENT A-B to each well and shake plate gently by hand for 10 seconds.
- 7. Incubate 45 minutes at R.T.
- 8. Read at 550-570 nm and calculate the results by plotting the standard curve on linear graph paper (See "Calculations" Section).

C. Quality Control

BLOOD SPOT CONTROLS containing low, intermediate and high levels of phenylalanine are included in this kit. They should be included in each assay run as unknowns in order to monitor the performance and reliability of the assay. Likewise, external blood spot controls (i.e., CDC, etc.) containing phenylalanine at several different levels should be routinely included in each run. Analysis of the results obtained should be done according to acceptance criteria established by the individual laboratory.





Revised 6 Jul 2005

IVD

PROTOCOL

WELLS	SAMPLE (3/16"- Blood Spot)	TC EX SO	CA TR. LN	Rotate & Incubate	Vacuum Manifold	NEI SO	UT. LN	REAG A-	BENT B	Incubate	Read Abs.
A1, B1	Reagent Blank	100	μL			50	μL	100	μL		
C1, D1	Standard 1										
E1, F1	Standard 2										
G1, H1	Standard 3				Transfer						
A2, B2	Standard 4			60	Extracts					45	550-570 nm
C2, D2	Standard 6			minutes	to					minutes	
E2, F2	Control I			at	Collection					at	
G2 H2	Control II			R.T.	Plate					R.T.	
A3, B3	Control III										
C3, D3	Unknown 1										
E3, F3	Unknown 2										
Etc.	Etc.	,	~			~	~		ļ		

CALCULATIONS

- A. Average the absorbance duplicates for all standards, controls and patients. Subtract the average reagent blank absorbance from each of the averages obtained above. This yields the *net absorbance (Abs.)*.
- B. Construct the standard curve by plotting the net absorbance (y-axis) versus the concentration of the phenylalanine standards (x-axis) using linear graph paper and a weighted linear curve fit. This yields the *standard curve*.
- C. Using the standard curve, determine the phenylalanine concentrations of each patient sample. Read patient samples (Abs. *Blank*) directly off curve as mg/dL phenylalanine in whole blood. A sample assay and calibration curve are provided below.

<u>۲</u>۲(۲ D



Revised 6 Jul 2005

|--|

SAMPLE	ABS	ABS-BL	CONC (mg/dL)	PHE (mg/dL)
Blank	0.121			
Std. 1	0.166	0.045	1.0	
Std. 2	0.242	0.121	4.0	
Std. 3	0.286	0.165	7.0	
Std. 4	0.431	0.310	15.0	
Std. 5	0.671	0.550	27.0	
Control I	0.206	0.085		2.9
Control II	0.318	0.197		8.7
Control III	0.515	0.394		18.9

The following is a sample standard curve and related information as illustrated from a computer assisted data reduction program.







Revised 6 Jul 2005

EXPECTED VALUES

NEONATES (24 - 48 hours): 0.6 to 3.6 mg/dL whole blood

Literature indicates that serum phenylalanine levels in normal, full term infants are 2.1 ± 0.5 mg/dL (4). The actual values range from 0.6 - 3.6 mg/dL. A study performed using the DRG PHEA-MW EA produced similar results of 2.0 ± 0.7 mg/dL.

As with any diagnostic test, differences in physiological ranges may be encountered from laboratory to laboratory due to patient demographics, laboratory techniques, and population sampling. These ranges should only be used as a guideline. We recommend each laboratory establish its own ranges using a statistically significant number of characterized patient specimens.

A. Normal Range Analysis

The DRG PHEA-MW EA was compared with a commercially available quantitative method for phenylalanine determination. Data from this study (n = 1046) are highly compatible with a Gaussian (normal) distribution, with an estimated mean of 2.04 mg/dL and an estimated standard deviation of 0.685 mg/dL. The 95% confidence interval for the mean is 2.00 to 2.09 mg/dL. The actual observed range of the assay in the database is 0.13 to 4.9 mg/dL. The 95th percentile corresponds to the value 3.23 mg/dL and the 99th percentile to the value 3.7 mg/dL. The median and the mode are both 2.03 mg/dL.



B. Cut-Off Derivation

The table below illustrates the statistically derived cut-off for phenylalanine levels in the DRG PHEA-MW EA and a comparison assay at the 95th, 99th and 99.9th percentile. These data are based on results generated at an established testing facility.





Revised 6 Jul 2005

	DRG PHEA-MW EA (mg/dL)	COMPARISON ASSAY (mg/dL)
95th percentile	3.2	2.3
99th percentile	3.7	2.7
99.9th percentile	4.2	3.1

C. Referral Rate Analysis

The table below demonstrates the variable presumptive positive rate observed in the DRG PHEA-MW EA utilizing a broad range of cut-off values for phenylalanine. The data below was generated at a state testing facility utilizing Inter-assay variation, or precision was calculated by evaluating the a population of normally distributed neonatal blood spot samples (n = 1044). Definition of a cut-off lower than 4.0 mg/dL may enhance screening security through re-testing of presumptive positive samples. National reference intervals, which also determine the need for re-testing, might exist.

DRG PHEA-MW EA Cut-Off (mg/dL)	Expected Presumptive Positive Rate (%)
3.0	8.33
3.1	6.42
3.2 (95% Normal Range)	5.36
3.3	3.93
3.4	2.87
3.5	2.01
3.6	1.44
3.7 (99% Normal Range)	0.96
3.8	0.67
3.9	0.57
4.0	0.48







Revised 6 Jul 2005

PERFORMANCE CHARACTERISTICS

A. Inter-Assay Variation

Inter-assay variation, or precision was calculated by evaluating the same tri-level samples in multiple assay runs.

	Level I n = 42	Level II n = 42	Level III n = 42
Mean	2.85	8.63	19.0
S.D.	0.39	0.77	1.46
% C.V.	13.7	8.9	7.7

B. Intra-Assay Variation

Intra-assay variation, or precision was calculated by assaying 10 replicates each of three tri-level samples.

	Level I n = 10	Level II n = 10	Level III n = 10
Mean	3.09	9.33	20.4
S.D.	0.38	0.45	1.18
% C.V.	12.3	4.8	5.8

C. Functional Sensitivity

Although lower doses are distinguishable from zero, it is recommended that, based on the functional sensitivity of the assay, all values which fall below the lowest standard (Standard 1), should be reported out as "less than the value of the lowest standard".

D. Parallelism (Linearity Of Dilutions)

Results shown represent actual values multiplied back by the indicated dilution factor.

SAMPLE	PHE (mg/dL)	1:2 PHE (mg/dL)	1:4 PHE (mg/dL)
1	27.8	25.8	25.8
2	35.4	29.9	31.1
3	32.6	29.6	30.8

E. Comparison Study

1. Patient Sample Correlation

Results from samples with values distributed throughout the quantitative range of this assay were compared with those obtained with an a commercially available method. The correlation





Revised 6 Jul 2005



coefficient was 0.969 (slope = 1.125, y-intercept = -0.494 mg/dL). A paired t-test did not show a difference in the values obtained by the two methods.

415 patients were selected from a statewide screening program for PKU. They were stratified over a broad range of measurements and were chosen from low, intermediate and presumptive positive sample pools, as well as from elevated and known positive PKU patients.

Two samples in this study gave poor correlation between the DRG and the comparison methods. However, both methods indicated that the samples were well above the usual clinical cut-off.



2. Relative Presumptive Postitive Analysis

The table below classifies the 1044 cases of the normal range study according to normal and presumptive positive results for the 95th percentile cut-off. Utilizing the statistically derived cutoff values, the presumptive positive rates are similar for both methods (4.2% and 5.36% respectively). Based on statistical analysis, the discordance observed between the two is due to random method variation.

No confirmed positives were encountered in this naturalistic population. We have shown, however, that inclusion of all elevated and known PKU positive samples available for reanalysis in both assays (n = 54), produced all positive results at the highest examined cut-off of 4.0 mg/dL.

Frequency	DRG Group < 2.3 mg/dL	Comparison Group > 2.3 mg/dL	Total
< 3.2 mg/dL	959	29	988 (94.64%)
> 3.2 mg/dL	42	14	56 (5.36%)
Total	n = 1001 (95.88%)	n = 43 (4.12%)	1044 (100%)









F. Recovery

Results shown illustrate the percent recoveries of samples to which known amounts of phenylalanine were added.

Sample Number	Expected PHE (mg/dL)	Observed PHE (mg/dL)	% PHE Recovered
1	4.25	4.33	102
2	8.29	7.72	93
3	14.0	13.1	94
4	6.17	5.54	90
5	10.2	9.96	98
6	15.9	15.8	99
7	12.2	12.1	100
8	17.8	17.7	99
9	21.9	25.2	115
			Average Recovery 98.9%

G. Interfering Substances

The following substances were added to a blood spot standard in order to determine their effect on the result. Ascorbic acid was found to interfere with the phenylalanine assay at superphysiologic levels. No interference was observed when ascorbic acid was tested at normal physiologic levels (0.2 - 2.0 mg/dL). Additionally, tetracycline was found to present a possible interference. Infants receiving tetracycline antibiotic therapy should be evaluated accordingly, especially if elevated phenylalanine results are observed. None of the other compounds interfered at the concentrations tested.





Revised 6 Jul 2005

	Amount Added (mg/dL)	Baseline PHE (mg/dL)	Observed PHE (mg/dL)
Amoxycillin	100	4.20	4.30
Ampicillin	100	4.20	4.20
Erythromycin	100	4.20	4.20
Gentamycin	100	4.20	4.50
Penicillin	100	4.20	4.50
Tetracycline	100	4.20	6.30
Bilirubin	20	4.10	4.20
Hemoglobin	100	4.20	4.20
Triglycerides	125	3.71	3.90
Cysteine	1.44	7.13	6.66
Cystine	1.44	7.13	6.73
Tryptophan	19.4	7.13	7.37
Tyrosine	17.9	7.13	7.44
Leucine	15.8	7.13	7.87
Methionine	3.7	7.13	6.65
Ascorbic Acid	11.0	7.13	10.8
Ascorbic Acid	1.1	0.0	0.19
Ascorbic Acid	1.44	0.0	0.52

IVD

ALTERNATE PROCEDURE (for manual transfer of extracts)

- 1. For the five (5) STANDARDS, three (3) CONTROLS, and variable unknowns: Punch one 3/16" (5 mm) blood spot or two 1/8" (3 mm) blood spot(s), into the appropriate wells of the plate. Be sure to leave an empty well for the reagent blank.
- 2. Add 200 µL of TCA Extraction Solution to each well.
- 3. Rotate mix for 60 minutes at room temperature.



Revised 6 Jul 2005

- 4. Using a single channel or multichannel pipet, transfer 100 μ L of each extract to the corresponding well of a new plate.
- 5. Add Neutralizing Solution and Reagent A-B as indicated in package insert.



IVD







VD

REFERENCES

- 1. Wolf, L.I., Inherited Metabolic Disorders: Errors of Phenylalanine and T y r o s i n e i n S e r u m . A d v . C l in. C h e m . e d . by Sobotka and Stewart. 1963 6. pp. 97-230.
- 2. Berry, H.K., Screening Newborns for Genetic Disease: The PKU Model. *Diagnostic Medicine* 7(1):50,1984.
- 3. Wong, P.K, O'Flynn, M.E., Inouye T., Micromethods for Measuring Phenylalanine and Tyrosine in Serum, Clin . Chem. 10(12):1098, 1964. (2,3,4,5).
- 4. Hsia, D.Y.Y., Berman, J.L, Slatis, H.M., Screening Newborn Infants for Phenylketonuria . *JAMA* 188: 203, 1964.
- 5. Jervis, G.A., Phenylpyruvic Oligophrenia Deficiency of Phenylalanine-Oxidizing System. *Proc*. Soc. Exp. Biol. Med 82:514, 1953.
- 6. Voss, J.C. and Weisman, H.A., The Phenylalanine Hydroxylase Contents of Livers of Various Vertebrates . *Comp Biochem, Physiol* 17:49, 1966.
- 7. Wallace, H.W., et. al. Studies on Conversion of Phenylalanine to Tyrosine in Phenylpyruvic Oligophrenia. *Proc. Soc Exp Biol Med* 94:632, 1957.
- 8. Hoffman, G.L. et. al., Dual Channel Continuous Flow System for Determination of Phenylalanine and Galactose: Application to Newborn Screening. *Clin Chem* 30(2): 287, 1984.
- 9. Scriver, C.R., and Clow, C.L., Phenylketonuria: Epitome of Human Biochemical Genetics Part 1. *N Engl J Med* 303:1336, 1980.
- 10. West, E.S. and Todd, W.R. Textbook of Biochemistry, 3rd ed. MacMillan, New York, 1964, pp. 1095-1098.
- 11. Kirkman, H.N., Carroll, C.L., Moore, E.G., et. al., Fifteen-Year Experience with Screening for Phenylke to nuria with an Automated Fluorometric Method. *Am J Hum Genet* 34:743-752, (1982).
- 12. Committee for the Study of Inborn Errors of Metabolism. National Research Council: Genetic Screening Programs, Principles and Research. Washington, D.C., National Academy of Sciences, 1975.
- 13. Reilly, P., Genetic Screening Legislation. Adv Hum Genet 5: 319, 1975.
- 14. Guthrie, R., Letters to the Journal: Blood Screening for Phenylketonuria. JAMA 178(8):863, 1961.
- 15. Guthrie, R., and Susi, A., A Simple Phenylalanine Method for Detecting Phenylketonuria in Large Populations of Newborn Infants. *Pediatrics* 32:338, 1963.
- 16. Scriver, C.R. and Clow, C.L., Phenylketonuria: Epitome of Human Biochemical Genetics Part II. *N Engl J Med* 303:1336 & 1394, 1980.
- 17. McCabe, E.R.B., McCabe, L., Mosher, G.A., et al., Newborn Screening for Phenylketonuria: PredictiveValidity as a Function of Age. *Pediatrics* 72(3):390, 1983.
- 18. American Academy of Pediatrics, Committee on Genetics: New Issues In NewbornScreening for Phenylketonuria and Congenital Hypothyroidism. Pediatrics 69(1):104, 1982.
- 19. Matalon, R., Michals, K., Lee, C.L., et. al., Screening for Biopterin Defects in Newborns with Phenylke to nuria and Other Hyperphenylalaninemias. *Annals of Clin Lab Sci* 12(5):411, 1982.
- 20. Jew, K., Kan, K., Koch, R., et al., Validity of Screening Early Collected Newborn Specimens for Phenylke to nuria Using a Fluorometric Method. *Screening* 3:1, 1994.
- 21. Madira, W.M., et.al., Clin Chem 38:2162,1992.
- 22. Dontanville, V.K., and Cunningham, G.C., Pediatrics 51:531, 1973.
- 23. Doherty, L.B., Rohr, F.J., Levy, H.L., Detection of Phenylketonuria in the Very Early Newborn Blood Specimen . *Pediatrics* 87(2):240, 1991.
- 24. Adam, B.W., et. al., Impact of differences in the filter paper matrix on outcomes of PKU testing. In Proceedings: 10th National Neonatal Screening Symposium, 0-10, pp. 29, 1994.
- 25. American Academy of Pediatrics, Committee on Genetics: Issues in newborn screening. *Pediatrics* 89(2):345,1992.