

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA**Introduction****Intended Use**

The **DRG β -HCG ELISA Kit** is an enzyme immunoassay for the measurement of total human chorionic gonadotropin (hCG and β -hCG) in serum.

Summary and Explanation

Chorionic Gonadotropin (hCG) is a glycoprotein hormone which is normally produced by the placenta during pregnancy. After conception, the hCG concentration increases rapidly to reach a peak near the end of the first trimester. High concentrations are observed throughout pregnancy. After delivery, hCG levels fall rapidly and become undetectable after a few days.

Structurally intact hCG molecules are composed of an alpha and a beta subunit. The alpha subunit is nearly identical to the alpha subunits of other glycoprotein hormones, such as Thyroid Stimulating Hormone (TSH), Luteinizing Hormone (LH), and Follicle Stimulating Hormone (FSH): The differences in the beta subunit of the respective hormones account for their biological specificity and immunochemical distinctiveness.

Monoclonal antibodies recognizing unique sites on the beta chain of the β -hCG/hCG molecule are essential for differentiation between hCG and LH, FSH and TSH.

Specific assays for beta-hCG permit the early detection of pregnancy.

In addition to the elevated hCG levels during pregnancy, high concentrations of β hCG/hCG may be associated with neoplasms of trophoblastic and nontrophoblastic origin such as hydatiform mole, chorionepithelioma, embryonal cell carcinoma, seminoma and many others.

PRINCIPLE of the test

The DRG β -HCG ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle.

The microtiter wells are coated with a monoclonal antibody [mouse] directed towards a unique antigenic site on a β -HCG molecule. An aliquot of Donor sample containing endogenous β -HCG and/or HCG is incubated in the coated well with enzyme conjugate, which is an anti- β -HCG antibody conjugated with horseradish peroxidase. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase is proportional to the concentration of β -HCG/HCG in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of β -HCG/HCG in the Donor sample.

Warnings and Precautions

1. This kit is for in vitro use only. For professional use only. In the United States, this kit is intended for Research Use Only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA

4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the Donor samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Material Safety Data Sheets. Material Safety Data Sheets for this product are available upon request directly from DRG.

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA

Reagents

Reagents provided

1. **Microtiterwells**, 12x8 (break apart) strips, 96 wells;
Wells coated with anti- β -HCG antibody (monoclonal).
2. **Standard (Standard 0-5)**, 6 vials (lyophilized), 1 mL;
Concentrations: 0, 5; 25; 50; 100; 200 mIU/mL
Conversion: 1pg/mL = 0,00916 mIU/mL
The standards are calibrated against WHO approved Reference material IRR β -HCG , code 75/537
See „Preparation of Reagents“;
contain 0.3% Proclin as a preservative.
3. **Sample Diluent**, 1 vial, 10 mL, ready to use,
contains 0.3% Proclin as a preservative.
4. **Enzyme Conjugate**, 1 vial, 11 mL, ready to use,
Anti- β -HCG antibody conjugated to horseradish peroxidase;
contains 0.3% Proclin as a preservative.
5. **Substrate Solution**, 1 vial, 14 mL, ready to use,
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.

Note: Additional *Sample Diluent* for sample dilution is available upon request.

Materials required but not provided

- A microtiter plate calibrated reader (450±10 nm) (e.g. the DRG Instruments Microtiter Plate Reader).
- Calibrated variable precision micropipettes.
- Distilled or deionized water
- Timer
- Semi logarithmic graph paper or software for data reduction

Storage Conditions

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again.

Reagent Preparation

Allow all reagents and required number of strips to reach room temperature prior to use.

Standards

Reconstitute the lyophilized contents of the standard vial with 1.0 mL Aqua dest.

Note: *The reconstituted standards are stable for 2 months at 2-8°C. For longer storage freeze at -20°C.*

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA

Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets (see chapter 13).

Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

SPECIMEN Collection and Preparation

Serum should be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette # 02.1388.001), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Donors receiving anticoagulant therapy may require increased clotting time.

Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 5 days at 2-8°C prior to assaying.

Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.

Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 μ L Serum + 90 μ L *Sample Diluent* (mix thoroughly)
- b) dilution 1:100: 10 μ L dilution a) 1:10 + 90 μ L *Sample Diluent* (mix thoroughly).

Note:

Samples with expected values greater 200 mIU/ml should be diluted with *Sample Diluent* before assaying.

Assay procedure

General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA

- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

Test Procedure (quantitative method)

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the holder.
2. Dispense **25 μ L** of each *Standard*, *Control* and samples with new disposable tips into appropriate wells.
3. Dispense **100 μ L** *Enzyme Conjugate* into each well.
4. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature.
6. Briskly shake out the contents of the wells.
Rinse the wells 5 times with distilled water (400 μ L per well). Strike the wells sharply on absorbent paper to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

7. Add **100 μ L** of *Substrate Solution* to each well.
8. Incubate for **15 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **50 μ L** of *Stop Solution* to each well.
10. Determine the absorbance (OD) of each well at **450 \pm 10 nm** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

Calculation of Results (quantitative)

1. Calculate the average absorbance values for each set of standards, controls and Donor samples.
2. Construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 200 mIU/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

Revised 28 July 2010 rm (Vers. 6.1)

RUO in the USA

EXAMPLE OF TYPICAL STANDARD CURVE

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 mIU/mL)	0.04
Standard 1 (5 mIU/mL)	0.15
Standard 2 (25 mIU/mL)	0.28
Standard 3 (50 mIU/mL)	0.53
Standard 4 (100 mIU/mL)	0.94
Standard 5 (200 mIU/mL)	1.50

Assay Procedure (qualitative method)

This procedure differentiates positive (pregnant) from negative samples by comparing the sample beta hCG levels with the *Standard 0* (0 mIU/mL) and *Standard 3* (50 mIU/mL).

Donor samples are run with the Standard 0 and the 50 mIU/mL Standard. The assay procedure is identical with the Quantitative Method, but step 9 and 10 is omitted.

Calculation of Results (qualitative)

For a qualitative analysis of the β -hCG level the color development of the specimen is compared with the color of the 0 mIU/mL and 50 mIU/mL standards.

If the blue color is less intense than the color of the 50 mIU/mL standard, the sample is considered as negative.

If the blue color is more intense than or equal to the color of the 50 mIU/mL standard the sample is considered as positive.

CAUTION:

1. For the detection of pregnancy in serum, a qualitative assay is used with a cut-off point of 50 mIU/mL. Negative or borderline results should be repeated on a fresh specimen obtained at least 48 hours after the first specimen.
2. It has been shown that immunological pregnancy tests may yield false results in cases of several medications and diseases such as rheumatoid arthritis. In such cases, the interpretation of the pregnancy test should be done carefully.

REFERENCES

1. Sancken, U., Bartels, I., Der sogenannte Triple-Test. Reproduktionsmedizin 15:276-284, 1999
2. Cole L.A. Clinical Chemistry 43:12, 2233-2243, 1997
3. Iles R.K., Javid M.K., Gunn L.K., Chard T. Clinical Chemistry 45:4, 532-538, 1999
4. Iles R.K., Javid M.K., Gunn L.K., Chard T. Clinical Chemistry 40:3, 484-485, 1994
5. Korhonen J., Alfthan H., Ylöstalo P., Veldhuis J., Stenman U.-H. Clinical Chemistry 43:11, 2155-2163, 1997