

**DRG<sup>®</sup> Ultrasensitive Insulin ELISA (EIA-2337)**

Revised 18 Nov. 2010 rm (Vers. 8.1)

*Please use only the valid version of the package insert provided with the kit.***1 INTENDED USE**

DRG Ultrasensitive Insulin ELISA provides a method for the quantitative determination of insulin in human serum or plasma.

**2 SUMMARY AND EXPLANATION OF THE TEST**

Insulin is the principal hormone responsible for the control of glucose metabolism. It is synthesized in the  $\beta$ -cells of the islets of Langerhans as the precursor, proinsulin, which is processed to form C-peptide and insulin. Both are secreted in equimolar amounts into the portal circulation. The mature insulin molecule comprises two polypeptide chains, the A chain and B chain (21 and 30 amino acids respectively). The two chains are linked together by two inter-chain disulphide bridges. There is also an intra-chain disulphide bridge in the A chain.

Secretion of insulin is mainly controlled by plasma glucose concentration, and the hormone has a number of important metabolic actions. Its principal function is to control the uptake and utilisation of glucose in peripheral tissues via the glucose transporter. This and other hypoglycaemic activities, such as the inhibition of hepatic gluconeogenesis and glycogenolysis are counteracted by the hyperglycaemic hormones including glucagon, epinephrine (adrenaline), growth hormone and cortisol.

Insulin concentrations are severely reduced in insulin-dependent diabetes mellitus (IDDM) and some other conditions such as hypopituitarism. Insulin levels are raised in non-insulin-dependent diabetes mellitus (NIDDM), obesity, insulinoma and some endocrine dysfunctions such as Cushing's syndrome and acromegaly.

**3 PRINCIPLE OF THE PROCEDURE**

DRG Ultrasensitive Insulin ELISA is a solid phase two-site enzyme immunoassay. It is based on the direct sandwich technique in which two monoclonal antibodies are directed against separate antigenic determinants on the insulin molecule. During incubation insulin in the sample reacts with peroxidase-conjugated anti-insulin antibodies and anti-insulin antibodies bound to microplate. A simple washing step removes unbound enzyme labelled antibody. The bound conjugate is detected by reaction with 3,3',5,5'-tetramethylbenzidine. The reaction is stopped by adding acid to give a colorimetric endpoint that is read spectrophotometrically.

**4 WARNINGS AND PRECAUTIONS**

- For in vitro diagnostic use.
- The contents of this kit and their residues must not be allowed to come into contact with ruminating animals or swine.
- The Stop solution in this kit contains 0.5 M H<sub>2</sub>SO<sub>4</sub>. Follow routine precautions for handling hazardous chemicals.
- All patient specimens should be handled as of capable of transmitting infections.

**5 MATERIAL REQUIRED BUT NOT PROVIDED**

- Pipettes for 25, 50, 200 and 1000  $\mu$ L (repeat pipettes preferred for addition of Enzyme Conjugate 1X solution, Substrate TMB and Stop Solution)

## DRG® Ultrasensitive Insulin ELISA (EIA-2337)



Revised 18 Nov. 2010 rm (Vers. 8.1)



- Beakers and cylinders for reagent preparation
- Redistilled water
- Micro plate reader with 450 nm filter
- Plate shaker (The recommended velocity is 700–900 cycles per minute, orbital movement).
- Microplate washing device

### 6 REAGENTS

Each Ultrasensitive Insulin ELISA kit contains reagents for 96 wells, sufficient for 42 samples and one calibration curve in duplicate.

For larger series of assays, use pooled reagents from packages bearing identical Lot numbers.

The expiry date for the complete kit is stated on the outer label.

The recommended storage temperature is 2–8°C.

<b>Coated Plate</b>	1 plate	96 wells	Ready for use
(mouse monoclonal anti-insulin)		8-wells strips	

For unused microplate strips, reseal the bag using adhesive tape, store at 2–8°C and use within 8 weeks.

<b>Calibrators 1, 2, 3, 4, 5</b>	5 vials	1000 µL	Ready for use
----------------------------------	---------	---------	---------------

Color coded yellow

Concentration indicated on vial label

<b>Calibrator 0</b>	1 vial	5 mL	Ready for use
---------------------	--------	------	---------------

Colour coded yellow

<b>Enzyme Conjugate 11X</b>	1 vial	1.2 mL	Preparation see below,
-----------------------------	--------	--------	------------------------

Peroxidase conjugated mouse monoclonal anti-insulin

<b>Enzyme Conjugate Buffer</b>	1 vial	12 mL	Ready for use
--------------------------------	--------	-------	---------------

Colour coded blue

<b>Wash Buffer 21X</b>	1 bottle	50 mL	
------------------------	----------	-------	--

Dilute with 1000 mL redistilled water to make Wash Buffer 1X solution

Storage after dilution: 2–8°C for 8 weeks.

<b>Substrate TMB</b>	1 vial	22 mL	Ready for use
----------------------	--------	-------	---------------

Colourless solution

*Note! Light sensitive!*

<b>Stop Solution</b>	1 vial	7 mL	Ready for use
----------------------	--------	------	---------------

0.5 M H<sub>2</sub>SO<sub>4</sub>

## DRG<sup>®</sup> Ultrasensitive Insulin ELISA (EIA-2337)



Revised 18 Nov. 2010 rm (Vers. 8.1)



### 6.1 Preparation of Enzyme Conjugate 1X solution

Prepare the needed volume of enzyme conjugate 1X solution by dilution of Enzyme Conjugate 11X, (1+10) in Enzyme Conjugate Buffer or according to the table below. When preparing enzyme conjugate 1X solution for the whole plate, pour all of the Enzyme Conjugate Buffer into the Enzyme Conjugate 11X vial.

Mix gently. Use within one day.

Number of strips	Enzyme Conjugate 11X	Enzyme Conjugate Buffer
12 strips	1 vial	1 vial
8 strips	700 µL	7.0 mL
6 strips	500 µL	5.0 mL
4 strips	350 µL	3.5 mL

## 7 SPECIMEN COLLECTION AND HANDLING

### Serum

Collect blood by venipuncture, allow to clot, and separate the serum by centrifugation.

Samples can be stored at 2–8°C up to 24 hours. For longer periods, store samples at –20°C. Avoid repeated freezing and thawing.

### Plasma

Collect blood by venipuncture into tubes containing heparin or EDTA as anticoagulant, and separate the plasma fraction.

Samples can be stored at 2–8°C up to 24 hours. For longer periods store samples at -20°C. Avoid repeated freezing and thawing.

### Preparation of samples

No dilution is normally required, however, samples containing >20mU/L should be diluted e.g. 1/10v/v with Calibrator 0.

## 8 TEST PROCEDURE

All reagents and samples must be brought to room temperature before use.

Prepare a calibration curve for each assay run.

1. Prepare enzyme conjugate 1X solution and wash buffer 1X solution.
2. Prepare sufficient microplate wells to accommodate Calibrators and samples in duplicate.
3. Pipette 25 µL each of Calibrators and samples into appropriate wells.
4. Add 100 µL of enzyme conjugate 1X solution to each well.
5. Incubate on a plate shaker (700-900 rpm) for 1 hour at room temperature (18–25°C)

**DRG<sup>®</sup> Ultrasensitive Insulin ELISA (EIA-2337)****Revised 18 Nov. 2010 rm (Vers. 8.1)**

6. Wash 6 times with 700  $\mu$ L per well using an automatic plate washer with overflow-wash function. Do not include soak step in washing procedure.  
Or manually,  
Discard the reaction volume by inverting the microplate over a sink. Add 350  $\mu$ L wash solution to each well. Discard the wash solution, tap firmly several times against absorbent paper to remove excess liquid. Repeat 5 times. Avoid prolonged soaking during washing procedure.
7. Add 200  $\mu$ L Substrate TMB into each well
8. Incubate for 30 minutes at room temperature (18–25°C)
9. Add 50  $\mu$ L Stop Solution to each well.  
Place plate on a shaker for approximately 5 seconds to ensure mixing.
10. Read optical density at 450 nm and calculate results. Read within 30 minutes

*Note! To prevent contamination between the conjugate and substrate, separate pipettes are recommended.*

**9 INTERNAL QUALITY CONTROL**

Commercial controls such as DRG Diabetes antigen Control (EIA-2338) and/or internal serum pools with low, intermediate and high insulin concentrations should routinely be assayed as samples, and results charted from day to day. It is good laboratory practice to record the following data for each assay: kit lot number; dilution and/or reconstitution dates of kit components, OD values for the blank, Calibrators and controls..



**Revised 18 Nov. 2010 rm (Vers. 8.1)**



## **10 CALCULATION OF RESULTS**

### **Computerized calculation**

The concentration of insulin is obtained by computerized data reduction of the absorbance for the Calibrators, except for Calibrator 0, versus the concentration using cubic spline regression.

### **Manual calculation**

1. Plot the absorbance values obtained for the Calibrators, except Calibrator 0, against the insulin concentration on a lin-log paper and construct a calibration curve.
2. Read the concentration of the samples from the calibration curve.

### **Example of results**

Mean Wells	Identity	A <sub>450</sub>	conc. mU/l
1A–B	Calibrator 0	0.071/0.073	
1C–D	Calibrator 1*	0.087/0.088	
1E–F	Calibrator 2*	0.186/0.188	
1G–H	Calibrator 3*	0.405/0.424	
2A–B	Calibrator 4*	1.172/1.241	
2C–D	Calibrator 5*	2.262/2.290	
2E–F	Unknown 1	0.466/0.470	3.5
2G–H	Unknown 2	1.987/1.926	17.0
3A–B	Unknown 3	0.327/0.317	2.2

\* Concentration indicated on vial label.

## **11 LIMITATIONS OF THE PROCEDURE**

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical findings have been evaluated.

Application of this test to individuals already undergoing insulin therapy is complicated by formation of anti-insulin antibodies that are capable of interfering in the assay.

Grossly lipemic, icteric or hemolysed samples do not interfere in the assay.

## **12 EXPECTED VALUES**

Good practice dictates that each laboratory establishes its own expected range of values. The following results may serve as a guide until the laboratory has gathered sufficient data of its own.

Fasting levels for 137 tested, apparently healthy individuals, yielded a mean of 9.2 mU/L.

a median of 6.9 mU/L and a range- corresponding to the central 95% of the observations – of 2–25 mU/L.



**Revised 18 Nov. 2010 rm (Vers. 8.1)**



## **13 PERFORMANCE CHARACTERISTICS**

### **13.1 Detection limit**

The detection limit is 0.07 mU/L calculated as two standard deviations above the Calibrator 0

### **13.2 Recovery**

Recovery upon addition is 97–107% (mean 102%)

### **13.3 Hook effect**

Samples with a concentration of up to at least 800 000 mU/L can be measured without giving falsely low results.

### **13.4 Precision**

Each sample was analysed in 6 replicates on six different occasions.

Sample	Mean value mU/L	Coefficient of variation		
		Within assay %	Between assay %	Total assay %
1	5.7	5.3	2.7	6.0
2	8.2	4.2	3.9	5.8
3	15.5	5.1	1.8	5.4

### **13.5 Specificity**

The following cross reactions have been found:

C-peptide	< 0.01%
Proinsulin	< 0.01%
Proinsulin des (31-32)	< 0.5%
Proinsulin split (32-33)	< 0.5%
Proinsulin des (64-65)	98%
Proinsulin split (65-66)	56%
Insulin lispro (Humalog®, Eli Lilly)	< 0.006%
Insulin aspart	< 0.006%
IGF-I	< 0.02%
IGF-II	< 0.02%
Rat insulin	0.7%
Mouse insulin	0.3%
Porcine insulin	374%
Sheep insulin	48%
Bovine insulin	31%

## **14 CALIBRATION**

DRG Ultrasensitive Insulin ELISA kit is calibrated against 1<sup>st</sup> International Reference Preparation 66/304 for human insulin.

**DRG<sup>®</sup> Ultrasensitive Insulin ELISA (EIA-2337)****Revised 18 Nov. 2010 rm (Vers. 8.1)****15 CONVERSION FACTOR**

1 µg/L = 29 mU/L; 1 mU/L = 6.0 pmol/L

**16 WARRANTY**

The performance data presented here was obtained using the procedure indicated. Any change or modification in the procedure not recommended by DRG may affect the results, in which event DRG disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use.

DRG and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

**17 REFERENCES/LITERATURE / LITERATUR**

1. Craig ME, Howard NJ, Silink M, Rawlinson WD. (2003) Reduced frequency of HLA DRB1\*03-DQB1\*02 in children with type 1 diabetes associated with enterovirus RNA. J Infect Dis 187:1562-1570
2. Gaines-Das, R.E. and Bristow, A.F. (1988) WHO International reference reagents for human proinsulin and human C-peptide. J Biol Stand 16: 179-186
3. Herkner H, Klein N, Joukhadar C, Lackner E, Langenberger H, Frossard M, Bieglmayer C, Wagner O, Roden M, Muller M. (2003) Transcapillary insulin transfer in human skeletal muscle. Eur J Clin Invest 33:141-146
4. Sjostrand M, Gudbjornsdottir S, Holmang A, Lonn L, Strindberg L, Lonnroth P. (2002) Delayed transcapillary transport of insulin to muscle interstitial fluid in obese subjects. Diabetes 51:2742-2748
5. Thule PM, Campbell AG, Kleinhenz DJ, Olson DE, Boutwell JJ, Sutliff RL, Hart CM. (2006) Hepatic Insulin Gene Therapy Prevents Deterioration of Vascular Function and Improves Adipocytokine Profile in STZ-Diabetic Rats. Am J Physiol Endocrinol Metab 290:144-122
6. Zhou J, Pineyro MA, Wang X, Doyle ME, Egan JM. (2002) Exendin-4 differentiation of a human pancreatic duct cell line into endocrine cells: involvement of PDX-1 and HNF3beta transcription factors. J Cell Physiol 192:304-314