

---

*This kit is intended for Research Use Only.*

*This kit is not intended for diagnostic purposes.*

### **Principle of the Test**

The fT4 test is a solid phase competitive enzyme immunoassay. Patient serum samples, standards, and Thyroxine-Enzyme Conjugate Working Reagent are added to wells coated with monoclonal T4 antibody. fT4 in the patient specimen and the T4 labeled conjugate compete for available binding sites on the antibody. After a 60 minutes incubation at room temperature, the wells are washed with water to remove unbound T4 conjugate. A solution of H<sub>2</sub>O<sub>2</sub>/TMB is then added and incubated for 20 minutes, resulting in the development of blue color. The color development is stopped with the addition of 3N HCl, and the absorbance is measured spectrophotometrically at 450 nm. The intensity of the color formed is proportional to the amount of enzyme present and is inversely related to the amount of unlabeled fT4 in the sample. By reference to a series of fT4 standards assayed in the same way, the concentration of fT4 in the unknown sample is quantified.

### **Reagents**

#### **Materials provided with the kit**

- Antibody **Coated Microplate** (1 plate, 96 wells)  
Microtiter wells coated with Anti-T4
- fT4-**Enzyme Conjugate** Reagent, ready to use (1 vial, 10.5 ml)
- Free T4 Reference **Standards** 0, 0.3, 0.95, 2.1, 3.6, and 7.0 ng/dl\*\*, 1 ml each  
\*\*Exact levels are given on the labels on a lot specific basis.
- **Color Reagent A** (1 bottle, 13 ml)
- **Color Reagent B** (1 bottle, 13 ml)
- **Stop Solution** (3N HCl) (1 bottle, 10 ml)

**Revised 9 Feb. 2011 rm (Vers. 15.1)****USA:RUO**

---

**Materials required but not provided**

- Pipette capable of delivering 50 µl volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.050 ml and 0.200 ml volumes with a precision of better than 1.5%.
- Microplate Reader with 450 nm wavelength absorbance capability.
- Test tubes for dilution of enzyme conjugate and for mixing Color Reagent A with Color Reagent B.
- Absorbent paper for blotting the microplate wells.
- Timer.
- Control Sera for quality control.

**Storage of Test Kit and Instrumentation**

Unopened test kits should be stored at 2°C – 8°C upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air.

Opened test kits will remain stable until the expiration date shown, provided it is stored as described above.

A microtiter plate reader with a bandwidth of 10 nm or less and an optical density range of 0-2 OD or greater at 450 nm wavelength is acceptable for use in absorbance measurement.

**Reagent Preparation****Working Substrate Solution – Prepare immediately before use**

To prepare H<sub>2</sub>O<sub>2</sub>/TMB solution, make a 1:1 mixing of Color Reagent A with Color Reagent B up to 1 hour before use. Mix gently to ensure complete mixing.

The prepared H<sub>2</sub>O<sub>2</sub>/TMB reagent should be made at least 15 minutes before use and is stable at room temperature in the dark for up to 3 hours. Discard excess after use.

**Assay Procedure**

Before proceeding with the assay, bring all reagents, serum references and controls to room temperature (18°C – 25°C).

1. Format the microplates' wells for each serum reference, control, and patient specimen to be assayed in duplicate.
2. Pipette 0.050 ml (50 µl) of the appropriate serum reference, control and specimen into the assigned well.
3. Add 0.100 ml (100 µl) of Free T4 Enzyme Conjugate Reagent to all wells.
4. Swirl the microplate gently for 20-30 seconds to mix.
5. Incubate 60 minutes at room temperature.
6. Remove the incubation mixture by emptying the plate content into a waste container. Rinse and empty the microtiter plate 5 times with distilled water. Strike the microtiter plate sharply onto absorbent paper or paper towels to remove all residual water droplets.
7. Add 0.200 ml (200 µl) of Working Substrate Solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.** Gently mix for 10 seconds.
8. Incubate at room temperature in the dark for 20 minutes.



Revised 9 Feb. 2011 rm (Vers. 15.1)

USA:RUO

9. Stop the reaction by adding 0.050 ml (50 µl) of 3N HCl (Stop Solution) to each well.
10. Gently mix for 30 seconds. **It is important to make sure that all the blue color changes to yellow color completely.**
11. Read absorbance at 450 nm with a microtiter well reader within 30 minutes.


### Calculation of Results

1. Calculate the mean absorbance value ( $A_{450}$ ) for each set of reference standards, controls and patient samples.
2. Construct a standard curve by plotting the mean absorbance obtained for each reference standard against its concentration in ng/dl on graph paper, with absorbance values on the vertical or Y axis, and concentrations on the horizontal or X axis.
3. Use the mean absorbance values for each specimen to determine the corresponding concentration of fT4 in ng/dl from the standard curve.

### Example of Standard Curve

Results of a typical standard run with optical density readings at 450 nm shown in the Y axis against fT4 concentrations shown in the X axis. This standard curve is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own data and standard curve in each experiment.

fT4 (ng/dl)	Absorbance (450 nm)
0	2.496
0.3	2.292
0.95	1.903
2.1	1.295
3.6	0.819
7.0	0.410



**References/Literature**

1. Tietz, N.W., Fundamentals of Clinical Chemistry, 2<sup>nd</sup> Ed., pg. 602, Saunders Press, Phila., 1976.
2. Horworth, P.J.N., Ward, R.L., J. Clin Pathol. 1972; 25:259-62.
3. Sati, C., Chatter, A.J., Watts, N. Fundamentals of Clinical Chemistry. Ed. Tietz, N.W. 3<sup>rd</sup> Ed., pg. 586. Saunders press Phila. 1987.
4. Lundberg, P.A., Jagenburg, R., Lindstedt, G., Nystrom, E., Clin. Chem. 1982, 28:1241.
5. Melmed, S., Geola, F.L., Reed, A.W., Pekary, A.E., Park, J., Hershmen, J. M., Clin Endocrin. Metabol. 1982, 54:300.
6. Ingbar, S.H., et al. J. Clin. Invest., 1965, 44:1679.
7. Selenkow, H.A., and Robin, N.I., J. Maine Med. Assoc. 1970, 61:199.
8. Oppenheimer, J.H., et al. J. Clin. Invest. 1962, 42: 1769.
9. Dick, M., Watson, F., Med J. Aust. 1980, 1:115.
10. Dussault, j. H., Turcotte, R., and Gieyda, H., Clin Endocrin. Metabol. 1976, 43: 232-285.
11. Tarnoky, A.L. Advan. Clin. Chem. 1981, 21:101-146.
12. Emrich, D., Schondube, H.,- Sehlen, S., and Schreivagel, I., Nuc. Compact, 1985, 16:392.
13. Procedures for Decontamination of Plumbing Systems Containing Copper and/or Lead azides, Dept. of H.E.W., N.I.O.S.H., Rockville, Maryland, 1976.

Version 04-02-03~rm

Version 01-13-11~rm