

**DRG® Leptin Receptor (human) (EIA-3527)**

Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA*Please use only the valid version of the package insert provided with the kit.**This kit is intended for Research Use Only. Not for use in diagnostic procedures.***INTENDED USE**

The Human Leptin Receptor ELISA is a sandwich enzyme immunoassay for measurement of human leptin receptor.

**Features**

- The total assay time is less than 2.5 hours.
- The kit measures leptin receptor in serum, plasma (EDTA, citrate, heparin) and tissue culture medium.
- Assay format is 96 wells.
- Quality Controls are human serum based.
- Standard is recombinant protein based.
- Components of the kit are provided ready to use, concentrated or lyophilized.

**STORAGE, EXPIRATION**

Store the complete kit at 2-8°C. Under these conditions, the kit is stable until the expiration date (see label on the box).

For stability of opened reagents see Chapter 9.

**TEST PRINCIPLE**

In this Human Leptin Receptor ELISA, standards, quality controls and samples are incubated in microplate wells pre-coated with monoclonal anti-human leptin receptor antibody. After 60 minutes incubation and washing, monoclonal anti-human leptin receptor antibody, conjugated with horseradish peroxidase (HRP) is added to the wells and incubated for 60 minutes with captured leptin receptor. Following another washing step, the remaining HRP conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured spectrophotometrically at 450 nm. The absorbance is proportional to the concentration of leptin receptor. A standard curve is constructed by plotting absorbance values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

**PRECAUTIONS**

- **For professional use only.**
- Wear gloves and laboratory coats when handling materials.

**DRG® Leptin Receptor (human) (EIA-3527)****Revised 3 Dec. 2010 rm (Vers. 4.1)****RUO in the USA**

- Do not drink, eat or smoke in the areas where immunodiagnostic materials are being handled.
- This kit contains components of human origin. These materials were found non-reactive for HBsAg, HCV antibody and for HIV 1/2 antigen and antibody. However, these materials should be handled as potentially infectious, as no test can guarantee the complete absence of infectious agents.
- This kit contains components of animal origin. However, these materials should be handled as potentially infectious.
- Avoid contact with the acidic Stop Solution and Substrate Solution, which contains hydrogen peroxide and tetramethylbenzidine (TMB). Wear gloves and eye and clothing protection when handling these reagents. Stop and/or Substrate Solutions may cause skin/eyes irritation. In case of contact with the Stop Solution and the Substrate Solution wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- The materials must not be pipetted by mouth.

**TECHNICAL HINTS**

- Reagents with different lot numbers should not be mixed.
- Use thoroughly clean glassware.
- Use deionized (distilled) water, stored in clean containers.
- Avoid any contamination among samples and reagents. For this purpose, disposable tips should be used for each sample and reagent.
- Substrate Solution should remain colourless until added to the plate. Keep Substrate Solution protected from light.
- Stop Solution should remain colourless until added to the plate. The colour developed in the wells will turn from blue to yellow immediately after the addition of the Stop Solution. Wells that are green in colour indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.
- Dispose of consumable materials and unused contents in accordance with applicable national regulatory requirements.

## DRG® Leptin Receptor (human) (EIA-3527)



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

### REAGENT SUPPLIED

Kit Components	State	Quantity
Antibody Coated Microtiter Strips	ready to use	96 wells
Conjugate Solution	ready to use	13 ml
Human Leptin Receptor Standards (2, 5, 10, 20, 50, 100 ng/ml)	liquid	6 x 0.35 ml
Quality Control High	lyophilized	1 vial
Quality Control Low	lyophilized	1 vial
Dilution Buffer	ready to use	13 ml
Wash Solution Concentrate (10x)	concentrated	100 ml
Substrate Solution	ready to use	13 ml
Stop Solution	ready to use	13 ml
Product Data Sheet + Certificate of Analysis		1 pc

### MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized (distilled) water
- Test tubes for diluting samples
- Glassware (graduated cylinder and bottle) for Wash Solution (Dilution Buffer)
- Precision pipettes to deliver 10-1000 µl with disposable tips
- Multichannel pipette to deliver 100 µl with disposable tips
- Absorbent material (e.g. paper towels) for blotting the microtiter plate after washing
- Vortex mixer
- Orbital microplate shaker capable of approximately 300 rpm
- Microplate washer (optional). [Manual washing is possible but not preferable.]
- Microplate reader with 450 ± 10 nm filter
- Software package facilitating data generation and analysis (optional)

### PREPARATION OF REAGENTS

All reagents need to be brought to room temperature prior to use.

Always prepare only the appropriate quantity of reagents for your test.

Do not use components after the expiration date marked on their label.

**DRG® Leptin Receptor (human) (EIA-3527)**

Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

Assay reagents supplied ready to use

**Antibody Coated Microtiter Strips**Stability and storage:

Return the unused strips to the provided aluminium zip-sealed bag with desiccant and seal carefully.  
Remaining Microtiter Strips are stable 3 month when stored at 2-8°C and protected from the moisture.

**Conjugate Solution****Dilution Buffer****Substrate Solution****Stop Solution**Stability and storage:

Opened reagents are stable 3 month when stored at 2-8°C.

Assay reagents supplied concentrated or lyophilized

**Human Leptin Receptor Standards**

Dilute each concentration of standard 3x with the Dilution Buffer prior to the assay,  
e.g. 50 µl of standard + 100 µl of Dilution Buffer for singlets,  
or preferably 100 µl of standard + 200 µl of Dilution Buffer for duplicates. Mix well (not to foam).

Stability and storage:

Opened standards are stable 3 month when stored at 2-8°C. **Do not store the diluted (3x) standards.**

**Quality Controls High, Low**

**Refer to the Certificate of Analysis for current Quality Control concentrations!!!**

**Reconstitute** each Quality Control (High and Low) with **350 µl of distilled (deionized) water** just prior to the assay. Let it dissolve at least 15 minutes with occasional gentle shaking (not to foam).

**Dilute reconstituted** Quality Controls **3x with Dilution Buffer**,

e.g. 50 µl of Quality Control + 100 µl of Dilution Buffer when assaying samples in singlets,  
or preferably 100 µl of Quality Control + 200 µl of Dilution Buffer for duplicates. **Mix well** (not to foam).

Stability and storage:

The reconstituted Quality Controls must be used immediately or aliquoted and frozen at -20°C for 3 month.  
Avoid repeated freeze/thaw cycles. **Do not store the diluted Quality Controls.**

**DRG® Leptin Receptor (human) (EIA-3527)****Revised 3 Dec. 2010 rm (Vers. 4.1)****RUO in the USA****Wash Solution Concentrate (10x)**

Dilute Wash Solution Concentrate (10x) ten fold in distilled water to prepare a 1x working solution.

Example: 100 ml of Wash Solution Concentrate (10x)+ 900 ml of distilled water for use of all 96-wells.

Stability and storage:

The diluted Wash Solution is stable 1 month when stored at 2-8°C. Opened Wash Solution Concentrate (10x) is stable 3 months when stored at 2-8°C.

**PREPARATION OF SAMPLES**

The kit measures leptin receptor in serum or plasma (EDTA, citrate, heparin) and tissue culture medium. Samples should be assayed immediately after collection or should be stored at -20°C. Mix thoroughly thawed samples just prior to the assay and avoid repeated freeze/thaw cycles, which may cause erroneous results. Avoid using hemolyzed or lipemic samples.

**Dilute samples 3x** with Dilution Buffer just prior to the assay

e.g. 50 µl of sample + 100 µl of Dilution Buffer for singlets,  
or preferably 100 µl of sample + 200 µl of Dilution Buffer for duplicates.

**Mix well** (not to foam). Vortex is recommended.

Stability and storage:

Samples should be stored at -20°, or preferably at -70°C for long-term storage. Avoid repeated freeze/ thaw cycles.

**Do not store the diluted samples.**

See Chapter 13 for stability of serum and plasma samples when stored at 2-8°C, effect of freezing/thawing and effect of sample matrix (serum/plasma) on the concentration of leptin receptor.

*Note: It is recommended to use a precision pipette and a careful technique to perform the dilution in order to get precise results.*

Ask for protocol at DRG if assaying other samples.

**ASSAY PROCEDURE**

1. Pipet **100 µl** of diluted Standards, Quality Controls, Dilution Buffer (=Blank) and samples, preferably in duplicates, into the appropriate wells. See Figure 1 for example of work sheet
2. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
3. Wash the wells **5-times** with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.

**DRG® Leptin Receptor (human) (EIA-3527)**

**Revised 3 Dec. 2010 rm (Vers. 4.1)**
**RUO in the USA**

4. Add **100 µl** of Conjugate Solution into each well.
5. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
6. Wash the wells **5-times** with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
7. Add **100 µl** of Substrate Solution into each well. Avoid exposing the microtiter plate to direct sunlight. Covering the plate with e.g. aluminium foil is recommended.
8. Incubate the plate for **10 minutes** at room temperature. The incubation time may be extended [up to 20 minutes] if the reaction temperature is below than 20°C. Do not shake the plate during the incubation.
9. Stop the colour development by adding **100 µl** of Stop Solution.
10. Determine the absorbance of each well using a microplate reader set to 450 nm, preferably with the reference wavelength set to 630 nm (acceptable range: 550 – 650 nm). Subtract readings at 630 nm (550 – 650 nm) from the reading at 450 nm.  
The absorbance should be read within 5 minutes following step 9.

**Note:**

*If the microplate reader is not capable of reading absorbance greater than the absorbance of the highest standard, perform a second reading at 405 nm. A new standard curve, constructed using the values measured at 405 nm, is used to determine leptin receptor concentration of off-scale samples. The readings at 405 nm should not replace the on-scale readings at 450 nm.*

**Note 2:**

*Manual washing: Aspirate wells and pipet 0.35 ml Wash Solution into each well. Aspirate wells and repeat twice. After final wash, invert and tap the plate strongly against paper towel. Make certain that Wash Solution has been removed entirely.*

	strip 1+2	strip 3+4	strip 5+6	strip 7+8	strip 9+10	strip 11+12
<b>A</b>	<b>Standard 100</b>	Blank	Sample 8	Sample 16	Sample 24	Sample 32
<b>B</b>	<b>Standard 50</b>	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
<b>C</b>	<b>Standard 20</b>	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
<b>D</b>	<b>Standard 10</b>	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
<b>E</b>	<b>Standard 5</b>	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
<b>F</b>	<b>Standard 2</b>	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
<b>G</b>	<b>QC High</b>	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
<b>H</b>	<b>QC Low</b>	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39

Figure 1: Example of a work sheet.

**DRG® Leptin Receptor (human) (EIA-3527)**



**Revised 3 Dec. 2010 rm (Vers. 4.1)**


**RUO** in the USA

**CALCULATIONS**

Most microplate readers perform automatic calculations of analyte concentration. The standard curve is constructed by plotting the mean absorbance at 450 nm (Y) of Standards against log of the known concentration (X) of Standards, using the four-parameter algorithm. Results are reported as concentration of leptin receptor (ng/ml) in samples.

Alternatively, the logit log function can be used to linearize the standard curve (i.e. logit of the mean absorbance (Y) is plotted against log of the known concentration (X) of Standards.

**Samples, Quality Controls and Standards are all diluted 3x prior to analysis, so there is not need to take this dilution factor into account.**



*Figure 2: Typical Standard Curve for Human Leptin Receptor ELISA.*

# DRG® Leptin Receptor (human) (EIA-3527)



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

## Stability of samples stored at 2-8°C

Samples should be stored at -20°C. However, no decline in concentration of leptin receptor was observed in serum and plasma samples after 10 days when stored at 2-8°C. To avoid microbial contamination, samples were treated with ε-aminocaproic acid and sodium azide, resulting in the final concentration of 0.03% and 0.1%, respectively.

Sample	Incubation Temp, Period	Serum (ng/ml)	Plasma (ng/ml)		
			EDTA	Citrate	Heparin
1	-20°C	49.75	41.48	38.82	35.69
	2-8°C, 1 day	47.13	43.03	41.38	41.19
	2-8°C, 10 days	45.04	44.28	46.02	40.49
2	-20°C	22.36	23.27	21.70	22.97
	2-8°C, 1 day	21.79	24.54	21.81	24.10
	2-8°C, 10 days	23.56	22.69	22.93	19.35
3	-20°C	33.28	33.34	35.09	34.46
	2-8°C, 1 day	35.80	33.49	30.82	32.29
	2-8°C, 10 days	35.15	33.96	31.77	35.45

## Effect of Freezing/Thawing

No decline was observed in concentration of human leptin receptor in serum and plasma samples after repeated (5x) freeze/thaw cycles. However it is recommended to avoid unnecessary repeated freezing/thawing of the samples.

Sample	Number of f/t cycles	Serum (ng/ml)	Plasma (ng/ml)		
			EDTA	Citrate	Heparin
1	1x	21.39	20.89	16.77	21.79
	3x	18.69	18.31	16.20	22.72
	5x	20.06	20.47	16.43	21.15
2	1x	26.89	24.71	23.33	27.14
	3x	26.77	26.91	22.29	26.58
	5x	25.24	24.80	20.93	24.36
3	1x	18.07	18.51	16.00	18.96
	3x	17.57	17.95	17.19	19.73
	5x	19.39	18.55	17.29	19.62



**DRG® Leptin Receptor (human) (EIA-3527)****Revised 3 Dec. 2010 rm (Vers. 4.1)****RUO in the USA****Reference ranges**

It is recommended that each laboratory include its own panel of control sample in the assay. Each laboratory should establish its own normal and pathological reference ranges for leptin receptor levels with the assay.

**DEFINITION OF THE CALIBRATOR**

The Standard used in this kit is recombinant human IgG-Fc-fragment - human OB-R dimeric chimera, which is different from the native soluble OB-R that is measured in human serum. Mature OB-R is a disulfide-linked homodimeric protein. As a result of glycosylation, recombinant human leptin receptor/Fc chimera migrates as a 205 kDa protein in SDS-PAGE and we used to employ the unit U/ml.

From the lot RD-877 we started using the unit ng/ml.

1 ng OB-R/ml = 1 U OB-R/ml (previously used).

It can be to recalculate previous results with factor 1.0; e.g. concentration of leptin receptor 24.50 U/ml in a sample measured in previous assays corresponds to 24.50 ng/ml of leptin receptor measured in this assay.

**METHOD COMPARISON**

The Human Leptin Receptor ELISA has not been compared to any other immunoassay.

**DRG® Leptin Receptor (human) (EIA-3527)**



**Revised 3 Dec. 2010 rm (Vers. 4.1)**

**RUO in the USA**

**TROUBLESHOOTING AND FAQs**

**Weak signal in all wells**

---

Possible explanations:

- Omission of a reagent or a step
- Improper preparation or storage of a reagent
- Assay performed before reagents were allowed to come to room temperature
- Improper wavelength when reading absorbance

**High signal and background in all wells**

---

Possible explanations:

- Improper or inadequate washing
- Overdeveloping; incubation time with Substrate Solution should be decreased before addition of Stop Solution
- Incubation temperature over 30°C

**High coefficient of variation (CV)**

---

Possible explanation:

- Improper or inadequate washing
- Improper mixing Standards, Quality Controls or samples

## DRG® Leptin Receptor (human) (EIA-3527)



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

### REFERENCES / Literature

- Lahlou N. et al.: Soluble leptin receptor in serum of subjects with complete resistance to leptin: Diabetes; 19:1347-1352 (2000)
- Auwerx J. and Staels B.: Leptin (Review article). The Lancet; 13:737 (1998)
- Chen H. et al.: Evidence that the diabetes gene encodes the leptin receptor: identification of a mutation in the leptin receptor gene in db/db mice. Cell; 84:491 495 (1996)
- Ciofii J. A. et al.: Novel B219/OB receptor isoforms: possible role of leptin in hematopoiesis and reproduction. Nature Med.; 2:585 589 (1996)
- Houseknecht KL and Portocarrero CP: Leptin and its receptors: Regulators of whole body energy homeostasis. Domestic animal endocrinology; 15 (6):457 475 (1998)
- Lee GH et al.: Abnormal splicing of the leptin receptor in diabetic mice. Nature; 379:632 635 (1996)
- Tartaglia LA et al.: Identification and Expresion Cloning of a Leptin Receptor, OB R. Cell; 83:1263 1271 (1995)
- Tartaglia LA: The leptin receptor. J. Biol. Chem.; 272:6093 6096 (1997)
- Haluzikova D, Dostalova I, Kavalkova P, Roubicek T, Mraz M, Papezova H, Haluzik M. Serum concentrations of adipocyte fatty acid binding protein in patients with anorexia nervosa. Physiol Res. 2009;58 (4):577-81
- Dostalova I, Kavalkova P, Haluzikova D, Lacinova Z, Mraz M, Papezova H, Haluzik M. Plasma concentrations of fibroblast growth factors 19 and 21 in patients with anorexia nervosa. J Clin Endocrinol Metab. (2008)
- Anderson PD, Mehta NN, Wolfe ML, Hinkle CC, Pruscino L, Comiskey LL, Tabita-Martinez J, Sellers KF, Rickels MR, Ahima RS, Reilly MP . Innate immunity modulates adipokines in humans. J Clin Endocrinol Metab . Jun;92(6):2272-9 (2007)
- Hoefle G, Saely CH, Risch L, Rein P, Koch L, Schmid F, Aczel S, Marte T, Langer P, Drexel H . Leptin, leptin soluble receptor and coronary atherosclerosis. Eur J Clin Invest . Aug;37(8):629-36 (2007)
- Proto C, Romualdi D, Cento RM, Romano C, Campagna G, Lanzone A . Free and total leptin serum levels and soluble leptin receptors levels in two models of genetic obesity: the Prader-Willi and the Down syndromes. Metabolism . Aug;56(8):1076-80 (2007)
- Martos-Moreno GA, Barrios V, Argente J . Normative data for adiponectin, resistin, interleukin 6, and leptin/receptor ratio in a healthy Spanish pediatric population: relationship with sex steroids. Eur J Endocrinol . Sep;155(3):429-34 (2006)
- Paul G, Schaffler A, Neumeier M, Furst A, Bataille F, Buechler C, Muller-Ladner U, Scholmerich J, Rogler G, Herfarth H . Profiling adipocytokine secretion from creeping fat in Crohn's disease. Inflamm Bowel Dis . Jun;12(6):471-7 (2006)

## DRG® Leptin Receptor (human) (EIA-3527)



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

- Presle N, Pottie P, Dumond H, Guillaume C, Lapique F, Pallu S, Mainard D, Netter P, Terlain B . Differential distribution of adipokines between serum and synovial fluid in patients with osteoarthritis. Contribution of joint tissues to their articular production. *Osteoarthritis Cartilage* . Jul;14(7):690-5 (2006)
- Stejskal D, Svoboda I . Stanovení koncentrace solubilního leptinového receptoru (Ob-Re) pomocí nové metody ELISA. *Chem Listy* . 809-812 (2006)
- Dostalova I, Kopsky V, Duskova J, Papezova H, Pacak K, Nedvidkova J . Leptin concentrations in the abdominal subcutaneous adipose tissue of patients with anorexia nervosa assessed by in vivo microdialysis. *Regul Pept* . May 15;128(1):63-8 (2005)
- Housova J, Wilczek H, Haluzik MM, Kremen J, Krizova J, Haluzik M . Adipocyte-derived hormones in heroin addicts: the influence of methadone maintenance treatment. *Physiol Res* . 54(1):73-8 (2005)
- Jarkovska Z, Hodkova M, Sazamova M, Rosicka M, Dusilova-Sulkova S, Marek J, Justova V, Lacinova Z, Haluzik M, Haas T, Krsek M . Plasma levels of active and total ghrelin in renal failure: a relationship with GH/IGF-I axis. *Growth Horm IGF Res* . Dec;15(6):369-76 (2005)
- Papadopoulos DP, Makris TK, Krespi PG, Poulakou M, Paizis IA, Hatzizacharias AN, Perrea D, Voteras VV . Human soluble leptin receptor number in healthy normotensive individuals with high normal blood pressure. *Am J Hypertens* . Jul;18(7):1001-4 (2005)
- Popruk S, Tungtrongchitr R, Pongpaew P, Phonrat B, Tungtongchitr A, Tribunyatkul S, Paksanont S, Vudhivai N, Schelp FP . Relationship between soluble leptin receptor, leptin, lipid profiles and anthropometric parameters in overweight and obese Thai subjects. *J Med Assoc Thai* . Feb;88(2):220-7 (2005)
- Szalay F, Folhoffer A, Horvath A, Csak T, Speer G, Nagy Z, Lakatos P, Horvath C, Habior A, Tornai I, Lakatos PL . Serum leptin, soluble leptin receptor, free leptin index and bone mineral density in patients with primary biliary cirrhosis. *Eur J Gastroenterol Hepatol* . Sep;17(9):923-8 (2005)
- Verhaeghe J, van Bree R, Lambin S, Caluwaerts S . Adipokine profile and C-reactive protein in pregnancy: effects of glucose challenge response versus body mass index. *J Soc Gynecol Investig* . Jul;12(5):330-4 (2005)
- Zemanova M, Petruzela L, Pazdrova G, Haluzik M, Novak F, Svobodnik A . Changes in leptin, leptin receptor and fatty acids levels during chemoradiotherapy for oesophageal cancer. *Cas Lek Cesk* . 144(12):811-6 (2005)
- Tsiotra PC, Pappa V, Koukourava A, Economopoulos T, Tsigos C, Raptis SA . Expression of leptin receptors in mononuclear cells from myelodysplastic syndromes and acute myeloid leukemias. *Acta Haematol* . 114(2):71-7 (2005)
- Housova J, Anderlova K, Krizova J, Haluzikova D, Kremen J, Kumstyrova T, Papezova H, Haluzik M . Serum adiponectin and resistin concentrations in patients with restrictive and binge/purge form of anorexia nervosa and bulimia nervosa. *J Clin Endocrinol Metab* . Dec 14 (2004)

**DRG® Leptin Receptor (human) (EIA-3527)****Revised 3 Dec. 2010 rm (Vers. 4.1)****RUO in the USA**

- Chatzantoni K, Papathanassopoulos P, Gourzoulidou E, Mouzaki A . Leptin and its soluble receptor in plasma of patients suffering from remitting-relapsing multiple sclerosis (MS) In vitro effects of leptin on type-1 and type-2 cytokine secretion by peripheral blood mononuclear cells, T-cells and monocytes of MS patients. *J Autoimmun* . Sep;23(2):169-77 (2004)
- Krsek M, Silha JV, Jezkova J, Hana V, Marek J, Weiss V, Stepan JJ, Murphy LJ . Adipokine levels in Cushing's syndrome; elevated resistin levels in female patients with Cushing's syndrome. *Clin Endocrinol (Oxf)* . Mar;60(3):350-7 (2004)
- Meier U, Gressner AM . Endocrine regulation of energy metabolism: review of pathobiochemical and clinical chemical aspects of leptin, ghrelin, adiponectin, and resistin. *Clin Chem* . Sep;50(9):1511-25 (2004)
- Nuamah MA, Yura SA, Sagawa NA, Itoh H, Miese H, Korita D, Kakui K, Takemura M, Ogawa Y, Nakao K, Fujii S . Significant increase in maternal plasma leptin concentration in induced delivery: a possible contribution of pro-inflammatory cytokines to placental leptin secretion. *Endocr J* . Apr;51(2):177-87 (2004)
- Wolfe BE, Jimerson DC, Orlova C, Mantzoros CS . Effect of dieting on plasma leptin, soluble leptin receptor, adiponectin and resistin levels in healthy volunteers. *Clin Endocrinol (Oxf)* . Sep;61(3):332-8 (2004)
- Buyan N, Ozkaya O, Bideci A, Soylemezoglu O, Cinaz P, Gonen S, Kalman S, Bakkaloglu S, Hasanoglu E. Leptin, soluble leptin receptor, and transforming growth factor-beta1 levels in minimal change nephrotic syndrome. *Pediatr Nephrol*. Jul;29(8) (2003)
- Gavrilu A, Peng CK, Chan JL, Mietus JE, Goldberger AL, Mantzoros CS. Diurnal and ultradian dynamics of serum adiponectin in healthy men: comparison with leptin, circulating soluble leptin receptor, and cortisol patterns. *J Clin Endocrinol Metab* . Jun;88(6):2838-43 (2003)
- Kalousova M, Sulkova S, Zima T, Deppisch R, Beck W, Bednarova V, Fortova M, Tesar V . Advanced glycation end products in hemodialyzed patients with diabetes mellitus correlate with leptin and leptin/body fat ratio. *Ren Fail* . Mar;25(2):277-86 (2003)
- Kado N, Kitawaki J, Koshiba H, Ishihara H, Kitaoka Y, Teramoto M, Honjo H . Relationships between the serum levels of soluble leptin receptor and free and bound leptin in non-pregnant women of reproductive age and women undergoing controlled ovarian hyperstimulation. *Hum Reprod* . Apr;18(4):715-20 (2003)
- Kalman S, Buyan N, Bayazit AK, Bideci A, Bakkaloglu S, ozkaya O, Soylemezoglu O. The role of leptin and soluble leptin receptor on the nutritional parameters in pediatric CPD patients. *Nephrology Dialysis Transplantation* . Volume 18 Supplement 4 June 2003
- Krizova J, Sulkova S, Bednarova V, Parizkova J, Kotrlíkova E, Haluzik M . Soluble leptin receptor levels in patients with chronic renal failure. *Physiol Res* . 52(3):347-351 (2003)
- Lenchik L, Register TC, Hsu FC, Lohman K, Nicklas BJ, Freedman BI, Langefeld CD, Carr JJ, Bowden DW. Adiponectin as a novel determinant of bone mineral density and visceral fat. *Bone* . Oct;33(4):646-51 (2003)

# DRG® Leptin Receptor (human) (EIA-3527)



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

- Nuamah MA, Sagawa N, Yura S, Mise H, Itoh H, Ogawa Y, Nakao K, Fujii S . Free-to-total leptin ratio in maternal plasma is constant throughout human pregnancy. *Endocr J* . Aug;50(4):421-8 (2003) Welt CK, Schneyer AL, Heist K, Mantzoros CS . Leptin and soluble leptin receptor in follicular fluid. *J Assist Reprod Genet* . Dec;20(12):495-501 (2003)
- Rosicka M, Krsek M, Matoulek M, Jarkovska Z, Marek J, Justova V V, Lacinova Z . Serum ghrelin levels in obese patients: the relationship to serum leptin levels and soluble leptin receptors levels. *Physiol Res* . Feb;52(1):61-66 (2003)
- Sandhofer A, Laimer M, Ebenbichler CF, Kaser S, Paulweber B, Patsch JR. Soluble leptin receptor and soluble receptor-bound fraction of leptin in the metabolic syndrome. *Obes Res* . Jun;11(6):760-8 (2003)
- Stejskal D, Ruzicka V, Adamovska S, Jurak Mann DR, Johnson AO, Gimpel T, Castracane VD . Changes in circulating leptin, leptin receptor, and gonadal hormones from infancy until advanced age in humans. *J Clin Endocrinol Metab* Jul . Jul;88(7):3339-45 (2003)
- Teplan V, Schuck O, Stollova M, Vitko S . Obesity and hyperhomocysteinaemia after kidney transplantation. *Nephrol Dial Transplant* . Jul;18 Suppl 5:v71-3 (2003)
- Yiannakouris N, Melistas L, Yannakoulia M, Mungal K, Mantzoros CS . The -2548G/A polymorphism in the human leptin gene promoter region is associated with plasma free leptin levels; interaction with adiposity and gender in healthy subjects. *Hormones* . 2(4):229-236 (2003)
- Yannakoulia M, Yiannakouris N, Bluher S, Matalas AL, Klimis-Zacas D, Mantzoros CS. Body fat mass and macronutrient intake in relation to circulating soluble leptin receptor, free leptin index, adiponectin, and resistin concentrations in healthy humans. *J Clin Endocrinol Metab* . Apr;88(4):1730-6 (2003)
- Chan JL, Bluher S, Yiannakouris N, Suchard MA, Kratzsch J, Mantzoros CS . Regulation of Circulating Soluble Leptin Receptor Levels By Gender, Adiposity, Sex Steroids, and Leptin. *Observational and Interventional Studies in Humans. Diabetes* . 51(7):2105-2112 (2002)
- Haluzik M, Markova M, Slaby JJ, Jiskra J, Krizova J, Hass T . The changes of serum leptin and soluble leptin receptor levels in patients undergoing mobilization of peripheral blood stem cells before autologous stem cells transplantation. *Endocr Res* . Aug;28(3):189-97 (2002)
- Krsek M, Rosicka M, Haluzik M, Svobodova J, Kotrlíkova E, Justova V, Lacinova Z, Jarkovska Z . Plasma ghrelin levels in patients with short bowel syndrome. *Endocr Res* . Feb-May;28(1-2):27-33 (2002)
- Lahlou N, Issad T, Lebouc Y, Carel JC, Camoin L, Roger M, Girard J. Mutations in the human leptin and leptin receptor genes as models of serum leptin receptor regulation. *Diabetes* . 51(6):1980-5 (2002)
- Laimer M, Ebenbichler CF, Kaser S, Sandhofer A, Weiss A, Nehoda H, Aigner F, Patsch JR . Weight Loss Increases Soluble Leptin Receptor Levels and the Soluble Receptor Bound Fraction of Leptin. *Obesity Research* . 10(7):597-601 (2002)
- Manzella D, Parillo M, Razzino T, Gnasso P, Buonanno S, Gargiulo A, Caputi M, Paolisso G . Soluble leptin receptor and insulin resistance as determinant of sleep apnea. *Int J Obes Relat Metab Disord* . 26(3):370-5 (2002)

**DRG® Leptin Receptor (human) (EIA-3527)****Revised 3 Dec. 2010 rm (Vers. 4.1)****RUO in the USA**

- Monteleone P, Fabrazzo M, Tortorella A, Fuschino A, Maj M. Opposite modifications in circulating leptin and soluble leptin receptor across the eating disorder spectrum. *Mol Psychiatry* . 7(6):641-6 (2002)
- Ogier V, Ziegler O, Mejean L, Nicolas JP, Stricker-Krongrad A. Obesity is associated with decreasing levels of the circulating soluble leptin receptor in humans. *Int J Obesity* . 26(4):496-503 (2002)
- Pecoits-Filho R, Nordfors L, Heimbürger O, Lindholm B, Anderstam B, Marchlewska A, Stenvinkel P. Soluble leptin receptors and serum leptin in end-stage renal disease: relationship with inflammation and body composition. *Eur J Clin Invest Nov* . 32(11):811-7 (2002)
- Shimizu H, Shimomura K, Negishi M, Masunaga M, Uehara Y, Sato N, Shimomura Y, Kasai K, Mori M . Circulating concentrations of soluble leptin receptor: influence of menstrual cycle and diet therapy. *Nutrition* . 18(4):309-12 (2002)
- van Dielen FM, van 't Veer C, Buurman WA, Greve JW . Leptin and soluble leptin receptor levels in obese and weight-losing individuals. *J Clin Endocrinol Metab* . 87(4):1708-16 (2002)
- Malincikova J, Stejskal D, Hrebicek J . Serum leptin and leptin receptors in healthy prepubertal children: relations to insulin resistance and lipid parameters, body mass index (BMI), tumor necrosis factor alpha (TNF alpha), heart fatty acid binding protein (hFABP), and IgG anticardiolipin. *Acta Univ Palacki Olomuc Fac Med* . 143:51-7 (2000)
- Stejskal D, Bartek J, Ruzicka V, Stejskal P, Jedelsky L, Oral I, Horalik D, Zurek M . Determination of leptin receptor in the serum and relations to laboratory and anthropological parameters in patients with atherosclerotic complications. *Acta Univ Palacki Olomouc Fac Med* . 142:89-92 (1999)


**DRG® Leptin Receptor (human) (EIA-3527)**



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

**Assay Procedure Summary**





**DRG® Leptin Receptor (human) (EIA-3527)**



Revised 3 Dec. 2010 rm (Vers. 4.1)

**RUO** in the USA

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H