

DRG® Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA

INTENDED USE:


For Determination of Anti-H. Pylori Specific Antibodies of the IgG, IgA or IgM type in Human Serum or Plasma by Microplate Enzyme Immunoassay

PRINCIPLE

A Sequential ELISA Method

The reagents required for the sequential ELISA assay include immobilized antigen, circulating autoantibody and enzyme-linked species-specific antibody. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenous added biotinylated H. Pylori antigen.

Upon mixing biotinylated antigen, and a serum containing the antibody, reaction results between the antigen and the antibody to form an immune-complex. The interaction is illustrated by the followed equation:



k_a = Rate Constant of Association
 k_{-a} = Rate Constant of Disassociation
 $\text{B}^{\text{tn}}\text{Ag}_{(\text{H.Pylori})}$ = Biotinylated Antigen (Constant Quantity)
 $h\text{-Ab}_{(X\text{-H.Pylori})}$ = Human Auto-Antibody (Variable Quantity)
 $h\text{-Ab}_{(X\text{-H.Pylori})} - \text{B}^{\text{tn}}\text{Ag}_{(\text{H.Pylori})}$ = Immune Complex (Variable Quantity)

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antigen. This interaction is illustrated below:

$h\text{-Ab}_{(X\text{-H.Pylori})} - \text{B}^{\text{tn}}\text{Ag}_{(\text{H.Pylori})} + \text{Streptavidin}_{\text{C.W.}} \Rightarrow \text{immobilized complex (IC)}$
 $\text{Streptavidin}_{\text{C.W.}}$ = Streptavidin immobilized on well
 Immobilized complex = sandwich complex bound to the solid surface

After the incubation time, the well is washed to separate the unbound components by aspiration and/or decantation. The enzyme linked species-specific antibody (anti-h-IgG, M or A) is then added to the microwells. This conjugates binds to the immune complex that formed.

$\text{IC}_{(h\text{-IgG, M or A})} + \text{ENZ}\text{Ab}_{(X\text{-h-IgG, M or A})} \Rightarrow \text{ENZ}\text{Ab}_{(X\text{-h-IgG, M or A})} - \text{IC}_{(h\text{-IgG, M or A})}$
 $\text{IC}_{(h\text{-IgG, M or A})}$ = Immobilized Immune complex (Variable Quantity)
 $\text{ENZ}\text{Ab}_{(X\text{-h-IgG, M or A})}$ = Enzyme-antibody Conjugate (Constant Quantity)
 $\text{ENZ}\text{Ab}_{(X\text{-h-IgG, M or A})} - \text{I.C.}_{(h\text{-IgG, M or A})}$ = Ag-Ab Complex (Variable)

The anti-h-IgG, IgM or IgA enzyme conjugate that binds to the immune complex in a second incubation is separated from unreacted material by a wash step. The enzyme activity in this fraction is directly proportional to the antibody concentration in the specimen. By utilizing several different serum references of known antibody activity, a reference curve can be generated from which the antibody activity of an unknown can be ascertained.

DRG® Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA**REAGENTS****Materials provided:****A. Anti-H.Pylori Calibrators -- 1ml/vial**

Five (5) vials of references for anti-H.Pylori at levels of 0(A), 10(B), 25(C), 50(D), and 100(E) U/ml* of the IgG, IgM or IgA type. Store at 2-8°C. A preservative has been added. *Manufacturers' Reference Value

B. H.Pylori Biotin Conjugate -- 13ml/vial

One (1) vial of biotinylated inactivated H.Pylori in a buffering matrix. A preservative has been added. Store at 2-8°C.

C. Enzyme-antigen Conjugate -- 13ml/vial

One (1) vial of anti-human IgG, IgM or IgA-horseradish peroxidases (HRP) conjugate in a buffering matrix. A preservative has been added. Store at 2-8°C.

D. Streptavidin Coated Microplate -- 96 wells

One 96-well microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

E. Serum Diluent Concentrate - 20ml

One (1) vial of serum diluent containing buffer salts and a dye. Store at 2-8°C.

F. Wash Solution Concentrate -- 20ml

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-30°C

G. Substrate A --7ml/vial

One (1) bottle containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

H. Substrate B -- 7ml/vial

One (1) bottle containing hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.

I. Stop Solution -- 8ml/vial

One (1) bottle containing a strong acid (1N HCl). Store at 2-30°C.

J. Product Instructions

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: Opened reagents are stable for sixty (60) days when stored at 2-8°C.

Note 3: Above reagents are for a single 96-well microplate.

Required But Not Provided:

1. Pipette capable of delivering 10µl & 50µl volumes with a precision of better than 1.5%.
2. Dispenser(s) for repetitive deliveries of 0.100ml and 0.300ml volumes with a precision of better than 1.5%.
3. Microplate washers or a squeeze bottle (optional).
4. Microplate Reader with 450nm and 620nm wavelength absorbance capability.
5. Absorbent Paper for blotting the microplate wells.
6. Plastic wrap or microplate cover for incubation steps.
7. Vacuum aspirator (optional) for wash steps.
8. Timer.
9. Quality control materials.

PRECAUTIONS

This kit is intended for Research Use Only.

Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good

DRG® Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA

laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

SPECIMEN COLLECTION AND PREPARATION

Collect sample(s) by venipuncture in ten (10) ml silicone evacuated tube(s) or evacuated tube(s) containing EDTA or heparin. The usual precautions in the collection of venipuncture samples should be observed. Separate the red blood cells by centrifugation use serum or plasma for the anti-H.Pylori procedure. Specimen(s) may be refrigerated at 2-8°C for a maximum period of 48 hours. If the specimen(s) cannot be assayed within 48 hours, the sample(s) may be stored at temperatures of -20°C for up to 30 days.

REAGENT PREPARATION:

1. Serum Diluent

Dilute the serum diluent concentrate to 200ml in a suitable container with distilled or deionized water. Store at 2-8°C.

2. Wash Buffer

Dilute contents of Wash Concentrate to 1000ml with distilled or deionized water in a suitable storage container. Store at room temperature 20-27°C for up to 60 days.

3. Working Substrate Solution

Pour the contents of the vial labeled Solution 'A' into the vial labeled Solution 'B'. Mix and store at 2-8°C. Use within 60 days. Or for longer periods of usage determine the amount of reagent needed and prepare by mixing equal portions of Substrate A and Substrate B in a suitable container. For example, add 1ml of A and 1ml of B per two (2) eight well strips (A slight excess of solution is made. Discard the unused portion).

Note: Do not use the working substrate if it looks blue.

TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum references and controls to room temperature (20-27°C).

Donor Sample Dilution (1/100)

Dispense 0.010ml (10µl) of each donor specimen into 1ml of serum diluent. Cover and vortex or mix thoroughly by inversion. Store at 2-8°C for up to forty-eight (48) hours.

1. Format the microplates' wells for each serum reference, control and donor specimen to be assayed in duplicate.
Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.
2. Pipette 0.025 ml (25µl) of the appropriate serum reference, control or diluted donor specimen into the assigned well for IgG determination. **For IgM or IgA, pipette 0.050ml (50µl) of the appropriate serum reference, control or diluted donor specimen into the assigned well.**
3. Add 0.100 ml (100µl) of H.Pylori Biotinylated Conjugate Solution.
4. Swirl the microplate gently for 20-30 seconds to mix and cover.
5. Incubate 60 minutes at room temperature.
6. Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
7. Add 300µl of wash buffer (see Reagent Preparation Section), decant (blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**

DRG® Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA

8. Add 0.100 ml (100µl) of Enzyme anti-h-IgG, IgM or IgA Conjugate Solution to all wells. **Always add reagents in the same order to minimize reaction time differences between wells.**
9. Swirl the microplate gently, cover and incubate for thirty (30) minutes at room temperature.
10. Repeat steps (6 & 7) as explained above.
11. Add 0.100 ml (100µl) of Working Substrate Solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**
12. Incubate at room temperature for fifteen (15) minutes.
13. Add 0.050ml (50µl) of stop solution to each well and swirl the microplate gently for 15-20 seconds to mix. **Always add reagents in the same order to minimize reaction time differences between wells.**
14. Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

Note: For re-assaying specimens with concentrations greater than 100 U/ml, dilute the sample an additional 1:5 or 1:10 using the original diluted material in the serum diluent. Multiply by the dilution factor to obtain the concentration of the specimen.

RESULTS

A reference curve is used to ascertain the concentration of anti-*H. Pylori* in unknown specimens.

1. Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
2. Plot the absorbance for each duplicate serum reference versus the corresponding anti-*H. Pylori* activity in U/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
3. Draw the best-fit curve through the plotted points.
4. To determine the level of anti-*H. Pylori* activity for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in U/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance 1.603 (intersects the dose response curve) at 64.0 U/ml anti-H. Pylori concentration (See Figure 1). *


DRG[®] Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA

EXAMPLE 1 (Typical results for IgG, M or A)

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (U/ml)
Cal A	A1	0.042	0.044	0
	B1	0.046		
Cal B	C1	0.424	0.406	10
	D1	0.388		
Cal C	E1	0.810	0.791	25
	F1	0.772		
Cal D	G1	1.351	1.312	50
	H1	1.273		
Cal E	A2	2.377	2.328	100
	B2	2.279		
Donor 1	C2	0.163	0.172	5.2
	D2	0.182		
Donor 2	A3	1.534	1.603	64.0
	B3	1.671		



*The data presented in Example 1 and Figure 1 are for illustration only and **should not** be used in lieu of a standard curve prepared with each assay.

DRG® Anti-H. Pylori IgA (EIA-4099)

Revised 21 Sept. 2010 rm (Vers. 1.1)

RUO in the USA

REFERENCES

1. Warren, J.R., *Lancet*, **1**, 1273 (1983).
2. Marshall, B., *Lancet*, **1**, 1273 (1983).
3. Strickland, R. G., and Mackay, I.R., *Amer. J. Diag. Dis* , **18** 426 (1973).
4. Morris, A. G., and Nicholson, G., *Amer. J. Gastroenterol*, **82**, 192 (1987).
5. Sethi, P., et. al., *Post Grad. Med. J.*, **63**, 543 (1987).
6. Marshall, B.J., et. al., *Med. J. Aust.* **42**, 436 (1985)
7. Steer, H: *J. Pathology.*, 146:355.1985
8. Strickland, R. and Mackay,I: *Amer. J. Diag. Dis.*,18:426,1973.
9. McKenna,D: *Gastroenterology.*,912:1528,1987
10. Blaser,M. (ed): *Campylobacter Pylori in Gastritis and Peptic Ulcer Disease*. New york, Ikagu-Shion. 1989.