

DRG® 1,25 (OH) Vitamin D (EIA-4193)



RUO in the USA

Revised 12 Nov. 2010 rm (Vers. 3.1)

Please use only the valid version of the package insert provided with the kit.

INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for determination of 1,25-dihydroxy vitamin D in serum and plasma.

This kit is intended for Research Use Only. Not for use in diagnostic procedures.

MATERIAL SUPPLIED

Content	Kit Components	Quantity
PLATE	One holder with precoated strips	12 x 8 wells
WASHBUF	ELISA wash buffer concentrate 10x	100 mL
ETHANOL	Ethanol , ready-to-use	1.5 mL
TRIS-HCL	Tris-HCL buffer , ready-to-use	30 mL
AB	Detection antibody , anti 1,25-(OH) ₂ vitamin D, ready-to-use	25 mL
STD	Standard incl. NSB, ready-to-use (for range see specification or label)	7 x 2.5 mL
CTRL	Controls , ready-to-use (for range see specification)	2 x 2.5 mL
CONJ	Conjugate , polyclonal peroxidase labeled antibody, ready-to-use	22 mL
SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	2 x 15 mL
STOP	ELISA stop solution , ready-to-use	15 mL

MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- 48 Chromabond columns (Catalog No.: Se2112)
- 48 Silica Cartridges (solid phase extraction cartridges, Catalog No.: Sb2221)
- Diisopropylether (p.A.) 99.0 %
- Isopropanol (p.A.) 99.9 %
- n-Hexan (p.A.) 98.3 %
- Methanol (p.A.) 99.9 %

DRG® 1,25 (OH) Vitamin D (EIA-4193)

**RUO** in the USA

Revised 12 Nov. 2010 rm (Vers. 3.1)

- 75 x 12 glass tubes (no plastic)
- Extraction rack (Catalog No.: K2221sv)
- Precision pipettors and disposable tips to deliver 10-1000 µL
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Vacuum centrifuge or nitrogen distributor
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label.
Prepare just the appropriate amount necessary for the assay.
The kit can be used up to 2 times within the expiry date stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 mL concentrate + 900 mL aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be redissolved at room temperature or at 37°C using a water bath before dilution of the buffer solutions.
The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label.
Diluted **buffer solution** could be stored in a closed flask at **2-8°C for no longer than one month.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

DRG® 1,25 (OH) Vitamin D (EIA-4193)



RUO in the USA

Revised 12 Nov. 2010 rm (Vers. 3.1)

SAMPLE PREPARATION**Serum/plasma samples**


Fresh collected blood should be centrifuged within one hour. Vitamin D is an inert substance. The serum samples can be stored at room temperature. However, serum storage at 2-8°C is recommended when the analysis is performed within 24 h after collection. Otherwise, the serum samples must be stored at -20°C until analyzed. Avoid repeated freeze-thaw cycles. Serum samples can be shipped at 4-8 °C (for example with Coolpacks) and remain stable for up to 3 days.

Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.

We recommend to apply 1000 µL sample per the cartridge.

If the sample volume is less than 1000 µL, load an appropriate amount of Tris-HCl buffer in the column, and then the sample (minimum 500 µL) for a total volume of 1000 µL.

To calculate the actual concentration, each result should be multiplied with the respective dilution factor.

Extraction of the samples

– The extraction unit consists of three parts, which were put on top of each other.

DRG® 1,25 (OH) Vitamin D (EIA-4193)**RUO** in the USA**Revised 12 Nov. 2010 rm (Vers. 3.1)**

- The upper part is used for the Chromabond columns (extraction rack I), the lower part for the silica cartridges (extraction rack II).
- During sample application and the entire washing procedure, the whole unit should be put into a container big enough to collect the extraction solvents (extraction rack I and extraction rack II). After the first extraction step (ether) remove the extraction rack I with the Chromabond columns.
- It is recommended to place the glass-tubes (extraction rack III) directly under the cartridges (extraction rack II) for the last elution step. The tubes can then be used directly for the next step of the assay.

ASSAY PROCEDURE**Principle of the test**

The assay utilizes of a competitive Enzyme-Immuno-Assay (EIA) technique with a selected monoclonal antibody recognizing 1,25-dihydroxy vitamin D.

Standards, NSB (non-specific binding), controls and samples which are assayed for 1,25-dihydroxy vitamin D are incubated after the extraction step with the detection antibody. The pre-incubated solution is then transferred to the microplate coated with 1,25-dihydroxy vitamin D. During this incubation step, 1,25-dihydroxy vitamin D in the sample and a fixed amount of 1,25-dihydroxy vitamin D bound to the microtiter well compete for the binding of the detection antibodies. Then a peroxidase-conjugated anti-mouse antibody is added into each microplate well and a complex of 1,25-dihydroxy vitamin D - detection antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 1,25-dihydroxy vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 1,25-dihydroxy vitamin D in the samples is determined from this curve.

Extraction

1. Apply **1000 µL** of standards, NSB, control and sample (plasma or serum) on the Chromabond columns and incubate for 10 minutes. For sample volumes less than 1000 µL wet the cartridges with Tris–HCl buffer, e.g. pipette 500 µL Tris–HCl buffer in the cartridge and 500µl sample
2. Extract vitamin D from the chromabond columns with **4 x 1 mL** diisopropylether (3 min for each elution). The eluate should drip from the chromabond column directly on an untreated and dry silica cartridge.
After the extraction the chromabond columns should be removed (extraction rack I)
3. Wash the silica cartridges (extraction rack II) with **5 x 2 mL Isopropanol/ Hexane (4/96 v/v)**
4. Wash the silica cartridges (extraction rack II) with **3 x 2 mL Isopropanol/ Hexane (6/94 v/v)**
5. Elute 1,25-dihydroxy vitamin D from the silica cartridges with **2 x 2 mL Isopropanol/Hexane (25/75 v/v)**. Note: the glass tubes (extraction rack III) should be placed directly under the silica cartridges
6. Evaporate the eluate under a nitrogen stream at 37 °C or in a vacuum centrifuge

DRG® 1,25 (OH) Vitamin D (EIA-4193)**RUO** in the USA**Revised 12 Nov. 2010 rm (Vers. 3.1)****Pre-incubation**

1. Add **20 µL** of ethanol into each glass tube. Immediately after adding ethanol gently vortex each tube to avoid any possible evaporation.
2. Add **450 µL** antibody solution into each glass tube. The antibody solution is viscous. Please pipette slowly and carefully. Mix thoroughly.
3. Cover glass tubes with a plastic film and incubate exactly for **1 hour** at room temperature.

Test procedure

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well
2. Mark the positions of STD/NSB/SAMPLE/CTRL (Standards/non-specific binding/Sample/Control) on a protocol sheet
3. Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label
4. Add 200 µL of STD/NSB/SAMPLE/CTRL in duplicate into respective well. All these solutions are viscous. Please pipette slowly and carefully. We recommend to wet the pipette tip before using it to transfer the preincubate.
5. Cover the plate tightly and incubate for 18-22 hours at **6 - 10 °C***
6. Discard the contents of each well. Wash 5 times by dispensing 250 µL of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
7. Add 200 µL CONJ (conjugate) into each well
8. Cover the plate tightly and incubate for exactly 1 hour at room temperature (18-26°C) on a horizontal mixer
9. Discard the contents of each well. Wash 5 times by dispensing 250 µL of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
10. Add 200 µL of SUB (substrate) into each well
11. Incubate for 20 - 30 minutes at room temperature (18-26°C) in the dark**
12. Add 50 µL of STOP (stop solution) into each well, mix thoroughly
13. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

**As with any competitive immunoassay, consistent incubation times and temperature are essential for accurate plate-to-plate comparisons. Fluctuations in overnight incubation can lead to increased inter-assay CV's. It is therefore recommended to use always the same incubation time, i.e. 20 hours.*

***The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.*

DRG® 1,25 (OH) Vitamin D (EIA-4193)**RUO** in the USA**Revised 12 Nov. 2010 rm (Vers. 3.1)****RESULTS**

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. *4-Parameter-algorithm*

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration.

When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. *Point-to-point-calculation*

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. *Spline-algorithm*

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e.g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

REGENERATION OF THE SILICA CARTRIDGES

The silica cartridges can be used up to 5 times if regenerated as follows. The regeneration can be performed directly after the extraction. The silica cartridges have to be dry before the next use.

- **2 x 2 mL** methanol
- **2 x 2 mL** n-hexan
- Dry the columns in the hood

PRECAUTIONS

- This kit is intended for Research Use Only. Not for use in diagnostic procedures.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.

DRG® 1,25 (OH) Vitamin D (EIA-4193)**RUO** in the USA**Revised 12 Nov. 2010 rm (Vers. 3.1)**

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for research use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. DRG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to DRG together with a written complaint.

REFERENCES / Literature

1. Wildermuth S. et al., 1993; Clinica chimica Acta, 220, 61
2. Schilling M. et al., 1987; Clinical Chemistry, 33, 187
3. Armbruster F.P. et al., 1990; Ärztl. Lab., 36, 75
4. Durham B.W. et al., 1995; Ann. Clin. Biochem., 32, 77
5. Hollis B.W., 1996; Calcif. Tissue Int., 58, 4
6. Hollis B.W., 1995; Clinical Chemistry, 41, 1313
7. Iqbal S.J. et al., 1996; Clinical Chemistry, 42, 112
8. Withold W. et al., 1995; Eur. J. Clin. Chem. Clin. Biochem., 33, 15