

DRG[®] HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



For Research Use Only

INTRODUCTION

Human Antibodies to mouse immunoglobulins are known as HAMA (Human Anti-Mouse Antibody) ¹⁻¹⁰.

The HAMA ELISA kit is an enzyme immunoassay that provides materials for the sensitive detection and quantitation of HAMA.

The HAMA kit is completed with all reagents ready to use including standards and a positive control. Test samples are diluted and quantitated in µg HAMA/mL of sample from the standard curve. Test Samples and an enzyme-linked mouse IgG are added to coated wells of a 96-well microplate. A TMB-substrate is added and the enzyme reaction stopped with the absorbance read at 450 nm. The assay can be completed in less than one hour.

MATERIAL PROVIDED

1. **Micro-well Strips:** 8x12 microplate strips coated with Mouse- IgG, 96 wells.
2. **Enzyme Conjugate** (11 mL): Mouse IgG conjugated to horseradish peroxidase
3. **Sample Diluent** or Zero Standard (50 mL).
4. Reference **Standard Set** (0.75 mL/each):
Calibrated to 25, 100, 400, and 1000 ng/mL in BSA-containing diluent.
5. HAMA **positive Control** (0.75 mL)
6. Concentrated **Wash Buffer (100x)** (10 mL)
7. **TMB Solution** (11 mL): Buffer solution containing hydrogen peroxide and TMB
8. **Stop Solution** (10 mL): 2N HCl.
9. **Well holder:** for securing wells.
10. Package Insert.

*** The reagents except TMB Solution contain 0.01% Thimerosal as a preservative***

MATERIALS REQUIRED BUT NOT PROVIDED

1. Micro-well reader
2. Disposal tips and pipettor for measuring 25L µL, 50 µL and 100 µL.

WARNING AND PRECAUTION

1. The HAMA ELISA is designed for in Vitro Research Use only.
2. The Components in the kits are intended for usage as an integral unit. The components from different lots should not be mixed, and not be used beyond expiration date.
3. The material should be used in a designated work area, the bench surface should be cleaned with detergent and the contaminated materials should be disposed properly.

DRG[®] HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



4. Some components have been tested using FDA-approved methods and has been found negative for antibody to human immuno-deficiency virus (HIV-I, HIV-II), antibody to Hepatitis C and Hepatitis B surface antigen (HBsAg). No known test method can offer total assurance that HIV-I, HIV-II, Hepatitis B & C virus or other infectious agents are absent. Handle these reagents as if they were potentially infectious. Information on handling human serum is provided in the CDC/NIH manual A Biosafety in Microbiological and Biomedical Laboratories@ (U.S.A. HHS publication No. (NIH 88-8395.)
5. Avoid microbial contamination of reagents when removing aliquots from the vials.

STORAGE AND STABILITY

1. Store the kits at 2-8°C in a refrigerator.
2. Keep micro-wells in a dry bag with desiccants.
3. The reagents are stable until expiration of the kit. TMB Solution should be colorless; if the solution turns blue, it must be replaced. Do not expose these reagents to strong light during storage or usage.

SPECIMEN AND COLLECTION

Collect blood by venipuncture, allow to clot and separate the serum by centrifugation at room temperature.

Serum is required for the HAMA ELISA, and do not add sodium azide as preservative.

Samples should not be stored at room temperature or 4 °C for more than 24 hours.

Serum samples are recommended to be frozen for longer storage.

Avoid repeated freezing and thawing of serum samples.

Mild hemolysis and lipemia have been shown not to interfere with the results. Specimens that are grossly lipemic, hemolyzed or contaminated may interfere and should not be used.

PREPARATION FOR ASSAY

1. Bring all reagents and samples to room temperature (20 °C – 25 °C) and shake gently before beginning the test. Have all reagents and samples ready before the start of the assay. Once the test is begun it must be performed without any interruption to get the most reliable and consistent results.
2. Use new disposable tips for each specimen.

ASSAY PROCEDURE

Note: Standards and Positive Control are ready to use, do not dilute.

Sample dilution: Add 10 µL of test sample to sample diluent in the tube and mix (1:101 dilution).

It is recommended that samples, standards and positive control be run in duplicate.

1. Secure the desired number of coated wells in the holder. Mark data sheet with sample identification.
2. Dispense 25 µL of references, controls or diluted serum samples (1:101) into the appropriate wells.
3. Dispense 100 µL of enzyme conjugate into wells.

DRG[®] HAMA (EIA-4195)**Revised 8 Jan. 2009 (Vers. 3.0)**

4. Incubate for 30 minutes at room temperature.
5. Remove incubation mixture and rinse the wells 5 times with diluted washing buffer.
6. Dispense 100 μ L of TMB solution into each well.
7. Incubate for 15 minutes at room temperature.
8. Stop reaction by adding 50 μ L of 2 N HCl to each well.
9. Blank with the substrate only well(s) and read absorbance at 450 nm for dual wavelength readers use 570 nm as the reference wavelength.

PROCEDURE NOTE

1. Wash the microwells and remove water thoroughly.
2. Pipet all reagents and samples into bottom of the well. Vortex-mixing or shaking is not required.
3. Absorbance is a function of the time and temperature of incubations. It is recommended to have reagents, samples and needed wells ready for ensure the equal elapsed time for each pipetting without interruption.
4. For the same reason run no more than 20 patient samples with a set of reference standards in duplicate for each assay.

QUALITY CONTROL

Good laboratory practices include the use of control specimens to ensure that all reagents and protocols are performing properly. The HAMA ELISA kit does include serum control.

CALCULATION OF RESULTS


1. Plot the concentration (X) of each reference standards against its absorbance (Y) on a 3 cycle log-log paper (see below) paper.
2. Obtain the HAMA value of patient by reference to the standard curve as follows: (These data are for demonstration purpose only and must not be used in place of data generated for each assay).

DRG[®] HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



Well No.	Description (ng/mL)	Absorbance (450 nm)	HAMA (ng/mL)
A1	0	0.028	
B1	(Blank)	0.027	
A2	25	0.128	
B2		0.134	
A3	100	0.377	
B3		0.382	
A4	400	1.441	
B4		1.438	
A5	1000	2.590	
B5		2.624	
A6	CONTROL	0.053	8.7 x 10 ¹
B6		0.052	8.5 x 10 ¹
A7	CONTROL	0.608	162.6 x 10 ¹
B7		0.611	163.5 x 10 ¹



Determine the level of HAMA (ng/mL) in the test sample by reading of the standard curve and multiply result by dilution factor (x101).

Samples with absorbances greater than the highest standard should be diluted further and retest to quantitate the HAMA. Record results in µg HAMA/mL of sample.

DRG® HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



PERFORMANCE CHARACTERISTICS

Accuracy

Recovery studies were performed by mixing equal volume of test samples negative for HAMA spiked with known concentrations of HAMA. The HAMA values were measured and percentage of recovery determined.

Initial Values (ng/mL)	Concentration Spiked (ng/mL)	Expected Values (ng/mL)	Observed Values (ng/mL)	Recovery (%)
0	22	11	11	100
0	30	15	14	93
0	52	26	23	89
0	60	30	28	93
0	90	45	42	93

Precision

Intra-assay and inter-assay coefficient of variation were evaluated at three different pooled serum samples.

Intra-assay	Pool A	Pool B	Pool C
N	18	18	18
Mean (ng/mL)	26.3	203.6	528.1
S.D. (ng/mL)	1.7	11.1	21.5
C.V. (%)	6.5	5.5	4.1

Inter-assay	Pool A	Pool B	Pool C
N	18	18	18
Mean (ng/mL)	24.1	203.5	523.6
S.D. (ng/mL)	2.0	12.4	22.2
C.V. (%)	8.3	6.1	4.2

Specificity

In studies of this HAMA ELISA Kit the presence of rheumatoid factor as well as other autoantibodies (anti-DNA and anti-ENA) were shown not to interfere with the accurate detection and quantitation of HAMA. The presence of anti-viral (CMV, HSV and Hepatitis A and B) and anti-Toxoplasma antibodies in either the acute or convalescent phase of infection were also shown not to interfere with the HAMA ELISA kit (see data in table on reference range).

NOTE: Circulating mouse immunoglobulin (Ig) in samples may interfere with the accurate detection and HAMA quantitation. In these samples the mouse Ig should be quantitated.

DRG® HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



REFERENCE RANGE

Various samples were tested for presence of HAMA and the results are listed below.

Groups	n	Mean Abs	Abs. Range
Group 1 Healthy Controls	50	0.081	0.043-0.119
Group 2 Rheumatic Disease patients	30	0.080	0.069-0.092
Group 3 Patients with Acute Infections	30	0.088	0.069-0.110
Group 4 Patients treated with Mouse Monoclonal Antibodies	12	0.559	0.010-2.500

REFERENCES





















1. Hammond EH, Wittwer CT, Greenwood J, Knape WA, Yowell RL, Menlove RL, Craven C, Renlund DG, Bristow MR, DeWitt CW and O'Connell JB. 1990. Relationship of OKT3 sensitization and vascular rejection in cardiac transplant patients receiving OKT3 rejection prophylaxis. Transplantation 50:776-782.
2. Jaffers Gj, Colvin RB, Cosimi AB, Giorgi JV, Goldstein G, Fuller TC, Kurnick JT, Lillehei C and Pursell PS. 1983. The human immune response to murine OKT3 monoclonal antibody. Transplant. Proc. 15:646-648
3. Miller Ra, Maloney DG, McKillop J. and Levy R. 1981. In vivo effects of murine hydridoma monoclonal antibody in a patient with T-cell leukemia. Blood 58:78-84.
4. Schoff RW, Foon Ka, Beatty SM, Oldham RK and Morgan AC. 1985. Human Anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res. 45:879-885.
5. Shawler DL, Bartholomew RM, Smith LM and Dillman RO. 1985. Human immune response to multiple injections of murine monoclonal IgG. J. Immunol. 135:1530-1535.
6. Sasakuma F, Shimisu T, Wada H, Morii T, Sasaki A and Ehara M. 2000. Human anti-murine antibodies interfere with CRP assays performed with commercial kits. Diabetes Res Clin Pract. 48 (2) 105-11
7. Oltrogge JB, Baum RP, Lema KN, Donnerstag B, and Hor G. 1997. How to overcome the disturbing effects of human anti-mouse antibodies (HAMA) on in vitro assays. Int. J. Biol. Markers. 12 (1) 15-7.
8. Hasholzner U, Stieber P, Meeier W., and Lamerz R. 1997. Value of HAMA-determination in clinical practice. Anticancer Res. 17(4B): 3055-8
9. Kimball JA, Norman DJ, Shield CF, Schroeder TJ, Lisi P, Garovoy M, O'Connell JB, Stuart F, McDiarmid SV, and Wall W. 1995. The OKT3 antibody response study: a multicentre study of human anti-mouse antibody (HAMA) production following OKT3 use in solid organ transplantation. Transpl. Immunol. 3 (3) 212-21.
10. Perosa F, Rizzi R, Pulpito V, and Dammacco F. 1995. Soluble CD4 antigen reactivity in intravenous immunoglobulin preparations: Is it specific? Clin Exp Immunol. 99(1) 16-20.

DRG® HAMA (EIA-4195)

Revised 8 Jan. 2009 (Vers. 3.0)



Symbols used with DRG Assays

Symbol	English	Deutsch	Français	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
Content	Content	Inhalt	Conditionnement	Contenido	Contenuto
Volume/No.	Volume / No.	Volumen/Anzahl	Volume/Quantité	Volumen/Número	Volume/Quantità
Symbol	Portugues	Dansk	Svenska	Ελληνικά	
	Consulte as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Εγχειρίδιο χρήστη	
	Conformidade com as normas europeias	Europeisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση	
	Diagnóstico in vitro	In vitro diagnostik	Diagnostik in vitro	in vitro διαγνωστικό	
					
	Catálogo n.º	Katalognummer	Katalog nummer	Αριθμός καταλόγου	
	No do lote	Lot nummer	Batch-nummer	Αριθμός Παρτίδος	
		Indeholder tilstrækkeligt til "n" test	Innehåller tillräckligt till "n" tester	Περιεχόμενο επαρκές για «n» εξετάσεις	
	Temperatura de conservação	Opbevarings-temperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης	
	Prazo de validade	Udløbsdato	Bäst före datum	Ημερομηνία λήξης	
	Fabricante	Producent	Tillverkare	Κατασκευαστής	
Distributed by					
Content	Conteúdo	Indhold	Innehåll	Περιεχόμενο	
Volume/No.	Volume/Número	Volumen/antal	Volym/antal	Όγκος/αριθ..	