

REVISED 8 MAR. 2011 RM (VERS. 2.1)

Please use only the valid version of the package insert provided with the kit.

This kit is intended for Research Use Only.

Not for use in diagnostic procedures.

1 INTENDED USE

Anti-ASGPR is used for determination of IgG antibodies to asialoglycoprotein receptor (ASGPR) in human serum or plasma.

McFarlane BM, McSorley CG, Vergani D, McFarlane IG, Williams R: Serum autoantibodies reacting with the hepatic asialoglycoprotein receptor protein (hepatic lectin) in acute and chronic liver disorders. J Hepatol 1986 3, 196-205

Treichel U, Poralla T., Hess G, Manns M, Meyer zum Buschenfelde KH: Autoantibodies to human asialoglycoprotein receptor in autoimmune-type chronic hepatitis. Hepatology 1990 11, 606-612

PRINCIPLE OF THE TEST 2

Anti-ASGPR is an enzyme immunoassay for determination of IgG antibodies to ASGPR.

The antibodies of the calibrators and diluted samples react with ASGPR immobilized on the solid phase of microtiter plates. ASGPR highly purified from rabbit liver and coated on the microtiter plate guarantees the specific binding of ASGPR IgG antibodies of the specimen under investigation. Following an incubation period of 60 min at 37 °C, unbound serum components are removed by a washing step.

The bound IgG antibodies react specifically with anti-human-IgG conjugated to horseradish peroxidase (HRP) within an incubation period of 30 min at 37 °C. Excessive conjugate is separated from the solid-phase immune complexes by the following washing step.

Horseradish peroxidase converts the colorless substrate solution of 3,3',5,5'-tetramethylbenzidine (TMB) added into a blue product. The enzyme reaction is stopped by dispensing an acidic solution (H_2SO_4) into the wells after 10 min at room temperature turning the solution from blue to yellow.

The optical density (OD) of the solution at 450 nm is directly proportional to the amount of specific antibodies bound.

3 SPECIMEN SAMPLES

3.1 **Specimen collection and storage**

Blood is taken by venipuncture. Serum is separated after clotting by centrifugation. Plasma can be used, too. Lipaemic, haemolytic and contaminated samples should not be used.

Repeated freezing and thawing should be avoided. If samples are to be used for several assays, initially aliquot samples and keep at -20 °C.







C E REVISED 8 MAR. 2011 RM (VERS. 2.1)



3.2 Preparation before use

Allow samples to reach room temperature prior to assay. Take care to agitate serum samples gently in order to ensure homogeneity.

<u>Note:</u> Samples have to be diluted 1 + 100 (v/v), e.g. $10 \mu L$ sample $+ 1.0 \mu L$ sample diluent (C), prior to assay.

The samples may be kept at 2 °C - 8 °C for up to two days. Long-term storage requires - 20 °C.









4 TEST COMPONENTS

for 96 determinations			
A 96 Ag	Microtiter plate , 12 breakable strips per 8 wells (total 96 individual wells) coated with ASGPR (rabbit)	1 vacuum sealed with desiccant	
B BUF WASH 10x	Concentrated wash buffer sufficient for 1000 mL solution	100 mL concentrate capped white	
C DIL	Sample diluent	100 mL ready for use capped black	
D Conj	Conjugate containing anti-human-IgG-(sheep) coupled with HRP	15 mL ready for use capped red	
E SOLN TMB	Substrate 3,3',5,5'-tetramethylbenzidine in citrate buffer containing hydrogen peroxide	15 mL ready for use capped blue	
F H2SO4 0.25M	Stop solution 0.25 M sulfuric acid	15 mL ready for use capped yellow	
P CONTROL	Positive control (diluted serum) +	1 mL ready for use	
	Cut-off control (diluted serum)	1 mL ready for use	
	Negative control (diluted serum) -	1 mL ready for use	



DRG[®] A CE REVISED 8 MAR. 2011 RM (VERS. 2.1)



RUO IN THE USA

4.1 Materials required

- micropipette 100 1000 μL
- micropipette 10 100 μL
- multi-channel pipette 50 200 μL
- rough for multi-channel pipette
- 8-channel wash comb with vacuum pump and waste bottle or microplate washer
- incubator (37 °C)
- microplate reader with optical filters for 450 nm and 620 nm or 690 nm
- graduated cylinders
- distilled or de-ionized water

4.2 Size and storage

Anti-ASGPR has been designed for 96 determinations.

The expiry date of each component is reported on its respective label that of the complete kit on the box labels. Upon receipt, all components of the Anti-ASGPR have to be kept at 2 °C - 8 °C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage.

4.3 Preparation before use

Allow all components to reach room temperature prior to use in the assay.

The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed microplate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

<u>Prepare a sufficient amount of wash solution</u> by diluting the concentrated wash buffer 10 times (1 + 9) with de-ionized or distilled water. For example, dilute 8 mL of the concentrate with 72 mL of distilled water per strip. The wash solution prepared is stable at 2 °C - 8 °C up to 30 days.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle.

Avoid exposure of the TMB substrate solution to light!





REVISED 8 MAR. 2011 RM (VERS. 2.1)



- ASSAY PROCEDURE 5
- Dilute sera with sample diluent (C) 1 + 100 (v/v), • e.g. 10 µL serum + 1.0 mL sample diluent (C).
- Avoid any time shift during pipetting of reagents and samples.
- Bring all reagents to room temperature before use. Mix gently without causing foam. 1.
- Dispense 2. 100 µL controls (P, Co, N) 100 µL diluted samples into the respective wells.
- Seal plate, incubate 60 min at 37 °C. 3.
- 4. Decant, then wash each well five times using 300 µL wash buffer (B).
- 5. Add $100 \ \mu L$ of conjugate (D) solution to each well.
- Seal plate, incubate **30 min** at 37 °C. 6.
- Decant, then wash each well five times using $300 \ \mu L$ wash buffer (B). 7.
- 8. Add 100 μ L of substrate (E) to each well.
- Incubate 10 min protected from light at room temperature. 9.
- 10. Add 100 μ L of stop solution (F) to each well and mix gently.
- 11. Read the optical density at 450 nm versus 620 or 690 nm within 30 min after adding the stop solution.

6 **DATA PROCESSING**

Results are interpreted by calculating the following ratio of OD:

ratio = OD_{sample}/OD_{cut-off control}

This calculation can be done by the integrated evaluation software of the microplate reader used, too.





INCUBATION SCHEME

Dilute s sample 10 µL serum + 1.0 mL sample diluent (C)				
1	Bring all reagents to room temperature (18-25°C)			
2	Pipette controls (P, Co, N) 1 + 100 prediluted sera	100 μL	100 μL	
3	Seal plate and incubate	60 min, 37 °C		
4	Wash	Decant, 5 x 300 µL (made of B)		
5	Pipette conjugate (D)	100 µL	100 µL	
6	Seal plate and incubate 30 min, 37 °C			
7	Vash Decant, 5 x 300 μL (made of B)			
8	Pipette substrate (E)	100 µL	100 µL	
9	Incubate protected from light	10 min, room temperature (18-25 °C)		
10	Pipette stop solution (F)	100 µL	100 µL	
11	Read at 450 nm against 620 (690) nm within 30 min.			

7 SAFETY PRECAUTIONS

• Follow the working instructions carefully. DRG and its authorized distributors shall not be liable for damages indirectly or consequentially brought about by changing or modifying the procedure indicated. The kit should be performed by trained technical staff only.









C E REVISED 8 MAR. 2011 RM (VERS. 2.1)

- The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for reconstituted reagents.
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept at 2 °C 8 °C before use in the original shipping container.
- Some of the reagents contain small amounts of Thimerosal (< 0.1 % w/v) and Kathon (1.0 % v/v) as preservatives. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Source materials derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and for HIV as well as HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed:
 - Do not smoke, eat or drink while handling kit material,
 - Always use protective gloves,
 - Never pipette material by mouth,
 - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.

Version 2011-02-15~rm



