# **Renin-EASIA**

***KAP1531***

Read entire protocol before use.

# RENIN-EASIA

## I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative determination of active Renin in human plasma .

## II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource Renin-EASIA Kit
- B. **Catalogue number :** KAP1531 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)67 88.99.99      Fax : +32 (0)67 88.99.96**

## III. CLINICAL BACKGROUND

Renin, a polypeptidic enzyme (MW~ 40000) (1) also known as angiotensinogenase, is a circulating protease secreted by juxtaglomerular cells in the juxtaglomerular apparatus of the kidneys in response to low blood volume or low body NaCl content.

Renin activates the renin-angiotensin system by cleaving angiotensinogen produced in the liver into angiotensin I (inactive) which is further converted into angiotensin II (active) in the vascular epithelium of the lung. Angiotensin II can cause vasoconstriction by stimulating the central nervous system, in addition it stimulates ADH (antidiuretic hormone) secretion and aldosterone secretion from the adrenal gland (6).

Regulation of blood pressure and renal glomerular filtration control (2) are the most important functions of renin -angiotensin system.

Plasmatic concentration of renin is influenced by concentration of circulating angiotensinogen and subsequently the concentration of angiotensin II. High plasmatic levels of angiotensin II reduce renin secretion (negative feed back).

Determination of renin plasma levels is useful in the diagnosis of hypertension and in the therapeutic follow up of hypertensive patients (3).

Plasmatic concentration of renin decreases in patients with hypertension due to a primary hyperaldosteronism (4), contrary to renovascular hypertension (5) where concentrations of renin and aldosterone are both elevated.

## IV. PRINCIPLES OF THE METHOD

The DAsource Renin-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplates. Calibrators and samples react with the capture monoclonal antibody (MAb 1) coated on microtiter well and with a monoclonal antibody (MAb 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAB 1 – human Renin – MAB 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the human Renin concentration. A calibration curve is plotted and Renin concentration in samples is determined by interpolation from the calibration curve.

## V. REAGENTS PROVIDED

| Reagents                                                                                                                                                         | 96 tests Kit           | Color Code | Reconstitution                                                                              |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------|------------|---------------------------------------------------------------------------------------------|
|  Microtiterplate with 96 anti Renin monoclonal antibodies coated breakable wells | 96 wells               | blue       | <b>Ready</b> for use                                                                        |
| <div>AbHRPCONC</div> Conjugate: HRP labelled anti-Renin (monoclonal antibodies) in TRIS buffer with bovine serum albumin and thymol                              | 1 vial<br>0.05 ml      | red        | <b>To be diluted</b> with conjugate buffer (see on the label for the exact dilution factor) |
| <div>CONJBUF</div> Conjugate buffer: Tris buffer with bovine serum albumin and thymol                                                                            | 1 vial<br>10.5 ml      | red        | <b>Ready</b> for use                                                                        |
| <div>CALN</div> Calibrators - N = 0 to 5 , in human serum and thymol. See exact value on vial labels.                                                            | 6 vials<br>lyophilized | yellow     | <b>Add</b> 2 ml distilled water                                                             |
| <div>CONTROLN</div> Controls - N = 2 in human serum with thymol.                                                                                                 | 2 vials<br>lyophilized | silver     | <b>Add</b> 2 ml distilled water                                                             |
| <div>WASHSOLNCONC</div> Wash solution (Tris-HCl)                                                                                                                 | 1 vial<br>10 ml        | brown      | <b>Dilute</b> 200 x with distilled water (use a magnetic stirrer).                          |
| <div>CHROMTMB</div> Chromogenic TMB Solution (Tetramethylbenzidine)                                                                                              | 1 vial<br>12 ml        | white      | <b>Ready</b> for use                                                                        |
| <div>STOP SOLN</div> Stop solution: HCl 1.0 N                                                                                                                    | 1 vial<br>12 ml        | black      | <b>Ready</b> for use                                                                        |

**Note:** 1. Use the zero calibrator for sample dilution.  
 2. 1 pg of the calibrator preparation is equivalent to 2.2 +/- 0.2 µIU of NIBSC 68/356.  
 Values obtained in pg/ml must be multiplied by 2.2 to obtain results in µIU/ml or mIU/l.

## VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 25 µl ,100 µl, 200 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Plastic tubes for dilution of samples
- Vortex mixer
- Magnetic stirrer

- Horizontal microtiterplate shaker capable of 400 rpm ± 100 rpm
- Washer for microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (monochromatic reading)
- Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XIA.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

## VII. REAGENT PREPARATION

- A. Calibrators:** Reconstitute the calibrators with 2 ml distilled water.  
**! For a complete solubilisation : after reconstitution, let the vials 30 min on a shaker, then vortex them.**
- B. Controls:** Reconstitute the controls with 2 ml distilled water .

**C Working Renin-HRP conjugate :** Prepare an adequate volume of conjugate solution by diluting the concentrated Renin-HRP conjugate with conjugate buffer. See on the label for the exact dilution factor. Use a vortex to homogenize. Extemporaneous preparation is recommended.

**D. Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

## VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, the calibrators and controls are unstable, use them immediately after reconstitution, freeze them immediately in aliquots and keep them at –20°C for maximum 6 weeks. They are stable after 1 freeze-thaw cycle .
- Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

## IX. SPECIMEN COLLECTION AND PREPARATION

- Samples must be EDTA plasma.
- If the test is not run within 4 h , plasma should be aliquoted and stored at –20°C.
- Avoid subsequent freeze-thaw cycles.
- Do not use haemolysed samples.

## X. PROCEDURE

### A. Handling notes

Do not use the kit or components beyond expiry date.  
 Do not mix materials from different kit lots.  
 Bring all the reagents to room temperature prior to use.  
 Thoroughly mix all reagents and samples by gentle agitation or swirling.  
 Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.  
 Use a clean plastic container to prepare the Wash Solution.  
 In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.  
 For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.  
 High precision pipettes or automated pipetting equipment will improve the precision.  
 Respect the incubation times.  
 Prepare a calibration curve for each run, do not use data from previous runs.  
 Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.  
 During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

### B. Procedure

- Select the required number of wells for the run. The unused wells should be resealed in the bag with a desiccant and stored at 2-8°C.
- Secure the wells into the holding frame.
- Pipette 200 µl of each calibrator, control and sample into the appropriate wells.
- Incubate for 1 hour at room temperature on a horizontal shaker set at 400 rpm ± 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times

7. Pipette 100µl working anti-Renin-HRP conjugate into all the wells.
8. Incubate for 1 hour at room temperature on a horizontal shaker set at 400 rpm  $\pm$  100 rpm.
9. Aspirate the liquid from each well.
10. Wash the plate 3 times .
11. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
12. Incubate the microtiterplate for 30 minutes at room temperature, avoid direct sunlight
13. Pipette 100 µl of Stop Solution into each well.
14. Read at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

## XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of Hu Renin (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

## XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

| Hu Renin-EASIA |           | OD units |
|----------------|-----------|----------|
| Calibrator     |           |          |
|                | 0 pg/ml   | 0.074    |
|                | 3.9 pg/ml | 0.121    |
|                | 11 pg/ml  | 0.2      |
|                | 55 pg/ml  | 0.68     |
|                | 91 pg/ml  | 1.068    |
|                | 270 pg/ml | 2.6      |

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection Limit

Twenty four zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.8 pg/ml.

### B. Specificity

The Cross-reactivity of Prorenin in this Renin EASIA assay was determined by adding various concentrations of Prorenin to a plasma matrix and by measuring the apparent Renin response. Prorenin cross-reactivity was found to be 0.2 %.

### C. Hemoglobin and bilirubin interferences

Hemoglobin at 7.5 mg/ml and bilirubin at 0.2 mg/ml were tested in Dia source Renin EASIA kit. Results are shown in the table below. No significant interference was found with bilirubin, but hemolytic samples should be avoided.

| sample   | Renin value (pg/ml) | Renin value + Hemoglobin (pg/ml) | Renin value + Bilirubin (pg/ml) |
|----------|---------------------|----------------------------------|---------------------------------|
| Plasma 1 | 12.8                | 7.5                              | 13.1                            |
| Plasma 2 | 62.7                | 41.0                             | 60.0                            |
| Plasma 3 | 10.8                | 7.1                              | 11.7                            |

### D Precision

| INTRA ASSAY |    |                      |        | INTER ASSAY |    |                      |        |
|-------------|----|----------------------|--------|-------------|----|----------------------|--------|
| Sample      | N  | <X> $\pm$ SD (pg/ml) | CV (%) | Sample      | N  | <X> $\pm$ SD (pg/ml) | CV (%) |
| A           | 22 | 16.7 $\pm$ 0.9       | 5.4    | A           | 10 | 16.8 $\pm$ 1.1       | 6.2    |
| B           | 24 | 67.3 $\pm$ 1.5       | 2.2    | B           | 10 | 64.6 $\pm$ 2.2       | 3.3    |

SD : Standard Deviation; CV: Coefficient of variation

### E Accuracy

## RECOVERY TEST

| Sample | Added Renin (pg/ml) | Recovered Renin (pg/ml) | Recovery (%) |
|--------|---------------------|-------------------------|--------------|
| Serum  | 37.1                | 38.0                    | 102.4        |
|        | 78.0                | 62.0                    | 79.5         |
|        | 131.2               | 109.1                   | 83.2         |
|        | 228.5               | 241.5                   | 105.7        |

## DILUTION TEST

| Sample | Dilution | Theoretical Concent. (pg/ml) | Measured Concent. (pg/ml) |
|--------|----------|------------------------------|---------------------------|
| 1      | 1/1      | -                            | 103.8                     |
|        | 1/2      | 53.6                         | 51.9                      |
|        | 1/4      | 24.2                         | 25.9                      |
|        | 1/8      | 11.1                         | 12.9                      |
| 2      | 1/1      | -                            | 105.8                     |
|        | 1/2      | 52.9                         | 51.6                      |
|        | 1/4      | 26.4                         | 26.6                      |
|        | 1/8      | 13.2                         | 13.8                      |
|        | 1/16     | 6.6                          | 5.5                       |

Samples were diluted with zero calibrator.

## XIV. INTERNAL QUALITY CONTROL

- § If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- § If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- § Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- § It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- § It is good practise to check visually the curve fit selected by the computer.

## XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Normal range (0.8 to 16.5 pg/ml) has been determined on 20 normal subjects with a mean value of 4.9 pg/ml.

## XVI. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains H2SO4. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## ***XVII. BIBLIOGRAPHY***

- (1) Primary structure of the human Renin gene : Hardman J.A; and al. ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

## ***XVIII. SUMMARY OF THE PROTOCOL***

|                                                                                                                                                                               | CALIBRATORS<br>(µl) | SAMPLE(S)<br>CONTROLS<br>(µl) |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|-------------------------------|
| Calibrators (0-6)<br>Controls, Samples                                                                                                                                        | 200<br>-            | -<br>200                      |
| Incubate for 1 hour at room temperature with continuous shaking at 400 rpm.<br>Aspirate the contents of each well.<br>Wash 3 times with 400 µl of Wash Solution and aspirate. |                     |                               |
| Anti Renin Conjugate                                                                                                                                                          | 100                 | 100                           |
| Incubate for 1 hour at room temperature with continuous shaking at 400 rpm.<br>Aspirate the contents of each well.<br>Wash 3 times with 400 µl of Wash Solution and aspirate. |                     |                               |
| Chromogenic Solution                                                                                                                                                          | 100                 | 100                           |
| Incubate for 30 min at room temperature.                                                                                                                                      |                     |                               |
| Stop Solution                                                                                                                                                                 | 100                 | 100                           |
| Read on a microtiterplate reader.<br><br>Record the absorbance of each well at 450 nm (versus 630 or 650 nm).                                                                 |                     |                               |

|                                     |                             |                           |
|-------------------------------------|-----------------------------|---------------------------|
| DIAsource Catalogue Nr :<br>KAP1531 | P.I. Number :<br>1700482/en | Revision nr :<br>110120/1 |
|-------------------------------------|-----------------------------|---------------------------|

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XVII. BIBLIOGRAPHIE

(1) Primary structure of the human Renin gene : Hardman J.A; and al. ; DNA (1984); 3(6):457-468.  
(2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372  
(3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .  
(4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .  
(5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180  
(6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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XVII. BIBLIOGRAFIA

(1) Primary structure of the human Renin gene : Hardman J.A; and al. ; DNA (1984): 3(6):457-468.  
(2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372  
(3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .  
(4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .  
(5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180  
(6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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|                 | <b><u>Used symbols</u></b>         |
|-----------------|------------------------------------|
|                 | Consult instructions for use       |
|                 | Storage temperature                |
|                 | Use by                             |
| <b>LOT</b>      | Batch code                         |
| <b>REF</b>      | Catalogue number                   |
| <b>CONTROL</b>  | Control                            |
| <b>IVD</b>      | In vitro diagnostic medical device |
|                 | Manufacturer                       |
|                 | Contains sufficient for <n> tests  |
| WASH SOLN CONC  | Wash solution concentrated         |
| CAL 0           | Zero calibrator                    |
| CAL N           | Calibrator #                       |
| CONTROL N       | Control #                          |
| Ag 1251         | Tracer                             |
| Ab 1251         | Tracer                             |
| Ag 1251 CONC    | Tracer concentrated                |
| Ab 1251 CONC    | Tracer concentrated                |
|                 | Tubes                              |
| INC BUF         | Incubation buffer                  |
| ACETONITRILE    | Acetonitrile                       |
| SERUM           | Serum                              |
| DIL SPE         | Specimen diluent                   |
| DIL BUF         | Dilution buffer                    |
| ANTISERUM       | Antiserum                          |
| IMMUNOADSORBENT | Immunoabsorbent                    |
| DIL CAL         | Calibrator diluent                 |
| REC SOLN        | Reconstitution solution            |
| PEG             | Polyethylene glycol                |
| EXTR SOLN       | Extraction solution                |
| ELU SOLN        | Elution solution                   |
| GEL             | Bond Elut Silica cartridges        |
| PRE SOLN        | Pre-treatment solution             |
| NEUTR SOLN      | Neutralization solution            |
| TRACEUR BUF     | Tracer buffer                      |
|                 | Microtiterplate                    |
| Ab HRP          | HRP Conjugate                      |
| Ag HRP          | HRP Conjugate                      |
| Ab HRP CONC     | HRP Conjugate concentrate          |
| Ag HRP CONC     | HRP Conjugate concentrate          |
| CONJ BUF        | Conjugate buffer                   |
| CHROM TMB CONC  | Chromogenic TMB concentrate        |
| CHROM TMB       | Chromogenic TMB solution           |
| SUB BUF         | Substrate buffer                   |
| STOP SOLN       | Stop solution                      |
| INC SER         | Incubation serum                   |
| BUF             | Buffer                             |
| Ab AP           | AP Conjugate                       |
| SUB PNPP        | Substrate PNPP                     |
| BIOT CONJ CONC  | Biotin conjugate concentrate       |
| AVID HRP CONC   | Avidine HRP concentrate            |
| ASS BUF         | Assay buffer                       |
| Ab BIOT         | Biotin conjugate                   |
| Ab              | Specific Antibody                  |
| SAV HRP CONC    | Streptavidin HRP concentrate       |
| NSB             | Non-specific binding               |
| 2nd Ab          | 2nd Antibody                       |
| ACID BUF        | Acidification Buffer               |
| DIST            | Distributor                        |

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