IGF-1-EASIA

KAP1581

Read entire protocol before use.

IGF1-EASIA

I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative measurement of human Insulin-like Growth Factor-I (IGF-I) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource IGF1-EASIA Kit
- B. **Catalogue number :** KAP1581 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

A. Biological activities

Insulin-like growth factor I (IGF-I) or Somatomedin-C (SM-C) is a basic 70 amino acid single chain polypeptide (MW : 7649 Da) similar to proinsulin (50% sequence homology), and to the other well-characterized member of the somatomedin family : IGF II (67AA, 70 % sequence homology with IGF-I). IGF-I is the most important factor, which mediates the growth promoting actions of growth hormone, a pituitary hormone with highly fluctuating blood levels due to pulsatile release. The blood concentration of IGF-I is more stable due to the binding to carrier proteins. The concentration of the predominant binding protein (MW 53000) as well as the production of IGF-I, are regulated by growth hormone. IGF-I is produced by the liver, and other tissues, and it has endocrine, paracrine and autocrine activities. It stimulates growth and regulates differentiation of various tissues, displays insulin-like activities and promotes cartilage growth. Although GH is the most important factor controlling IGF-I secretion and concentration, other factors are also determinant: the age (with a peak at adolescence), the sex, the nutritional status, and other hormones (oestrogen, thyroxine, prolactin, ...). Specific trophic stimuli mainly control IGF-I secretion in the local microenvironment of a particular organ (paracrine activities), while blood IGF-I concentration is the most important variable for balanced systemic growth (endocrine activities).

B. Clinical applications

· **Growth retardation:** Growth retardation may be due to several causes, among which deficient GH production (hypopituitarism), which is associated with low IGF-I blood levels. Because of the difficulties to get interpretable results from GH measurements (by dynamic multiple or stimulation tests), the determination of the stable IGF-I concentration in plasma is often considered as a simple screening test to evaluate "GH impregnation" of the patient before deciding more extensive investigations. In several clinical situations with impaired growth, low IGF-I levels may be observed despite normal or high GH production (i.e. malnutrition, chronic diseases states, some genetic dwarfs like Pygmies, ...). Interestingly, children with discrete GH neuro-secretory dysfunction may display low IGF-I values despite normal GH levels by conventional testing. The results of IGF-I assay must be interpreted cautiously by considering the normal variations of IGF-I during childhood and adolescence (see Rosenfeld et al).

· **Acromegaly:** IGF-I levels are elevated in acromegaly (excess production of GH) and may serve as an indicator of disease severity. Results are more readily interpreted because the normal values are more easily defined in adults. IGF-I measurements are also useful to monitor treatment.

· **Research:** The IGF-I EASIA kit is an invaluable tool to study the modifications of this growth factor during physiologic (i.e. pregnancy) or pathologic (i.e. diabetes) situations, and the local regulation of IGF-I production in relation to its paracrine and autocrine activities (wound healing, organ regeneration, neoplastic growth, foetal development, gonadal regulation, etc).

IV. PRINCIPLES OF THE METHOD

The DIAsource IGF1-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on breakable microtiterplates. In the present kit, DIAsource has introduced a pre-treatment step in order to improve the clinical performance of the assay. It is well established that the binding proteins interfere with the EASIA assay for IGF1. The pre-treatment step used by DIAsource is the acid-ethanol procedure of Daughaday et al.(8) A fixed amount of IGF1-labelled with horseradish peroxidase (HRP), compete with unlabelled IGF1 present in the calibrators ,controls and samples for a limited number of binding sites on a specific antibody. After 1 hour incubation at room temperature, the microtiterplate is washed to stop the competition reaction. . The Chromogenic solution (TMB) is added and incubated for 30 min. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is inversely proportional to the IGF1 concentration. A calibration curve is plotted and IGF1 concentration in samples is determined by interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution			
 Microtiterplate with 96 anti IGF1 coated breakable wells	96 wells	blue	Ready for use			
<table><tr><td>Ag</td><td>HRP</td><td>CONC</td></tr></table> Conjugate: HRP labelled IGF1 in phosphate buffer with caseine and thymol	Ag	HRP	CONC	1 vial 0.2 ml	red	Dilute 100 x with conjugate buffer
Ag	HRP	CONC				
<table><tr><td>CONJ</td><td>BUF</td></tr></table> Conjugate buffer: Phosphate buffer with caseine and thymol	CONJ	BUF	1 vial 11,5 ml	red	Ready for use	
CONJ	BUF					
<table><tr><td>PRE</td><td>SOLN</td></tr></table> Pre-treatment Solution containing ethanol	PRE	SOLN	1 vial 20 ml	black	Ready for use	
PRE	SOLN					
<table><tr><td>REC</td><td>SOLN</td></tr></table> Reconstitution Solution containing ethanol	REC	SOLN	1 vial 10 ml	blue	Ready for use	
REC	SOLN					
<table><tr><td>NEUTR</td><td>SOLN</td></tr></table> Neutralization Solution containing phosphate buffer with bovine casein and azide (<0.1%)	NEUTR	SOLN	1 vial 30 ml	green	Ready for use	
NEUTR	SOLN					
<table><tr><td>CAL</td><td>0</td></tr></table> Zero calibrator in phosphate buffer with bovine casein	CAL	0	1 vial lyophilized	yellow	Add 1.0 ml reconstitution solution	
CAL	0					
<table><tr><td>CAL</td><td>N</td></tr></table> Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein .	CAL	N	5 vials lyophilized	yellow	Add 1.0 ml reconstitution solution	
CAL	N					
<table><tr><td>CONTROL</td><td>N</td></tr></table> Controls - N = 1 or 2 in human serum with thymol	CONTROL	N	2 vials lyophilized	silver	Add 0.5 ml distilled water	
CONTROL	N					
<table><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table> Wash Solution (Tris-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table><tr><td>CHROM</td><td>TMB</td></tr></table> Chromogenic TMB Solution (Tetramethylbenzidine)	CHROM	TMB	1 vial 12 ml	brown	Ready for use	
CHROM	TMB					
<table><tr><td>STOP</td><td>SOLN</td></tr></table> Stop solution: HCl 1.0 N	STOP	SOLN	1 vial 12 ml	white	Ready for use	
STOP	SOLN					

Note: 1. Use the zero calibrator for sample dilutions.
2. 1 ng of the calibrator preparation is equivalent to 1 ng of the NIBSC 1st IRR 87/518.

VI. SUPPLIES NOT PROVIDED

- The following material is required but not provided in the kit:
1. High quality distilled water
 2. Pipettes for delivery of: 50 µl, 100 µl, 400 µl, 500 µl, 600 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
 3. Optional: microcentrifuge tubes in polypropylene for pre-treatment of samples (for example:Eppendorf 3810)
 4. Polypropylene tubes for neutralization step (for example:Falcon 352063)
 5. Vortex mixer
 6. Magnetic stirrer
 7. Microcentrifuge or Centrifuge
 8. Washer for microtiterplates
 9. Microtiterplate reader capable of reading at 450 nm and 650 nm (monochromatic reading)

VII. REAGENT PREPARATION

- A. Calibrators:** Reconstitute the zero calibrator and other calibrators with 1.0ml reconstitution solution.
B. Controls : Reconstitute the controls with 0.5 ml distilled water.
C Working IGF1-HRP conjugate : Prepare an adequate volume of conjugate solution by adding for example :20 µl of the 100 x concentrated IGF1-HRP conjugate to 2 ml of conjugate buffer. Use a vortex to homogenize. Extemporaneous preparation is recommended.
D. Working Wash solution : Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- § Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
§ Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
§ After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing
§ The concentrated Wash Solution is stable at room temperature until expiration date.
§ Freshly prepared Working Wash solution should be used on the same day..
§ The Working IGF1-HRP conjugate is stable for 4 hours at room temperature, avoid direct sunlight.
§ Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- § Serum must be kept at 2 - 8°C.
§ If the test is not run within 24 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
§ Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
§ Do not use haemolysed samples.

X. PROCEDURE

- A. Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.
During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

B.1. Pre-treatment step for microfuge tubes
! Do not extract calibrators

- 1. Label one microcentrifuge tube(for extraction) and one polypropylene tube (for neutralization) for each sample and control.
- 2. Dispense 100 µl of each sample and control into microfuge tube.
- 3. Add 400 µl of pre-treatment solution into this tube.
- 4. Close the tubes ,vortex and incubate 30 minutes at room temperature.
- 5. Centrifuge for 2 minutes at 10000 rpm.
- 6. Take 100 µl of the supernatant and transfer it into the polypropylene labelled tube.
- 7. Add 600 µl of the neutralisation solution to this tube.
- 8. Vortex each tube.

B.2. Pre-treatment step for other polypropylene tubes
! Do not extract calibrators

- 1. Label two polypropylene tubes(one fro extraction and the other for neutralization) for each sample and control.
- 2. Dispense 100 µl of each sample and control into extraction tube.
- 3. Add 400 µl of pre-treatment solution into this tube.
- 4. Close the tubes ,vortex and incubate 30 minutes at room temperature.
- 5. Centrifuge for 15 minutes at > or = 3000 rpm.
- 6. Take 100 µl of the supernatant and transfer it into the polypropylene labelled tube.
- 7. Add 600 µl of the neutralisation solution to this tube.
- 8. Vortex each tube.

C. Procedure

- 1. Select the required number of wells for the run. The unused wells should be resealed in the bag with a desiccant and stored at 2-8°C.
- 2. Secure the wells into the holding frame.
- 3. Pipette 50 µl of each Calibrator, extracted control and extracted sample into the appropriate wells.
- 4. Pipette 100 µl of IGF1-HRP conjugate solution into all the wells.
- 5. Incubate for 1 hour at room temperature .
- 6. Aspirate the liquid from each well.
- 7. Wash the plate 3 times by:
§ dispensing 0.4 ml of Wash Solution into each well
§ aspirating the content of each well
- 8. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
- 9. Incubate the microtiterplate for 30 minutes at room temperature, avoid direct sunlight
- 10. Pipette 100 µl of Stop Solution into each well.
- 11. Read the absorbences at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

XI. CALCULATION OF RESULTS

- 1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
- 2. Calculate the mean of duplicate determinations.
- 3. Calculate for each calibrator, control and sample:

$$B/B0(\%) = \frac{OD \text{ (Calibrator, Control or Sample)}}{OD \text{ (Zero Calibrator)}} \times 100$$

- 4. Using either linear-linear of semi-logarithmic graph paper, plot the (B/B0(%)) values for each calibrator point as a function of the IGF1 concentration of each calibrator point. Reject obvious outliers.
- 5. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

By interpolation of the sample (B/B0 (%)) values, determine the IGF1 concentrations of the samples from the calibration curve

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

IGF1-EASIA		OD units
Calibrator		
	0 ng/ml	2.33
	75.2 ng/ml	1.71
	190 ng/ml	1.18
	454 ng/ml	0.76
	1253 ng/ml	0.38
	3132 ng/ml	0.23

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Sixteen zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average OD at zero binding, was 1.1 ng/ml

B. Specificity

The percentages of cross-reaction estimated by comparison of the concentration yielding a 50% inhibition are respectively:

Compound	Cross-Reactivity (%)
IGF-I	100.0
IGF-II	0.7
Insuline	ND
GH	ND

Note: this table shows the cross-reactivity for the anti IGF-I.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A B	10	198.2±11.2	5.6	A B	9	193.3±22.27	11.5
	10	511..2 ± 21.5	4.2		9	466.15±27.2	5.2

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added IGF1 (ng/ml)	Recovered IGF1 (ng/ml)	Recovery (%)
Serum 1	750	897	120
	375	400	106
Serum 2	750	808	107
	375	329	88

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum A	1/1	-	773
	1/2	386	407
	1/4	193	239
	1/8	97	105
Serum B	1/1		603
	1/2	301	318
	1/4	150.5	176
	1/8	75.2	96

Samples were diluted after extraction with zero calibrator

Conversion factor:

From ng/ml to nmol/L: x 0.1307
From nmol/L to ng/ml: x 7.649

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

TIME DELAY		
	0 min	30 min
S1	209.14	191.7
S2	312.5	270.8

XIV. INTERNAL QUALITY CONTROL

- § If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- § If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- § Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- § It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- § It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Age Group	MALES (ng/ml)			FEMALES (ng/ml)		
	Mean	Range	N	Mean	Range	N
0 - 2 years	81	19-241	15	70	15-217	11
3 - 5 years	62	15-168	12	78	35-184	13
6 - 8 years	106	46-208	17	121	45-262	13
9 - 11 years	146	32-376	21	215	91-611	10
12 - 14 years	203	75-373	18	281	102-486	14
15 - 17 years	382	63-735	21	327	125-624	20
18 - 20 years	318	78-584	10	203	39 -402	9
21 – 30 years	157	85-240	11	191	86-338	10
31 - 40 years	114	54-190	13	129	75-237	12
41 -50 years	151	79-327	20	126	59-210	15
51 - 60 years	147	51-297	16	120	31-227	15
> 60 years	101	24-186	22	99	53-208	21

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.

2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.

3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D’ERCOLE A.J. (1977)
Estimation of spmatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.

4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.

5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.

6. ALBERTSSON-WIKLAND K. and HALL K. (1987)

Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.

7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.

8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.

9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.

10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.

11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.

12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.

13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

XVIII.**SUMMARY OF THE PROTOCOL**

	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
<u>PRE-TREATMENT</u>		
Samples, controls	-	100
Pre-treatment solution	-	400
Incubation Centrifugation	30 minutes RT 2 minutes at 10000 rpm (microfuge tubes) or 15 minutes at > or = 3000 rpm (polypropylene tubes)	
Supernatant	-	100
Neutralization Solution	-	600
Shaking	Vortex	
<i>INCUBATION</i>		
Calibrators (0 to 5)	50	-
Pre-treated Samples, controls	-	50
Diluted Conjugate	100	100
Incubate for 1 hour at room temperature. Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 30 min at room temperature		
Stop Solution	100	100
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (versus 630 or 650 nm)		

DIAsource Catalogue Nr: KAP1581	P.I. Number : 1700772 / en	Revision nr : 090505/1
------------------------------------	-------------------------------	---------------------------

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XVII. BIBLIOGRAFIA

- DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
- FROESCH E.R. and ZAPP J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
- FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
- ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
- RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
- ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
- WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.
- DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
- BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
- SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
- GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
- TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
- KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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