TNF- α - EASIA

KAP1751

Read entire protocol before use.

TNF- α -EASIA

I. INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human Tumor Necrosis Factor α (TNF- α) in serum.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource TNF- α -EASIA Kit
- B. Catalogue number :** KAP1751 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

A. Biological activities

Human Tumor Necrosis Factor Alpha (TNF- α) also named cachectin, is a 157 A.A. unglycosylated polypeptide cytokine mainly produced by activated macrophages (monocytes). Lipopolysaccharide (LPS), the cell-wall component of gram-negative bacteria (endotoxin), is a potent stimulus for TNF- α production by macrophages and TNF- α is an important mediator of the well-known in vivo effects of LPS such as tumour hemorrhagic necrosis, fever, shock and activation of neutrophils. The various biological activities of TNF- α may be classified as :

- *Antitumoral and growth regulatory activities* : TNF- α displays a selective toxicity for tumor and virus-infected cells. Conversely, it is angiogenic and stimulates the growth of cultured fibroblasts.
- *Immunomodulatory and proinflammatory activities* : TNF- α activates macrophages, neutrophils and eosinophils, as well as endothelial cells (which display procoagulant activity). It regulates the production of antibodies by B cells and stimulates cytotoxic T cells. It induces the production of several other inflammatory mediators such as IL-1, IL-6, colony stimulating factors, prostaglandins, platelet-activating factor (PAF), collagenases, etc.
- *Metabolic activities* : TNF- α strongly inhibits lipoprotein lipase and adipocyte gene expression.

B. Clinical application

TNF- α has a major pathogenic role : in cachexia associated with chronic infectious or cancerous diseases ; in septic shock where the neutralization of TNF- α protects against the associated acute lethality ; in graft rejection and graft-versus-host disease ; and in parasitic infections where TNF- α may provide some protection but also favours more severe forms of the disease (e.g. the cerebral form of malaria). TNF- α often in combination with other cytokines, has also been involved in several autoimmune diseases and even in the pathogenesis of arteriosclerosis. Abnormal high levels of serum TNF- α have been described in septic shock, graft rejection, parasitic infections, cancer, post hemofiltrations, during in vivo cytokine (IL-2) therapy, etc. Besides an insight into pathogenesis, these determinations might provide an aid in diagnosis (e.g. in graft rejection) and have prognostic value (e.g. in systemic infections).

IV. PRINCIPLES OF THE METHOD

The DIAsource TNF- α -EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. The assay uses monoclonal antibodies (MABs) directed against distinct epitopes of TNF- α . Calibrators and samples react with the capture monoclonal antibody (MAB 1) coated on microtiter well and with a monoclonal antibody (MAB 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAB 1 – human TNF- α – MAB 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the TNF- α concentration.

A calibration curve is plotted and TNF- α concentration in samples is determined by interpolation from the calibration curve. The use of the EASIA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution			
 Microtiterplate with 96 anti TNF- α (monoclonal antibodies) coated wells	96 wells	blue	Ready for use			
<table border="1"><tr><td>Ab</td><td>HRP</td></tr></table> Conjugate: HRP labelled anti-TNF- α (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin and thymol	Ab	HRP	1 vial 0.75 ml	red	Add conjugate buffer (see section VII)	
Ab	HRP					
<table border="1"><tr><td>CAL</td><td>0</td></tr></table> Zero calibrator in human serum, benzamadin and thymol	CAL	0	2 vials lyophil.	yellow	Add distilled water (see on the label for the exact volume)	
CAL	0					
<table border="1"><tr><td>CAL</td><td>N</td></tr></table> Calibrator N = 1 to 5 (see exact values on vial labels) in human serum, benzamadin and thymol	CAL	N	5 vials lyophil.	yellow	Add 2 ml distilled water	
CAL	N					
<table border="1"><tr><td>CONJ</td><td>BUF</td></tr></table> Conjugate buffer: TRIS-Maleate buffer with bovine serum albumin, EDTA and thymol	CONJ	BUF	1 vial 6 ml	red	Ready for use	
CONJ	BUF					
<table border="1"><tr><td>INC</td><td>BUF</td></tr></table> Incubation buffer: TRIS-Maleate buffer with bovine serum albumin, EDTA and thymol	INC	BUF	1 vial 6 ml	black	Ready for use	
INC	BUF					
<table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table> Wash Solution (Tris-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1"><tr><td>CONTROL</td><td>N</td></tr></table> Controls - N = 1 or 2 in human serum and thymol	CONTROL	N	2 vials lyophil.	silver	Add 2 ml distilled water	
CONTROL	N					
<table border="1"><tr><td>CHROM</td><td>TMB</td><td>CONC</td></tr></table> Chromogen TMB (Tetramethylbenzidine) in Dimethylformamide	CHROM	TMB	CONC	1 vial 1 ml	green	Dilute 0.2 ml into 1 vial of substrate buffer
CHROM	TMB	CONC				
<table border="1"><tr><td>SUB</td><td>BUF</td></tr></table> Substrate buffer: H ₂ O ₂ in acetate / citrate buffer	SUB	BUF	3 vials 21 ml	white	Ready for use	
SUB	BUF					
<table border="1"><tr><td>STOP</td><td>SOLN</td></tr></table> Stopping solution: H ₂ SO ₄ 1.8N	STOP	SOLN	1 vial 6 ml	black	Ready for use	
STOP	SOLN					

Note: 1. Use the zero calibrator for sample dilutions.
2. 1 pg of the calibrator preparation is equivalent to 40 mIU of the NIBSC IS 87/650.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 μ l, 200 μ l, 1 ml and 10 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Horizontal microtiterplate shaker capable of 700 rpm $\pm$ 100 rpm
- Washer for Microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)
- Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XI.A.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

VII. REAGENT PREPARATION

- A. **Calibrators** : Reconstitute the zero calibrator to the volume specified on the vial label with distilled water and the other calibrators with 2 ml distilled water.
- B. **Controls** : Reconstitute the controls with 2 ml distilled water.
- C. **Conjugate Solution** : following the number of wells to be used, dilute the concentrated conjugate with the conjugate buffer in a clean glass vial : see below table for the volumes to pipette. Extemporaneous preparation is recommended. Diluted conjugate is stable for max. 1 week at 2-8°C.

TABLE CONJUGATE DILUTION			
Number of wells	Concentrated conjugate	Conjugate buffer	Working volume
8	50 μ l	500 μ l	550 μ l
16	100 μ l	1000 μ l	1100 μ l
24	150 μ l	1500 μ l	1650 μ l
32	200 μ l	2000 μ l	2200 μ l
48	300 μ l	3000 μ l	3300 μ l
96	600 μ l	6000 μ l	6600 μ l

- D. **Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.
- E. **Revelation Solution**: pipette 0.2 ml of the chromogen TMB into one of the vials of substrate buffer (H₂O₂ in acetate/citrate buffer). Extemporaneous preparation is recommended.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- § Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- § Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- § After reconstitution, calibrators and controls are stable for 4 days at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 2 months. Avoid successive freeze thaw cycles.
- § The concentrated Wash Solution is stable at room temperature until expiration date.
- § Freshly prepared Working Wash solution should be used on the same day.
- § After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- § The freshly prepared revelation solution is stable, before use, for maximum 15 minutes at room temperature and must be discarded afterwards.
- § Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- § Serum must be removed as soon as possible from the clot of red cells after clotting and centrifugation, and kept at 4°C. If the samples are not used immediately, they must be kept at -20°C for maximum 2 months, and at -70°C for longer storage (maximum one year).
- § Avoid subsequent freeze thaw cycles.
- § Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.

- § Sampling conditions can affect values, therefore, strict precautions have to be taken during sampling to avoid impurities contained in sampling materials that would stimulate TNF- α production by blood cells and thus falsely increase serum TNF- α values.
- § Collection tubes must be pyrogen-free.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Revelation Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
The Revelation Solution should be colourless. If a blue colour develops within a few minutes after preparation, this indicates that the reagent is unusable, and must be discarded.
Dispense the Revelation Solution within 15 minutes following the washing of the microtiterplate.
During incubation with Revelation Solution, avoid direct sunlight on the microtiterplate.

B. Procedure

- Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
- Secure the strips into the holding frame.
- Pipette 50 μ l of incubation buffer into all the wells
- Pipette 200 μ l of each Calibrator, Control and Sample into the appropriate wells.
- Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm $\pm$ 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times by:
 - § Dispensing 0.4 ml of Wash Solution into each well
 - § Aspirating the content of each well
- Pipette 100 μ l of zero calibrator into all the wells
- Pipette 50 μ l of anti- TNF- α -HRP conjugate into all the wells.
- Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm $\pm$ 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times by:
 - § Dispensing 0.4 ml of Wash Solution into each well
 - § Aspirating the content of each well
- Pipette 200 μ l of the freshly prepared revelation solution into each well within 15 minutes following the washing step.
- Incubate the microtiterplate for 30 minutes at room temperature on a horizontal shaker set at 700 rpm $\pm$ 100 rpm, avoid direct sunlight.
- Pipette 50 μ l of Stop solution into each well.
- Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 3 hours and calculate the results as described in section XI.

XI. CALCULATION OF RESULTS

A. Polychromatic Reading:

- In this case, the ELISA-AID™ software will do the data processing.
- The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
- A second reading is performed at 490 nm against the same reference filter.
- The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
- The principle of polychromatic data processing is as follows:
 - § X_i = OD at 450 nm
 - § Y_i = OD at 490 nm
 - § Using a standard unweighted linear regression, the parameters A & B are calculated : $Y = A \cdot X + B$
 - § If $X_i < 3$ OD units, then X calculated = X_i

- § If $X_i > 3$ OD units, then X calculated = $(Y_i - B)/A$
- § A 4 parameter logistic curve fitting is used to build up the calibration curve.
- § The TNF- α concentration in samples is determined by interpolation on the calibration curve.

B. Bichromatic Reading

- Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
- Calculate the mean of duplicate determinations.
- On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of TNF- α (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

TNF- α -EASIA		OD units Polychromatic model
Calibrator	0 pg/ml	0.045
	6.8 pg/ml	0.120
	18 pg/ml	0.259
	52 pg/ml	0.619
	176 pg/ml	1.435
	518 pg/ml	3.237

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.7 pg/ml.

B. Specificity

No significant cross-reaction was observed in presence of 50 ng of IL-1 α , IL-1 β , IL-2, IL-3, IL-4, IL-6, IL-7, IL-8, IL-10, TNF- β , IFN- α , IFN- β , IFN- γ , TGF- β , GM-CSF, OSM, MIP-1 α , MIP-1 β , LIF, MCP-1, G-CSF and RANTES. This TNF- α assay is specific for human natural and recombinant TNF- α .

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> $\pm$ SD (pg/ml)	CV (%)	Serum	N	<X> $\pm$ SD (pg/ml)	CV (%)
A	20	91 $\pm$ 6	6.6	A	24	122 $\pm$ 5	4.5
B	20	526 $\pm$ 33	6.3	B	24	431 $\pm$ 14	3.3

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added TNF- α (pg/ml)	Recovered TNF- α (pg/ml)	Recovery (%)
Serum 1	0	6.2	-
	38.4	43.3	97
	83.9	90.0	100
	188.3	192.5	99
	408.2	376.2	91
Serum 2	0	3.8	-
	38.4	45.5	108
	83.9	91.2	104
	188.3	162.2	84
	408.2	379.2	92

DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
Serum 1	1	-	436.5
	2	218.3	212.4
	4	109.1	104.8
	8	54.6	59.5
	16	27.3	31.7
Serum 2	1	-	420.2
	2	210.1	211.2
	4	105.0	98
	8	52.5	58.3
	16	26.3	30.7

Samples were diluted with zero calibrator.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

	T0	30 min	45 min
SC1	202	183	222
SC2	506	520	565

XIV. INTERNAL QUALITY CONTROL

- § If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- § If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- § Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices
- § It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- § It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

For guidance, the results of 30 serum samples from apparently healthy persons with low CRP levels, ranged between 4.6 and 12.4 pg/ml.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains H₂SO₄, the chromogen contains TMB in Dimethylformamide, Substrate buffer contains H₂O₂. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

- BEUTLER B., CERAMI A. (1987)
Cachectin : more than a tumor necrosis factor.
N. Engl. J. Med., 316 ; 379-385.
- TRACEY K.J., FONG Y., HESSE D.G., MANOGUE K.R., LEE A.T., KUO G.C., LOWRY S.F. and CERAMI A. (1987)
Anti cachectin/TNF monoclonal antibodies prevent septic shock during lethal bacteraemia.
Nature, 330 : 662-664.
- PIGUET P.F., GRAU G.E., ALLET B. and VASSALLI P. (1987)
Tumor necrosis factor/cachectin is an effector of skin and gut lesions of the acute phase of graft-versus-host disease.
J. Exp. Med., 166 ; 1280-1289.
- AUKRUST P., LIABAKK N-B., MÜLLER F., LIEN E., ESPEVIK T. and FROLAND S.S. (1994)
Serum Levels of Tumor Necrosis Factor- α (TNF α) and Soluble TNF Receptors in Human Immunodeficiency Virus Type 1 Infection - Correlations to Clinical Immunologic, and Virologic Parameters.
J. Inf. Dis., 169:420-424.
- WAAGE A., HALSTENSEN A. and ESPEVIK T. (1987)
Association between tumor necrosis factor in serum and fatal outcome in patients with meningococcal disease.
Lancet, 1 ; 355-357.
- LEROUX-ROELS G., OFFNER F., PHILIPPE J. and VERMEULEN A. (1988)
Influence of Blood-Collecting Systems on Concentrations of Tumor Necrosis Factor in Serum and Plasma.
Clin. Chem., 34 ; 2373-2374.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (μ l)	SAMPLE(S) CONTROLS (μ l)
Incubation buffer Calibrators (0-5) Samples, Controls	50 200 -	50 - 200
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 μ l of Wash Solution and aspirate.		
Zero Calibrator Anti-TNF- α -HRP conjugate	100 50	100 50
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 μ l of Wash Solution and aspirate.		
Revelation Solution	200	200
Incubate for 30 min at room temperature with continuous shaking at 700 rpm.		
Stop Solution	50	50
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

DIAsource Catalogue Nr : KAP1751	P.I. Number : 1700571/en	Revision nr : 090508/1
-------------------------------------	-----------------------------	---------------------------

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- BEUTLER B., CERAMI A. (1987)
Cachectin : more than a tumor necrosis factor.
N. Engl. J. Med., 316 ; 379-385.
- TRACEY K.J., FONG Y., HESSE D.G., MANOGUE K.R., LEE A.T., KUO G.C., LOWRY S.F. and CERAMI A. (1987)
Anti cachectin/TNF monoclonal antibodies prevent septic shock during lethal bacteraemia. Nature, 330 ; 662-664.
- PIGUET P.F., GRAU G.E., ALLET B. and VASSALLI P. (1987)
Tumor necrosis factor/cachectin is an effector of skin and gut lesions of the acute phase of graft-versus-host disease.
J. Exp. Med., 166 ; 1280-1289.
- AUKRUST P., LIABAKK N-B., MÜLLER F., LIEN E., ESPEVIK T. and FROLAND S.S. (1994)
Serum Levels of Tumor Necrosis Factor-α (TNFα) and Soluble TNF Receptors in Human Immunodeficiency Virus Type 1 Infection - Correlations to Clinical Immunologic, and Virologic Parameters.
J. Inf. Dis., 169:420-424.
- WAAGE A., HALSTENSEN A. and ESPEVIK T. (1987)
Association between tumor necrosis factor in serum and fatal outcome in patients with meningococcal disease. Lancet, 1 ; 355-357.
- LEROUX-ROELS G., OFFNER F., PHILIPPE J. and VERMEULEN A. (1988)
Influence of Blood-Collecting Systems on Concentrations of Tumor Necrosis Factor in Serum and Plasma. Clin. Chem., 34 ; 2373-2374.

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