sTNF-RII-EASIA

KAP1771

Read entire protocol before use.

sTNF-RII-EASIA

I. INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human soluble Tumor Necrosis Factor Receptor N°. 2 (sTNF-RII) in serum.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource sTNF-RII-EASIA Kit
- B. Catalogue number :** KAP1771 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

A. Biological activities

sTNF-RII, also named TBPII or p75, is one of the two receptors of tumor necrosis factors that are present at the surface of many cells. Different processes modulate their presence. IL-2 and the activation of T-lymphocytes increase the presence of both TNF-RI and RII. An increase is also noticed during the maturation of the macrophages or in presence of protein kinases activators. On the opposite, TNF-Rs decrease in the presence of H₂O₂, epinephrine, Insulin and somatostatin. The two receptors of TNF are able to bind TNF- α and TNF- β . TNF-RII is the form present on monocytes. The gene present in the chromosome 1 codes for a protein of 461 aa. The MW is about 75 kDa ; that suggests an important glycosylation. According to Tartagliata et al, the TNF-RII could be at the origin of the thymocytes and lymphocytes T cytotoxic proliferation. The proliferation of circulating human mononuclear cells under the influence of PHA involves the participation of the two receptors. Soluble forms of the receptors are shedded from the cell membrane and are present in urine, plasma or culture supernatants.

B. Clinical application

The presence of the soluble receptors for TNF has been proved in the serum of cancer patients, chronic renal deficiency and in the broncho-alveolar lavage of patients suffering from ARDS. The soluble receptors for TNF are also putative markers of disease progression in HIV infection and, compared to controls, sTNF-RII is already higher in asymptomatic carriers. sTNF-R correlates also with parasitemia and disease severity in human malaria. These forms bind perfectly the TNF and, in high concentration, inhibit the biological activity of TNF. In some conditions, these soluble forms are able to protect TNF and increase its half live.

IV. PRINCIPLES OF THE METHOD

The DAsource sTNF-RII-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. The assay uses monoclonal antibodies (MAbs) directed against distinct epitopes of sTNF-RII. Calibrators and samples react with the capture monoclonal antibody (MAb 1) coated on microtiter well and with a monoclonal antibody (MAb 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAb 1 – human sTNF-RII – MAb 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the sTNF-RII concentration.

A calibration curve is plotted and sTNF-RII concentration in samples is determined by interpolation from the calibration curve. The use of the EASIA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution			
 Microtiterplate with 96 anti sTNF-RII (monoclonal antibodies)	96 wells	blue	Ready for use			
<table border="1"><tr><td>Ab</td><td>HRP</td></tr></table> Conjugate: HRP labelled anti-sTNF-RII (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin and thymol	Ab	HRP	1 vial 21 ml	red	Ready for use	
Ab	HRP					
<table border="1"><tr><td>CAL</td><td>0</td></tr></table> Zero Calibrator: bovine serum with benzamidin and thymol	CAL	0	1 vial lyophil.	yellow	Add distilled water (see on the label for the exact volume)	
CAL	0					
<table border="1"><tr><td>CAL</td><td>N</td></tr></table> Calibrator N = 1 to 5 (see exact values on vial labels) in bovine serum with benzamidin and thymol	CAL	N	5 vials lyophil.	yellow	Add 0.5 ml distilled water	
CAL	N					
<table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table> Wash Solution (Tris-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1"><tr><td>CONTROL</td><td>N</td></tr></table> Controls - N = 1 or 2 in human serum with benzamidin and thymol	CONTROL	N	2 vials lyophil.	silver	Add 0.5 ml distilled water	
CONTROL	N					
<table border="1"><tr><td>CHROM</td><td>TMB</td></tr></table> Chromogenic TMB Solution	CHROM	TMB	1 vial 25 ml	white	Ready for use	
CHROM	TMB					
<table border="1"><tr><td>STOP</td><td>SOLN</td></tr></table> Stop Solution: HCl 2N	STOP	SOLN	1 vial 25 ml	white	Ready for use	
STOP	SOLN					

Note: 1. Use zero calibrator for sample dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 200 µl, 1 ml and 10 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Horizontal microtiterplate shaker capable of 700 rpm ± 100 rpm
- Washer for Microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)

- Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XI.A.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

VII. REAGENT PREPARATION

- Calibrators:** Reconstitute the zero calibrator to the volume specified on the vial label with distilled water and the other calibrators with 0.5 ml distilled water.
- Controls:** Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators and controls are stable for 4 days at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 2 months. Avoid successive freeze thaw cycles.
- The concentrated Wash Solution is stable at room temperature until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum must be removed as soon as possible from the clot of red cells after clotting and centrifugation, and kept at 4°C. If the samples are not used immediately, they must be kept at -20°C for maximum 2 months, and at -70°C for longer storage (maximum one year).
- Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- Sampling conditions can affect values, therefore, strict precautions have to be taken during sampling to avoid impurities contained in sampling materials that would stimulate sTNF-RII production by blood cells and thus falsely increase plasma sTNF-RII values.
- Collection tubes must be pyrogen-free.

X. PROCEDURE

- Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.
During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.
- Procedure**
 - Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
 - Secure the strips into the holding frame.
 - Pipette 50 µl of each Calibrator, Control and Sample into the appropriate wells.

- Pipette 200 µl of anti-sTNF-RII-HRP conjugate into all the wells.
- Incubate for 1 hour at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times by:
 - § Dispensing 0.4 ml of Wash Solution into each well
 - § Aspirating the content of each well
- Pipette 50 µl of the Chromogenic Solution into each well within 15 minutes following the washing step.
- Incubate the microtiterplate for 15 minutes at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm, avoid direct sunlight.
- Pipette 200 µl of Stop Solution into each well.
- Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 3 hours and calculate the results as described in section XI.

XI. CALCULATION OF RESULTS

- A. Polychromatic Reading:**
- In this case, the ELISA-AID™ software will do the data processing.
 - The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
 - A second reading is performed at 490 nm against the same reference filter.
 - The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
 - The principle of polychromatic data processing is as follows:
 - § $X_i = OD$ at 450 nm
 - § $Y_i = OD$ at 490 nm
 - § Using a standard unweighted linear regression, the parameters A & B are calculated : $Y = A \cdot X + B$
 - § If $X_i < 3$ OD units, then X calculated = X_i
 - § If $X_i > 3$ OD units, then X calculated = $(Y_i - B)/A$
 - § A 4-parameter logistic curve fitting is used to build up the calibration curve.
 - § The sTNF-RII concentration in samples is determined by interpolation on the calibration curve.

- B. Bichromatic Reading**
- Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
 - Calculate the mean of duplicate determinations.
On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of sTNF-RII (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
 - Read the concentration for each control and sample by interpolation on the calibration curve.
 - Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

sTNF-RII-EASIA		OD units
		Polychromatic model
Calibrator	0 ng/ml	0.028
	1.9 ng/ml	0.141
	3.9 ng/ml	0.287
	17 ng/ml	0.926
	48 ng/ml	1.840
	142 ng/ml	3.407

XIII. PERFORMANCE AND LIMITATIONS

- A. Detection Limit**
- Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.04 ng/ml.

B. Specificity

Cross-reactivity and interference were analysed by the addition of different analytes to sTNF-RII samples and measuring the apparent sTNF-RII concentration.

Analyte	sTNF-RII 2.86 ng/ml (Ratio [*])	sTNF-RII 8.85 ng/ml (Ratio [*])	sTNF-RII 28.76 ng/ml (Ratio [*])
+ sTNF-RI (p55) (1500 ng/ml)	3.0 ng/ml (104.9 %)	9.08 ng/ml (102.6 %)	29.48 ng/ml (102.5 %)
Analyte	sTNF-RII 0.80 ng/ml (Ratio [*])	sTNF-RII 4.36 ng/ml (Ratio [*])	sTNF-RII 27.66 ng/ml (Ratio [*])
+ TNF-α (500 ng/ml)	0.75 ng/ml (93.8 %)	3.83 ng/ml (87.8 %)	23.63 ng/ml (85.4 %)

Ratio^{*} : $\frac{\text{sTNF-RII measured in the presence of analyte} \times 100}{\text{sTNF-RII added in the absence of analyte}}$

This demonstrates that the sTNF-RII EASIA does not cross react with sTNF-RI and that TNF-α does not interfere with the assay.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A	20	3.03 ± 0.06	2.1	A	20	5.2 ± 0.7	13.5
B	9	18.5 ± 0.7	3.7	B	20	16.2 ± 0.8	4.9

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added sTNF-RII (ng/ml)	Recovered sTNF-RII (ng/ml)	Recovery (%)
Serum	0	4.2	-
	3.9	8.4	107.7
	15.4	20.7	106.5
	34.9	39.1	99.7

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum	1/1	-	39.1
	1/2	19.5	19.3
	1/4	9.8	9.6
	1/8	4.9	5.1
	1/16	2.4	2.6

Samples were diluted with zero calibrator.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

TIME DELAY			
	T0	15 min	30 min
S1	4.7	5.3	5.1
S2	3.6	3.9	3.9
S3	4.3	4.3	3.9
S4	4.6	4.9	5.4

F. Hook effect

A sample spiked with sTNF-RII up to 3000 ng/ml gives higher OD's than the last calibrator point.

XIV. INTERNAL QUALITY CONTROL

- § If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- § If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- § Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices
- § It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- § It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

For guidance, the results of 36 serum samples from apparently healthy persons with low CRP levels, ranged between 0 and 60.2 ng/ml. 35 of the 36 samples obtained values below 20 ng/ml).

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. DIETRICH M. and PORZSOLT F. (1992)
Soluble tumor necrosis factor receptors correlate with parasitemia and disease severity in human malaria.
J. Infec. Dis., 166 : 930-934.

2. GODFRIED M.H., VAN DER POLL T., JANSEN J., ROMIJN J.A., EEFTINCK-SCHATTENKERK J.K.M., ENDERT E., VAN DEVENTER S.J.H. and SAUERWEIN H.P. (1993)
Soluble receptors for tumor necrosis factor : a putative marker of disease progression in HIV infection.
AIDS, 7 : 33-36.
3. RIMANIOL A.C., BOUSSIN F., HERBELIN A., DE GROOTE D., DORMONT D., BACH J.F., DESCAMPS-LATSCHA B. and ZAVALA F. (1994)
Induction of tumor necrosis factor soluble receptor (TNF-sR75) release by HIV adsorption on cultured human monocytes.
Eur. J. Immunol., 24 : 2055-2060.
4. HERBELIN A., CHATENAUD L., ROUX-LOMBARD P., DE GROOTE D., LEGENDRE C., DAYER J.M., DESCAMPS-LATSCHA B., KREIS H. and BACH J.F. (1995)
In vivo soluble tumor necrosis factor receptor release in OKT3-treated patients.
Transplantation, 59 : 1470-1475.
5. DESCAMPS-LATSCHA B., HERBELIN A., NGUYEN A.T., ROUX-LOMBARD P., ZINGRAFF J., MOYNOT A., VERGER C., DAHMANE D., DE GROOTE D., JUNGERS P. and DAYER J.M. (1995)
Balance between IL-1 β , TNF- α and their specific inhibitors in Chronic Renal Failure and maintenance dialysis : relationship with activation markers of T-cells, B-cells and monocytes.
The J. of Immunol., 154 : 882-892.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (μ l)	SAMPLE(S) CONTROLS (μ l)
Calibrators (0-5) Samples, Controls Anti-sTNF-RII -HRP conjugate	50 - 200	- 50 200
Incubate for 1 hour at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 μ l of Wash Solution and aspirate.		
Chromogenic Solution	50	50
Incubate for 15 min at room temperature with continuous shaking at 700 rpm.		
Stop Solution	200	200
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

DIAsource Catalogue Nr : KAP1771	P.I. Number : 1700601/en	Revision nr : 090506/1
-------------------------------------	-----------------------------	---------------------------

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2. GODFRIED M.H., VAN DER POLL T., JANSEN J., ROMIJN J.A., EEFTINCK-SCHATTENKERK J.K.M., ENDERT E., VAN DEVENTER S.J.H. and SAUERWEIN H.P. (1993)
Soluble receptors for tumor necrosis factor : a putative marker of disease progression in HIV infection.
AIDS, 7 : 33-36.
3. RIMANIOLO A.C., BOUSSIN F., HERBELIN A., DE GROOTE D., DORMONT D., BACH J.F., DESCAMPS-LATSCHA B. and ZAVALA F. (1994)
Induction of tumor necrosis factor soluble receptor (TNF-sR75) release by HIV adsorption on cultured human monocytes.
Eur. J. Immunol., 24 : 2055-2060.
4. HERBELIN A., CHATENAUD L., ROUX-LOMBARD P., DE GROOTE D., LEGENDRE C., DAYER J.M., DESCAMPS-LATSCHA B., KREIS H. and BACH J.F. (1995)
In vivo soluble tumor necrosis factor receptor release in OKT3-treated patients.
Transplantation, 59 : 1470-1475.
5. DESCAMPS-LATSCHA B., HERBELIN A., NGUYEN A.T., ROUX-LOMBARD P., ZINGRAFF J., MOYNOT A., VERGER C., DAHMANE D., DE GROOTE D., JUNGERS P. and DAYER J.M. (1995)
Balance between IL-1 β , TNF- α and their specific inhibitors in Chronic Renal Failure and maintenance dialysis : relationship with activation markers of T-cells, B-cells and monocytes.
The J. of Immunol., 154 : 882-892.

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- GODFRIED M.H., VAN DER POLL T., JANSEN J., ROMIJN J.A., EEFTINCK-SCHATTENKERK J.K.M., ENDERT E., VAN DEVENTER S.J.H. and SAUERWEIN H.P. (1993)
Soluble receptors for tumor necrosis factor : a putative marker of disease progression in HIV infection.
AIDS, 7 : 33-36.
- RIMANIOL A.C., BOUSSIN F., HERBELIN A., DE GROOTE D., DORMONT D., BACH J.F., DESCAMPS-LATSCHA B. and ZAVALA F. (1994)
Induction of tumor necrosis factor soluble receptor (TNF-sR75) release by HIV adsorption on cultured human monocytes.
Eur. J. Immunol., 24 : 2055-2060.
- HERBELIN A., CHATENAUD L., ROUX-LOMBARD P., DE GROOTE D., LEGENDRE C., DAYER J.M., DESCAMPS-LATSCHA B., KREIS H. and BACH J.F. (1995)
In vivo soluble tumor necrosis factor receptor release in OKT3-treated patients.
Transplantation, 59 : 1470-1475.
- DESCAMPS-LATSCHA B., HERBELIN A., NGUYEN A.T., ROUX-LOMBARD P., ZINGRAFF J., MOYNOT A., VERGER C., DAHMANE D., DE GROOTE D., JUNGERS P. and DAYER J.M. (1995)
Balance between IL-1β, TNF-α and their specific inhibitors in Chronic Renal Failure and maintenance dialysis : relationship with activation markers of T-cells, B-cells and monocytes.
The J. of Immunol., 154 : 882-892.

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