AFP-IRMA

KIP0081 - KIP0084

Read entire protocol before use.

AFP-IRMA

I. INTENDED USE

Immunoradiometric assay for the *in vitro* quantitative measurement of human Alpha – Fetoprotein (AFP) in serum and plasma.

This kit is not intended to be used for the risk evaluation of trisomy 21.

II. GENERAL INFORMATION

A. **Proprietary name :** DIAsource AFP-IRMA

B. **Catalog number :** KIP0081 : 96 tests
KIP0084 : 4 x 96 tests

C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

A. Introduction

α -Fetoprotein (AFP) is a 70.000 Da MW oncofetal protein synthesized by liverparenchymal cells, yolk sac and gastrointestinal tract of human fetus. The peak of AFP concentration occurs between weeks 12 and 15 of gestation. After birth AFP concentration in plasma rapidly decreases to less than 5 IU/ml. AFP levels are elevated in the following clinical situation:

B. Clinical applications

The main clinical applications of measurements of AFP are found in the monitoring of cancer following treatment. However, AFP measurement may also be of clinical interest in monitoring of pregnancy when applied to serum or amniotic fluid.

- **Cancer**
 - . Hepatocellular carcinoma
 - . Teratocarcinomas and embryonal cell carcinoma of testis and ovaries
 - . Yolk sac tumor
 - . Other cancers (less than 5 %).
- **Viral diseases**
 - . Acute hepatitis (usually < 100 IU/ml)
 - . Chronic active hepatitis (usually < 100 IU/ml).

IV. PRINCIPLES OF THE METHOD

DIAsource AFP-IRMA is a two-step immunoradiometric assay based on coated tube separation. Mabs 1 - the capture antibodies - are attached to the lower and inner surface of the plastic tube. Calibrators or samples are added to the tubes. After incubation, washing removes the occasional excess of antigen. Mab 2 - the ¹²⁵I-labelled-antibody is added. After incubation and washing, the remaining radioactivity bound to the tube reflects the antigen concentration. Because of the occasional extremely high concentration of AFP in serum (> 1 mg/ml) in case of cancer, the present DIAsource AFP-IRMA has been developed on a two steps procedure to avoid high dose hook effects.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	4 x 96 Test Kit	Colour Code	Reconstitution
 Tubes coated with anti AFP (monoclonal antibodies)	2 x 48	8 x 48	dark red	Ready for use
Ab^{125I} TRACER: ¹²⁵ Iodine labelled anti - AFP (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide (0.1%) and an inert red dye	1 vial 10.5 ml 350 kBq	4 vials 10.5 ml 4x350 kBq	red	Ready for use
CAL0 Zero calibrator (0 IU/ml) in bovine serum with azide (0.5%) and merthiolate	1 vial 10 ml	2 vials 10 ml	yellow	Ready for use
CALN Calibrators AFP N = 1 to 5 (see exact values on vial labels) in bovine serum with thymol	5 vials lyophil.	2 x 5 vials lyophil.	yellow	Add 0.5 ml distilled water
WASH SOLN CONC Wash solution (TRIS-HCl)	1 vial 10 ml	4 vials 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
CONTROLN Controls - N = 1 or 2 in human serum with merthiolate	2 vials lyophil.	2 x 2 vials lyophil.	silver	Add 0.5 ml distilled water

Note : Use the zero calibrator for sample dilutions.

1 IU of the calibrator preparation is equivalent to 1 IU of the 1st IRP 72/225

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 100 µl, 200 µl, 500 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. 5 ml automatic syringe (Cornwall type) for washing
6. Aspiration system (optional)
7. Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators** : Reconstitute the calibrators 1-5 with 0.5 ml distilled water.
- Controls** : Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months.
- Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
 - If the test is not run within 24 hrs, storage in aliquots at -20°C is recommended.
 - Avoid successive freezing and thawing.
 - Serum or plasma (EDTA or Heparin) provide similar results.
- $$y(\text{serum}) = 1.05 \times (\text{hep. plasma}) - 0.2 \quad r = 0.98 \quad n = 31$$
- $$y(\text{serum}) = 1.05 \times (\text{EDTA plasma}) - 0.5 \quad r = 0.99 \quad n = 32$$

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes
2. Briefly vortex calibrators, controls and samples and dispense 50µl of each into respective tubes.
3. Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
4. Incubate for 1 hour at room temperature.
5. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
6. Wash tubes twice with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
7. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
8. Dispense 100 µl of anti-AFP-¹²⁵I tracer into each tube, including the uncoated tubes for total counts.
9. Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
10. Incubate for 1 hour at room temperature.
11. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
12. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
13. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
14. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of AFP (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

AFP-IRMA		cpm	B/T (%)
Total count		133479	
Calibrator	0.0 IU/ml	268	0.2
	5.3 IU/ml	2821	2.11
	15.0 IU/ml	6076	4.55
	45.0 IU/ml	19122	14.33
	155.0 IU/ml	44001	32.96
	440.0 IU/ml	66466	49.79

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.5 IU/ml.

B. Specificity

The AFP-IRMA does not cross-react with structurally related substances. Human serum albumin (HSA) at a concentration of 60 ng/ml was added to a low (5 IU/ml) and a high (50 IU/ml) AFP Calibrator, the apparent AFP response was measured.

	Observed AFP value IU/ml
Cal 5 IU/ml + HSA (60 ng/ml)	4.8
Cal 50 IU/ml + HSA (60 ng/ml)	49.0

Note : this table shows the cross-reactivity for the AFP-IRMA.

C. Precision

INTRA-ASSAY				INTER-ASSAY			
Serum	N	<X> ± SD (IU/ml)	CV (%)	Serum	N	<X> ± SD (IU/ml)	CV (%)
A	10	4.2 ± 0.2	4.8%	E	22	7.8 ± 0.5	5.8%
B	10	11.4 ± 0.2	1.8%	F	22	19.8 ± 1.2	6.0%
C	10	26.2 ± 0.6	2.3%				
D	10	357 ± 19.6	5.5%				

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST			
Sample	Dilution	Theoretical Concent. (IU/ml)	Measured Concent. (IU/ml)
Serum 1	1/1	-	444
	1/3	148.0	158
	1/6	74.0	86
	1/10	44.4	42
	1/50	8.9	9.8
Serum 2	1/1	-	87.4
	1/2	43.7	39.6
	1/4	21.9	21.8
	1/10	8.7	9.0
	1/20	4.4	4.2

Samples were diluted with the zero calibrator.

RECOVERY TEST			
Sample	added AFP (IU/ml)	Recovered AFP (IU/ml)	Recovered (%)
Serum 1	7.5	9.6	128%
	22.5	23.2	103%
	80	82.2	103%
Serum 2	2.5	2.4	96%
	22.5	21.4	95%
	225	218	97%

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to coated tubes.

TIME DELAY			
Serum (ng/ml)	0'	20'	30'
Serum 1	12.7	11.0	10.1
Serum 2	31.3	29.4	27.3

F. Hook effect

A sample spiked with an AFP concentration up to 1460000 IU/ml gives higher counts than the last calibrator point.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophelic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.


XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XVI. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

A. Healthy non-pregnant individuals

B. Pregnancy (*)

Maternal serum


The distribution of serum values in normal pregnancy. The upper normal limits for serum samples were obtained by expressing the data as the median multiplied by 2.5.

Week of gestation	Lower values 2.5 Percentile (IU AFP/ml)	Median (IU AFP/ml)	Higher values 97.5 Percentile (IU AFP/ml)	2.5 Median (IU AFP/ml)
15	11.1	21	59	52.5
16	15.1	28	73	70
17	19.0	36	84	90
18	23	42	93	105
19	26	50	104	125
20	30	58	113	145

(*) Studies performed for DIAsource ImmunoAssays S.A., Belgium based on samples from Belgium.

C. Oncology

AFP concentrations were measured in serum during disease- monitoring in three individual cases of cancer of different cell type origin.

XVII. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.
This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.
This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.
All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.
Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.
The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.
All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.
Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.
Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. CUCKLE, H.S., WALD, N.J. 1984.
Maternal serum alpha-fetoprotein measurement : a screening test for Down Syndrome.
Lancet, i:926
2. HUNTER, W.M. BUDD, P.S. 1981.
Immunoradiometric versusradioimmunoassay a comparison using alpha-fetoprotein as the model analyte.
J. of Immunol. Meth.45:255

3. KOHN, J. WEAVER, P.C. 1974.
Serum alpha-fetoprotein in hepatocellular carcinoma.
Lancet, i:334
4. KOHN, J., ORR, A.H. MC ELWAIN, T.J. et al.
Serum alpha-fetoprotein in patients with testicular tumours.
Lancet; II:433
5. NORGAARD-PENDERSEN, B. 1976.
Human alpha-fetoprotein - A review of recent methodological and clinical studies.
Scand. J. Immunol., Suppl. 4:7.
6. U.K. Collaborative study, First report. 1977.
Maternal serum alpha-fetoprotein measurement in antenatal screening for anencephaly and spina bifidein early pregnancy.
Lancet, i:1323.
7. JEVY, L.B. et al. 1995
Hepatocellular carcinoma in a pregnant woman detected by routine screening of maternal alpha-fetoprotein.
Am. J. Obstet. Gynecology, 172(1pt1):219-220.
8. VERLOES, A. et al. 1995.
A prenatal trisomy 21 screening program using AFP, hCG and free estriol assays on maternal dried blood.
Am. J. Obstet. Gynecology, 172(1pt1):167-179.
9. CHRISTMOS, J.T. et al. 1994.
The effect of fetomaternal bleeding of the risk of adverse pregnancy outcome in patients with elevated second trimester maternal serum alpha-fetoprotein levels.
Am. J. Obstet. Gynecology, 171(2):315-320.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Calibrators (0 to 5) Samples, controls	- -	0.05 -	- 0.05
Incubation	1 hour at room temperature		
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate carefully	
Tracer	0.1	0.1	0.1
Incubation	1 hour at room temperature		
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate carefully	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP0081 – KIP0084	P.I. Number : 1700503/en	Revision nr : 101028/1
---	-----------------------------	---------------------------

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XVIII. BIBLIOGRAPHIE

- CUCKLE, H.S., WALD, N.J. 1984.
Maternal serum alpha-fetoprotein measurement : a screening test for Down Syndrome.
Lancet, i:926
- HUNTER, W.M. BUDD, P.S. 1981.
Immunoradiometric versus radioimmunoassay a comparison using alpha-fetoprotein as the model analyte.
J. of Immunol. Meth.45:255
- KOHN, J. WEAVER, P.C. 1974.
Serum alpha-fetoprotein in hepatocellular carcinoma.
Lancet, i:334
- KOHN, J., ORR, A.H. MC ELWAIN, T.J. et al.
Serum alpha-fetoprotein in patients with testicular tumours.
Lancet; II:433
- NORGAARD-PENDERSEN, B. 1976.
Human alpha-fetoprotein - A review of recent methodological and clinical studies.
Scand. J. Immunol., Suppl. 4:7.
- U.K. Collaborative study, First report. 1977.
Maternal serum alpha-fetoprotein measurement in antenatal screening for anencephaly and spina bifida in early pregnancy.
Lancet, i:1323.
- JEVY, L.B. et al. 1995
Hepatocellular carcinoma in a pregnant woman detected by routine screening of maternal alpha-fetoprotein.
Am. J. Obstet. Gynecology, 172(1pt1):219-220.
- VERLOES, A. et al. 1995.
A prenatal trisomy 21 screening program using AFP, hCG and free estriol assays on maternal dried blood.
Am. J. Obstet. Gynecology, 172(1pt1):167-179.
- CHRISTMOS, J.T. et al. 1994.
The effect of fetomaternal bleeding of the risk of adverse pregnancy outcome in patients with elevated second trimester maternal serum alpha-fetoprotein levels.
Am. J. Obstet. Gynecology, 171(2):315-320.

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Data wydania : 2010-10-28

	<u>Used symbols</u>
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
U	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer

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