25OH-VIT.D3-RIA-CT

KIP1961

Read entire protocol before use.

25OH-VIT.D3-RIA-CT

I. INTENDED USE

Radioimmunoassay for the in vitro quantitative measurement of 25OH-Vit.D3 in human serum and heparin plasma.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource 25OH-Vit.D3-Ria-CT Kit
- B. Catalog number :** KIP1961 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)67 88.99.99 Fax : +32 (0)67 88.99.96

III. CLINICAL BACKGROUND

A. *Physiological function of 25OH-Vit.D₃*

25-Hydroxyvitamin D₃ (25OHD₃) is the trivial name of 9, 10-secocholesta-5, 7, 10(19)-triene-3 β , 25-diol. This secosteroid is produced in the liver by 25-hydroxylation of cholecalciferol or Vitamin D₃. 25OHD₃ is a precursor for other Vitamin D metabolites and has only a limited biological activity in itself. The most active derivative is 1 α ,25-Hydroxyvitamin D₃, produced in the kidney (or placenta) by 1 α -hydroxylation of 25OHD₃. This hormonally regulated steroid stimulates the intestinal absorption of both calcium and phosphorus. It also stimulates bone resorption and mineralisation thereby preventing the development of rickets, osteoporosis or osteomalacia. This Vitamin D hormone might also be active in other tissues responsible for calcium transport (placenta, kidney, mammary gland, ...) and endocrine glands (beta-cells, parathyroid glands, ...). 25-hydroxyvitamin D₃ is a main precursor for the metabolites.

B. *Regulatory mechanism*

The production of 25-hydroxyvitamin D mainly occurs in the liver although other tissues (intestine, kidney) might perform the same hydroxylation. Although there might be some feedback inhibition of 25-hydroxyvitamin D on its own production, a higher substrate availability (Vitamin D₃) results in higher 25-hydroxyvitamin D production and higher circulating 25-hydroxyvitamin D concentrations in blood.

C. *Clinical applications*

This assay is of importance for the diagnosis of Vitamin D deficiency, insufficiency or intoxication.

IV. PRINCIPLES OF THE METHOD

At first calibrators, controls and samples (serum or heparin plasma) are extracted with acetonitrile.

A fixed amount of ^{125}I labelled 25OH Vitamin D₃ competes with the 25OH Vitamin D₃ from either extracted samples, controls or calibrators for a fixed amount of specific antibody sites immobilized to the lower and inner surface of plastic tubes.

After 2 hours incubation at room temperature, an aspiration step stops the competition reaction. The tubes are then washed with 2 ml washing solution and counted in a gamma counter.

V. REAGENTS PROVIDED

Reagents	Quantity	Colour Code	Reconstitution			
 Tubes coated with anti 25OH-Vitamin.D ₃	2 x 48	pink	Ready for use			
<table border="1"><tr><td>Ag</td><td>¹²⁵I</td></tr></table> ¹²⁵ I 25OH-Vitamin.D ₃	Ag	¹²⁵ I	1 vial lyophil. 160 kBq	red	Reconstitute extemporaneously by adding 6 ml of the tracer buffer	
Ag	¹²⁵ I					
<table border="1"><tr><td>CAL</td><td>0</td></tr></table> Calibrator 0: horse serum / phosphate buffer, with gentamycin	CAL	0	1 vial lyophil.	yellow	Add 0.5 ml distilled water	
CAL	0					
<table border="1"><tr><td>CAL</td><td>N</td></tr></table> Calibrators 1-5 in horse serum / phosphate buffer, with gentamycin (see exact values on vial labels)	CAL	N	5 vials lyophil.	yellow	Add 0.5 ml distilled water	
CAL	N					
<table border="1"><tr><td>INC</td><td>BUF</td></tr></table> Incubation Buffer with sodium azide (< 0.1%)	INC	BUF	1 vial 45 ml	black	Ready for use	
INC	BUF					
<table border="1"><tr><td>TRACER</td><td>BUF</td></tr></table> Ethanol solution with sodium azide (< 0.1%)	TRACER	BUF	1 vial 7 ml	red	Ready for use	
TRACER	BUF					
<table border="1"><tr><td>ACETONITRILE</td></tr></table> Acetonitrile	ACETONITRILE	1 vial 25 ml	black	Ready for use		
ACETONITRILE						
<table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table> Wash solution	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1"><tr><td>CONTROL</td><td>N</td></tr></table> Controls 1 and 2 in human serum with thymol	CONTROL	N	2 vials lyophil.	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note : Use Calibrator 0 for dilution of samples with values above the highest calibrator before extraction step.
No international reference material is available.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 µl, 100 µl, 400 µl, 500 µl.
3. Vortex mixer
4. Magnetic stirrer
5. Glass tubes (12x75 mm) for extraction step
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration and washing device
8. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).
10. Centrifuge operating at 800-1500 g

VII. REAGENT PREPARATION

- A. Calibrators :** Reconstitute the calibrators with 0.5 ml distilled water.
- B. Controls :** Reconstitute the controls with 0.5 ml distilled water.
- C. ^{125}I 25OH-Vit.D₃ :** Reconstitute with 6 ml of the tracer buffer.
- D. Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- The reconstituted tracer has to be frozen after first use. Then, it is stable until the expiry date.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum and heparin plasma samples must be kept at 2-8°C.
- This kit is suitable for serum and heparin plasma samples. A correlation has been established between 23 serum and heparin plasma samples from same patients: Plasma = 0.95 Serum + 1.25, R = 0.89.
- If the test is not run within 48 h., storage at -20°C is recommended
- After thawing serum and heparin plasma samples should be mixed (Vortex), then centrifuged.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

Extraction step :

1. Label glass tubes (12x75 mm) for extraction: 6 calibrators, 2 controls and up to 40 samples in duplicate.
2. Dispense 100 µl of each calibrator, control or sample in the respective tubes.
3. Add 0.5 ml acetonitrile to each tube.
4. Mix for 7 seconds with a vortex.
5. Centrifuge for 5 minutes at room temperature (at 800-1500 g).

Incubation step :

1. Label coated tubes in duplicate for each calibrator, sample, control. For the determination of total counts, label 2 normal tubes.
2. Add 100 µl of the supernatant obtained after the extraction step in the corresponding tubes. Pipette tips have to be saturated with corresponding supernatant before the addition in the tube.
3. Dispense 400 µl Incubation Buffer in each tube, except those for total counts.
4. Add 50 µl Tracer in each tube, including total counts.
5. Incubate for 2 hours, under stirring (300-700 RPM), at room temperature.
6. Aspirate the content of each tube (except total counts).
7. Wash tubes twice with 2 ml Wash Solution and aspirate. Avoid foaming during the addition of the wash solution.
8. After the washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
9. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- 1. Calculate the mean of duplicate determinations.
- 2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

B/B0 x 100 = (Counts (Calibrator or sample) / Counts (Zero Calibrator)) x 100

- 3. Using a 3 cycle semi-logarithmic or logit-log graph paper plot the (B/B0 x 100) values for each calibrator point as a function of the 25OH.D3 concentration of each calibrator point, reject obvious outliers.
- 4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.
- 5. By interpolation of the sample (B/B0 x 100) values, determine the 25OH.D3 concentrations of the samples from the reference curve.
- 6. For each assay, the percentage of total tracer bound in the absence of unlabelled 25OH.D3 (B0/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

25OH-Vit.D3		cpm	B/Bo x 100
Total count		40301	-
Calibrator	0 ng/ml	23803	100 %
	3.0 ng/ml	20171	85 %
	9.0 ng/ml	16237	68 %
	20.0 ng/ml	12156	51 %
	65.0 ng/ml	7090	30 %
	130.0 ng/ml	5330	22 %

Note : 1 ng/ml = 2.5 pmol/ml

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent 25OH-Vitamin.D3 concentration two standard deviations below the average counts at zero binding, was 1.2 ng/ml.

B. Specificity

The percentage of cross reaction estimated by comparison of the concentration yielding a 50 % inhibition are respectively :

Compound	Cross-Reactivity (%)
25OH-Vitamin D3	100
25OH-Vitamin D2	<0.3
1,25(OH)2-Vitamin.D3	89*
1,25(OH)2-Vitamin.D2	<0.01
Vitamin D3	<0.03
Vitamin D2	<0.03
24,25(OH)2-Vitamin.D3	<0.8

* As 1,25(OH)2-Vit.D3 concentrations are practically 1000 times lower than 25-OH-Vit.D3, this cross-reactivity is insignificant and does not interfere in this 25-OH-Vit.D3 assay.

The assay performance is not affected by hemolysis (5 g/L haemoglobin tested) or bilirubinemia (0.25 g/L bilirubin tested).

C. Precision

INTRA-ASSAY				INTER-ASSAY			
Sample	N	<X> ± SD (ng/ml)	C.V. (%)	Sample	N	<X> ± SD (ng/ml)	C.V. (%)
A	49	19.1 ± 1.4	7.3	A	9	6.7 ± 0.5	7.2
B	43	10.4 ± 0.9	8.8	B	9	13.1 ± 1.0	7.3

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST	
Added 25OH-Vit.D3 (ng/ml)	Recovery (%)
0	100
15	109
40	97
60	106

DILUTION TEST			
Sample dilution	Theoretical concent. (ng/ml)	Measured concent. (ng/ml)	Recovery (%)
1/1	42.5	42.5	100
1/2	21.3	20.9	98.4
1/4	10.6	11.2	105.3
1/8	5.3	5.7	107.0
1/16	2.7	2.6	98.8
1/32	1.3	1.3	98.0

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 15 and 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY			
	0' (ng/ml)	15' (ng/ml)	30' (ng/ml)
Sample 1	9.5	9.2	10.1
Sample 2	36.8	33.9	36.2

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplo results of the samples should rely on Good Laboratory Practises

XV. EXPECTED VALUES

Dietary intake, race, season and age are known to affect the normal levels of 25OH.Vit.D3.

The expected values given hereafter should not be considered as absolute. DIAsource evaluated serum collected from 39 male and 40 female, healthy according to Ca, PTH and albumin values, from Western Europe. The age of the volunteers fell within the range of 17 - 58 years. Samples were collected during the months of December 2010 and January 2011. The mean for the population (n = 79) was 12.5 ng/mL, ranging from 4.1 to 28.7 ng/ml (based on 2.5 % to 97.5 % percentiles).

Each laboratory should establish its own range based on their local population.

Recent literature has suggested the following ranges for the classification of 25 OH Vitamin D status: Deficiency: 0-10 ng/mL; Insufficiency: 10-30 ng/mL; Sufficiency: 30 to 150 ng/mL; Toxicity: >150 ng/mL.

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves. All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

XVII. BIBLIOGRAPHY

1. ZERWEKH J.E. (2008)
Blood biomarkers of Vitamin D status.
Am. J. Clin. Nutr., 87(suppl):1087S-91S.
2. HOLICK M.F. (2006)
Resurrection of Vitamin D deficiency and rickets.
J. Clin. Invest., 116:2062-2072.
3. HEANEY R.P. (2000)
Vitamin D: how much do we need, and how much is too much.
Osteoporos. Int., 11:553-555.
4. DAWSON-HUGHES B., HEANEY R.P., HOLICK M.F., LIPS P., MEUNIER P.J. (1997)
Prevalence of Vitamin D insufficiency in an adult normal population.
Osteoporos. Int., 7:439-443.
5. BISCHOFF-FERRARI H.A., GIOVANNUCCI E., WILLETT W.C., DIETRICH T., DAWSON-HUGHES B. (2006)
Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes.
Am. J. Clin. Nutr., 84:18-28.
6. HOLICK M.F. (2007)
Vitamin D deficiency.
N. Engl. J. Med., 357:266-281.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
<u>EXTRACTION</u>			
Calibrators	-	0.1	-
Samples/Controls	-	-	0.1
Acetonitrile	-	0.5	0.5
Vortex Centrifugation	Vortex for 7 seconds 5 minutes at 800-1500 g		
<u>INCUBATION</u>			
Supernatant of extraction	-	0.1	0.1
Incubation Buffer	-	0.4	0.4
Tracer	0.05	0.05	0.05
Incubation	2 hours at room temperature under stirring (300-700 RPM)		
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate	
Wash solution		2.0	
Separation		aspirate carefully	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP1961	P.I. Number : 1700543/en	Revision nr : 110328/1
-------------------------------------	-----------------------------	---------------------------

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