# **25OH-Vitamin D total- RIA-CT**

***KIP1971***

Read entire protocol before use.

## 25OH Vitamin D total -RIA-CT

### I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of 25-hydroxyvitamin D2 and D3 (25-OH-D2 and 25-OH-D3) in serum or heparin plasma.

### II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource 25OH Vitamin D total -RIA-CT Kit
- B. **Catalog number :** KIP 1971 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue de l'Industrie, 8, B-1400 Nivelles, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)67 88.99.99      Fax : +32 (0)67 88.99.96**

### III. CLINICAL BACKGROUND

Vitamin D is the generic term used to designate Vitamin D2 or ergocalciferol and Vitamin D3 or cholecalciferol.

Humans naturally produce Vitamin D3 when the skin is exposed to ultraviolet sun rays.

In the liver mainly, Vitamin D3 is metabolised into 25-Hydroxyvitamin D3 (25 OH D3) which is the main form of Vitamin D circulating in the body.

25 OH D3 is a precursor for other Vitamin D metabolites and has also a limited activity by itself.

The most active derivative is 1,25-Hydroxyvitamin D3, produced in the kidney (or placenta) by 1  $\alpha$  -hydroxylation of 25OHD3.

25OHVitamin D stimulates the intestinal absorption of both calcium and phosphorus and also bone resorption and mineralisation. 25OH Vitamin D might also be active in other tissues responsible for calcium transport (placenta, kidney, mammary gland...) and endocrine gland (parathyroid glands, beta cells...).

Vitamin D3 and Vitamin D2 are also available by ingestion through food or dietary supplementation.

As Vitamin D2 is metabolised in a similar way to vitamin D3, both contribute to the overall Vitamin D status of an individual.

It is the reason why it is very important to measure both forms of 25 OH Vitamin D equally for a correct diagnosis of Vitamin D deficiency, insufficiency or intoxication.

Vitamin D deficiency is an important risk factor for rickets, osteomalacia, senile osteoporosis, cancer and abnormal pregnancy outcomes.

The measurement of both 25 OH Vitamin D forms is also required to determine the cause of abnormal serum calcium concentrations in patients.

Vitamin D intoxication has been shown to cause kidney and tissue damages.

IV. PRINCIPLES OF THE METHOD

At first calibrators, controls and samples (serum or heparin plasma) are pre-treated to release 25OH Vitamin D<sub>2</sub> and 25OH Vitamin D<sub>3</sub> from Vitamin D Binding Protein (DBP). A fixed amount of <sup>125</sup>I labelled 25OH Vitamin D competes with the 25OH Vitamin D<sub>2</sub> and 25OH Vitamin D<sub>3</sub> from treated samples, controls or calibrators, for a fixed amount of a specific monoclonal antibody sites immobilized to the lower and inner surface of plastic tubes. After 2 hours incubation at room temperature on a tube shaker, an aspiration step terminates the competition reaction. The tubes are then washed twice and aspirated again. A calibration curve is plotted and the total 25 OH Vitamin D (D<sub>2</sub> and D<sub>3</sub>) concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

| Reagents                                                                                                                                                           | 96 Test Kit      | Colour Code      | Reconstitution                   |                 |                                           |                                                                   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|------------------|----------------------------------|-----------------|-------------------------------------------|-------------------------------------------------------------------|
|  Tubes coated with Mab anti 25OH Vit D2 and 25OH Vit D3                             | 2 x 48           | green            | <b>Ready</b> for use             |                 |                                           |                                                                   |
| <table border="1"><tr><td>Ag</td><td><sup>125</sup>I</td></tr></table><br><sup>125</sup> Iodine labelled 25OH Vit D (HPLC grade).                                  | Ag               | <sup>125</sup> I | 1 vial<br>160 kBq<br>lyophilised | red             | <b>Add</b> 6 ml of Tracer Buffer          |                                                                   |
| Ag                                                                                                                                                                 | <sup>125</sup> I |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>CAL</td><td>0</td></tr></table><br>Calibrator 0: horse serum and phosphate buffer with gentamycin.                                       | CAL              | 0                | 1 vial<br>lyophilised            | yellow          | <b>Add</b> 0.5 ml distilled water         |                                                                   |
| CAL                                                                                                                                                                | 0                |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>CAL</td><td>N</td></tr></table><br>Calibrators 1-5 in horse serum and phosphate buffer with gentamycin (see exact values on vial labels) | CAL              | N                | 5 vials<br>lyophilised           | yellow          | <b>Add</b> 0.5 ml distilled water         |                                                                   |
| CAL                                                                                                                                                                | N                |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table><br>Wash solution (TRIS-HCl)                                                             | WASH             | SOLN             | CONC                             | 1 vial<br>10 ml | brown                                     | <b>Dilute</b> 70 x with distilled water (use a magnetic stirrer). |
| WASH                                                                                                                                                               | SOLN             | CONC             |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>CONTROL</td><td>N</td></tr></table><br>Controls N = 2 in human serum with thymol                                                         | CONTROL          | N                | 2 vials<br>lyophilised           | silver          | <b>Add</b> 0.5 ml distilled water         |                                                                   |
| CONTROL                                                                                                                                                            | N                |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>TRACER</td><td>BUF</td></tr></table><br>Tracer Buffer with casein and gentamycin                                                         | TRACER           | BUF              | 1 vial<br>7 ml                   | red             | <b>Ready</b> for use                      |                                                                   |
| TRACER                                                                                                                                                             | BUF              |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>PRE</td><td>SOLN</td></tr></table><br>Pre-treatment Solution with Proclin 300 (0.1 %)                                                    | PRE              | SOLN             | 1 vial<br>lyophilised            | green           | <b>Add</b> 5 ml distilled distilled water |                                                                   |
| PRE                                                                                                                                                                | SOLN             |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>PRE</td><td>BUF</td></tr></table><br>Pre-treatment Buffer with sodium azide (<0.1%)                                                      | PRE              | BUF              | 1 vial<br>20 ml                  | blue            | <b>Ready</b> for use                      |                                                                   |
| PRE                                                                                                                                                                | BUF              |                  |                                  |                 |                                           |                                                                   |
| <table border="1"><tr><td>DIL</td><td>SOLN</td></tr></table><br>Dilution Solution with sodium azide (<0.1%)                                                        | DIL              | SOLN             | 1 vial<br>100 ml                 | black           | <b>Ready</b> for use                      |                                                                   |
| DIL                                                                                                                                                                | SOLN             |                  |                                  |                 |                                           |                                                                   |

**Note :** Use Calibrator 0 for dilution of samples with values above the highest calibrator before pre-treatment step.  
No international reference material is available.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 50 µl ,100 µl , 200µl , 500 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)

- Vortex mixer
- Magnetic stirrer
- Tube shaker (300 to 700 rpm)
- Glass tubes (12x 75 mm) or polypropylene tubes (for example: Falcon 352063) for pre-treatment step.
- A water bath at 37°C
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system
- Any gamma counter capable of measuring <sup>125</sup>I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators :** Reconstitute the calibrators with 0.5 ml distilled water
- Controls:** Reconstitute the controls with 0.5 ml distilled water.
- Tracer:** Reconstitute the lyophilized tracer with 6 ml of tracer buffer.
- Pre-treatment Solution :**Reconstitute extemporaneously the lyophilized pre-treatment solution with 5 ml distilled water.
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at -20°C (with one thawing).
- After its first use, pre-treatment solution is stable 2 months at -20°C(with one thawing)
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Don't use EDTA plasma.**
- This kit is suitable for serum and heparin plasma samples. A correlation has been established between 23 serum and heparin plasma samples from same patients: Plasma = 0.997 Serum, R = 0.91.
- Serum and heparin samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, samples storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.

X. PROCEDURE

- Handling notes**  
Do not use the kit or components beyond expiry date.  
Do not mix materials from different kit lots.  
Bring all the reagents to room temperature prior to use.  
Thoroughly mix all reagents and samples by gentle agitation or swirling.  
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.  
Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

B1 Pre-treatment :

- Label one glass tube or one polypropylene tube for each calibrator, control and sample.
- Briefly vortex calibrators, controls and samples and dispense 100 µl of each into the respective tubes.
- Dispense 50 µl of pre-treatment solution into each tube.
- Dispense 200 µl of pre-treatment buffer into each tube.
- Briefly vortex tubes, cover them and incubate 30 minutes at 37°C in a water bath.  
**! Be careful incubation volume must be completely immersed during this step.**  
**! Respect incubation time: 30 - 35 minutes is necessary.**  
**! Respect incubation temperature: 37 - 40 °C is necessary.**
- Dispense 1 ml of Dilution Solution into each tube and vortex briefly.

B2 Incubation :

- Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes.
- Dispense 500µl of pre-treated calibrator or control or sample

- Dispense 50 µl of <sup>125</sup>Iodine labelled 25OH Vitamin D into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand to liberate any trapped air bubbles.
- Incubate for 2 hours at room temperature on a tube shaker (300 to 700 rpm)
- Aspirate the content of each tube (except total counts). Be sure that the tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate. Avoid foaming during the addition of the Working Wash solution.
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate .
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

**XI. CALCULATION OF RESULTS**

- Calculate the mean of duplicate determinations.
- Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

- Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B<sub>0</sub>(%)) values for each calibrator point as a function of 25OH vitamin D concentration of each calibrator point. Reject obvious outliers.
- Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
- By interpolation of the sample (B/B<sub>0</sub> (%)) values, determine the total 25OH vitamin D concentrations of the samples from the calibration curve.
- For each assay, the percentage of total tracer bound in the absence of unlabelled 25 OH vitamin D (B<sub>0</sub>/T) must be checked.

**XII. TYPICAL DATA**

The following data are for illustration only and should never be used instead of the real time calibration curve.

| 25OH Vitamin D total |           | cpm   | B/Bo (%) |
|----------------------|-----------|-------|----------|
| Total count          |           | 45756 |          |
| Calibrator           | 0.0 ng/ml | 21958 | 100.0    |
|                      | 2 ng/ml   | 17916 | 81.6     |
|                      | 6 ng/ml   | 14170 | 64.5     |
|                      | 15 ng/ml  | 10710 | 48.8     |
|                      | 55 ng/ml  | 4949  | 22.5     |
|                      | 120 ng/ml | 3161  | 14.4     |

**Note :** 1 ng/ml = 2.5 pmol/ml

**XIII. PERFORMANCE AND LIMITATIONS**

**A. Detection Limit**

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent 25OH-Vitamin D total concentration two standard deviations below the average counts at zero binding, was 0.4ng/ml.

**B. Specificity**

The percentage of cross reaction estimated by comparison of the concentration yielding a 50 % inhibition are respectively :

| Compound                                       | Cross-Reactivity (%) |
|------------------------------------------------|----------------------|
| 25OH-Vitamin D <sub>3</sub>                    | 100                  |
| 25OH-Vitamin D <sub>2</sub>                    | 95                   |
| 1,25(OH) <sub>2</sub> -Vitamin.D <sub>3</sub>  | 4                    |
| 1,25(OH) <sub>2</sub> -Vitamin.D <sub>2</sub>  | 4                    |
| Vitamin D3                                     | <0.3                 |
| Vitamin D2                                     | <0.3                 |
| 24,25(OH) <sub>2</sub> -Vitamin.D <sub>3</sub> | ≥100                 |

The assay performance is not affected by hemolysis (5 g/L hemoglobin tested), bilirubinemia (0.25 g/L bilirubin tested) or triglycerides (5 g/L tested).

**C. Precision**

| INTRA-ASSAY |    |                  |          | INTER-ASSAY |    |                  |          |
|-------------|----|------------------|----------|-------------|----|------------------|----------|
| Sample      | N  | <X> ± SD (ng/ml) | C.V. (%) | Sample      | N  | <X> ± SD (ng/ml) | C.V. (%) |
| A           | 23 | 14.7 ± 0.8       | 5.4      | A           | 16 | 13.8 ± 1.3       | 9.4      |
| B           | 23 | 42.4 ± 2         | 4.7      | B           | 16 | 38.5 ± 3.5       | 9.1      |

SD : Standard Deviation; CV: Coefficient of variation

**D. Accuracy**

| RECOVERY TEST                         |              |
|---------------------------------------|--------------|
| Added 25OH-Vit.D <sub>3</sub> (ng/ml) | Recovery (%) |
| 0                                     | 100          |
| 25                                    | 105          |
| 50                                    | 105          |
| 100                                   | 91           |
| Added 25OH-Vit.D <sub>2</sub> (ng/ml) | Recovery (%) |
| 0                                     | 100          |
| 25                                    | 101          |
| 50                                    | 101          |
| 100                                   | 96           |

| DILUTION TEST   |                              |                           |              |
|-----------------|------------------------------|---------------------------|--------------|
| Sample dilution | Theoretical concent. (ng/ml) | Measured concent. (ng/ml) | Recovery (%) |
| 1/1             | 102.8                        | 102.8                     | 100          |
| 1/2             | 51.4                         | 60.1                      | 117          |
| 1/4             | 25.7                         | 25.6                      | 100          |
| 1/8             | 12.9                         | 12.9                      | 100          |
| 1/1             | 80.3                         | 80.3                      | 100          |
| 1/2             | 40.2                         | 42.0                      | 105          |
| 1/4             | 20.1                         | 24.7                      | 123          |
| 1/8             | 10.0                         | 10.9                      | 108          |

**E. Time delay between last calibrator and sample dispensing**

As shown hereafter, assay results remain accurate even when a sample is dispensed 10 and 20 minutes after the calibrator has been added to the coated tubes.

| TIME DELAY |            |             |             |
|------------|------------|-------------|-------------|
|            | 0' (ng/ml) | 10' (ng/ml) | 20' (ng/ml) |
| Sample 1   | 13.5       | 13.6        | 11.7        |
| Sample 2   | 39.5       | 39.3        | 37.4        |

#### XIV. EXPECTED VALUES

Dietary intake, race, season and age are known to affect the normal levels of 25OH Vitamin D3.

The expected values given hereafter should not be considered as absolute. DIAsource evaluated serum collected from 39 male and 40 female, healthy according to Ca, PTH and albumin values, from Western Europe. The age of the volunteers fell within the range of 17 - 58 years. Samples were collected during the months of December 2010 and January 2011. The mean for the population (n = 79) was 12.9 ng/mL, ranging from 6.6 to 24.5 ng/ml (based on 2.5 % to 97.5 % percentiles).

Each laboratory should establish its own range based on their local population.

Recent literature has suggested the following ranges for the classification of 25 OH Vitamin D status: Deficiency: 0-10 ng/mL; Insufficiency: 10-30 ng/mL; Sufficiency: 30 to 150 ng/mL; Toxicity: >150 ng/mL.

#### XV. PRECAUTIONS AND WARNINGS

##### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

#### XIV. PRECAUTIONS AND WARNINGS

1. ZERWEKH J.E. (2008)  
**Blood biomarkers of Vitamin D status.**  
Am. J. Clin. Nutr., 87(suppl):1087S-91S.
2. HOLICK M.F. (2006)  
**Resurrection of Vitamin D deficiency and rickets.**  
J. Clin. Invest., 116:2062-2072.
3. HEANEY R.P. (2000)  
**Vitamin D: how much do we need, and how much is too much.**  
Osteoporos. Int., 11:553-555.
4. DAWSON-HUGHES B., HEANEY R.P., HOLICK M.F., LIPS P., MEUNIER P.J. (1997)  
**Prevalence of Vitamin D insufficiency in an adult normal population.**  
Osteoporos. Int., 7:439-443.
5. BISCHOFF-FERRARI H.A., GIOVANNUCCI E., WILLETT W.C., DIETRICH T., DAWSON-HUGHES B. (2006)  
**Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes.**  
Am. J. Clin. Nutr., 84:18-28.
6. HOLICK M.F.(2004)  
**Sunlight and vitamin D for bone health and prevention of autoimmune diseases, cancers and cardiovascular disease.**  
Am. J. Clin. Nutr., 80:1678S-1688S.

7. HEANEY R.P. (2010)  
**Defining deficiency of vitamin D .**  
Clinical Laboratory International , October 2010, vol.34 : 16-19.
8. HOLICK M.F. (2007)  
**Vitamin D deficiency.**  
N. Engl. J. Med., 357:266-281.
9. TAHA N. M. , VIETH R.(2010)  
**The problem of an optimal target level for 25-Hydroxyvitamin D, the test for vitamin D nutritional status .**  
Clinical Laboratory International , November 2010, vol.34 : 28-30

#### XV. SUMMARY OF THE PROTOCOL

|                                                                        | TOTAL<br>COUNTS<br>µl                                                                   | CALIBRATORS<br>µl              | SAMPLE (S)<br>CONTROLS<br>µl |
|------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------|------------------------------|
| <b><u>PRE-TREATMENT</u></b><br>(in glass tubes or polypropylene tubes) |                                                                                         |                                |                              |
| Samples, controls, calibrators                                         | -                                                                                       | 100                            | 100                          |
| Pre-treatment solution                                                 | -                                                                                       | 50                             | 50                           |
| Pre-treatment buffer                                                   | -                                                                                       | 200                            | 200                          |
| Incubation:                                                            | Vortex<br>30-35 minutes at 37-40°C in a water bath (incubation volume must be immersed) |                                |                              |
| Dilution Solution                                                      | -                                                                                       | 1000                           | 1000                         |
|                                                                        | Vortex                                                                                  |                                |                              |
| <b><u>INCUBATION</u></b><br>(in coated tubes)                          |                                                                                         |                                |                              |
| Pre-treated Calibrators, samples, controls                             | -                                                                                       | 500                            | 500                          |
| Tracer                                                                 | 50                                                                                      | 50                             | 50                           |
| Incubation                                                             | 120 minutes at RT on a shaker (300 to 700 rpm)                                          |                                |                              |
| Separation<br>Working Wash solution                                    | -                                                                                       | Aspirate<br>2.0 ml             |                              |
| Separation<br>Working Wash solution                                    | -                                                                                       | Aspirate<br>2.0 ml<br>Aspirate |                              |
| Counting                                                               | Count tubes for 60 seconds                                                              |                                |                              |

|                                      |                             |                           |
|--------------------------------------|-----------------------------|---------------------------|
| DIAsource Catalogue Nr :<br>KIP 1971 | P.I. Number :<br>1700528/en | Revision nr :<br>110324/1 |
|--------------------------------------|-----------------------------|---------------------------|

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- ZERWEKH J.E. (2008)  
**Blood biomarkers of Vitamin D status.**  
Am. J. Clin. Nutr., 87(suppl):1087S-91S.
- HOLICK M.F. (2006)  
**Resurrection of Vitamin D deficiency and rickets.**  
J. Clin. Invest., 116:2062-2072.
- HEANEY R.P. (2000)  
**Vitamin D: how much do we need, and how much is too much.**  
Osteoporos. Int., 11:553-555.

- DAWSON-HUGHES B., HEANEY R.P., HOLICK M.F., LIPS P., MEUNIER P.J. (1997)  
**Prevalence of Vitamin D insufficiency in an adult normal population.**  
Osteoporos. Int., 7:439-443.
- BISCHOFF-FERRARI H.A., GIOVANNUCCI E., WILLETT W.C., DIETRICH T., DAWSON-HUGHES B. (2006)  
**Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes.**  
Am. J. Clin. Nutr., 84:18-28.
- HOLICK M.F.(2004)  
**Sunlight and vitamin D for bone health and prevention of autoimmune diseases, cancers and cardiovascular disease.**  
Am. J. Clin. Nutr., 80:16788S-1688S.
- HEANEY R.P. (2010)  
**Defining deficiency of vitamin D .**  
Clinical Laboratory International , October 2010, vol.34 : 16-19.
- HOLICK M.F. (2007)  
**Vitamin D deficiency.**  
N. Engl. J. Med., 357:266-281.
- TAHA N. M. , VIETH R.(2010)  
**The problem of an optimal target level for 25-Hydroxyvitamin D, the test for vitamin D nutritional status .**  
Clinical Laboratory International , November 2010, vol.34 : 28-30

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|                 | <b><u>Used symbols</u></b>         |
|-----------------|------------------------------------|
|                 | Consult instructions for use       |
|                 | Storage temperature                |
|                 | Use by                             |
| <b>LOT</b>      | Batch code                         |
| <b>REF</b>      | Catalogue number                   |
| <b>CONTROL</b>  | Control                            |
| <b>IVD</b>      | In vitro diagnostic medical device |
|                 | Manufacturer                       |
|                 | Contains sufficient for <n> tests  |
| WASH SOLN CONC  | Wash solution concentrated         |
| CAL 0           | Zero calibrator                    |
| CAL N           | Calibrator #                       |
| CONTROL N       | Control #                          |
| Ag 1251         | Tracer                             |
| Ab 1251         | Tracer                             |
| Ag 1251 CONC    | Tracer concentrated                |
| Ab 1251 CONC    | Tracer concentrated                |
|                 | Tubes                              |
| INC BUF         | Incubation buffer                  |
| ACETONITRILE    | Acetonitrile                       |
| SERUM           | Serum                              |
| DIL SPE         | Specimen diluent                   |
| DIL BUF         | Dilution buffer                    |
| ANTISERUM       | Antiserum                          |
| IMMUNOADSORBENT | Immunoabsorbent                    |
| DIL CAL         | Calibrator diluent                 |
| REC SOLN        | Reconstitution solution            |
| PEG             | Polyethylene glycol                |
| EXTR SOLN       | Extraction solution                |
| ELU SOLN        | Elution solution                   |
| GEL             | Bond Elut Silica cartridges        |
| PRE SOLN        | Pre-treatment solution             |
| NEUTR SOLN      | Neutralization solution            |
| TRACEUR BUF     | Tracer buffer                      |
|                 | Microtiterplate                    |
| Ab HRP          | HRP Conjugate                      |
| Ag HRP          | HRP Conjugate                      |
| Ab HRP CONC     | HRP Conjugate concentrate          |
| Ag HRP CONC     | HRP Conjugate concentrate          |
| CONJ BUF        | Conjugate buffer                   |
| CHROM TMB CONC  | Chromogenic TMB concentrate        |
| CHROM TMB       | Chromogenic TMB solution           |
| SUB BUF         | Substrate buffer                   |
| STOP SOLN       | Stop solution                      |
| INC SER         | Incubation serum                   |
| BUF             | Buffer                             |
| Ab AP           | AP Conjugate                       |
| SUB PNPP        | Substrate PNPP                     |
| BIOT CONJ CONC  | Biotin conjugate concentrate       |
| AVID HRP CONC   | Avidine HRP concentrate            |
| ASS BUF         | Assay buffer                       |
| Ab BIOT         | Biotin conjugate                   |
| Ab              | Specific Antibody                  |
| SAV HRP CONC    | Streptavidin HRP concentrate       |
| NSB             | Non-specific binding               |
| 2nd Ab          | 2nd Antibody                       |
| ACID BUF        | Acidification Buffer               |
| DIST            | Distributor                        |

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