**NSE-IRMA**

***KIP2471***

Read entire protocol before use.

## NSE-IRMA

### ***I. INTENDED USE***

Immunoradiometric assay kit for the in vitro quantitative measurement of human Neuron Specific Enolase (NSE) in serum.

### ***II. GENERAL INFORMATION***

- A. Proprietary name :** DIAsource NSE-IRMA Kit
- B. Catalog number :** KIP2471 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue de l'Industrie 8, B-1400 Nivelles, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)67 88.99.99      Fax : +32 (0)67 88.99.96**

### ***III. CLINICAL BACKGROUND***

#### **A. Biological activities**

The glycolytic enzyme enolase (2-phospho-D-glycerate hydrolase, EC 4.2.1.11) exists as several dimeric isoenzymes ( $\alpha\alpha$ ,  $\alpha\beta$ ,  $\alpha\gamma$ ,  $\beta\beta$  and  $\gamma\gamma$ ). The  $\gamma\gamma$  and  $\alpha\gamma$  enolase isoenzymes are also known as neuron-specific enolase (NSE) as they are produced in central and peripheral neurons and malignant tumors of neuroectodermal origin.

Elevated serum concentrations of NSE are detected in patients with medullary carcinoma of the thyroid, pancreatic islet cell tumor, pheochromocytoma carcinoid, neuroblastoma and small-cell lung carcinoma (SCLC).

#### **B. Clinical applications**

The diagnosis of lung cancer generally requires medical imaging, endoscopy, intra-operative findings and histology.

An important indication for tumor marker determinations in lung cancer is in assessing the efficacy of therapy and post-operative follow-up care.

In patients with SCLC, serum levels of NSE reflect the response to chemotherapy.

#### IV. PRINCIPLES OF THE METHOD

The DIAsource NSE-IRMA is an immunoradiometric assay based on coated-tube separation. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mab1. Addition of Mab2, the signal antibody labelled with  $^{125}\text{I}$ , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

#### V. REAGENTS PROVIDED

| Reagents                                                                                                                                                                                                    | 96 tests Kit     | Colour Code      | Reconstitution              |                 |                                   |                                                                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|------------------|-----------------------------|-----------------|-----------------------------------|------------------------------------------------------------------|
|  Tubes coated with anti NSE (monoclonal antibodies)                                                                        | 2 x 48           | orange           | <b>Ready</b> for use        |                 |                                   |                                                                  |
| <table border="1"><tr><td>Ab</td><td><sup>125</sup>I</td></tr></table><br>Anti-NSE- <sup>125</sup> I (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide (<0.1%) and inert red dye | Ab               | <sup>125</sup> I | 1 vial<br>5.5 ml<br>700 kBq | red             | <b>Ready</b> for use              |                                                                  |
| Ab                                                                                                                                                                                                          | <sup>125</sup> I |                  |                             |                 |                                   |                                                                  |
| <table border="1"><tr><td>CAL</td><td>O</td></tr></table><br>Calibrator 0 in bovin serum<br>With thymol                                                                                                     | CAL              | O                | 1 vial<br>lyophil.          | yellow          | <b>Add</b> 3ml distilled water    |                                                                  |
| CAL                                                                                                                                                                                                         | O                |                  |                             |                 |                                   |                                                                  |
| <table border="1"><tr><td>CAL</td><td>N</td></tr></table><br>Calibrators 1-5 in bovin serum with thymol (see exact value on vial labels)                                                                    | CAL              | N                | 5 vials<br>lyophil.         | yellow          | <b>Add</b> 0.5 ml distilled water |                                                                  |
| CAL                                                                                                                                                                                                         | N                |                  |                             |                 |                                   |                                                                  |
| <table border="1"><tr><td>WASH</td><td>SOLN</td><td>CONC</td></tr></table><br>Wash solution (TRIS-HCl)                                                                                                      | WASH             | SOLN             | CONC                        | 1 vial<br>10 ml | brown                             | <b>Dilute</b> 70x with distilled water (use a magnetic stirrer). |
| WASH                                                                                                                                                                                                        | SOLN             | CONC             |                             |                 |                                   |                                                                  |
| <table border="1"><tr><td>CONTROL</td><td>N</td></tr></table><br>Controls 1 and 2 in human serum and thymol                                                                                                 | CONTROL          | N                | 2 vials<br>lyophil.         | silver          | <b>Add</b> 0.5 ml distilled water |                                                                  |
| CONTROL                                                                                                                                                                                                     | N                |                  |                             |                 |                                   |                                                                  |

**Note:** Use the zero calibrator for sera dilutions.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50  $\mu\text{l}$ , 500  $\mu\text{l}$  and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Pipette for delivery of 5 to 10 ml of distilled water
4. Vortex mixer
5. Magnetic stirrer
6. Tube shaker (400rpm)
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- A. Calibrators** : Reconstitute the calibrator 0 with 3 ml distilled water and the calibrators 1-5 with 0.5 ml.
- B. Controls** : Reconstitute the controls with 0.5 ml distilled water.
- C. Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- After reconstitution of the calibrators and controls, aliquots should be made and kept at -20°C for maximum 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

**NSE is also found in erythrocytes, plasma cells and platelets, it may be released into serum if separation from red cells does not occur within 60 minutes of venipuncture.**

- Serum must be kept at 2-8°C.
- If the test is not run within 24 h., storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- Do not use plasma samples.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. Thoroughly mix all reagents and samples by gentle agitation or swirling.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, samples, controls and dispense 50  $\mu\text{l}$  of each into the respective tubes.
3. Dispense 50  $\mu\text{l}$  of anti-NSE- $^{125}\text{I}$  tracer into each tube, including the uncoated tubes for total counts.
4. Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
5. Incubate for 2 hours at room temperature on a tube shaker (400 rpm).
6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
7. Wash tubes with 2 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
8. Aspirate (or decant) the content of each tube (except total counts).
9. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
10. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
11. Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of NSE (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

| NSE-IRMA    |             | cpm    | B/T (%) |
|-------------|-------------|--------|---------|
| Total count |             | 275069 | 100     |
| Calibrator  | 0.0 ng/ml   | 399    | 0.15    |
|             | 2.5 ng/ml   | 2669   | 0.97    |
|             | 7.8 ng/ml   | 6710   | 2.44    |
|             | 26.0 ng/ml  | 20604  | 7.49    |
|             | 78.0 ng/ml  | 57102  | 20.76   |
|             | 270.0 ng/ml | 144748 | 52.62   |

Since no international reference material is available for NSE antigen, BSE NSE-IRMA calibrator values are assigned against a set of in-house reference standards.

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.19 ng/ml.

B. Precision

| INTRA ASSAY |           |                  |        | INTER ASSAY |           |                  |        |
|-------------|-----------|------------------|--------|-------------|-----------|------------------|--------|
| Serum       | Replicate | <X> ± SD (ng/ml) | CV (%) | Serum       | Replicate | <X> ± SD (ng/ml) | CV (%) |
| A           | 20        | 12.9 ± 0.5       | 3.6    | A           | 20        | 12.7 ± 0.5       | 3.8    |
| B           | 20        | 104.4 ± 1.3      | 1.2    | B           | 20        | 101.4 ± 3.8      | 3.8    |

SD : Standard Deviation; CV: Coefficient of variation

C. Accuracy

| RECOVERY TEST     |                       |              |
|-------------------|-----------------------|--------------|
| Added NSE (ng/ml) | Recovered NSE (ng/ml) | Recovery (%) |
| 20.0              | 20.0                  | 100.0        |
| 39.0              | 38.7                  | 99.0         |
| 79.0              | 76.0                  | 96.2         |
| 156.4             | 149.4                 | 95.5         |

| DILUTION TEST |          |                              |                           |
|---------------|----------|------------------------------|---------------------------|
| Sample        | Dilution | Theoretical Concent. (ng/ml) | Measured Concent. (ng/ml) |
| 1             | 1/1      |                              | 164.7                     |
|               | 1/2      | 82.4                         | 79.2                      |
|               | 1/4      | 41.2                         | 39.4                      |
|               | 1/8      | 20.6                         | 20.3                      |
|               | 1/16     | 10.3                         | 10.2                      |
|               | 1/32     | 5.1                          | 5.1                       |
|               | 1/64     | 2.6                          | 2.5                       |
|               | 1/128    | 1.3                          | 1.2                       |
| 2             | 1/1      |                              | 240.5                     |
|               | 1/2      | 120.3                        | 120.4                     |
|               | 1/4      | 60.1                         | 61.3                      |
|               | 1/8      | 30.1                         | 29.9                      |
|               | 1/16     | 15.0                         | 14.7                      |
|               | 1/32     | 7.5                          | 6.6                       |
|               | 1/64     | 3.8                          | 3.1                       |
|               | 1/128    | 1.9                          | 1.4                       |

Samples were diluted with zero calibrator .

D. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

| TIME DELAY  |      |      |      |      |
|-------------|------|------|------|------|
|             | 0'   | 10'  | 20'  | 30'  |
| S 1 (ng/ml) | 11.5 | 11.8 | 11.4 | 11.9 |
| S 2 (ng/ml) | 45.1 | 44.7 | 43.6 | 43.1 |

E. Hook-effect

Samples containing 15000 ng/ml NSE give a result higher than the last calibration point.

F. Specificity

The monoclonal antibodies used are specific for the g-subunit of enolase. No measurable cross-reactions with other enolase have been observed.

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XV. REFERENCE INTERVALS

Among 332 apparently healthy individuals, 98% of the results were below 12.5 ng/ml.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.  
This kit contains <sup>125</sup>I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.  
This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.  
All radioactive handling should be executed in a designated area. away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.  
Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.  
All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.  
Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.  
Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## XVII. BIBLIOGRAPHY

1. LAMERZ R.  
**NSE (neuron-specific enolase).**  
Thomas L, ed. Clinical laboratory diagnostics. Verlag Frankfurt : TH-Books, 1998; 979-981.
2. LIIPPO KK, TERHO T.  
**Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.**  
Acta Oncol 1991; 30:321-324.
3. OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTÉ G., NIKSSON K.  
**Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.**  
J. Pediatr Surg 1982; 17:381-385.
4. TAPIA FL., POLAK JM., BARBOSA AJ., et al.  
**Neuron-specific enolase is produced by neuroendocrine tumours.**  
Lancet 1981 ; 1 :808-811.
5. FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.  
**Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.**  
Cancer 1989; 63:1143-1149.
6. BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.  
**Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.**  
Clin Cancer Res 2000; 6:597-601.
7. EBERT W., HOPPE M., MULEY T., DRINGS P.  
**Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.**  
Anticancer RES 1997; 17:2875-2878.

## XVIII. SUMMARY OF THE PROTOCOL

|                                                  | TOTAL<br>COUNTS<br>ml                               | CALIBRATORS<br>ml           | SAMPLE(S)<br>CONTROLS<br>ml |
|--------------------------------------------------|-----------------------------------------------------|-----------------------------|-----------------------------|
| Calibrators (0-5)<br>Samples, Controls<br>Tracer | -<br>-<br>0.05                                      | 0.05<br>-<br>0.05           | -<br>0.05<br>0.05           |
| Incubation                                       | 2 hours at room temperature with shaking at 400 rpm |                             |                             |
| Separation<br>Working Wash<br>solution           | -<br>-                                              | Aspirate (or decant)<br>2.0 |                             |
| Separation<br>Working Wash<br>solution           | -<br>-                                              | Aspirate (or decant)<br>2.0 |                             |
| Separation                                       | -                                                   | Aspirate (or decant)        |                             |
| Counting                                         | Count tubes for 60 seconds                          |                             |                             |

|                                     |                             |                               |
|-------------------------------------|-----------------------------|-------------------------------|
| DIAsource Catalogue Nr :<br>KIP2471 | P.I. Number :<br>1700863/en | Revision Number :<br>090429/1 |
|-------------------------------------|-----------------------------|-------------------------------|

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2. LIIPPO KK, TERHO T.  
**Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.**  
Acta Oncol 1991; 30:321-324.

3. OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTÉ G., NIKSSON K.  
**Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.**  
J. Pediatr Surg 1982; 17:381-385.

4. TAPIA FL., POLAK JM., BARBOSA AJ., et al.  
**Neuron-specific enolase is produced by neuroendocrine tumours.**  
Lancet 1981 ; 1 :808-811.

5. FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.  
**Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.**  
Cancer 1989; 63:1143-1149.

6. BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.  
**Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.**  
Clin Cancer Res 2000; 6:597-601.

7. EBERT W., HOPPE M., MULEY T., DRINGS P.  
**Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.**  
Anticancer RES 1997; 17:2875-2878.

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2.

LIIPPO KK, TERHO T.  
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.  
Acta Oncol 1991; 30:321-324.

3. OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTÉ G., NIKSSON K.  
**Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.**  
J. Pediatr Surg 1982; 17:381-385.
  
4. TAPIA FL., POLAK JM., BARBOSA AJ., et al.  
**Neuron-specific enolase is produced by neuroendocrine tumours.**  
Lancet 1981 ; 1 :808-811.
  
5. FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.  
**Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.**  
Cancer 1989; 63:1143-1149.
  
6. BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.  
**Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.**  
Clin Cancer Res 2000; 6:597-601.
  
7. EBERT W., HOPPE M., MULEY T., DRINGS P.  
**Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.**  
Anticancer RES 1997; 17:2875-2878.

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