# Acetylcholine Receptor

*KIPC512B*

# Acetylcholine Receptor, Binding Ria

en

For the determination of anti-acetylcholine receptor antibody in human serum.

## KIPC512B

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

### INTENDED USE

The DIAsource AChRab Assay kit is for the determination of anti-acetylcholine receptor antibody in human serum.

### SUMMARY AND EXPLANATION

Myasthenia gravis (MG) is a skeletal muscle disorder characterized by muscular weakness. In such cases, muscular weakness is due to anti-acetylcholine receptor (AChR) antibodies. Anti-AChR antibodies are present in approximately 90% of patients with MG. Anti-AChR antibodies could be: binding antibodies (multitudes of wide populations of antibodies directed to hydrophilic domains of receptors), blocking antibodies (preventing binding of acetylcholine to receptors) and modulating antibodies (accelerating endocytosis resulting in loss of receptors).

### PRINCIPLE OF THE TEST

Anti-AChR antibodies are detectable by a radiobinding that is followed by precipitation of the antibodies. Patient specimens and references are incubated with detergent-solubilized fetal and adult AChRs labeled with  $^{125}\text{I}$ -alpha-bungarotoxin. The resulting bound complexes of labeled AChRs and autoantibodies are then immunoprecipitated with anti-human IgG. After centrifugation, the supernatant is aspirated and the pellet containing labeled AChR/autoantibody-bound complexes is counted in a gamma counter. Counts are directly proportional to the amount of autoantibodies present.

### MATERIALS PROVIDED

|     |   |
|-----|---|
| CAL | N |
|-----|---|

Calibrators : 6 vials.N = 0 to 5.  
Ready to use  
120 Tests : 0.20 ml

|    |                  |
|----|------------------|
| Ag | $^{125}\text{I}$ |
|----|------------------|

$^{125}\text{I}$ AChR 4 vials. 120 Tests.  
Lyophilized. 4 x 76 kBq.  
To be reconstituted with 3.2 ml of  $^{125}\text{I}$ -AChR Diluent

|        |     |
|--------|-----|
| TRACER | BUF |
|--------|-----|

$^{125}\text{I}$ -AChR Diluent  
120 Tests : 15 ml

|         |   |
|---------|---|
| Control | H |
|---------|---|

Positive Control  
120 Tests : 0.20 ml

|         |   |
|---------|---|
| Control | L |
|---------|---|

Negative Control  
120 Tests : 0.20 ml

|    |
|----|
| Ab |
|----|

Goat Anti-Human IgG  
120 Tests : 25 ml

|         |
|---------|
| Diluent |
|---------|

Normal Human Serum : 2 vials  
120 Tests : 1 ml

|      |      |
|------|------|
| WASH | SOLN |
|------|------|

Wash Buffer : 2 vials  
120 Tests : 125 ml

### MATERIALS NOT PROVIDED

1. Distilled or deionized water
2. Precision pipettes
3. Disposable pipette tips
4. Absorbance paper or paper towel
5. Graph paper
6. Tubes > 2 ml
7. Gamma Counter

### STORAGE AND STABILITY

1. Store the kit at 2 - 8° C.
2. Keep microwells sealed in a dry bag with desiccants.
3. The reagents are stable until expiration of the kit.
4. Do not expose test reagents to heat, sun, or strong light.

### WARNINGS AND PRECAUTIONS

For in vitro diagnostic use only.

This kit contains  $^{125}\text{I}$  (half-life: 60 days), emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations.

1. Potential biohazardous materials:  
The calibrator and controls contain human source components which have been tested and found non-reactive for hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents. However, as there is no test method that can offer complete assurance that HIV, Hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984
2. Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
3. The components in this kit are intended for use as an integral unit. The components of different lots should not be mixed.
4. It is recommended that serum samples be run in duplicate.
5. Optimal results will be obtained by strict adherence to this protocol. Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from this may yield invalid data.
6. All radioactive handling should be executed in a designated area, away from regular passage.
7. Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory.
8. Use protective clothing and disposable gloves.

### SPECIMEN COLLECTION HANDLING

1. Collect blood specimens and separate the serum immediately.
2. Specimens may be stored refrigerated at (2-8°C) for 5 days. If storage time exceeds 5 days, store frozen at (-20° C) for up to one month.
3. Avoid multiple freeze-thaw cycles.
4. Prior to assay, frozen sera should be completely thawed and mixed well.
5. Do not use grossly lipemic specimens.

### ASSAY PROCEDURE

Prior to assay, allow reagents to stand at room temperature.

1. Up to 30 minutes before use, reconstitute the lyophilized  $^{125}\text{I}$  Acetylcholine Receptor with 3200  $\mu\text{L}$  of  $^{125}\text{I}$ -AChR diluent. (Use within 2 weeks of reconstitution; store at 2-8°C) Bring all components and patient specimens to room temperature and mix well before use.
2. Pipette 20  $\mu\text{L}$  of each Calibrator in duplicate into their respective tubes
3. Pipette 20  $\mu\text{L}$  of each control in duplicate into their respective tubes
4. Pipette 20  $\mu\text{L}$  of the patient specimen(s) into their respective tubes.
5. Pipette 100  $\mu\text{L}$  of the diluted  $^{125}\text{I}$  Acetylcholine Receptor into all of the tubes and Vortex. Incubate all tubes at room temperature (22-28°C) for two (2) hours.
6. Pipette 200  $\mu\text{L}$  of anti-human IgG into all tubes and Vortex gently. Incubate at room temperature (22-28°C) for 30 minutes.
7. Add 1000  $\mu\text{L}$  of Wash Buffer to each tube and Vortex gently. Centrifuge all tubes for 30 minutes at 3500 rpm at a temperature of 2-8°C.
8. Carefully decant or aspirate the supernatants without disturbing the pellets.
9. Add 1000  $\mu\text{L}$  of Wash Buffer to each tube and Vortex gently. Centrifuge all tubes for 15 minutes at 3500 rpm at a temperature of 2-8°C.
10. Carefully decant or aspirate the supernatants without disturbing the pellets
11. Pipette 100  $\mu\text{L}$  of the reconstituted  $^{125}\text{I}$  Acetylcholine Receptor into each of two tubes for counting Total Activity.
12. Determine the CPM for each tube, including the two Total Activity tubes, for one (1) minute in a gamma counter set for measuring  $^{125}\text{I}$ .

## CALCULATIONS

For each run, a calibration curve must be calculated from the mean CPM of each of the six calibrators.

1. Calculate the ratio of the counts per minute (CPM) of calibrators (1-6), control serums and samples to the CPM of calibrator 1 ( $B/B_0$ ).
2. Fit the curve with Log(dose)-Log( $B/B_{max}$ ) model or four parameter model, plot the  $B/B_0$  of each calibrator on the ordinate (Y-axis) against the calibrator concentration on the abscissa (X-axis).
3. Read the concentrations (ng/ml) of AChR in the samples from the calibration curve.

### Example of a Calibration Curve

| Calibrator | Anti-AChR (nmol/L) | Mean CPM |
|------------|--------------------|----------|
| 1          | 0                  | 518      |
| 2          | 0,20               | 1853     |
| 3          | 0,50               | 2817     |
| 4          | 1,20               | 5067     |
| 5          | 2,50               | 7483     |
| 6          | 7,50               | 11159    |

### Typical Interpretation:

< 0.25 nmol/L: Negative  
 0.25-0.40 nmol/L: Equivocal  
 >0.4 nmol/L: Positive  
 > 2.0 nmol/L: Dilute sample and multiply to verify value

## EXPECTED VALUE

Serum samples from apparently healthy volunteers showed values between 0.00 and 0.35 nmol/L. (mean 0.093, median 0.080, S.D. 0.092, 95% interval 0.0 – 0.30 nmol/L).

It is recommended that each laboratory establish its own range of normal values.

## PERFORMANCE CHARACTERISTICS

### 1. Precision

Intra-Assay (Within-run) precision was determined with 20 replicates of four different serum samples:

| Serum Sample | Replicates | Mean (nmol/L) | SD (nmol/L) | CV (%) |
|--------------|------------|---------------|-------------|--------|
| 1            | 20         | 0.41          | 0.03        | 2.4    |
| 2            | 20         | 0.98          | 0.04        | 4.8    |
| 3            | 20         | 3.20          | 0.05        | 6.2    |
| 4            | 20         | 7.50          | 0.19        | 2.1    |

Inter-Assay (Between-run) precision was determined with 20 daily assays of two different serum samples:

| Sample Serum | Days | Mean (nmol/L) | SD (nmol/L) | CV (%) |
|--------------|------|---------------|-------------|--------|
| 1            | 20   | 0.98          | 0.15        | 15     |
| 2            | 20   | 7.50          | 0.12        | 1.6    |

### 2. Sensitivity

The minimum detectable concentration (calculated from the 95% confidence limit of 20 determinations of the 0 nmol/L calibrator (A1)) was 0,01nmol/L.

## 3. Parallelism

Four high serum samples with different acetylcholine receptor antibody levels were diluted in normal human serum.

| Sample | Dilution  | Measured (nmol/L) | Recovery (%) |
|--------|-----------|-------------------|--------------|
| 1      | undiluted | 7,69              | 100          |
|        | 1:2       | 3,98              | 104          |
|        | 1:4       | 1,57              | 82           |
|        | 1:8       | 0,75              | 78           |
|        | 1:16      | 0,39              | 81           |
|        | 1:32      | 0,27              | 112          |
| 2      | undiluted | 2,82              | 100          |
|        | 1:2       | 1,56              | 111          |
|        | 1:4       | 0,74              | 105          |
|        | 1:8       | 0,46              | 130          |
|        | 1:16      | 0,19              | 108          |
|        | undiluted | 1,75              | 100          |
| 3      | undiluted | 1,75              | 100          |
|        | 1:2       | 0,84              | 96           |
|        | 1:4       | 0,44              | 101          |
|        | 1:8       | 0,25              | 114          |
|        | 1:16      | 0,11              | 101          |
|        | undiluted | 1,37              | 100          |
| 4      | undiluted | 1,37              | 100          |
|        | 1:2       | 0,65              | 95           |
|        | 1:4       | 0,32              | 93           |
|        | 1:8       | 0,15              | 88           |

## 4. Specificity

Nine patients' serum samples containing non-anti-Acetylcholine Receptor autoantibodies (anti -sm-mp, -ro(ss-b), la(ss-b), -dsdna, -rf, and ana) were tested. No cross reactivities were found with any of these antibodies.

## REFERENCES

Vincent A, Newsom-Davis J. Acetylcholine receptor antibody as a diagnostic test for myasthenia gravis: results in 153 validated cases and 1967 diagnostic assays. Journal of Neurology, Neurosurgery, and Psychiatry 48:1246-1252 (1985).

Lindstrom JM, Seybold ME, Lennon VA, Whittingham S, and Duane DD. Antibody to acetylcholine receptor in myasthenia gravis. Neurology 26:1054-1059 (1976).

MacLennan C, Beeson D, Buijs A-M, Vincent A and Newsom-Davis J. Acetylcholine receptor expression in human extraocular muscles and their susceptibility to myasthenia gravis. Annals of Neurology 41:423-431 (1997).

Beeson D, Jacobson L, Newsom-Davis J and Vincent A. A transfected human muscle cell line expressing the adult subtype of the human muscle acetylcholine receptor for diagnostic assays in myasthenia gravis Neurology 47:1552-1555 (1996).

Revision date : 2011-01-25

# Acetylcholine Receptor, Binding Ria

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MacLenan C, Beeson D, Buijs A-M, Vincent A and Newsom-Davis J.Acetylcholine receptor expression in human extraocular muscles and their susceptibility to myasthenia gravis. Annals of Neurology 41:423-431 (1997).

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Date de révision: 2011-01-25

|                 | <b><u>Used symbols</u></b>         |
|-----------------|------------------------------------|
|                 | Consult instructions for use       |
|                 | Storage temperature                |
|                 | Use by                             |
| <b>LOT</b>      | Batch code                         |
| <b>REF</b>      | Catalogue number                   |
| <b>CONTROL</b>  | Control                            |
| <b>IVD</b>      | In vitro diagnostic medical device |
|                 | Manufacturer                       |
|                 | Contains sufficient for <n> tests  |
| WASH SOLN CONC  | Wash solution concentrated         |
| CAL 0           | Zero calibrator                    |
| CAL N           | Calibrator #                       |
| CONTROL N       | Control #                          |
| Ag 1251         | Tracer                             |
| Ab 1251         | Tracer                             |
| Ag 1251 CONC    | Tracer concentrated                |
| Ab 1251 CONC    | Tracer concentrated                |
|                 | Tubes                              |
| INC BUF         | Incubation buffer                  |
| ACETONITRILE    | Acetonitrile                       |
| SERUM           | Serum                              |
| DIL SPE         | Specimen diluent                   |
| DIL BUF         | Dilution buffer                    |
| ANTISERUM       | Antiserum                          |
| IMMUNOADSORBENT | Immunoabsorbent                    |
| DIL CAL         | Calibrator diluent                 |
| REC SOLN        | Reconstitution solution            |
| PEG             | Polyethylene glycol                |
| EXTR SOLN       | Extraction solution                |
| ELU SOLN        | Elution solution                   |
| GEL             | Bond Elut Silica cartridges        |
| PRE SOLN        | Pre-treatment solution             |
| NEUTR SOLN      | Neutralization solution            |
| TRACEUR BUF     | Tracer buffer                      |
|                 | Microtiterplate                    |
| Ab HRP          | HRP Conjugate                      |
| Ag HRP          | HRP Conjugate                      |
| Ab HRP CONC     | HRP Conjugate concentrate          |
| Ag HRP CONC     | HRP Conjugate concentrate          |
| CONJ BUF        | Conjugate buffer                   |
| CHROM TMB CONC  | Chromogenic TMB concentrate        |
| CHROM TMB       | Chromogenic TMB solution           |
| SUB BUF         | Substrate buffer                   |
| STOP SOLN       | Stop solution                      |
| INC SER         | Incubation serum                   |
| BUF             | Buffer                             |
| Ab AP           | AP Conjugate                       |
| SUB PNPP        | Substrate PNPP                     |
| BIOT CONJ CONC  | Biotin conjugate concentrate       |
| AVID HRP CONC   | Avidine HRP concentrate            |
| ASS BUF         | Assay buffer                       |
| Ab BIOT         | Biotin conjugate                   |
| Ab              | Specific Antibody                  |
| SAV HRP CONC    | Streptavidin HRP concentrate       |
| NSB             | Non-specific binding               |
| 2nd Ab          | 2nd Antibody                       |
| ACID BUF        | Acidification Buffer               |
| DIST            | Distributor                        |

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