VASOPRESSIN-RIA

KIPERB319

VASOPRESSIN-RIA

en

KIPERB319

IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

1 INTENDED USE

The DIAsource vasopressin-ria kit contains reagents and instructions for the quantitative measurement of vasopressin in plasma or urine. After solid phase extraction (SPE) or ethanol extraction the plasma vasopressin concentrations are measured by radioimmunoassay (RIA). Urine vasopressin concentrations can be measured directly.
For professional use within a laboratory.

2 CLINICAL APPLICATION

Vasopressin, or Antidiuretic Hormone (ADH) is a cyclic nanopeptide with a molecular weight of 1083. Its structure is very similar to that of oxytocin, differing in only two amino acids. Endogenous ADH has antidiuretic and pressor activity, both approaching 400 units per mg, with an antidiuretic-to-vasopressin ratio of 1, and a biphasic plasma half-life of 2.5 and 14.5 minutes. ADH is synthesized in the hypothalamic supraoptic nucleus and paraventricular nucleus of primates and transported via axonal flow to the posterior pituitary for storage and eventual release. The clinical application of a vasopressin radioimmunoassay is in diabetes insipidus, psychogenic water intoxication, hyponatraemia, stress conditions, ADH as a neurotransmitter and hypertension studies. ADH values can be influenced by cigarettes, tea, coffee, alcohol and some drugs.

3 PRINCIPLE OF TEST

After solid phase extraction (SPE) or ethanol extraction of the plasma samples, vasopressin is assayed by a competitive radioimmunoassay. Urinary vasopressin can be measured directly. This assay uses a rabbit anti-vasopressin antiserum and a radioiodinated vasopressin [¹²⁵I] tracer. Bound and free phases are separated by a second antibody bound to solid phase particles, followed by a centrifugation step. The radioactivity in the bound fractions is measured and a typical calibration curve can be generated. The values of the extracted samples are corrected for extraction recovery.

4 PRECAUTIONS

Materials derived from human blood and used in the preparation of this kit were tested and found negative for hepatitis B surface antigen (HBsAg), antibodies to HCV and for antibodies to HIV-1 and HIV-2. However, handle all components as a possible source of infection.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

The radioactive material included may be received, acquired, possessed and used only by physicians, clinical laboratories or hospitals for in-vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulation of each country.

Adherence to the basic rules of radiation safety should provide adequate protection.

- Do not eat, drink, smoke or apply cosmetics where radioactive materials are used.
- Do not pipette radioactive solutions by mouth.
- Avoid direct contact with all radioactive materials by using protective articles such as lab coats and disposable gloves.
- All radiological work should be done in a designated area.
- Radioactive materials should be stored in original containers in a designated area.
- Laboratory equipment and glassware, which are subject to contamination, should be segregated to prevent cross-contamination of different radioisotopes.
- Any radioactive spills should be taken care of immediately in accordance with established procedures.
- All radioactive materials must be disposed of in accordance with the prevailing regulations and guidelines of the agencies jurisdiction over the laboratory.

The reagents in this kit contain sodium azide (0.05%). Contact with copper or lead drain pipes may result in the cumulative formation of highly explosive azide deposits. On disposal of the reagents in the sewerage, always flush with copious amounts of water, which prevents metallic azide formation. Plumbing suspected of being contaminated with these explosive deposits should be rinsed thoroughly with 10% sodium hydroxide solution.

5 SPECIMEN COLLECTION

Careful standardization of the patient preparation and sampling conditions is recommended.

5.1 Vasopressin in plasma

- Draw blood from fasting patient into a chilled tube, containing EDTA or Heparin.
- Centrifuge at 4° C to separate the plasma.
- Freeze the sample in plastic tubes at -20° C until assayed.

NOTE: Vasopressin (ADH) in plasma is stable at -20° C only for 4 weeks, or stable up to 3 months after addition of 500 KIU Trasylol (Bayer) per mL blood; after extraction.

5.2 Vasopressin in urine

Vasopressin can be determined directly, in unextracted human urine.

- Collect 24 hours urine sample.
- Register urine volume.
- Measure urinary osmolarity.
- If the sample is not assayed immediately, keep an aliquot at -20° C.
- Measure urine samples undiluted and in dilutions of 1:2, 1:4 or higher.

6 MATERIALS AND EQUIPMENT REQUIRED

- Pipettes (100 µL, 200 µL, 300 µL, 1.00 mL, 2.00 mL, 5.00 mL)
- Repeating dispensers (100 µL, 200 µL)
- Measuring cylinder 25 mL
- Polystyrene RIA tubes (12 x 75 mm)
- Ethanol absolute (99%)
- Vortex
- Centrifuge
- Icebath
- Vac-concentrator
- Nitrogen gas
- Polystyrene or glass tubes for extraction (16 x 100 mm)
- Sep-pak C18 (Waters Ass Inc. art. 051910)
- Acetic acid 4%
- Methanol absolute (99%)
- 1N HCl
- Gamma counter

7 QUALITY CONTROL

Controls should be carried out in each assay run. Two controls are included in the kit, the value (without extraction procedure) is indicated on the Quality Control sheet and on the label of the vials. Use also controls as recommended by the control plasma manufacturer and in accordance with reference laboratories practice to monitor the accuracy and precision of reagents and techniques. Each laboratory should establish its own extraction recovery under their own experimental conditions.

8 SHELF LIFE AND STORAGE

This kit is stable until the stated expiry date if stored as specified.

Upon receipt of the kit, all reagents should be stored at 2-8° C.

The reconstituted reagents should be stored according to table on section 10.

The reconstituted reagents are stable according to table on section 10, but no longer than to the expiry date.

9 COMPONENTS

Contents of the Kit

- | | | | | | |
|----|------|---|---------|--|---|
| 1. | REAG | A | Ab | 1 vial, lyophilized anti-vasopressin (Rabbit) for 100 tubes.
Colour : yellow. | |
| 2. | REAG | B | Aq | ¹²⁵ I | 1 vial, 28 KBq or 0.75 µCi. Lyophilized. Specific activity : 62-77 MBq/nmol
(1700-2100 µCi/nmol).
Colour : blue. |
| 3. | REAG | C | DASP | 1 vial, goat anti-rabbit IgG's bound to solid phase.
11 mL suspension. | |
| 4. | REAG | D | ASSAY | BUF | 2 vials, 0.05 M phosphate buffer with 0.25 % HSA, 0.25 % EDTA disodium salt,
0.05 % NaN ₃ , 500 KIU Trasylol/mL. 2 x 50 mL (liquid) |
| 5. | REAG | E | CAL | 1 vial, 5.00 mL vasopressin calibrator, 60 pmol/L. Lyophilized in assay buffer. | |
| 6. | REAG | F | Control | 1 vial, 2.00 mL lyophilized vasopressin control, low level. | |
| 7. | REAG | G | Control | 1 vial, 2.00 mL lyophilized vasopressin control, high level. | |

10 PREPARATION OF REAGENTS

Item	Reconstitute each vial with		Stable at	Special remarks
Anti-vasopressin (Reagent A)	22 mL distilled water	Mix gently	-20° C for at least 3 months after recon- stitution	
¹²⁵ I-vasopressin (Reagent B)	25 mL distilled water	Mix gently	-20° C until expiry date	
Double antibody solid phase (Reagent C)	Ready for use. The separation reagent should be placed on a magnetic stirrer for 10 minutes at room temp		2-8° C until expiry date	It is possible to pipette the reagent with a repeating dispenser
Assay buffer (Reagent D)	Ready for use		2-8° C until expiry date	
Vasopressin calibrator 60 pmol/L (Reagent E)	5.00 mL distilled water	Mix gently	-20° C for at least 3 months after recon- stitution	Refer to table for calibration curve preparation
Vasopressin low control (Reagent F)	2.00 mL distilled water	Mix gently	-20° C for at least 3 months after recon- stitution	The value of the control is found on the label of the vial and in the QC sheet (without extraction)
Vasopressin high control (Reagent G)	2.00 mL distilled water	Mix gently	-20° C for at least 3 months after recon- stitution	The value of the control is found on the label of the vial and in the QC sheet (without extraction)

PREPARE ALL REAGENTS 15 MINUTES BEFORE USE!

11 SAMPLE PREPARATION

Before proceeding in the RIA procedure two different sample preparation methods can be used:

A1 : Sep-pak C18 extraction

A2 : Ethanol extraction

A1: Sep-pak C18 extraction procedure

Column: Sep-pak C18 cartridge (Waters Ass. Inc. art. 051910)

Procedure: DO NOT EXTRACT CALIBRATORS AND CONTROLS.

1. Wash the column with 10 mL distilled water, 5 mL methanol and 10 mL distilled water, respectively.
2. Acidify 1.0 mL plasma sample with 150 µL 1N HCl.
3. Bring this acidified sample into the column.
4. Wash the column with 20 mL 4% acetic acid.
5. Elute with 4 mL methanol.
6. Dry the methanol under a stream of nitrogen or air.
7. Reconstitute the residue with 1.0 mL of assay buffer.
8. Proceed the RIA procedure immediately or store the extracted samples at – 20°C.

To estimate recovery, add an aliquot (200 µL) ¹²⁵I-vasopressin tracer to a random plasma sample and submit the recovery sample for the same extraction procedure.

Recovery calculation

1. Prepare a recovery estimation tube (R).
 - Pipette 1.0 mL of a random plasma sample into the recovery tube (R). The sample used for this recovery assay should have a protein matrix similar to the samples being tested.
 - Add 200 µL ¹²⁵I-vasopressin tracer into tube R and mix.
 - Extract this sample along with samples in the above procedure.
2. Prepare 2 Total Recovery tubes (TR).
 - Pipette 200 µL ¹²⁵I-vasopressin tracer into two TR tubes.
 - Add 100 µL assay buffer and mix.
 - Cap and set aside these tubes to be counted for recovery calculation.
3. Reconstitute the dried recovery sample (R) by adding 1.0 mL assay buffer and vortex thoroughly.

4. Pipette 300 µL of the reconstituted recovery sample tube (R) into two assay tubes.
5. Count the total recovery (TR) and recovery (R) tubes for at least two minutes in a gamma counter.
Calculate % recovery by dividing the cpm in the recovery tubes (R) by cpm in the total recovery tubes (TR) and multiply by 3.33:

$$\% \text{ Recovery} : \frac{\text{cpm recovery tube (R)}}{\text{cpm total recovery tube}} \times 3.33 \times 100\%$$

Extraction recoveries should score values between 60-80%.

A2: Ethanol extraction procedure

DO NOT EXTRACT CALIBRATORS AND CONTROLS.

1. Label one extraction tube (E) for each patient sample. Label one additional tube (R) in order to estimate the extraction recovery.
2. Place the extraction tubes and ethanol on ice.
3. Pipette 0.8 mL of each sample into the appropriately labelled extraction tubes (E).
4. Prepare a recovery estimation tube (R).
 - Pipette 0.8 mL of a random plasma sample into the recovery tube (R). The sample used for this recovery assay should have a protein matrix similar to the samples being tested.
 - Add 200 µL ¹²⁵I-vasopressin tracer into tube R and mix.
 - Extract this sample along with samples in step 6.
5. Prepare 2 Total Recovery tubes (TR).
 - Pipette 200 µL ¹²⁵I-vasopressin tracer into two TR tubes.
 - Add 100 µL assay buffer and mix.
 - Cap and set aside these tubes to be counted for recovery calculation.
6. Add 4 mL chilled ethanol to each sample (tubes E) and Recovery Tube (R).
7. Mix and vortex for 2 minutes.
8. Centrifuge all extraction tubes (E and R) at 2000 g. for 15 min. at 2-8° C.
9. Decant supernatant from each extraction tube into previous prepared clean, appropriately labelled 16 x 100 mm tubes.
10. Evaporate the supernatants under a stream of nitrogen to dryness (at max. 37° C), or evaporate by using a Vac-concentrator.
11. Reconstitute the dried samples by adding 0.8 mL assay buffer and vortex thoroughly.
12. Proceed RIA procedure immediately or store the extracted samples at -20° C up to two weeks before using in the assay.
13. Reconstitute the dried recovery sample (R) by adding 0.8 mL assay buffer and vortex thoroughly.
14. Pipette 300 µL of the reconstituted recovery sample tube (R) into two assay tubes.
15. Count the total recovery (TR) and recovery (R) tubes for at least two minutes in a gamma counter.

Recovery calculation

Calculate % recovery by dividing the cpm in the recovery tubes (R) by cpm in the total recovery tubes (TR) and multiply by 2.67:

$$\% \text{ Recovery} : \frac{\text{cpm recovery tube (R)}}{\text{cpm total recovery tube (TR)}} \times 2.67 \times 100\%$$

Extraction recoveries should score values between 40-50%.

12 ASSAY PROTOCOL

All tests (calibrators, controls and samples) should be performed in duplicate.

A. Preparation of calibrator solutions

Dilution	Reagent E (=calibrator a)	Vasopressin concentration 60 pmol/L
1000 µL of Reagent E + 1000 µL assay buffer vortex	Calibrator b	30 pmol/L
1000 µL of calibrator b + 1000 µL assay buffer Vortex	Calibrator c	15 pmol/L
1000 µL of calibrator c + 1000 µL assay buffer Vortex	Calibrator d	7.5 pmol/L
1000 µL of calibrator d + 1000 µL assay buffer Vortex	Calibrator e	3.8 pmol/L
1000 µL of calibrator e + 1000 µL assay buffer vortex	Calibrator f	1.9 pmol/L
	assay buffer	0 pmol/L

B. Assay Procedure

- After preparation of the calibrator solutions pipette 300 µL of each calibrator, controls, each extract from plasma, or diluted urine into the correspondingly labelled tubes.
- Add 300 µL assay buffer (Reagent D) to the max binding (0 pmol/L binding) tubes and 500µL assay buffer to the blank tubes (NSB).
- Add 200 µL vasopressin antiserum (Reagent A) to all tubes, except blank (NSB) and total counts (TC) tubes.
- Vortex all tubes and incubate at 4° C for 18-24 hours.
- Add 200 µL ¹²⁵I-vasopressin (Reagent B) to all tubes.
- Vortex all tubes and incubate at 4° C for 18-24 hours.
- While stirring continuously add 100 µL double antibody solid phase (Reagent C) to all tubes, except TC tubes.
- Vortex and incubate 30-60 minutes at 4° C.
- Centrifuge all tubes for 15 min. at 1700 g at 4° C.
- Decant or aspirate supernatant.
- Count residue for 2-4 min.

	Assay buffer (D)	Calibrator or sample or contr.	Anti- Vasopressin (A)		¹²⁵ I-Vasopressin (B)		Double antibody solid phase (C)		
Total Counts (TC)	-	-	-	Vortex	200 µL	Vortex	-	Vortex	Centrifuge
Blank (NSB)	500 µL	-	-	and	200 µL	and	100 µL	and	for 15
St. 0 pmol/L	300 µL	-	200 µL	incubate	200 µL	incubate	100 µL	incubate	minutes at
St. 1.9 pmol/L	-	300 µL	200 µL	for 18-24	200 µL	for 18-24	100 µL	for 30-60	1700 g at
St. 3.8 pmol/L	-	300 µL	200 µL	hrs	200 µL	hrs	100 µL	minutes	4° C.
St. 7.5 pmol/L	-	300 µL	200 µL	at +4° C.	200 µL	at +4° C	100 µL	at	Aspirate or
St. 15 pmol/L	-	300 µL	200 µL		200 µL		100 µL	4° C	decant the
St. 30 pmol/L	-	300 µL	200 µL		200 µL		100 µL		supernatant.
St. 60 pmol/L	-	300 µL	200 µL		200 µL		100 µL		Count the
Controls	-	300 µL	200 µL		200 µL		100 µL		residue for 2-4
Unknown sample	-	300 µL	200 µL		200 µL		100 µL		minutes.

C. Calculation

- Subtract the average count rate (cpm) of the NSB from the average count rate (cpm) of the replicates of calibrators, controls and patient samples.
- A calibration curve can be generated by plotting cpm, % B/Bo or % B/T of precipitated bound fraction against the concentration of the vasopressin calibrators.
- To obtain the vasopressin concentration in the extracted patient samples and controls, their cpm, % B/Bo or % B/T of precipitated bound fractions are interpolated now from the generated calibration curve.
- The calibration curve can also be constructed by computer methods. For automated data reduction, both logit/log and Spline methods can be used.
- Correct plasma and urine values for % extraction recovery.

Urine: calculate the 24 hour vasopressin excretion:


Vasopressin pmol/L x dilution x 24 hours urine volume in L.

This calculation provides vasopressin concentration in pmol/24 hours.

D. Calibration Curve Data

	Average cpm	Corrected cpm	% B/Bo	Results (pmol/L)
Total counts	11107			
NSB	555			
Calibrator 0 pmol/L	4770	4215	100	
Calibrator f 1.9 pmol/L	4340	3785	89.8	
Calibrator e 3.8 pmol/L	3739	3184	75.5	
Calibrator d 7.5 pmol/L	3135	2580	61.2	
Calibrator c 15 pmol/L	2199	1644	39.0	
Calibrator b 30 pmol/L	1490	935	22.2	
Calibrator a 60 pmol/L	1142	587	13.9	
Control low	3741	3186	75.6	4.0
Control high	1769	1214	28.8	21.1

13 EXAMPLE OF CALIBRATION CURVE

14 ASSAY CHARACTERISTICS

14.1 Precision									
Within-run					Between-run				
	n	mean pmol/L	SD	% c.v.		n	mean pmol/L	SD	% c.v.
sample A	18	4.17	0.27	6.5	sample A	6	4.38	0.26	6.0
sample B	16	20.2	1.00	4.9	sample B	6	21.4	1.48	6.9

14.2 Calibration
This assay is calibrated against the first international WHO calibrator 77/501

14.3 Recovery			
Two different samples are spiked with different amounts of vasopressin calibrator			
Sample	Expected conc. (pmol/L)	Observed conc. (pmol/L)	% Recovery
A1	9.2	9.8	106
A2	14.3	14.4	101
B1	9.6	10.0	104
B2	15.3	15.1	98

14.4 Specificity

Vasopressin antiserum is raised in rabbits.

The following cross reactivities were measured at 50% binding (B/Bo)

<u>Peptide</u>	<u>% Crossreactivity</u>
Arg ⁸ Vasopressin	100 %
Oxytocin	<0.1 %
Lys ⁸ -Vasopressin	<0.1 %
Desmopressin	<0.1 %
Arg ⁸ Vasotocin	80 %

14.5 Sensitivity

The sensitivity judged as 3 standard deviations change from zero calibrator is 0.5 pmol/L

14.6 Normal Range

Each laboratory should establish its own normal range of expected values.

Plasma: up to 13 pmol/L

Urine: 57 ± 22 pmol/24 hours urine

15 LITERATURE

Sakamoto, M.I. (1987). Atrial Natriuretic Peptide and Vasopressin in Human Plasma Peptides 9, 187-191.

Durr, J.A. (1987). Diabetes Insipidus in Pregnancy Associated with abnormally high circulating Vasopressinase activity. New Engl. J. Med. 316, 1070-1074.

Vokes, T.J. (1988). Disorders of antidiuretic hormone. Endocr. Metab. Clin. North. Am. 17(2) 281.

Miller, M. (1970). Potentiation of Vasopressin Action by Chlorpropamide in Vivo. Endocrinology 86, 1024-1027.

Beardwell, C.G. (1971). Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Endocr. 33, 254-260.

Roberson, Gary L. (1973). Development and Clinical Application of a New method for the Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Invest. 52, 2340-2352.

Uhlich, E. (1975). Radioimmunoassay of Arginine Vasopressine in Human Plasma. Hormon. Metab. res 7, 501-507.

Wagner, H. (1977). Improved Method and its Clinical Application of a Radioimmunoassay of Arginine Vasopressin in Human Serum. Horm. Metab. Res. 9, 223-227.

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Pullan, P.T. (1979). Plasma Vasopressin and Human Neurophysins in Physiological and Pathological States Associated with Changes in Vasopressin Secretion. J. of Clin. Endocrin. metab. 49, 580-578.

Rowe, J.W. (1980). Evidence in Man that Cigarette Smoking induces Vasopressin Release via an Airway-Specific mechanism. J. Clin. Endocrin. Metab. 51, 170-171.

Zerbe, R.A. (1981). Comparison of Plasma Vasopressin Measurements with a Standard Indirect Test in the Differential Diagnosis of Polyuria. New England. J. Med. 305.

Zipser, D. (1981). Dual Effects of Antidiuretic Hormone of Urinary Prostaglandin E2 Excretion in Man. J. Clin. Endocr. Metabolism 65, 522-526.

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VASOPRESSIN-RIA

de

KIPERB319

IN-VITRO-DIAGNOSTIKUM

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

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The DIAsource vasopressin-ria kit contains reagents and instructions for the quantitative measurement of vasopressin in plasma or urine. After solid phase extraction (SPE) or ethanol extraction the plasma vasopressin concentrations are measured by radioimmunoassay (RIA). Urine vasopressin concentrations can be measured directly.

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Urin: 57 ± 22 pmol/24 h Urin

15 LITERATUR

1. Sakamoto, M.I. (1987). Atrial Natriuretic Peptide and Vasopressin in Human Plasma Peptides 9, 187-191.
2. Durr, J.A. (1987). Diabetes Insipidus in Pregnancy Associated with abnormally high circulating Vasopressinase activity. New Engl. J. Med. 316, 1070-1074.
3. Vokes, T.J. (1988). Disorders of antidiuretic hormone. Endocr. Metab. Clin. North. Am. 17(2) 281.
4. Miller, M. (1970). Potentiation of Vasopressin Action by Chlorpropamide in Vivo. Endocrinology 86, 1024-1027.
5. Beardwell, C.G. (1971). Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Endocr. 33, 254-260.
6. Roberson, Gary L. (1973). Development and Clinical Application of a New method for the Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Invest. 52, 2340-2352.
7. Uhlich, E. (1975). Radioimmunoassay of Arginine Vasopressin in Human Plasma. Hormon. Metab. res 7, 501-507.
8. Wagner, H. (1977). Improved Method and its Clinical Application of a Radioimmunoassay of Arginine Vasopressin in Human Serum. Horm. Metab. Res. 9, 223-227.
9. Von Zur Mühlen, A. (1977). Untersuchung zur Stimulation der Vasopressin Sekretion bei Gesunden und Patienten mit Diabetes Insipidus. Schweiz med. Wschr. 107, 1097-1100.
10. Pullan, P.T. (1979). Plasma Vasopressin and Human Neurophysins in Physiological and Pathological States Associated with Changes in Vasopressin Secretion. J. of Clin. Endocrin. metab. 49, 580-578.
11. Rowe, J.W. (1980). Evidence in Man that Cigarette Smoking induces Vasopressin Release via an Airway-Specific mechanism. J. Clin. Endocrin. Metab. 51, 170-171.
12. Zerbe, R.A. (1981). Comparison of Plasma Vasopressin Measurements with a Standard Indirect Test in the Differential Diagnosis of Polyuria. New England. J. Med. 305.
13. Zipser, D. (1981). Dual Effects of Antidiuretic Hormone of Urinary Prostaglandin E2 Excretion in Man. J. Clin. Endocr. Metabolism 65, 522-526.

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Urine: 57 ± 22 pmol/24 ore

15 BIBLIOGRAFIA

1. Sakamoto, M.I. (1987). Atrial Natriuretic Peptide and Vasopressin in Human Plasma Peptides 9, 187-191.
2. Durr, J.A. (1987). Diabetes Insipidus in Pregnancy Associated with abnormally high circulating Vasopressinase activity. New Engl. J. Med. 316, 1070-1074.
3. Vokes, T.J. (1988). Disorders of antidiuretic hormone. Endocr. Metab. Clin. North. Am. 17(2) 281.
4. Miller, M. (1970). Potentiation of Vasopressin Action by Chlorpropamide in Vivo. Endocrinology 86, 1024-1027.
5. Beardwell, C.G. (1971). Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Endocr. 33, 254-260.
6. Roberson, Gary L. (1973). Development and Clinical Application of a New method for the Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Invest. 52, 2340-2352.
7. Uhlich, E. (1975). Radioimmunoassay of Arginine Vasopressine in Human Plasma. Hormon. Metab. res 7, 501-507.
8. Wagner, H. (1977). Improved Method and its Clinical Application of a Radioimmunoassay of Arginine Vasopressin in Human Serum. Horm. Metab. Res. 9, 223-227.
9. Von Zur Mühlen, A. (1977). Untersuchung zur Stimulation der Vasopressin Sekretion bei Gesunden und Patienten mit Diabetes Insipidus. Schweiz med. Wschr. 107, 1097-1100.
10. Pullan, P.T. (1979). Plasma Vasopressin and Human Neurophysins in Physiological and Pathological States Associated with Changes in Vasopressin Secretion. J. of Clin. Endocrin. metab. 49, 580-578.
11. Rowe, J.W. (1980). Evidence in Man that Cigarette Smoking induces Vasopressin Release via an Airway-Specific mechanism. J. Clin. Endocrin. Metab. 51, 170-171.
12. Zerbe, R.A. (1981). Comparison of Plasma Vasopressin Measurements with a Standard Indirect Test in the Differential Diagnosis of Polyuria. New England. J. Med. 305.
13. Zipser, D. (1981). Dual Effects of Antidiuretic Hormone of Urinary Prostaglandin E2 Excretion in Man. J. Clin. Endocr. Metabolism 65, 522-526.

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Urin: 57 ± 22 pmol/dygnsurin

15 REFERENSER

1. Sakamoto, M.I. (1987). Atrial Natriuretic Peptide and Vasopressin in Human Plasma Peptides 9, 187-191.
2. Durr, J.A. (1987). Diabetes Insipidus in Pregnancy Associated with abnormally high circulating Vasopressinase activity. New Engl. J. Med. 316, 1070-1074.
3. Vokes, T.J. (1988). Disorders of antidiuretic hormone. Endocr. Metb. Clin. North. Am. 17(2) 281.
4. Miller, M. (1970). Potentiation of Vasopressin Action by Chlorpropamide in Vivo. Endocrinology 86, 1024-1027.
5. Beardwell, C.G. (1971). Radioimmunoassay of Arginine Vasopressin in Human Plasma J. Clin. Endocr. 33, 254-260.
6. Roberson, Gary L. (1973). Development and Clinical Application of a New method for the Radioimmunoassay of Arginine Vasopressin in Human Plasma. J. Clin. Invest. 52, 2340-2352.
7. Uhlich, E. (1975). Radioimmunoassay of Arginine Vasopressine in Human Plasma. Hormon. Metab. res 7, 501-507.
8. Wagner, H. (1977). Improved Method and its Clinical Application of a Radioimmunoassay of Arginine Vasopressin in Human Serum. Horm. Metab. Res. 9, 223-227.
9. Von Zur Mühlen, A. (1977). Untersuchung zur Stimulation der Vasopressin Sekretion bei Gesunden und Patienten mit Diabetes Insipidus. Schweiz med. Wschr. 107, 1097-1100.
10. Pullan, P.T. (1979). Plasma Vasopressin and Human Neurophysins in Physiological and Pathological States Associated with Changes in Vasopressin Secretion. J. of Clin. Endocrin. metab. 49, 580-578.
11. Rowe, J.W. (1980). Evidence in Man that Cigarette Smoking induces Vasopressin Release via an Airway-Specific mechanism. J. Clin. Endocrin. Metab. 51, 170-171.
12. Zerbe, R.A. (1981). Comparison of Plasma Vasopressin Measurements with a Standard Indirect Test in the Differential Diagnosis of Polyuria. New England. J. Med. 305.
13. Zipser, D. (1981). Dual Effects of Antidiuretic Hormone of Urinary Prostaglandin E₂ Excretion in Man. J. Clin. Endocr. Metabolism 65, 522-526.

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