DHT-RIA

KIPI9900

DHT-RIA

Radioimmunoassay for the Quantitative Determination of 5 α -Dihydrotestosterone in Human Serum or Plasma

KIP19900

IN VITRO DIAGNOSTIC USE

en

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

1. INTENDED USE:

For IN VITRO determination of serum or plasma DHT levels.


5 α -Dihydrotestosterone (DHT; 17 β -Hydroxy-5 α -androstan-3-one), the most potent naturally-occurring androgen, is produced from testosterone through the action of cholesteryl 5 α -reductase [1]. The concentrations of 5 α -reductase are highest in certain peripheral tissues, including genital skin and hair follicles, and is localized intracellularly in apparent association with the nuclear membrane. DHT exerts its biological action by intracellular binding to the androgen receptor; this complex is then transferred to the nucleus where DNA-binding occurs with resultant effects on DNA transcription [2]. Most of the residual DHT undergoes intracellular metabolism to 3 α -androstenediol and 3 α -androstenediol glucuronide. Only a small proportion of DHT escapes into the peripheral circulation, where it is present primarily complexed to sex-hormone binding globulin [3].

Testosterone causes virilization of the Wolffian ducts during fetal life, while DHT is responsible for the development of the male external genitalia and prostate, and is primarily responsible for the physical changes which occur during male sexual maturation. An autosomal-recessive genetic deficiency of 5 α -reductase, sometimes called male pseudohermaphroditism or pseudovaginal perineoscrotal hypospadias, leads to inadequate differentiation of DHT-dependent peripheral tissues. Male infants with this disorder have ambiguous genitalia and are often raised as females, although significant virilization may occur later in life presumably due to the natural increase in testosterone levels [2].

Measurement of DHT concentrations can be complicated by antibody cross-reactivity to testosterone. The DHT Double Antibody radioimmunoassay utilizes a sample oxidation / extraction procedure to remove most of the testosterone, coupled with a specific immunoassay for DHT.

2. PRINCIPLE OF THE METHOD:

The 5 α dihydrotestosterone (DHT) DA RIA obeys the law of mass action according to the following equation :

The procedure follows the basic principle of radioimmunoassay where there is competition between the radioactive and non-radioactive antigen for a fixed number of antibody binding sites. The amount of (I-125)-labeled analyte bound to the antibody is inversely proportional to the concentration of unlabeled analyte present. The separation of free and bound antigen is achieved by using a double antibody system.

3. MATERIAL PROVIDED AND STORAGE:

Stored at 2 - 8°C, the material can be used up to the expiration date printed on each label.

3.1. CAL

Stock solution in ethyl alcohol, containing 100ng/ml.

Preparation of the curve.

Concentration (pg/ml)	Calibrator (ml)	Buffer (ml)	
2500	Stock : 0.1	3.9	= Cal 7
1000	Stock : 0.1	9.9	= Cal 6
500	Cal 6 : 1.0	1.0	= Cal 5
200	Cal 6 : 0.2	0.8	= Cal 4
100	Cal 6 : 0.1	0.9	= Cal 3
50	Cal 6 : 0.1	1.9	= Cal 2
25	Cal 2 : 0.5	0.5	= Cal 1

The calibrator 0 pg/ml is the diluting buffer.

Store at 2-8°C for up to 3 weeks.

3.2. Ag 125I

DHT ¹²⁵I (yellow)

1 vial, 10.5 ml, containing < 5 μ Ci (< 185 kBq) of ¹²⁵I labelled DHT in protein-based buffer containing sodium azide (< 0.1%) as a preservative. Store at 2-8°C.

3.3. ANTISERUM

Ab DHT

2 vials, 5.2 ml, containing lyophilised anti-DHT rabbit serum (blue coloured) and sodium azide (<0.1%) as a preservative. Reconstitute by adding 5.2 ml of buffer. Store at 2-8°C.

3.4. PEG

1 vial, 105 ml, containing sheep anti-rabbit gamma globulin in a buffer, polyethylene glycol as a precipitating aid and sodium azide (<0.1 %) as a preservative. Store at 2-8°C.

3.5. KMnO₄

1 vial, 5 ml, containing KMnO₄ solution. Store at 2-8°C.

3.6. BUF

1 vial, 100 ml, Phosphate Buffered Saline 0.1M containing sodium azide (<0.1%) as a preservative. Store at 2-8°C.

3.7. CONTROL N

N = 2

2 vials horse serum containing preservative (NaN₃ < 0.1%). The control sera are to be assayed along with the patient samples. The ranges for the control sera are printed on the vial labels.

Before use, reconstitute the contents of the controls with 2 ml of distilled water. Mix gently to avoid foaming. Wait at least 15 minutes after reconstitution before dispensing. After reconstitution, the solutions are stable for 3 weeks at 2-8°C. or for a longer period if stored at -20°C.

All reagents and samples must be allowed to reach room temperature and mixed thoroughly by gentle inversion before use.

4. MATERIAL REQUIRED BUT NOT PROVIDED :

- Refrigerator (2-8°C).
- 12 x 75 mm plastic and glass test tubes.
- Test tube rack for 12 x 75 mm tubes.
- Precision pipette to deliver 100 μ l.
- Precision repeating pipette to deliver 100 μ l, 200 μ l, 500 μ l and 1 ml.
- Controlled-temperature water bath, 37 $\pm$ 2°C.
- Centrifuge (preferably refrigerated) capable of 1500 x g.
- Absorbent material for blotting tubes.
- Gamma counter.
- Semi-log (linear-log) graph paper.
- Vortex mixer.
- Diethyl ether HPLC grade
- Hexane HPLC grade
- Ethanol HPLC grade

5. METHODOLOGY

5.1. Collection and handling of blood samples:

The blood sample may be collected into a dry tube or one containing an anticoagulant. If heparin is used, only the minimum required should be added to avoid clotting.

After separation from the red blood cells, plasma or serum samples may be assayed immediately, within 24 hours if stored at 2 - 8°C, or later, after period up to several months if stored at -20°C. Repeating freezing and thawing must be avoided.

5.2. Preparation of the samples:

1. Extract 300 μ l of serum or plasma and controls by adding 3 ml of Diethyl Ether in **glass** tubes (duplicate).
2. Vortex or shake for 2 minutes.
3. Centrifuge. The recommend speed is 1000 g for 5 minutes.
4. Freeze aqueous layer in a freezer.
5. Transfer solvent into glass tubes and evaporate to dryness.
6. Add 0.5 ml of buffer and vortex for 30 seconds.
7. Add 50 μ l of the oxidation solution (KMnO₄)
8. Incubate at 37°C for 20 minutes.
9. Extract DHT by adding 3 ml of a mixture of Hexane:Ethanol (98:2). Vortex or shake for 2 minutes. Centrifuge. The recommend speed is 1000 g for 5 minutes.
10. Transfer 2 ml of the organic layer into glass tubes and evaporate to dryness.

The tubes are now ready for the RIA determination.

5.3. Assay procedure:

Reagents stored at 2° - 8° C. must be brought at room temperature prior to use. Do not mix reagents of different lots. Label the tubes for T (Total Counts), calibrators, samples and control sera.

Perform the assay in duplicate. Calibrators, controls and samples must be assayed at the same time.

1. Label and arrange tubes in duplicate for total counts, non-specific binding, calibrators, controls and unknowns.
2. Add 200 μ l of the Buffer into extracted tubes and incubate for 20 minutes at room temperature. Vortex now and again.
3. Add 200 μ l of the prepared calibrators into glass tubes.

4. Add 100 µl of the antibody anti-DHT to all tubes except Total counts and NSB.
5. Add 100 µl of DHT ¹²⁵I Reagent to each tube.
6. Vortex, cover and incubate all tubes at 37°C for 1 hour.
7. Add 1 ml of Precipitating Reagent to all tubes except the Total Count tubes and immediately vortex all tubes. THE PRECIPITATING REAGENT MUST BE MIXED THOROUGHLY BEFORE USE.
8. Incubate at room temperature for 10 to 15 minutes.
9. Centrifuge all tubes except the Total Count tubes for 15 minutes at 1500 x g.
10. Aspirate or decant and discard the supernatant by inverting the tubes. Blot the tubes to remove any droplets adhering to the rim before returning the tubes to the upright position.
11. Count the pellet remaining in the tubes in a gamma counter for one minute.
12. Construct a calibration curve for each test run as described below.

5.4. Data processing:

Determine the mean count rate for each set of duplicate tubes.
Calculate the ratio B/B0 as follows:

$$B/B0 \% = [\text{Cal or Sample cpm} / B0 (\text{Cal 0}) \text{cpm}] \times 100$$

Draw the calibration curve on semi logarithmic paper by plotting the ratio B/B0 % (linear scale) obtained for each calibrator versus its respective concentration expressed in pg/ml (logarithmic scale). DHT concentrations in samples may be read directly from the calibration curve.

If a computer is used to calculate the results, the data can be fitted to the appropriate equation: 4 PL weighted.

5.5. Example of a typical assay:

	Contents (pg/ml)	Net average bound counts (cpm)	B/Bo (%)	DHT (ng/ml)
Total counts	-	36935	-	-
NSB	-	888	2.4	-
Cal 0	0	27000	100	-
Cal 1	25	24747	91.4	-
Cal 2	50	23559	86.8	-
Cal 3	100	19816	72.5	-
Cal 4	200	15742	56.9	-
Cal 5	500	9844	34.3	-
Cal 6	1000	6277	20.6	-
Cal 7	2500	3319	9.3	-
Sample 1	-	7750	27	0.6

Example of a typical assay, do not use for calculations

6. PERFORMANCE CHARACTERISTICS:

6.1. Specificity

Steroid	% Cross-reactivity
DHT	100.00
Testosterone	0.19
Estrilol	0.03
Estradiol	0.005
Progesterone	0.004
Free DHEA	N.D.

6.2. Minimum detectable concentration of DHT:

The minimum detectable concentration has been assayed at ± 7 pg/tube (± 20 pg/ml) and corresponds to the concentration given by two standard deviations below the mean cpm of 20 replicate determinations of the zero calibrators.

6.3. Recovery test:

When sera of known DHT contents have their DHT supplemented by addition of DHT in equal volumes (1/1), a satisfactory correlation between theoretical and assayed DHT is obtained.

Added DHT (pg/ml)	0	50	100	200	500	1000
Theoretical	-	180.5	205.5	255.6	405.6	655.6
Assayed DHT (pg/ml)	311.1	185.1	210.5	260.4	429.5	672.8
% recovery	-	102.5	102.4	101.9	105.9	102.6

6.4. Dilution test:

The dilution test indicates that there is immunological identity between the DHT present in the sample and the DHT used to calibrate the calibration curve.

Dilution Factor	1	1/2	1/4	1/8	1/16
Assayed DHT (pg/ml)	737	373.4	192.5	88	52.1
Expected DHT (pg/ml)	-	369	184	92	46

6.5. Reproducibility:

	Within assay variation 10 replicates		Between assay variation 12 replications	
	Mean (pg/ml)	%CV	Mean (pg/ml)	%CV
Pool 1	159.3	6	60.5	18.6
Pool 2	411.6	4.8	201.4	12.3
Pool 3	700.1	4.8	653.9	8.7

7. LIMITATION OF THE PROCEDURE

- 7.1. The results obtained from this or any other diagnostic kit should be used and interpreted only in the context of an overall clinical picture.
- 7.2. Do not use lipemic, haemolyzed, icteric or turbid specimens.

8. EXPECTED VALUES

It is recommended that each laboratory establishes its own reference values.

		DHT (pg/ml)
Males	Pre-puberty	50
	Adult	250 - 1000
Females	Follicular phase	50 - 200
	Luteal phase	100 - 300
	Menopausal	< 100

9. WARNING AND PRECAUTION

For IN VITRO DIAGNOSTIC use only

CAUTION: Radioactive material

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

WARNING: Sodium azide

Some components contain sodium azide as preservative agent (NaN₃ < 0.1%). Dispose of the reagents by flushing with large amount of water through the plumbing system.

WARNING: Potentially infectious material

Handle all components (and all patient samples) as if capable of transmitting viral diseases such as hepatitis B and C and the acquired immunodeficiency syndrome (AIDS).

Source material derived from human body fluids or organs and used in the preparation of this kit were tested and found negative for hepatitis B and C surface antigen by immunoassay. However, no known test can guarantee that such material does not contain the causative agent of viral hepatitis.

Likewise, all human materials used in the preparation of this kit were screened for the presence of the HIV-1 and 2 antibodies and found negative by enzyme-immunoassay. However, absence of these antibodies cannot guarantee the absence of the viral agent responsible for the acquired immunodeficiency syndrome.

WARNING: Ethanol

The calibrator stock solution contains ethanol, which is highly flammable.

10. BIBLIOGRAPHY

1. Miller WL: Molecular biology of steroid hormone synthesis. *Endocrin Rev* 9:295-318, 1988.
2. Pang S, Riddick L: Hirsutism. IN Lifshitz F: *Pediatric Endocrinology, A Clinical Guide, second edition*. Marcel Dekker, New York, 1990, pp. 259-291.
3. Wilson JD: Syndromes of androgen resistance. *Biol Reprod* 46:168-173, 1992.
4. Yalow R, Berson S: Introduction and general considerations. IN: Odell WD, Doughday WH (eds): *Principles of Competitive Protein Binding Assays*. J. B. Lippincott Co., Philadelphia, 1971, pp. 1-19.
5. Burtis CA, Ashwood ER: *Tietz Textbook of Clinical Chemistry, 2nd edition*. W.B. Saunders Company, Philadelphia, 1994, p.1821.

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