

CE

# Leptin-Ria-CT

## KIPMR44

DIAsource ImmunoAssays S.A. - Rue de l'Industrie, 8 - B-1400 Nivelles - Belgium



CE

Human-Leptin-RIA-CT

Radioimmunoassay with Coated Tubes for the Quantitative Detection of

human Leptin in serum and plasma

KIPMR44

### IN VITRO DIAGNOSTIC USE

Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

### **TECHNICAL FEATURES+APPLICATIONS**

- Measures total Leptin concentration in serum, EDTA/Heparin plasma
- Calibrated against the WHO International Standard: NIBSC Code 97/594
- High affinity antibodies, no influence of Leptin Bindingprotein
- ready-to-use calibrators with 0-64 ng/ml are
- Small serum volumina sufficient, measurement in capillar blood possible
- Easy handling

### INTRODUCTION

Leptin, the product of the ob gene (1,2), is a recently discovered single-chain proteohormone with a molecular weight of 16 kD, which is thought to play a key role in the regulation of body weight. Its amino acid sequence exhibits no major homologies with other proteins (1). Leptin is almost exclusively produced by differentiated adipocytes (3-5). It acts on the central nervous system, in particular the hypothalamus, thereby suppressing food intake and stimulating energy expenditure (2,6-9). Leptin receptors - alternatively spliced forms exist that differ in length - belong to the cytokine class I receptor family (10-12). They are found ubiquitously in the body (10, 11, 13, 14) indicating a general role of leptin, which is currently not fully understood. A circulating form of the leptin receptor exists which acts as one of several leptin binding proteins (15). Besides its metabolic effects, leptin was shown to have a strong influence on a number of endocrine axes. In male mice, leptin prevented the starvation-induced delay in ovulation (16). Ob/ob mice, which are leptin deficient due to an ob gene mutation, are infertile. This defect could be corrected by administration of leptin, but not through weight loss due to fasting (17), suggesting that leptin is pivotal for reproductive functions.

All these actions may, at least in part, be explained by the suppressive effect of leptin on neuropeptide Y (NPY) expression and secretion by neurons in the arcuate nucleus (6,18,19). NPY is a strong stimulator of appetite (20,21) and is known to be involved in the regulation of various pituitary hormones, e.g. suppression of GH through stimulation of somatostatin (22,23), suppression of gonadotropins (23) or stimulation of the pituitary adrenal axis (21).

The most important variable that determines circulating leptin levels is body fat mass (24-26). Obviously, under conditions of regular eating cycles, leptin reflects the proportion of adipose tissue (27) showing an exponential relationship (37). This constitutive synthesis of leptin is modulated by a number of non-hormonal and hormonal variables. Stimulators in both rodents and humans are overfeeding (28,29), insulin (3,5,30-33) and glucocorticoids (5,34-36). Suppression has been shown for fasting (27), cAMP and beta-3-adrenoceptor agonists (35). From these findings it becomes clear that leptin is an integral component of various metabolic and endocrine feedback loops (38).

For clinical purposes, it is important to note that serum leptin levels show a moderate circadian variation with a peak during the night at about 2 a.m. (37). The leptin values at this time are about 30 to 100 % higher than the levels measured in the morning or early afternoon. This variation together with the influence of food intake needs to be taken into account, when blood samples are collected.

Under fairly standardized conditions, i.e. normal eating cycles and blood sampling in the morning or early afternoon, a single leptin measurement is informative.

For the appropriate interpretation of measured leptin levels, reference ranges are required. Because body fat mass is the major confounding variable, these ranges should be referred to measures of the percentage body fat such as body mass index (BMI) or percent body fat determined by, e.g., bioelectric impedance assessment (BIA). Leptin levels are higher in females than in males (38,39) and an age dependence was shown in children and adolescents (40). Therefore, reference ranges referring to measures of body fat should be stratified according to gender and pubertal development.

Leptin levels are high in most obese patients suggesting the presence of leptin insensitivity (20,26,37,38,41,42). In a small percentage of patients, however, leptin levels have been found inappropriately low with respect to their fat mass. It remains for future studies to prove that these patients represent a new pathophysiologic entity: leptin deficiency. Since leptin has also been shown to be of great importance for reproductive functions, possible new pathophysiologic mechanisms may be discovered relating infertility to insufficient leptin production.

The discovery of leptin has released an avalanche of research activities seeking to understand the regulation and actions of this new hormone. Most importantly, it has provided a key to better understand the physiology of body weight regulation and to unveil possible pathophysiologic mechanisms in both obesity and eating disorders. Further, it may provide new insights into certain causes of infertility.

### INTENDED USE

The widespread importance makes leptin an interesting parameter for physicians dealing with **metabolic syndrome**, **obesity**, **cachexia** and other **metabolic disturbances**, as **diabetologists**, **endocrinologists**, **gynaecologists**, **andrologists**, and **psychiatrists** treating patients with **eating disorders**.

This radio immunoassay kit is suited for measuring human leptin in serum or plasma, and conditioned adipocyte culture media for scientific and diagnostic purposes.

Measuring leptin in anorectic or cachectic patients, young children or in specimen other than serum, such as urine, cerebrospinal fluid, and certain cell culture media, is also possible with this kit.

Under conditions of normal eating cycles, measurement in a single blood sample collected in the morning or early afternoon is sufficient.

The comparison with BMI-related reference ranges may be useful to detect conditions of relative **leptin deficiency** as a possible cause of obesity or provide an indication for **leptin resistance** respectively.

Due to its high correlation with body fat mass leptin measurements under standardized conditions may be used as a simple and inexpensive test for determination of **body fat**.

### PRINCIPLE

For the Radioimmunoassay for the determination of human Leptin a polyclonal rabbit-antibody of high specificity is used. Leptin is measured quantitatively and binding protein does not influence the test results. Tubes were coated with streptavidin and separation is facilitated by biotinylated anti-rabbit-antibody.

Calibrators are prepared of recombinant human leptin in concentrations between 1, 2, 4, 8, 16, 32 and 64 ng/ml.

### Calibration of the Assay

The Assays is calibrated against the International Standard for human leptin. The standard preparation of the WHO with code 97/594 (y) is available from the NIBSC (z). One ampoule of the preparation, reconstituted in 1 ml Assay Buffer, will be quantified with the nominal content of 5000 ng human leptin.

### PERFORMANCE CHARACTERISTICS

### Sensitivity

The analytical sensitivity of the assay yields 0.1 ng/ml (2x SD of zero calibrators)

### Intra-Assay-Variation

	Number of determinations	Mean value [ng/ml]	Standard Deviation [ng/ml]	VC%
Sample 1	6	19.67	0.87	4.4
Sample 2	6	6.95	0.33	4.8

### Inter-Assay-Variation

	Mean value (ng/ml)	Standard deviation	VC%
Sample 1	8.975	0.48	5.3
Sample 2	6.3	0.32	5.1
Sample 3	20.9	1.05	5.0

### Linearity

The Leptin RIA Leptin KIPMR44 can be used for diluted samples. The linearity of serum dilutions is excellent (Fig.9).

### Recovery

Serum spiking experiments with recombinant human leptin yielded a recovery of 97% (± 2%).

### SPECIMEN COLLECTION, PREPARATION, AND STORAGE

**Serum** as well as EDTA/Heparin **plasma samples** are suitable (significant deviation of hLeptin levels in corresponding Serum or -Plasma samples were not found). Conditioned adipocyte culture medium was found to be suitable. An external sample preparation prior to assay is not required. Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen at -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although Leptin levels were found to be unaffected by few cycles (5x) in our experiments.

The high sensitivity of the assay allows measurement of hLeptin in small sample volumes. Because of the wide effective range of this RIA kit a preparative sample dilution is generally not necessary. For most of the determinations (serum or plasma samples, and no extreme values expected) **the use of undiluted samples 25 µl per tube**, **should be appropriate.** In case Leptin levels of more than 64 ng/ml are expected, eg. obese patients (BMI>35), the sample should and can be diluted, e.g. 1:10.

The hLeptin concentrations may be completely different in body fluids of human origin other than serum or cell culture supernatants.

### REAGENTS PROVIDED

ASS	BUF	Assay buffer
		(1 bottle, 125 ml; ready for use)
Ab	ВЮТ	Capture-antibody: anti-rabbit IgG, biotin-conjugated. Reconstitute in 7 ml Assay buffer. (1 bottle, 7 ml; lyophilized)
Ab		Specific antibody: rabbit-anti-human leptin antibody. Reconstitute in 7 ml Assay buffer. (1 bottle, 7 ml; lyophilized)
Ag	1251	Tracer: <sup>125</sup> I-Ieptin; < 2.25 μCi or < 83 kBq respectively. Reconstitute in 13 ml Assay buffer. (1 bottle, 13 ml; lyophilized, red coloured)
CAL	N	Calibrators 1-7. (7 vials, 750 $\mu l$ each; ready for use). Contain recombinant Leptin. Use 25 $\mu l$ / tube. See exact values on the vial label



Control 1

(1 vial, 0.75 ml; human serum; lyophilized). Reconstitute with 750 µl Assay buffer. See exact values on the vial label

Tubes (125 streptavidin-coated tubes)

### MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes (50, 100 and 500 µl) Micropipettes and multichannel pipettes with disposable plastic tips Vortex-mixer

Device to aspirate the fluid from the tubes (recommended because of the potential danger of radioactivity and infection by human samples) Shaking Device : 350 rpm

Gamma Counter

### **REAGENT PREPARATION**

In conducting the assay, follow strictly the test protocol. Room temperature incubation means: Incubation at 20 - 25°C.

Reagents with different lot numbers should not be mixed. All reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

Control 1, Capture Antibody, Specific Antibody and Tracer have to be reconstituted in Assay Buffer. It is recommended to keep the reconstituted reagents at room temperature for 30 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The shelf life of the components after opening is not affected, if used appropriately. Reconstituted Components (Control 1, Capture Antibody, Specific Antibody and Tracer) should be stored at -20°C (or below). Repeated freeze-thaw cycles have to be avoided.

Before use, all kit components should be brought to room temperature. Precipitates, possible in buffers, should be dissolved before use through mixing and warming.

### WARNINGS AND PRECAUTIONS

### For in-vitro diagnostic use only. For professional use only.

Possession and use of the kit is subject to the regulations of the national nuclear regulatory authorities.

Reagents with different lot numbers should not be mixed.

Reagents contain Sodium-Azide as preservative, however, highly diluted (0.02%). Avoid any skin contact.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to room temperature at 20 - 25°C. Precipitates in buffers should be dissolved before use by thorough mixing and warming. Temperature WILL affect the the assay. However, Values for the patient samples will not be affected.

Caution: This kit contains material of human and/or animal origin. Source human serum for the Control 1 provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibody. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

**Radioactivity** - Before ordering or using radioactive materials, it is necessary to take the appropriate actions to ensure compliance with national regulations governing their use. Local rules in each establishment, which define actions and behaviour in the radioactivity working areas, should also be adhered to. The advice given here does not replace any local rules, instructions or training in the establishment, or advice from the radiation protection advisers. It is important to follow the code of good laboratory practice in addition to the specific precautions relating to the radionuclide I-125 used.

lodine-125 has a radioactive half-life T1/2 of 60 days and emits 35.5 keV gamma radiation, 27 – 32 keV x-rays and no beta radiation. Shielding is effective done by lead, first half value layer is 0.02 mm lead, reduction to 10 % is made by 0.2 mm.

To reduce the radiation dose time spent handling radioactivity should be minimized (plan ahead), and distance from source of radiation should be maximized (doubling the distance from the source quarters the radiation dose).

Formation of aerosols, e.g. by improper opening and mixing of vials or pipetting of solutions which may cause minute droplets of radioactivity become airborne, is a hazard and should be avoided.

Solutions containing iodine should not be made acidic, because this might lead to the formation of volatile elemental iodine.

As some iodo-compounds can penetrate rubber gloves, it is advisable to wear two pairs, or polyethylene gloves over rubber.

For cleaning of contaminated areas or equipment, the lodine-125 should be rendered chemically stable by using alkaline sodium thiosulphate solution together with paper or cellulose tissue.

### General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

### The handling of radioactive and potentially infectious material must comply with the following guidelines:

The material should be stored and used in a special designated area.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Avoid direct contact with these materials by wearing laboratory coats and disposable gloves.

Spilled material must be wiped off immediately. Clean contaminated areas and equipment with a suitable detergent.

Unused radioactive material and radioactive waste should be disposed according to the recommendations of the national regulatory authorities.

### ASSAY PROCEDURE

NOTES: All determinations (Calibrators, Control 1 and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Calibrators, Control 1 and the samples should be pipette as fast as possible.

### Flow Chart of Assay Protocol

#	Tube	Assay Buffer, Calibrators, Control 1, Samples µl	Capture-antibody µl	Specific antibody µI	Tracer µI
1,2	тс				100
3,4	B <sub>0</sub>	Assay Buffer :25	50	50	100
5-18	Calibrators	Cal 1-7: 25	50	50	100
19,20	Control 1	Control 1: 25	50	50	100
21,22	Sample 1	25	50	50	100
23,24	Sample 2	25	50	50	100
etc.					

1) Labelling of the assay tubes should be done in the following order:

1, 2	total counts ( <b>TC</b> ),
3, 4	zero calibrator ( <b>B</b> <sub>0</sub> ),
5 - 18	duplicates of calibrators (1-7),
19, 20	duplicates of control (1),
21. 22 etc.	duplicates of samples.

- 2 ) Add **25** μl of **Assay Buffer** to tubes 3 and 4. Add **25** μl of **Calibrators** to tubes 5 - 18: 5, 6 calibrator **1** 
  - 7, 8 calibrator **2**, etc.
- 3) Add 25  $\mu$ I of control 1 to tubes 19 and 20.
- 4) Add 25 µl of the first sample to tubes 21 and 22,etc.
- 5) Add **50** µl capture-antibody beginning with tube 3.
- 6) Add 50 µl specific antibody beginning with tube 3.
- 7) Add 100 µl tracer to all tubes.
- 8) Remove tubes 1 and 2 (total counts) or seal with a stopper.
- 9) Shake the tubes on a shaking device (350 rpm) overnight (at least 15 h) at room temperature. (steps 9 to 12 may not be performed with tubes 1 & 2)
- 10) Aspirate the liquid (except tubes 1 and 2 !) completely.

Take care that the coating of the tubes remains intact. Depending on laboratory equipment and common laboratory practice, aspiration of the liquid can be replaced by careful decantation.

- 11) Add **500** µl of **assay buffer** to the tubes (except tubes 1 and 2 !).
- 12) Aspirate the liquid (see step 10).
- 13) Count the radioactivity of **all** tubes.

### CALCULATION OF RESULTS

### Establishing the Calibration Curve

Calibrator	B <sub>0</sub>	1	2	3	4	5	6	7
ng/ml	0	1	2	4	8	16	32	64

1. Calculate the average counts of each pair of tubes. This equals to the binding value B.

2. The mean of the zero standard (tubes 3 and 4) equals  $B_0$ .

3. Calculate the percent bound (%  $B/B_0$ ) by dividing the B values by  $B_0$ : % $B/B_0 = B/B_0 x100\%$ .

- 4. Plot % B/B<sub>0</sub> versus the calibrator concentrations on a semi-logarithmic or logit-log paper respectively. For convenience, it is recommended to use computer assisted data reduction programs.
- Calculate the 'percentage bound of the zero calibrator': 5 Average of tubes 3, 4 (B<sub>0</sub>) divided by average of tubes 1, 2 (total counts) multiplied by 100%. It should be  $\%B_0/TC > 20\%$ .

### Evaluation of sample concentrations:

Read the concentration value (abscissa) corresponding to the %B/B<sub>0</sub> of the sample as in the example given below:

average o average o	counts of zero calibrator $(B_0)$ : counts of sample tube:	8648 cpm 4283 cpm
%B/B_	= sample counts / $B_0$ x 100%	

= 4283 / 8648 x 100% = 0.495 x 100% = 49.5%

For a 49.5% value on the y-axis (ordinate) a value of 7.28 ng/ml on the x-axis (abscissa) may be obtained. In case of dilution of the samples, multiplication of this value determined graphically or by the aid of a computer program with the dilution factor gives the leptin concentration of the sample.

### Concentration of the control:

The control 1 should fit with the labelled concentration range. The measured Leptin concentration are only valid, if the measured value of the control is in the labelled concentration range. Otherwise process analysis is required and depending on the results the measured values are accepted or not.

### EXPECTED NORMAL VALUES

Serum leptin levels are mainly determined by body fat mass with low levels in lean individuals and high levels in obese subjects. In addition, there is a clear gender difference with higher levels in females at a given percentage body fat. Further, leptin levels are influenced by pubertal development. Any attempt, therefore, to give ranges of expected leptin levels must account for these relationships.

Various methods for the estimation of body fat are available such as calculation of body mass index (weight (kg) divided by the square of height (m)) (BMI), bioelectric impedance assessment (BIA) or total body dual energy x-ray absorptiometry (DXA). Although the accuracy of BMI with respect to reflecting true fat mass is inferior to other more sophisticated methods such as BIA or DXA, BMI provides a number of advantages:

It is independent of the regression models applied.

- 1) 2) It is easy to determine, only weight and height measurements are required.
- 3) It is retrospectively mostly available.

4) It is the most precise measure during short-term changes of fat mass, e.g. during fasting.

Therefore, the following expectation ranges of serum leptin levels were referred to BMI as the major confounding independent variable and were stratified according to gender and pubertal development (45; see figures 1-8 and tables 4 - 11). After the age of 20 years, no significant age dependence was observed. These gender and age adjusted expectation ranges may be used to compare a measured leptin level at a given BMI with normal subjects to detect pathologic deviations.

The best-fit regression lines for the various subgroups are exponential curves of the form leptin =  $a \cdot e^{(b \cdot BMI)}$ . The 5th and 95th percentiles are given by the following equations:

leptin =  $a \cdot e^{(b \cdot BMI - c)}$  and leptin =  $a \cdot e^{(b \cdot BMI + c)}$ 

In a semi-logarithmic plot (y-axis = log leptin), these curves give straight lines. The values for a, b and c are given in table 3 according to gender and pubertal stage and also for adults. Using these values, the expectation ranges of leptin levels can be easily extended to lower or higher BMI ranges if required.

Example:

The 50th percentile for boys at Tanner stages 3 and 4 is given by the following curve: leptin =  $0.0181 \cdot e^{(0,2067 \cdot BMI)}$ 

	(0.0007 DMIL 4.4040)
The 5th percentile is given by:	$leptin = 0.0181 \cdot \mathbf{e}^{(0,2067 \cdot BMI - 1,1919)}$
	, (0.0007
and the 95th percentile is given by:	$leptin = 0.0181 \cdot \mathbf{e}^{(0,2067 \cdot BWI + 1,1919)}$

In a semi-logarithmic plot, these lines are parallel with an equal distance to the 50th percentile.

### Calculation of standard deviation scores (SDS; Z-scores)

A convenient method to detect any deviation of a measured leptin level from the corresponding reference range is to calculate its standard deviation score by relating the leptin level at the patient's BMI to the average leptin value of the corresponding sex and age group and expressing its deviation by the x-fold standard deviation. This method may be considered as normalization to the normal reference cohort. Thus, the leptin values can be adjusted for BMI, gender and pubertal stage/age (i.e., the influence of gender, age and BMI are removed) and may be pooled for further analysis.

Accounting for the logarithmic distribution of leptin levels, the leptin SDS can be calculated by the following equation:

leptin SDS = (In(leptin) - In(a) - b·BMI) / d

In this equation, In represents the natural logarithm (referring to the basis e). The constants a, b and d are given in table 3 according to gender and pubertal stage/age.

### Example:

A boy at Tanner stage 3, BMI = 25 kg/m<sup>2</sup>, measured leptin concentration = 5 ng/ml. leptin SDS = (ln(5) - ln (0.0181) - 0.2067.25) / 0.6850 = (1.6094 - (-4.0118) - 5.1675) / 0.6850 = 0.66

### Estimation of optimal dilution of samples

Because serum leptin levels vary widely over several orders of magnitude, depending mainly on the body fat mass, adequate dilution might be a prerequisite for precise measurement. Therefore, samples should be diluted such, that the leptin concentration is near this value in order to take maximum advantage of the precision of the kit. The reference ranges provide a useful tool for a good estimation of the expected leptin value according to BMI, gender and age.

### Example:

Adult woman, BMI = 45 kg/m<sup>2</sup> (130 kgs,1.70 m height). From the reference range for adult women, the average leptin level at a BMI of 45 kg/m<sup>2</sup> is approximately 224 ng/mI. The optimal dilution would be 1:30.

### APPENDIX

**Table 3:** Constants a, b, c and d for calculation of leptin reference ranges and leptin SDS based on BMI. Groups of normal healthy individuals were stratified according to gender and pubertal stage/age. TS= Tanner stage, n= number of subjects, a,b,c and d = constants as defined in the text.

Groups	n	а	b	с	d
Males					
TS 1&2	136	0,0146	0,2706	0,8821	0,5379
TS 3&4	50	0,0181	0,2067	1,1919	0,6850
TS 5	112	0,0316	0,1462	1,0821	0,6558
Adults	380	0,0130	0,2200	1,1053	0,6740
Females					
TS 1&2	136	0,0422	0,2499	0,7849	0,4786
TS 3&4	43	0,0543	0,2357	0,5745	0,3379
TS 5	157	0,2550	0,1508	0,7053	0,4301
Adults	587	0,3042	0,1467	0,8548	0,5212

	Pe	rcentile	(µg/L)		
BMI (kg/m²)	1	5	50	95	99
11	0.22	0.30	0.66	1.45	1.99
12	0.28	0.39	0.85	1.86	2.56
13	0.36	0.50	1.09	2.38	3.29
14	0.46	0.64	1.40	3.06	4.22
15	0.60	0.82	1.79	3.93	5.42
16	0.76	1.05	2.30	5.04	6.96
17	0.98	1.35	2.95	6.47	8.93
18	1.25	1.73	3.79	8.31	11.5
19	1.61	2.22	4.87	10.7	14.7
20	2.07	2.85	6.25	13.7	18.9
21	2.65	3.66	8.03	17.6	24.3
22	3.41	4.70	10.3	22.6	31.2
23	4.37	6.03	13.2	29.0	40.0
24	5.62	7.75	17.0	37.2	51.4
25	7.21	9.95	21.8	47.8	65.9
26	9.26	12.8	28.0	61.4	84.7
27	11.9	16.4	35.9	78.8	109.0
28	15.3	21.1	46.1	101.0	140.0
29	19.6	27.0	59.2	130.0	
30	15.2	34.7	76.1		
31	32.3	44.6	97.7		
32	41.5	57.2	125.0		
33	53.2	73.4			
34	68.4	94.3			
35	87.8	121.0			
36	113				
37	145				

Table 4: Girls Tanner stages 1 and 2



Figure 1: Reference ranges of human serum levels referring to BMI: GirlsTanner stage 1 & 2 (see text for details)

### Table 5: Boys Tanner stages 1 and 2

	Percer	ntile / Perze	entile (µg/L)	)	
BMI (kg/m²)	1	5	50	95	99
11	0.08	0.12	0.29	0.69	0.99
12	0.01	0.16	0.38	0.91	1.30
13	0.14	0.20	0.49	1.19	1.71
14	0.19	0.26	0.65	1.56	2.24
15	0.24	0.35	0.85	2.04	2.93
16	0.32	0.46	1.11	2.68	3.84
17	0.41	0.60	1.45	3.51	5.04
18	0.55	0.79	1.90	4.60	6.60
19	0.72	1.03	2.50	6.03	8.66
20	0.94	1.35	3.27	7.90	11.3
21	1.24	1.77	4.29	10.4	14.9
22	1.62	2.33	5.62	13.6	19.5
23	2.12	3.05	7.37	17.8	25.5
24	2.78	3.99	9.66	23.3	33.5
25	3.65	5.24	12.7	30.6	43.9
26	7.78	6.87	16.9	40.1	57.5
27	6.27	9.0	21.7	52.5	75.4
28	8.22	11.8	28.5	68.9	98.8
29	10.7	15.5	37.4	90.3	129.0
30	14.1	20.3	48.9	118.0	
31	18.5	26.6	64.2		
32	24.3	34.8	84.1		
33	31.8	45.6	110.0		
34	41.7	59.8	144.0		
35	54.6	78.4			
36	71.6	102.0			
37	93.9	134.0			
38	123.0				





Percentile (µg/L)							
BMI (kg/m²)	1	5	50	95	99		
11	0.32	0.41	0.73	1.29	1.63		
12	0.41	0.52	0.92	1.63	2.06		
13	0.52	0.66	1.16	2.07	2.61		
14	0.65	0.83	1.47	2.61	3.31		
15	0.83	1.05	1.87	3.31	4.19		
16	1.05	1.33	2.36	4.19	5.30		
17	1.33	1.68	2.99	5.30	6.71		
18	1.68	2.13	3.78	6.71	8.49		
19	2.13	2.69	4.79	8.5	10.8		
20	2.69	3.41	6.06	10.7	13.6		
21	3.41	4.31	7.67	13.61	17.2		
22	4.32	5.46	9.71	17.2	21.8		
23	5.46	6.91	12.3	21.8	27.6		
24	6.91	8.75	15.6	27.6	34.9		
25	8.75	11.1	19.7	34.9	44.2		
26	11.1	14.0	24.9	44.2	56.0		
27	14.0	17.7	31.6	56.0	70.9		
28	17.8	22.5	39.9	70.9	89.7		
29	22.5	28.4	50.5	89.7	114.0		
30	28.4	36.0	63.9	114.0	144.0		
31	36.0	45.6	80.9	144.0			
32	45.6	57.7	80.2	144.0			
33	57.7	73.0	102.0				
34	73.0	92.4	130.0				
35	92.4	117.0					
36	117.0	148.0					
37	148.0						

Table 6: Girls Tanner stages 3 and 4.



**Figure 3**: Reference ranges of human serum levels referring to BMI: Girls Tanner stage 3 & 4 (see text for details).

Table 7 : Boys Tanner stage 3 & 4

Percentile (µg/L)					
BMI (kg/m²)	1	5	50	95	99
11	0.03	0.05	0.18	0.58	0.94
12	0.04	0.07	0.22	0.71	1.16
13	0.49	0.08	0.27	0.88	1.43
14	0.06	0.10	0.33	1.08	1.75
15	0.07	0.12	0.40	1.32	2.16
16	0.09	0.15	0.49	1.63	2.65
17	0.11	0.18	0.61	2.00	3.26
18	0.14	0.23	0.75	2.46	4.01
19	0.17	0.28	0.92	3.03	4.93
20	0.21	0.34	1.13	3.72	6.06
21	0.26	0.42	1.39	4.58	7.46
22	0.32	0.52	1.71	5.63	9.17
23	0.39	0.64	2.10	6.92	11.3
24	0.48	0.78	2.58	8.51	13.9
25	0.59	0.96	3.18	10.5	17.0
26	0.73	1.19	3.91	12.9	21.0
27	0.89	1.46	4.80	15.8	25.8
28	1.10	1.79	5.90	19.4	31.7
29	1.35	2.20	7.26	23.9	39.0
30	1.66	2.71	8.93	29.4	48.0
31	2.05	3.33	11.0	36.2	58.9
32	2.51	4.09	13.5	44.5	72.4
33	3.09	5.04	16.6	54.7	89.1
34	3.80	6.20	20.4	67.2	109.0
35	4.68	7.62	25.1	82.6	134.0
36	5.75	9.37	30.9	101.0	
37	7.07	11.5	37.9	124.0	
38	8.7	14.2	46.7		
39	10.7	17.4	57.4		
40	13.1	21.4	70.5		



**Figure 4**: Reference ranges of human serum levels referring to BMI: Boys Tanner stage 3 & 4 (see text for details).

Table 8 : Girls Tanner stage 5.

Percentile (µg/L)					
BMI (kg/m²)	1	5	50	95	99
11	0.50	0.66	1.34	2.71	3.62
12	0.58	0.77	1.56	3.15	4.21
13	0.67	0.89	1.81	3.67	4.89
14	0.78	1.04	2.11	4.26	5.69
15	0.91	1.21	2.45	4.96	6.62
16	1.05	1.41	2.85	5.76	7.70
17	1.22	1.64	3.31	6.70	8.95
18	1.42	1.90	3.85	7.79	10.4
19	1.66	2.21	4.48	9.06	12.1
20	1.93	2.57	5.20	10.5	14.1
21	2.24	2.99	6.05	12.3	16.4
22	2.60	3.48	7.03	14.2	19.0
23	3.03	4.04	8.18	16.6	22.1
24	3.52	4.70	9.51	19.3	25.7
25	4.09	5.46	11.0	22.4	29.9
26	4.76	6.35	12.9	26.0	34.8
27	5.53	7.39	15.0	30.3	40.4
28	6.43	8.59	17.39	35.2	47.0
29	7.48	9.99	20.2	40.9	54.7
30	8.70	11.6	23.5	47.6	63.5
31	10.1	13.5	27.3	55.3	73.9
32	11.8	15.7	31.8	64.4	85.9
33	13.7	18.3	37.0	74.9	99.9
34	15.9	21.2	43.0	87.0	116.0
35	18.5	24.7	50.0	101.0	135.0
36	21.5	28.7	58.1	118.0	
37	25.0	33.4	67.6	137.0	
38	29.1	38.8	78.6		
39	33.8	45.1	91.4		
40	39.4	52.5	106.0		



Figure 5: Reference ranges of human serum levels referring to BMI: Girls Tanner stage 5 (see text for details)

Table 9 : Boys Tanner stage 5

	Percentile (µg/L)				
BMI (kg/m²)	1	5	50	95	99
11	0.03	0.05	0.16	0.47	0.73
12	0.04	0.06	0.18	0.54	0.84
13	0.05	0.07	0.21	0.62	0.97
14	0.05	0.08	0.24	0.72	1.12
15	0.06	0.10	0.28	0.84	1.30
16	0.07	0.11	0.33	0.97	1.51
17	0.08	0.13	0.38	1.12	1.74
18	0.1	0.15	0.44	1.3	2.02
19	0.11	0.17	0.51	1.50	2.34
20	0.13	0.2	0.59	1.74	2.7
21	0.15	0.23	0.68	2.01	3.13
22	0.17	0.27	0.79	2.33	3.62
23	0.20	0.31	0.91	2.69	4.19
24	0.23	0.36	1.05	3.12	4.85
25	0.27	0.41	1.22	3.61	5.62
26	0.31	0.48	1.41	4.17	6.5
27	0.36	0.55	1.63	4.83	7.52
28	0.41	0.64	1.89	5.59	8.71
29	0.48	0.74	2.19	6.47	10.1
30	0.55	0.86	2.54	7.49	11.7
31	0.64	1.00	2.94	8.67	13.5
32	0.74	1.15	3.4	10.0	15.6
33	0.86	1.33	3.94	11.6	18.1
34	0.99	1.54	4.55	13.4	20.9
35	1.15	1.79	5.27	15.6	24.2
36	1.33	2.07	6.10	18.0	28.1
37	1.54	2.39	7.06	20.8	32.5
38	1.78	2.77	8.17	24.1	37.6
39	2.06	3.21	9.46	27.9	43.5
40	2.38	3.71	10.9	32.3	50.3



**Figure 6:** Reference ranges of human serum levels referring to BMI: Boys Tanner stage 5 (see text for details).

Percentile (µg/L)					
BMI (kg/m²)	1	5	50	95	99
11	0.46	0.65	1.53	3.59	5.10
12	0.53	0.75	1.77	4.16	5.90
13	0.61	0.87	2.05	4.82	6.83
14	0.71	1.01	2.37	5.58	7.91
15	7.82	1.17	2.75	6.46	9.17
16	0.95	1.35	3.18	7.48	10.61
17	1.10	1.57	3.68	8.66	12.3
18	1.28	1.81	4.27	10.0	14.2
19	1.48	2.10	4.94	11.6	16.5
20	1.71	2.43	5.72	13.4	19.1
21	1.99	2.82	6.62	15.6	22.1
22	2.30	3.26	7.67	18.0	25.6
23	2.66	3.78	8.88	20.9	29.3
24	3.08	4.38	10.3	24.2	34.3
25	3.57	5.07	11.9	28.0	39.7
26	4.13	5.87	13.8	32.4	46.0
27	4.79	6.79	16.0	37.5	53.3
28	5.54	7.87	18.5	43.5	61.7
29	6.42	9.11	21.4	50.4	71.5
30	7.43	10.6	24.8	58.3	82.8
31	8.61	12.2	28.7	67.5	95.8
32	9.97	14.1	33.3	78.2	111.0
33	11.5	16.4	38.5	90.5	129.0
34	13.4	19.0	44.6	105.0	149.0
35	15.5	22.0	51.6	121.0	
36	17.9	25.4	59.8	141.0	
37	20.8	29.5	69.3		
38	24.0	34.1	80.2		
39	27.8	39.5	92.9		
40	32.2	45.7	108.0		

Table 10: Adult women

# Adult Women (20 - 80 yrs)

**Figure 7**: Reference ranges of human serum levels referring to BMI: Adult women (see text for details)

### Table 11 Adult men

Percentile (µg/L)					
BMI (kg/m²)	1	5	50	95	99
11	0.03	0.05	0.15	0.44	0.69
12	0.04	0.06	0.18	0.55	0.87
13	0.05	0.08	0.23	0.69	1.08
14	0.06	0.09	0.28	0.85	1.34
15	0.07	0.12	0.35	1.06	1.67
16	0.09	0.15	0.44	1.33	2.09
17	0.12	0.18	0.55	1.65	2.60
18	0.14	0.23	0.68	2.06	3.24
19	0.18	0.28	0.85	2.57	4.04
20	0.22	0.35	1.06	3.20	5.03
21	0.23	0.44	1.32	3.98	6.27
22	0.35	0.54	1.64	4.97	7.81
23	0.43	0.78	2.05	6.19	9.73
24	0.54	0.85	2.55	7.71	12.1
25	0.67	1.05	3.18	9.61	15.1
26	0.83	1.31	3.96	12.0	18.8
27	1.04	1.64	4.94	14.9	23.5
28	1.30	2.04	6.15	18.6	29.2
29	1.61	2.54	7.67	23.2	36.4
30	2.01	3.16	9.56	28.9	45.4
31	2.51	3.94	11.9	36.0	56.6
32	3.12	4.91	14.8	44.9	70.5
33	3.89	6.12	18.5	55.8	87.8
34	4.85	7.63	23.0	69.6	109.0
35	6.04	9.51	28.7	86.7	136.0
36	7.53	11.8	35.8	108.0	
37	9.38	14.8	44.6	135.0	
38	11.7	18.4	55.5		
39	14.6	22.9	69.2		
40	18.2	28.6	86.2		



**Figure 8**: Reference ranges of human serum levels referring to BMI: Adult men (see text for details).



Figure 9: Serial dilution of leptin calibrator (filled circle) and different sera with varying concentrations of leptin. The leptin level of the undiluted calibrator was 16 ng/ml. Means of duplicate measurements are shown

### REFERENCES

- 1. Zhang Y, Proenca R, Maffei M, Barone M, Leopold L, Friedman JM. 1994 Positional cloning of the mouse obese gene and its human homologue. Nature. 372:425-432.
- 2. Halaas JL, Gajiwala KS, Maffei M, et al. 1995 Weight-reducing effects of the plasma protein encoded by the obese gene. Science. 269:543-546.
- MacDougald OA, Hwang CS, Fan H, Lane MD. 1995 Regulated expression of the obese gene product (leptin) in white adipose tissue and 3T3-L1 adipocytes. Proc Natl Acad Sci U S A. 92:9034-9037.
- 4. Rentsch J, Chiesi M. 1996 Regulation of ob gene mRNA levels in cultured adipocytes. FEBS Lett. 379:55-59.
- 5. Wabitsch M, Jensen PB, Blum WF, et al. 1996 Insulin and cortisol promote leptin production in cultured human fat cells. Diabetes. 45:1435-1438.
- 6. Stephens TW, Basinski M, Bristow PK, et al. 1995 The role of neuropeptide Y in the antiobesity action of the obese gene product. Nature. 377:530-532.
- 7. Campfield LA, Smith FJ, Guisez Y, Devos R, Burn P. 1995 Recombinant mouse OB protein: evidence for a peripheral signal linking adiposity and central neural networks. Science. 269:546-549.
- 8. Pelleymounter MA, Cullen MJ, Baker MB, et al. 1995 Effects of the obese gene product on body weight regulation in ob/ob mice. Science. 269:540-543.
- 9. Levin N, Nelson C, Gurney A, Vandlen R, de-Sauvage F. 1996 Decreased food intake does not completely account for adiposity reduction after ob protein infusion. Proc Natl Acad Sci U S A. 93:1726-1730.

10. Tartaglia LA, Dembski M, Weng X, et al. 1995 Identification and expression cloning of a leptin receptor, OB-R. Cell. 83:1263-1271.

- 11. Chen H, Charlat O, Tartaglia LA, et al. 1996 Evidence that the diabetes gene encodes the leptin receptor: identification of a mutation in the leptin receptor gene in db/db mice. Cell. 84:491-495.
- 12. Lee GH, Proenca R, Montez JM, et al. 1996 Abnormal splicing of the leptin receptor in diabetic mice. Nature. 379:632-635.
- 13. Lynn RB, Cao GY, Considine RV, Hyde TM, Caro JF. 1996 Autoradiographic localization of leptin binding in the choroid plexus of ob/ob and db/db mice. Biochem Biophys Res Commun. 219:884-889.
- 14. Considine RV, Considine EL, Williams CJ, Hyde TM, Caro JF. 1996 The hypothalamic leptin receptor in humans: identification of incidental sequence polymorphisms and absence of the db/db mouse and fa/fa rat mutations. Diabetes. 45:992-994.
- 15. Sinha MK, Opentanova I, Ohannesian JP, et al. 1996 Evidence of free and bound leptin in human circulation. J Clin Invest. 98:1277-1282.
- 16. Ahima RS, Prabakaran D, Mantzoros C, et al. 1996 Role of leptin in the neuroendocrine response to fasting. Nature. 382:250-252.
- 17. Chehab FF, Lim ME, Lu R. 1996 Correction of the sterility defect in homozygous obese female mice by treatment with the human recombinant leptin. Nat Genet. 12:318-320.
- Schwartz MW, Baskin DG, Bukowski TR, et al. 1996 Specificity of leptin action on elevated blood glucose levels and hypothalamic neuropeptide Y gene expression in ob/ob mice. Diabetes. 45:531-535.
- 19.Schwartz MW, Seeley RJ, Campfield LA, Burn P, Baskin DG. 1996 Identification of targets of leptin action in rat hypothalamus. J Clin Invest. 98:1101-1106.
- 20. Campfield LA, Smith FJ, Burn P. 1996 The OB protein (leptin) pathway a link between adipose tissue mass and central neural networks. Horm Metab Res. 28:619-632.
- 21. Rohner-Jeanrenaud F, Cusin I, Sainsbury A, Zakrzewska KE, Jeanrenaud B. 1996 The loop system between neuropeptide Y and leptin in normal and obese rodents. Horm Metab Res. 28:642-648.
- 22. Chan YY, Steiner RA, Clifton DK. 1996 Regulation of hypothalamic neuropeptide-Y neurons by growth hormone in the rat. Endocrinol. 137:1319-1325.
- 23. Pierroz DD, Catzeflis C, Aebi AC, Rivier JE, Aubert ML. 1996 Chronic administration of neuropeptide Y into the lateral ventricle inhibits both the pituitary-testicular axis and growth hormone and insulin-like growth factor I secretion in intact adult male rats. Endocrinol. 137:3-12.
- 24. Frederich RC, Hamann A, Anderson S, Lollmann B, Lowell BB, Flier JS. 1995 Leptin levels reflect body lipid content in mice: evidence for diet-induced resistance to leptin action. Nat Med. 1:1311-1314.
- 25. Maffei M, Halaas J, Ravussin E, et al. 1995 Leptin levels in human and rodent: measurement of plasma leptin and ob RNA in obese and weight-reduced subjects. Nat Med. 1:1155-1161.
- 26. Considine RV, Sinha MK, Heiman ML, et al. 1996 Serum immunoreactive-leptin concentrations in normal-weight and obese humans. N Engl J Med. 334:292-295.
- 27. Kolaczynski JW, Considine RV, Ohannesian J, et al. 1996 Responses of leptin to short-term fasting and refeeding in humans: A link with keto-genesis but not ketones themselves. Diabetes. 45:1511-1515.
- 28. Harris RB, Ramsay TG, Smith SR, Bruch RC. 1996 Early and late stimulation of ob mRNA expression in meal-fed and overfed rats. J Clin Invest. 97:2020-2026.
- 29. Kolaczynski JW, Ohannesian J, Considine RV, Marco C, Caro JF. 1996 Response of leptin to short-term and prolonged overfeeding in humans. J Clin Endocrinol Metab. 91:4162-4165.
- 30. Saladin R, De-Vos P, Guerre-Millo M, et al. 1995 Transient increase in obese gene expression after food intake or insulin administration. Nature. 377:527-529.
- 31. Cusin I, Sainsbury A, Doyle P, Rohner-Jeanrenaud F, Jeanrenaud B. 1995 The ob gene and insulin. A relationship leading to clues to the understanding of obesity. Diabetes. 44:1467-1470.
- 32. Kolaczynski JW, Nyce MR, Considine RV, et al. 1996 Acute and chronic effects of insulin on leptin production in humans: studies in vivo and in vitro. Diabetes. 45:699-701.
- 33. Malström R, Taskinen M-R, Karonen S-L, Yki-Järvinen H. 1996 Insulin increases plasma leptin concentrations in normal subjects and patients with NIDDM. Diabetologia. 39:993-996.
- 34. De-Vos P, Saladin R, Auwerx J, Staels B. 1995 Induction of ob gene expression by corticosteroids is accompanied by body weight loss and reduced food intake. J Biol Chem. 270:15958-15961.
- 35. Slieker LJ, Sloop KW, Surface PL, et al. 1996 Regulation of expression of ob mRNA and protein by glucocorticoids and cAMP. J Biol Chem. 271:5301-5304.
- 36. Miell JP, Englaro P, Blum WF. 1996 Dexamethasone induces an acute and sustained rise in circulating leptin levels in normal human subjects. Horm Metab Res. 28:704-707.
- 37. Sinha MK, Ohannesian JP, Heiman ML, et al. 1996 Nocturnal rise of leptin in lean, obese, and non-insulin- dependent diabetes mellitus subjects. J Clin Invest. 97:1344-1347.
- Havel PJ, Kasim-Karakas S, Dubuc GR, Mueller W, Phinney SD. 1996 Gender difference in plasma leptin concentrations. Nature Med. 2:949-950.

- 39. Rosenbaum M, Nicolson M, Hirsch J, et al. 1996 Effects of gender, body composition, and menopause on plasma concentrations of leptin. J Clin Endocrinol Metab. 81:3424-3427.
- 40. Hassink SG, Sheslow DV, de Lancey E, Opentanova I, Considine RV, Caro JF 1996 Serum leptin in children with obesity: relationship to gender and development. Pediatrics. 98:201-203.
- 41. Scholz GH, Englaro P, Thiele I, et al. 1996 Dissociation of serum leptin concentration and body fat content during long term dietary intervention in obese individuals. Horm Metab Res. 28:718-723.
- 42. Caro JF, Kolaczynski JW, Nyce MR, et al. 1996 Decreased cerebrospinal-fluid/serum leptin ratio in obesity: a possible mechanism for leptin resistance. Lancet. 348:159-161.
- 43. Robinson CJ, Gaines-Das R,Woollacott D, et al. 2001 The first international standard for human leptin and the first international standard for mouse leptin: comparison of candidate preparations by in vitro bioassays and immunoassays. J Molecular Endocrinol. 27: 69-76.
- 44. Adresse NIBSC: Blanche Lane, South Mimms, Potters Bar, Hertfordshire EN6 3QG, Great Britain.
- 45. Blum WF, Juul A; Reference ranges of serum leptin,s.318-326. In:Leptin -the voice of adipose tissue, Blum WF et al, eds., Johann Ambrosius Verlag, Heidelberg, 1997.
- 46. Holtkamp K, Blum WF, et al., 2006, Elevated Physical Activity and Low Leptin Levels Co-Occur in Patients with Anorexia Nervosa, The Journal of Clinical Endocrinology Metabolism 88(11) 5169-5174
- 47. Guebra-Egziabher, Bernhard J., et al. 2005, Adiponectin in chronic kidney disease is related more to metabolic distrubances than to decline in renal function, Nephrol Dial Transplant 20: 129-134

### SUMMARY – DIAsource Hu LEPTIN RIA

Reagent Preparation:	Reconstitution:	
Capture Antibody	in 7 ml Assay Buffer	
Specific Antibody	in 7 ml Assay Buffer	
Tracer	in 13 ml Assay Buffer	
Control 1	in 750 µl Assay Buffer	
Dia ama (Camuna Camuniana ang karupadi un dikuta d		

**Proposed Assay Procedure for Double Detemination** 

Addition of Reagent [µl]							
Nr. of Tubes	Contents of Tubes	Assay Buffer Calibrators Control 1 Samples	Capture Antibody	Specific Antibody	Tracer		
1,2	TC	_	-	-	100		
3,4	B <sub>0</sub>	Ass Buf: 25	50	50	100		
5-18	Calibrators	CAL 1-7: 25	50	50	100		
19,20	Control 1	25	50	50	100		
21,22	Sample 1	25	50	50	100		
23,24 (etc)	Sample 2 (etc)	25	50	50	100		

 Tubes Nr.:1,2 remove until counting the activity

 Incubate, over night, at least. 15 hours at RT, 350 rpm

 Aspirate the liquid completely. Take care that the coating of the tubes remains intact.

 Add 500 µl of Assay Buffer to the tubes.

 Aspirate the liquid completely (see above).

 Count the radioactivity of all tubes.

Revision Date : 2011-02-04

	<u>Used symbols</u>		
Ţ.	Consult instructions for use		
1	Storage temperature		
Ω	Use by		
LOT	Batch code		
REF	Catalogue number		
CONTROL	Control		
IVD	In vitro diagnostic medical device		
	Manufacturer		
<u> </u>	Contains sufficient for <n> tests</n>		
WASH SOLN CONC	Wash solution concentrated		
CAL 0	Zero calibrator		
CAL N	Calibrator #		
CONTROL N	Control #		
Ασ 1251	Tracer		
Ab 1251	Tracer		
	Tracer concentrated		
Ab 1251 CONC			
Π	Tubes		
	Acetonitrile		
SERUM	Serum		
DIL SPE	Specimen diluent		
DIL BUF	Dilution buffer		
ANTISERUM	Antiserum		
IMMUNOADSORBENT	Immunoadsorbent		
DIL CAL	Calibrator diluent		
REC SOLN	Reconstitution solution		
PEG	Polyethylene glycol		
EXTR SOLN	Extraction solution		
ELU SOLN	Elution solution		
GEL	Bond Elut Silica cartridges		
PRE SOLN	Pre-treatment solution		
NEUTR SOLN	Neutralization solution		
TRACEUR BUF	Tracer buffer		
w	Microtiterplate		
Ab HRP	HRP Conjugate		
Ag HRP	HRP Conjugate		
Ab HRP CONC	HRP Conjugate concentrate		
Ag HRP CONC	HRP Conjugate concentrate		
CONJ BUF	Conjugate buffer		
CHROM TMB CONC	Chromogenic TMB concentrate		
CHROM TMB	Chromogenic TMB solution		
SUB BUF	Substrate buffer		
STOP SOLN	Stop solution		
INC SER	Incubation serum		
BUF	Buffer		
Ab AP	AP Conjugate		
SUB PNPP	Substrate PNPP		
BIOT CONJ CONC	Biotin conjugate concentrate		
AVID HRP CONC	Avidine HRP concentrate		
ASS BUF	Assay buffer		
Ab BIOT	Biotin conjugate		
Ab	Specific Antibody		
SAV HRP CONC	Streptavidin HRP concentrate		
NSB	Non-specific binding		
2nd Ab	2nd Antibody		
ACID BUF	Acidification Buffer		
DIST	Distributor		