Fecal Rotavirus Antigen Test Kit

RAPEPKT917

Fecal Rotavirus Antigen Test Kit

Rapid Immunochromatographic Test Device for the detection of Rotavirus Antigen in Feces – Instructions for Fecal Sample Collection


RAPEPKT917

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96


1

1. Take the sampling tube and unscrew the sampling lid, keeping the sampling tube in a **vertical position** to prevent loss of solution.


2

1. Hold the sampling lid by the **Thumb Grip**.
2. Use the tip of the sampling lid to collect a small amount of fecal sample at two or more sites. Only take the fecal sample that sticks to the sampling lid tip (never intentionally place any separate piece of fecal sample into the tube).

3

1. Insert and screw the sampling lid back into the sampling tube in a **vertical position**. Do not spill any solution from the tube.
2. Tightly seal the lid with the tube.

Fecal Rotavirus Antigen Test Kit

en

Rapid Immunochromatographic Test Device for the detection of Rotavirus Antigen in Feces – Instructions for Test Procedures


RAPEPKT917

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96


A

1. Shake to dissolve the stool into solution.
2. Turn the sampling tube **upside down vertically**.
3. Remove the test strip from foil pouch.

B

1. Insert and screw the test strip **in a vertical position** into the sampling tube by breaking the bottom seal of the sampling tube.
2. Allow the solution to flow into the bottom space of test strip, keeping the device **in a vertical position**.
3. You may soon see a red fluid moving across the white area of the test strip. Read test result after 5 minutes.

Revision date : 2010-03-12

Fecal Rotavirus Antigen Test Kit

en

Rapid Immunochromatographic Test Device for the detection of Rotavirus Antigen in Feces

RAPEPKT917

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue de l'Industrie 8, B-1400 Nivelles, Belgium - Tel: +32 67 88 99 99 - Fax : +32 67 88 99 96

INTENDED USE

This rotavirus antigen test kit is intended for the direct qualitative detection of the presence of rotavirus antigen in patient fecal samples. The test might be used as an aid for the diagnosis of a rotavirus infection resulting in acute gastroenteritis.

SUMMARY OF PHYSIOLOGY

Rotaviruses are the main cause of acute gastroenteritis and diarrhea, especially in children under the age of two years. Rotaviruses have been identified in almost 40% of the feces of children with gastroenteritis. Rotavirus is the cause of up to 50% of the hospitalized cases of diarrhea in infant and young children. If not treated, the infection may result in severe dehydration and disorders of body electrolyte balance. Therefore, it can be mortal in risk populations such as children, the elderly or immunosuppressed individuals. Rotavirus is transmitted by oral-fecal contact with an incubation period of 1-3 days. Characteristic symptoms include vomiting, hydrodiarrhoea for 3 to 8 days, high temperature and stomach pains. A large amount of rotavirus particles is shed during infection. Diagnosis of gastroenteritis with rotavirus infection can be established based on the detection of the virus particles by electron microscopy or the virus antigen by specific immunoassay methods.

ASSAY PRINCIPLE

The Rotavirus Rapid Test Strip employs dye-conjugated monoclonal antibodies against antigen VP6 of group A rotavirus, and solid-phase specific rotavirus antibodies. In this test the specimen is first treated with an extraction solution to extract rotavirus antigens from the feces. Following extraction, the only step required is to screw the rotavirus test strip into the sample collection device. As the sample extraction flows through the chamber and reach the test strip, the colored particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the colored particles. Different colored lines will be visible, depending on the virus content of the sample. These lines, after 5 minutes of incubation at room temperature, are used to interpret the result.

REAGENTS: Preparation and Storage

This test kit must be stored at 2 – 8°C upon receipt. For the expiration date of the kit refer to the label on the kit box. All components are stable until this expiration date.

Prior to use allow all reagents to come to room temperature if the kit is stored at refrigerated condition.

- | | | |
|---|--------------|---|
| § | TUBE | 30 Fecal specimen collection device: containing sampling tube, sampling lid and pre-added extraction solution (1.1 ml Tris Buffer containing 0.1% bovine serum albumin) in the sampling tube. This device should be stored at 2 to 8°C. Do not freeze. |
| § | STRIP | 30 Test strip: one dipstick for the Fecal Rotavirus Antigen test is assembled in a transparent housing and sealed in a foil pouch with desiccant. It should remain in its original sealed pouch until ready for use. The test strip should be stored at 2 to 30°C. Do not freeze. |
| § | | 30 Disposable spoons for sample collection. |
| § | | Instruction for use. |

SAFETY PRECAUTIONS

The reagents are for professional use only. Source material from which reagents of bovine serum was derived in the contiguous 48 United States. It was obtained only from donor health animals maintained under veterinary supervision and found free of contagious diseases. Wear gloves while performing this test and handle these reagents and patient samples as if they are infectious. Do not get in eyes, on skin, or on clothing. Do not ingest. On contact, flush with copious amounts of water for at least 15 minutes. When the assay procedure is completed, dispose of specimens (biohazard materials) carefully after autoclaving for at least one hour. Alternatively, treat with a 0.5 or 1% solution of sodium hypochlorite for one hour before disposal. Use Good Laboratory Practices.

MATERIALS REQUIRED BUT NOT PROVIDED

Positive Control

SPECIMEN COLLECTION

1. Stool specimens can be collected at any time of the day.
2. Collect a random sample of feces in a clean, dry cup or toilet paper.
3. Unscrew the sampling lid and keep the sampling tube in a vertical position to prevent the loss of any extraction solution (figure 1).
4. Insert and twist the tip of the sampling lid into the stool specimen at two or more different sites (figure 2).
5. Collect fecal sample that is stuck to the surface of the sampling lid. Do not intentionally collect any separate and large pieces of fecal sample into the tube.
6. Replace the sampling lid into the tube and secure tightly (figure 3).
7. The specimen is ready for testing, transportation or storage. It can be stored at 2-8°C for up to 14 days and at room temperature for up to 5 days.

Note: Two specimens from three consecutive bowel movements are recommended from American Cancer Society. Specimen should not be collected during digital rectal examination.

TEST PROCEDURE

1. Bring the sealed foil pouch test strips and collected specimens to room temperature.
2. Shake the sampling tube vigorously to ensure a good liquid suspension.
3. Position the sampling tube upside down vertically and let it settle for about 1 minute.
4. Remove the test strip from the sealed foil pouch (figure A).
5. Screw the test strip **in a vertical position** into the sampling tube by **breaking** the bottom seal of the sampling tube. Secure tightly! (Figure B)
6. Allow the solution to flow into the bottom space of the test strip and keeping the device **in a vertical position**.
7. Read test result at 5 minutes. Do not interpret test result after 10 minutes.

INTERPRETATION OF RESULTS

- **Positive:**
If two red/pink colored bands are visible within 10 minutes, the test result is positive and valid.
- **Negative:**
If the test area has no colored band and the control area displays a red/pink colored band, the test result is negative.
- **Invalid:**
If a colored band does not appear in the control area regardless whether there is any band in the test area, the test result is invalid.

Note: The control line is located at the upper half of the test membrane area and the test line is located at the lower half of it.

QUALITY CONTROL

Good laboratory practices recommend the use of appropriate controls. There are two types of controls for the Rotavirus antigen test, the internal procedural control and external controls.

1. **Internal procedural control:** Each fecal rotavirus antigen test contains a build in procedural control. A colored control line will appear if the test has been performed correctly and the reagents are reactive. This control does not ensure that the test line antibody is accurately detecting the presence or absence of viral antigen in the test fecal sample.
2. **External controls:** It is recommended to use external positive controls. The external positive controls are not provided with this kit. . External controls are used to assure that the test line antibody is reactive. However, external controls will not detect an error in the performance of a particular test with a patient sample It is recommended that the external control be tested once per kit.

You should always follow local, state, federal guidelines for running quality control.

LIMITATION OF THE PROCEDURE

1. The test should be used only for the detection of rotavirus antigen in fecal samples.
2. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result
3. Two hundred samples were evaluated to assure the correct performance of the test. The correlation of the results with other techniques (ELISA) was satisfactory. However, interferences in the performance of the tests should not be excluded.
4. No cross-reactions with other viruses or substances were observed during the evaluation of the test. A negative result does not totally exclude a possible rotavirus infection. The significance of the results must be evaluated in relation to the patient's clinical symptoms.
5. As with all diagnostic tests, the definitive clinical diagnosis must not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. Fecal rotavirus antigen test is designed for the aid of clinical diagnosis and should not replace other diagnostic procedures.

PERFORMANCE

The sensitivity and specificity of this rotavirus antigen test device are studied with 206 clinical samples and compared with a rotavirus antigen ELISA test. The results are shown below.

Strip \ ELISA	Positive	Negative	Total
Positive	68	2	70
Negative	2	134	136
Total	70	136	206

Specificity: 98.5 % (134/136)

Sensitivity: 97.1 % (68/70)

Accuracy: 98.1% (202/206)

Inter-series and intra-series accuracy: 100 %

Interference: Cross reactivity has been evaluated and found to be negative compared to positive specimens of *Cryptosporidium parvum*, to the Adenovirus group and to the 40/41 strain of Adenovirus.

LITERATURES

1. Set-up of a new rapid immunochromatographic diagnostic test for a Rotavirus detection. *D. Van Beers , M. DE Foor , R. Viehoff , D. Col , M. Venuti and T. Leclipteux.* Progress in Clinical Virology III , Bologne , Septembre 1997.
2. Detection of rotavirus in faecal specimens with a monoclonal antibody enzyme-linked immunosorbent assay : comparison with polyclonal antibody enzyme-immunoassays and a latex agglutination test. *Sneyers et al.* Comp. Immun. Microbiol. Infect. Dis. , vol 12 , n°4 , pp 95-104 , 1989
3. Comparison of Three Rapid Immunoassays for the Detection of Rotavirus Antigen in Stool Samples *I. Van der Donck et al.* ESCV Winter Meeting 1999, Rotterdam, the Netherlands
4. Evaluacion de tres Metodos de Deteccion de Rotavirus en Heces *I. Wilhelmi et al.* 6th Congreso Nacional de Virologia, Madrid, 26th Oct. 99

 Consult instructions for use	 Manufacturer
 Storage temperature	 Contains sufficient for n tests
 Use by	 In vitro diagnostic medical device
 Batch code	 Test Strip
 Catalogue number	 Collection device

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