ZenTech

RIAZEN Renin Plasma Activity

Deutsch (de)

English (en)

Español (es)

Français (fr)

Italiano (it)

Polski (pl)

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R-EX-125 125 Determinations

1. INTENDED USE

Radioimmunoassay for the measurement of Plasmatic Renin Activity in human Plasma.

2. EXPLANATION OF THE TEST

Renin-angiotensin-aldosterone system is a hormonal "cascade" system, whose first compound is represented by angiotensinogen. Renin, a proteolytic enzyme, transforms angiotensinogen into angiotensin I a polipeptide which has not biological activity. Angiotensin I is the substrate on which converting enzyme acts, producing biologically active angiotensin II. The main biological actions of angiotensin II consist of vasoconstriction, regulation of synthesis and release of aldosterone which modulates electrolytes balance. Thus the system is involved in regulation of blood pressure with several mechanisms. It is possible to evaluate system activity by measurement of plasmatic renin activity, because either angiotensin I or renin assay have methodological problems (low values for the first one, presence of structurally similar enzymes for the second one). Renin is synthesized by iuxta-glomerular apparatus of kidney and it has a molecular weight about 40.000 Dalton. The measurement of its activity is a valid tool in order to determine prognosis and therapy for patients affected by blood hypertension. Nevertheless it is important to know several factors which may alter the values such as: sodium intake, upright or supine position and drugs. Plasmatic renin activity increases with the reduction of sodium intake and it decreases changing the upright position into supine position. Moreover some drugs such as diuretics or oral contraceptives can produce an increase while other ones such as L-Dopa, Reserpine can lower renin activity.

Hyperreninism can be due to several situations:

- Hyperreninism with blood hypertension: reno-vascular hypertension, Wilms' tumor, Robertson-Kihara syndrome.
- Hyperreninism without blood hypertension: sodium depletion, nephrotic syndromes, Bartter's syndrome.

Hyporeninism can be due to several situations:

- Hyporeninism with blood hypertension: Conn's syndrome, primary pseudo-hyperaldosteronism, Biglieri's syndrome, some cases of essential hypertension.
- Hyporeninism without blood hypertension: hypoaldosteronism, hydro-saline retention, reduced activity of sympathetic nervous system.

In some situations it can be useful to evaluate changes of plasmatic renin activity after dynamic tests. One of the more useful test is furosemide stimulation.

3. PRINCIPLE OF THE METHOD

The kit has been designed for determination of plasmatic renin activity by radioimmunoassay of angiotensin I. Before radioimmunoassay, plasma must be incubated so that renin generates angiotensin I in standardized conditions. Enzymatic inhibitor avoids the enzymatic degradation of angiotensin I during generation step.

Procedure foresees that the user can choose pH generation.

The kit gives the possibility to performe the plasma generation at two different pH:

- Generation Buffer pH 6.0, utilizing the vial supplied with set reagents
- Generation Buffer pH 7.4, utilizing the vial provided only by request.

pH 7.4 generation is the same process involved in human physiology, blocking the activity of the other enzymes of the plasma.

pH 6.0 generation produces twice the amount as much pH 7.4, so that pH 6.0 method is more reliable for low renin activity samples.

Angiotensin generation is a linear function in the time, so the same results are reached changing generation time. Nevertheless the procedure foresees 90 minutes incubation at pH 6.0 and 3 hours incubation at pH 7.4 in order to standardize the method and to guarantee a good precision. After generation radioimmunoassay is performed. Antigen, whether from standards or from samples competes

Verdünnungs-faktor	gemessene Angiotensin I-Konz. (ng/ml)	erwartete Angiotensin I-Konz. (ng/ml)	Wiederfindung (%)
-	19,74	-	-
2	9,77	9,87	99,0
4	5,11	4,94	103,5
8	2,48	2,47	100,5
16	1,15	1,23	93,2
32	0,60	0,62	97,3
-	20,18	-	-
2	9,66	10,09	95,7
4	5,05	5,04	100,1
8	2,51	2,52	99,5
16	1,14	1,26	90,1
32	0,62	0,63	97,7

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with radioactive tracer (¹²⁵I-Angiotensin) for the binding sites of antibody. After incubation the amount of tracer bound to antibody will be inversely proportional to the amount of antigen present in standards or in samples. B/F separation is performed by addition of double antibody coupled to magnetic particles. The application of magnetic field allows the sedimentation of immunocomplexes, avoiding the centrifugation step.

The main steps of RIA assay are:

- incubation of reaction mixture for 18-20 hours at 4°C;
- B/F separation by adding double antibody coupled to magnetic particles;
- Decantation of supernatant;
- Counting of radioactivity in Gamma-counter.

4. REAGENTS

Do not interchange reagent from different kit lots. Lot numbers of reagents should be as stated on the "Certificate of Analysis". Do not use kit components beyond their expiration date.

The Renin kit contains sufficient reagents for 125 tubes. On receipt, store the kit at 2-8°C, until the expiration date on the label.

1. **Renin Tracer (Red), 1 vial:** Containing ¹²⁵I labelled Angiotensin I in phosphate buffer with bovine serum albumin and sodium azide 0.05% . 13.0 ml per vial. Maximum Radioactivity: 85 kBq. Expiry date is stated on vial label. The reagent is ready to use. Store at 2-8°C.
2. **Renin Calibrators, 8 vials:** Containing: 0- 0.2- 0.5- 1.0- 2.0- 5.0- 10.0- 25.0 ng/ml of lyophilized Angiotensin I in phosphate buffer with bovine serum albumin and sodium azide 0.2% w/v. Expiry date is stated on vial label. Reconstitute the contents of "0" standard with 2 ml distilled water and the other standards with 1 ml distilled water. Mix gently until the complete solubilization of freeze-dried residue. Store at 2-8°C for 1 day or at -20°C for longer length of time.
3. **Renin Antiserum, (color code green), 1 vial:** Containing lyophilized rabbit anti-Angiotensin I in phosphate buffer with rabbit gammaglobulin, bovine serum albumin, and sodium azide 0.05% w/v. Expiry date is stated on vial label. Reconstitute the contents of the vial with 12.5 ml distilled water and mix gently until the complete solubilization of freeze-dried residue. Store at 2-8°C for 3 days or at -20°C for longer length of time.
4. **Renin Separation Reagent, 1 vial:** Containing goat anti-rabbit gammaglobulin coupled to magnetic particles in Tris buffer with bovine serum albumin, and sodium azide <0.1 % w/v. 125 ml. Expiry date is stated on vial label. The solution is ready for use. Store at 2-8°C. Do not freeze.
5. **Renin Generation, Buffer pH 6.0, 1 vial:** Containing 6.0 ml phosphate buffer with yellow dye. Expiry date is stated on vial label. The suspension is ready for use. Store at 2-8°C. Do not freeze.
6. **Renin Generation Buffer, pH 7.4, 1 vial (Provided by request):** Containing 3.0 ml phosphate buffer with dyes (orange color). Expiry date is stated on vial label. The solution is ready for use. Store at 2-8°C. Do not freeze.
7. **Renin Enzyme Inhibitor, 1 vial:** Containing Phenilmethylsulfonyl Fluoride (PMSF) in ethanol. 0.5 ml. Store at 2-8°C. If the solution crystallizes after 4°C storage, let it stand at 37°C in stoppered vial. Warning: Avoid to keep the open vial besides the use. The solution is ready for use.

5. WARNING AND PRECAUTIONS FOR USERS

For In Vitro diagnostic use.

Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.

Radioactive Material - Not for Internal or External Use in Humans or Animals.

This radioactive material may be received, acquired, possessed, and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation there from to human

beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations of each country.

Physical characteristics of ¹²⁵I : see end of instructions

1. **Safety Precautions:** The following precautions should be observed in handling radioactive material:
 - o Store radioactive materials in a designated area.
 - o Do not eat, drink, smoke or apply cosmetics where radioactive materials are being handled.
 - o Do not pipette by mouth.
 - o Wear gloves when handling radioactive materials and wash hands thoroughly afterwards.
 - o Cover working area with disposable absorbent paper.
 - o Wipe up all spills immediately and thoroughly and dispose of the contaminated materials as radioactive waste.
 - o Dispose of the liquid radioactive waste into the sanitary sewage system if permitted by the local regulations.
2. **Chemical Hazard - Sodium Azide (NaN₃) Warning:** Some of the reagents in this kit contain sodium azide as a preservative. For all such reagents, the concentration of sodium azide is <0.2%. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Dispose of all non radioactive reagents by flushing with large amounts of water through the plumbing system.
Risk Phrases
R 21/22 Harmful in contact with skin and if swallowed.
Safety Phrases
S-26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S-28.1 After contact with skin, wash immediately with plenty of water.
S-46 If swallowed seek medical advice immediately and show the container or label.
3. **WARNING: Potential Biohazardous Material:** This kit contains some reagents made with human serum or plasma. The serum or plasma used has been tested by and FDA-approved method and found to be non-reactive for HIV-1/2 Antibodies, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 3rd Edition, 1993.

6. SPECIMEN COLLECTION AND STORAGE

Draw specimen at room temperature. **Do not chill EDTA-Na₂ (2 mg/ml blood) collection tube or store the collection tube on ice.** Centrifuge the tube in a **non-refrigerated** centrifuge, and separate the plasma from the cells immediately after centrifugation, then aliquot and freeze immediately. Do not refrigerate. Cryoactivation, that is conversion of prorenin to renin, occurs when samples are chilled or exposed to temperatures of 6°C or lower for extended periods of time and when samples are chilled but remain in a liquid state (i.e not frozen) (12, 13). Before analysis samples should be thawed rapidly to room temperature by standing them in lukewarm water. Patients taking diuretics, antihypertensives, vasodilators, oral contraceptives, and licorice should discontinue use of these substances for two to four weeks before the test. It should be noted that renin is increased in pregnancy and in diets with reduced salt intake. Also, since renin is affected by body position, as well as by diurnal (daily) variation, blood samples should be drawn in the morning, and the position of the patient (sitting or lying down) should be noted. Eventual dilutions can be performed after Angiotensin I generation, using "0" standard supplied with the kit.

Warning:

- Do not use heparin as anticoagulant because it causes interference in Angiotensin I generation.
- Do not use hemolyzed samples.

7. ASSAY PROCEDURE

1. Materials Provided

The Renin Kit (R-EX-125) is sufficient for 125 tubes and contains the following reagents:

Reagents	Quantity
Renin Antiserum	1
Renin Standards	8
Renin Tracer	1
Renin Separation Reagent	1
Renin Generation Buffer pH 6.0	1
Renin Enzyme Inhibitor	1

2. Reagents Not Supplied with kit provided by request

Generation Buffer pH 7.4 - Code: R-EX-LC.

3. Materials and Equipment Required but not Provided

- Distilled water;
- Plastic test tubes (about 1 x 7);
- Water bath at 37°C;
- Graduates pipettes;
- Ice bath;
- Automatic micropipettes with disposable tips (0.01 - 0.1 ml);
- 4°C Refrigerator;
- Vortex Mixer;
- Permanent rack: 50 tubes disposable rack to be slid into Magnetic Separator;
- Magnetic Separator: Magnetic plate for B/F separation;
- Gamma-counter.

4. Test Procedure

4.1 Angiotensin I Generation

The user can choose the generation pH by using generation buffer pH 7.4 (orange solution - provided only by request) or generation buffer pH 6.0 (yellow solution - supplied with the kit).

In both the cases follow this procedure, add to plastic tubes, following the indicated order:

- 1.0 ml sample;
- 0.01 ml enzymatic inhibitor (PMSF);
- 0.1 ml buffer;
- Mix on vortex;

Prepare for each sample 4 tubes and put them in 2 Rack-:

- 2 tubes for 37°C incubation (marked with 37);
- 2 tubes for ice-bath incubation (marked with 4).

These tubes will be used for radioimmunoassay.

- pipette 0.1 ml sample treated with inhibitor in each tube. Put the 4 series in ice-bath and 37 series in water-bath at 37°C.
- incubation time is 3 hours if pH 7.4 generation buffer is used or 90 minutes if pH 6.0 generation buffer is used.
- at the end of incubation put 37 series tubes in ice-bath where 4 series tubes were kept. Then RIA is performed.

4.2 Radioimmunoassay Procedure

Put the following groups of tubes into ice-bath, together with the two series samples tubes:

- 2 tubes for total radioactivity;
- 2 tubes for NSB;
- 2 tubes for B0;
- 2 tubes for each standard concentration;
- Pipette 0.1 ml standard "0" into B0 tubes;
- Pipette 0.2 ml standard "0" into NSB tubes;
- Pipette 0.1 ml each standard into conveniently marked tubes;
- Add 0.1 ml Renin Tracer into all tubes;
- Add 0.1 ml Renin Antiserum into all tubes except into total radioactivity;
- Mix on vortex and incubate for 18-20 hours in ice-bath put in a refrigerator (4°C);
- At the end of incubation take the ice-bath out of the refrigerator and always keeping the tubes in ice-bath add 1 ml Renin Separation Reagent except to total radioactivity. During dispensing, the magnetic second antibody should be occasionally mixed to ensure the uniformity.
- Mix the whole rack on vortex; remove the rack from ice-bath and incubate for 10 minutes at room temperature. Remove total tubes from the rack and slide the rack into Magnetic Separator; control that the bottoms of tubes are in contact with the surface of separator.
- Allow the magnetic suspension to sediment for 10 minutes.

- Decant the tubes by inversion of magnetic separator; place the inverted separator on absorbent paper in order to remove the drops of liquid adhering to the sides of the tubes.
- Count all tubes, also total radioactivity, for 1 at least minute in Gamma-counter.

4.3 Assay Procedure Scheme - (The volumes are expressed in ml)

Reagents	Tubes	Sample "4"- "37"	B0	Standard	NSB	Total Radioactivity
Sample	0.1 - 0.1	-	-	-	-	-
"0" Standard	- -	0.1	-	-	0.2	-
Standards	- -	-	-	0.1	-	-
Tracer	0.1 - 0.1	0.1	0.1	0.1	0.1	0.1
Antiserum	0.1 - 0.1	0.1	0.1	0.1	-	-
Mix on vortex and incubate for 18-20 hours at 4°C in ice-bath to be placed in refrigerator. At the end of incubation, take the ice-bath with the tubes out of the refrigerator.						
Separation Reagent	1.0 1.0	1.0	1.0	1.0	1.0	-
Mix the rack on vortex and incubate for 10 minutes at room temperature. Remove total tubes. Slide the rack into Magnetic Separator and wait for 10 minutes. Count all tubes, also total radioactivity, for at least 1 minute in Gamma-counter.						

8. CALCULATIONS AND RESULTS

NSB mean counts may be subtracted from all mean counts except from total radioactivity.

B0 mean counts are used to calculate the % binding in absence of cold antigen (maximum binding).

$$\frac{\text{B0 mean counts}}{\text{Total radioactivity mean counts}} \times 100 = \% \text{ binding}$$

% relative binding of standard and samples is calculated with the following formula:

$$\frac{\text{Mean counts (standard, sample)}}{\text{B0 mean counts}} \times 100 = \% \text{ relative binding}$$

Draw the dose-response curve by plotting the % relative binding of each standard (y axis) against the relative concentration (x axis) using semilogarithmic paper.

By interpolating on the standard curve the % relative binding of each sample, one obtains the sample antigen concentration in ng/ml.

The plasmatic renin activity (PRA) is calculated as ng/ml/hour with the following formula:

$$\text{PRA} = \frac{(\text{ng/ml "37"} - \text{ng/ml "4"}) \times 1.11^*}{\text{Hours}^{**}} \text{ ng/ml/hour}$$

*1.11: Dilution factor of the sample after addition of PMSF and generation buffer.

**Hours: Incubation time for angiotensin I generation.

For data computerization we suggest to refer to bibliography (11, 12). The user can choose, according to assay characteristic, the calculation method and the graphic expression which allows the best data treatment.

1. Example of Calculation - Representative Standard Curve

Tubes	Mean cpm	B/T (%)	B/BO (%)	Concentration (ng/ml)
Total radioactivity	26114	-	-	-
NSB	274	1.05	-	-
B0	21217	81.25	-	-
St. 0.2 ng/ml	18458	-	86.8	-
St. 0.50 ng/ml	15002	-	70.3	-
St. 1.0 ng/ml	10906	-	50.8	-
St. 2.0 ng/ml	7060	-	32.4	-
St. 5.0 ng/ml	3551	-	15.6	-
St. 10.0 ng/ml	2077	-	8.6	-
St. 25.0 ng/ml	1064	-	3.8	-
Sample 1 ("4")	17053	-	80.12	0.31
Sample 1 ("37")	2103	-	8.73	9.84

$$\text{Sample 1} = \frac{(9.84 - 0.31) \times 1.11}{1.5} = 7.05 \text{ ng/ml/hour.}$$

INTERNAL QUALITY CONTROL

- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

9. LIMITATIONS

1. To obtain reliable results it is necessary to follow carefully the instruction reported in this package insert.
2. Do not use the reagents beyond the expiry date stated on vial labels.
3. Do not mix reagents from different lots.
4. After reconstitution, store the reagents at 4°C for 3 days in the stoppered vials. After reconstitution and storage at 4°C, the reagents can be used only once. For longer than 3 days storage, freeze the reagents at -20°C. After freezing, allow the reagents to reach room temperature gradually, avoiding to use water-bath at 37°C. Avoid to thaw the reagents more than once.
5. Samples contaminated with exogenous radioactivity may give inaccurate results.
6. Since angiotensin I generation is a linear function in the time, for samples supposed at low renin activity (less than 0.15 ng/ml/hour) it is possible to prolong generation time to increase assay sensitivity. In this way it is possible to read on standard curve a value which could be out of standard curve (see generation curve on section 8).
7. Plasmatic renin activity is closely correlated with sodium intake of the patient. It is necessary to evaluate both sodium urinary excretion and plasmatic renin activity to obtain a complete clinical state.

10. EXPECTED VALUES

Sample	Plasmatic renin activity ng/ml/hour	
	pH 6.0	pH 7.4
Erect	0.98 - 4.18	0.15 - 2.12
Supine	0.51 - 2.64	0.12 - 1.59

The range has been calculated on patients who have a urinary sodium excretion of about 100 mEq/24 hours. We suggest that each laboratory establishes its own normal range.

11. PERFORMANCE CHARACTERISTICS

1. Accuracy

Recovery studies were performed by adding Angiotensin I to pooled plasma samples. Angiotensin I values were determined before and after addition, and the % recovery of added Angiotensin calculated.

Sample	Angiotensin I Added (ng/ml)	Angiotensin I Measured (ng/ml)	Recovery %
1	0	1.4	-
	10	10.29	88.9
	5	6.17	95.3
	2	3.45	102.3
	1	2.43	103.0
2	0.5	1.94	107.0
	0	1.97	-
	10	12.22	102.5
	5	6.42	89.0
	2	3.85	94.3
3	1	2.97	100.7
	0.5	2.46	98.8
	0	3.21	-
	10	13.90	106.9
	5	8.25	100.8
	2	5.37	108.0
	1	4.30	109.0
	0.5	3.67	92.0

2. Dilution

In a dilution study, samples were diluted in Zero standard. Proceed to the dilution after angiotensin generation. The % recovery determined are shown below:

Sample	Dilution Factor	Angiotensin I Measured (ng/ml)	Recovery %
1	-	14.51	-
	2	7.38	101.7
	4	3.92	108.1
	8	1.83	101.1
	16	0.83	91.7
2	32	0.39	85.6
	-	19.74	-
	2	9.77	99.0
	4	5.11	103.5
	8	2.48	100.5
3	16	1.15	93.2
	32	0.6	97.3
	-	20.18	-
	2	9.66	95.7
	4	5.05	100.1
	8	2.51	99.5
	16	1.14	90.1
	32	0.62	97.7

3. Precision

Intra-Assay coefficient of variation (C.V.) was evaluated on serum pools measured 20 time in the same assay. The within-assay variability is shown below:

Sample	Mean (ng/ml)	S.D. (ng/ml)	C.V. %
1	2.86	0.15	5.1
2	4.66	0.28	6.04
3	6.41	0.22	3.39

Inter-Assay coefficient of variation (C.V.) was evaluated at several points along the standard curve in 10 different assays. The between assay variability is shown below:

Sample	Mean (ng/ml)	S.D. (ng/ml)	C.V. %
1	2.96	0.15	5.15
2	4.61	0.19	4.18
3	6.33	0.24	3.82

4. Sensitivity

The detection limit of the assay, defined as the concentration of Angiotensin I equivalent to the mean counts of 20 replicates of the Zero Standard less two standard deviations, is typically 0.033 ng/ml.

5. Specificity

The specificity of the Renin assay system was assessed by measuring the apparent response of the assay to potentially cross-reactive analytes. The percentage of interference was calculated according to Abraham's method ($x/y \times 100$) where x and y are respectively the weight of Angiotensin I and interfering compound that causes a decrease of 50% binding capacity.

Cross-Reactant	Cross-Reactivity
Angiotensin I	100.00 %
Angiotensin II	0.0023 %
Angiotensin III	0.02 %
Angiotensin Fragment 1-13	0.32 %
Angiotensin II Pentapeptide	0.015 %

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1. Kosimen and Pakarinen - Correlation between plasma renin activity, Angiotensin and Aldosterone - J.C.E.M. vol. 47, n. 3, pag. 665, 1978.
2. Brooks C.S. et al. - Renin activity in plasma of patients with normal renin and low renin essential hypertension - J.C.E.M., 44, 322, 1977.
3. Weinberger M.H. - The effect of posture and saline loading on plasma Renin activity and Aldosterone concentration in pregnant, non-pregnant and estrogen-treated women - J.C.E.M. 44, 69, 1977.
4. Speckart P. et al. - The effect of sodium restriction and prostaglandin inhibition on the renin-angiotensin system in man. - J.C.E.M. 44, 832, 1977.
5. Swales D. and Thurston H. - Plasma renin and Angiotensin II measurement in hypertensive and normal subjects: correlation of basal and stimulated states - J.C.E.M. 45, 159, 1978.
6. Re R.N. Sancho J. et al. - The characterization of low renin hypertension by plasma renin activity and plasma aldosterone concentration J.C.E.M. 46, 189, 1978.
7. Pallas K.G. et al. - The effect of conjugated estrogens on the renin-angiotensin system - J.C.E.M. 44, 1061, 1977.
8. F. Mantero et al. - Plasma renin activity and urinary aldosterone in acromegaly - J. Endocrinol. Invest 2, 13, 1979.
9. Hsueh W.A. et al. - Changes in active and inactive renin throughout pregnancy - J.C.E.M. 54, 1010, 1982.
10. Giovannelli G. et al. - Plasma renin activity (PRA) in the perinatal agefrom «Problems in pediatric endocrinology» - Ed. La Cauza-Root, Academic Press 1980.
11. Rollieri E. - Descrizione parametrica della curva dose-risposta nel metodo radioimmunologico e nelle tecniche Affini. Giorn. It. Clin. 2(1), 1-13-1977.
12. Rollieri E. et al. - Automatic treatment of Radioimmunoassay data - The Journal of Nuclear Biology and Medicine. Vol. 17, n. 3, pag. 128-141, 1973.
13. Sealey I.E. Plasma renin activity and plasma pro-renin assays. Clinical Chemistry Vol. 37 (10 Pt 2), pag. 1811-9/1991-10.
14. Plouin PF; Chatellier g.; Guyene TT; Vincent N; Corvol P. Recent advance in the clinical study of renin system. Reference values and conditions of validity. La Presse medicale: Vol. 18 (18); pag. 917-21/1989-05-06.

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