

**Revised 6 Sept. 2006 (Vers. 1.0)****NOT FOR SALE IN THE USA**

## **1 INTRODUCTION**

The **DRG® CA 15-3 Enzyme Immunoassay Kit** provides materials for the quantitative determination of the cancer associated antigen CA 15-3 in serum and plasma.

**This assay is intended for in vitro use only.**

Breast cancer is one of the most common malignancies among women. Metastatic disease may be present at the time of initial diagnosis and can occur at any time following primary therapy. Up to 70% of patients with metastases will respond to systemic treatment with cytotoxic drugs or endocrine therapy; therefore, early detection of recurrence is important to patient management.(3)

In patients previously treated for stage II or stage III breast cancer, early detection of recurrence cannot be readily accomplished by routine clinical or diagnostic studies alone. The use of a circulating serum tumor marker assay, such as DRG® CA 15-3 ELISA, can be useful in the identification of these patients.

The DRG® CA 15-3 ELISA utilizes two monoclonal antibodies (115D8 and DF3) which react with DF3 reactive determinants, expressed by human breast carcinoma cells. The 115D8 antibody was raised against antigens of human milkfat globule membranes. (6-9) The DF3 antibody was prepared against a membrane-enriched fraction of a human breast carcinoma.(10,11) The results of research studies indicate that the concentration of the DF3-reactive determinants is frequently elevated in the serum of patients with breast cancer as well as with other malignancies, such as lung cancer, and in some non-malignant disorders.(12-15)

CA 15-3 assay values were not elevated in the sera of the majority of normal individuals or those with nonmalignant conditions.(12-15)

**NOTE:** CA 15-3 values determined with different assays and from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the assay used. Assay values obtained with different assay methods cannot be used interchangeably.

## **2 PRINCIPLE OF THE TEST**

The DRG® CA 15-3 ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle. The microtiter wells are coated with a monoclonal antibody directed towards a unique antigenic site on a CA 15-3 molecule.

An aliquot of patient sample containing endogenous CA 15-3 is incubated in the coated well with enzyme conjugate, which is an anti-CA 15-3 monoclonal antibody conjugated with horseradish peroxidase. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase is proportional to the concentration of CA 15-3 in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of CA 15-3 in the patient sample.

## **3 PRECAUTIONS**

- This kit is for in vitro use only.
- For information on hazardous substances included in the kit please refer to Material Safety Data Sheets.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Avoid contact with Stop Solution containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.

**Revised 6 Sept. 2006 (Vers. 1.0)****NOT FOR SALE IN THE USA**

- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even if the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according the national biohazard safety guideline or regulation.
- Safety Data Sheets for this product are available upon request directly from DRG International.  
The Safety Data Sheets fit the demands of: EU-Guideline 91/155 EC.



Revised 6 Sept. 2006 (Vers. 1.0)

NOT FOR SALE IN THE USA

## 4 KIT COMPONENTS

### 4.1 Contents of the Kit

1. **Microtiterwells**, 12x8 (break apart) strips, 96 wells  
Wells coated with anti-CA 15-3 monoclonal antibody
2. **Zero Standard**, 1 vial, 3 mL, ready to use  
contains 0.010% methylisothiazolone, 0.015% bromonitrodioxane as preservative.
3. **Standard (Standard 1-4)**, 4 vials, 0.5 mL, ready to use  
Concentrations 25, 50, 100, 200 U/mL;  
contain 0.010% methylisothiazolone, 0.015% bromonitrodioxane as preservative.
4. **Control**, 1 vial (lyoph.), 0.5 mL  
see „Reagent Preparation“  
Control values and ranges please refer to vial label or QC-Datasheet.  
Contain 0.3% Proclin 300 as a preservative.
5. **Assay Buffer**, 1 vial, 30 mL, ready to use:  
contains 0.010% methylisothiazolone, 0.015% bromonitrodioxane as preservative.
6. **Enzyme Conjugate**, 1 vial, 14 mL, ready to use  
Anti-CA 15-3 antibody conjugated to horseradish Peroxidase  
contains 0.01% methylisothiazolone, 0.01% bromonitrodioxane, 10 ppm Proclin 300 as preservative.
7. **Substrate Solution**, 1 vial, 14 mL, ready to use  
TMB.
8. **Stop Solution**, 1 vial, 14 mL, ready to use  
contains 0.5M H<sub>2</sub>SO<sub>4</sub>.  
Avoid contact with the stop solution. It may cause skin irritations and burns.
9. **Wash Solution**, 1 vial, 30 mL (40X concentrated);  
see „Preparation of Reagents“.

#### 4.1.1 Equipment and material required but not provided

- A microtiterplate calibrated reader (450±10 nm) (e.g. the DRG International Microtiterplate Reader).
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Aqua dest.

### 4.2 Storage and stability of the Kit

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. All opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foilbag has been opened, care should be taken to close it tightly again.

### 4.3 Preparation of Reagents

Allow all reagents and required number of strips to reach room temperature prior to use.

#### Control

**Revised 6 Sept. 2006 (Vers. 1.0)****NOT FOR SALE IN THE USA**

Reconstitute the lyophilized content with 0.5 mL Aqua dest. and let stand for 10 minutes in minimum. Mix the control several times before use.

*Note: The reconstituted control should be apportioned and stored at -20°C.*

**Wash Solution**

Dilute 30 mL of concentrated Wash Solution with 1170 mL deionized water to a final volume of 1200 mL.

*The diluted Wash Solution is stable for 2 weeks at room temperature.*

**4.4 Disposal of the Kit**

The disposal of the kit must be made according to the national official regulations. Special information for this product is given in the Material Safety Data Sheets (see chapter 13).

**4.5 Damaged Test Kits**

In case of any severe damage of the test kit or components, DRG® have to be informed written, latest one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

**5 SPECIMEN**

Serum or plasma (EDTA-, and Heparin plasma) can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

**5.1 Specimen Collection****Serum:**

Collect blood by venipuncture (e.g Sarstedt Monovette # 02.1388.001), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

**Plasma:**

Whole blood should be collected into centrifuge tubes containing anti coagulant and centrifuged immediately after collection.

(E.g for EDTA plasma Sarstedt Monovette – red cap - # 02.166.001; for Heparin plasma Sarstedt Monovette – orange cap - # 02.165.001.)

**5.2 Specimen Storage**

Specimens should be capped and may be stored for up to 48 hours at 2-8°C prior to assaying.

Specimen held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.



Revised 6 Sept. 2006 (Vers. 1.0)

NOT FOR SALE IN THE USA

### 5.3 Specimen Dilution

If in an initial assay, a serum specimen is found to contain more than the highest standard, the specimens can be diluted 10-fold or 100 fold with *Zero Standard* and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL Serum + 90 µL *Zero Standard* (mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Zero Standard* (mix thoroughly).

## 6 TEST PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipet tips for each standard, control of sample in order to avoid cross-contamination
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents be ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.



Revised 6 Sept. 2006 (Vers. 1.0)

NOT FOR SALE IN THE USA

## 6.2 Assay Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiterwells in the holder.
2. Dispense **10 µL** of each Standard, controls and samples with new disposable tips into appropriate wells of the first strip.
3. Dispense **250 µL Assay Buffer** into each well of the first strip.  
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Proceed like described in steps 2 and 3 for the following strips 2 – 12 until all secured strips are filled with sample and Assay Buffer.  
Please note: It is important to do the pipetting of Standard, controls, samples and Assay Buffer strip by strip, because incubation of the sample in absence of Assay Buffer may lead to increased CA 15-3 values.
5. Incubate for **60 minutes** at room temperature without covering the plate.
6. Briskly shake out the contents of the wells.  
Rinse the wells 3 times with diluted Wash Solution (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Dispense **100 µL Enzyme Conjugate** into each well.
8. Incubate for **30 minutes** at room temperature without covering the plate.
9. Briskly shake out the contents of the wells.  
Rinse the wells 3 times with diluted Wash Solution (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.
10. Add **100 µL of Substrate Solution** to each well.
11. Incubate for **30 minutes** at room temperature.
12. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
13. Read the OD at **450±10 nm** with a microtiterplate reader **within 10 minutes** after adding the *Stop Solution*.

## 6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.


**Revised 6 Sept. 2006 (Vers. 1.0)**
**NOT FOR SALE IN THE USA**

Below is listed a typical example of a standard curve with the CA 15-3 ELISA.

<b>Standard</b>	<b>Optical Units (450 nm)</b>
Zero Standard (0 U/mL)	0.03
Standard 1 (25 U/mL)	0.83
Standard 2 (50 U/mL)	1.24
Standard 3 (100 U/mL)	1.60
Standard 4 (200 U/mL)	2.05

## 7 EXPECTED VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DRG® CA 15-3 ELISA the following values are observed:

<b>Population</b>	<b>5% Percentile</b>	<b>95% Percentile</b>
females and males	8.3 U/mL	33.8 U/mL

The results correspond to published values.

## 8 ASSAY CHARACTERISTICS

### 8.1 Assay Dynamic Range

The range of the assay is between 0 – 200 U/mL.

### 8.2 Specificity of Antibodies (Cross Reactivity)

Cross-reactivities of the assay are not known.

### 8.3 Analytical Sensitivity

The analytical sensitivity was calculated from the mean plus two standard deviations of twenty (20) replicate analyses of the *Zero Standard* and was found to be 0.905 U/mL.

### 8.4 Precision

#### 8.4.1 Intra Assay Variation

The within assay variability is shown below:

<b>Sample</b>	<b>n</b>	<b>Mean (U/mL)</b>	<b>CV (%)</b>
1	20	18.07	3.78
2	20	28.59	5.53


**Revised 6 Sept. 2006 (Vers. 1.0)**
**NOT FOR SALE IN THE USA**

## 8.4.2 Inter Assay Variation

The between assay variability is shown below:

Sample	n	Mean (U/mL)	CV (%)
1	12	19.26	7.86
2	12	29.85	5,20

## 8.5 Accuracy

### 8.5.1 Quality Control

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG® directly.

### 8.5.2 Recovery

Samples have been spiked by adding CA 15-3 solutions with known concentrations in a 1:1 ratio.

The expected values were calculated by addition of half of the values determined for the undiluted samples and half of the values of the known solutions. The % Recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100.

Sample	Added Concentration 1:1 (v/v) (U/mL)	Measured Conc. (U/mL)	Expected Conc. (U/mL)	Recovery (%)
1	0	19.98		
	25	19.39	22.49	86.2
	50	30.12	34.99	86.1
	100	51.46	59.99	85.8
	200	124.04	109.99	112.8
2	0	35.29		
	25	34.40	30.15	114.1
	50	47.73	42.65	111.9
	100	74.54	67.65	110.2
	200	129.17	117.65	109.8





Revised 6 Sept. 2006 (Vers. 1.0)

**NOT FOR SALE IN THE USA**

### 8.5.3 Linearity

Sample	Dilution	Mean Conc. (U/mL)	Recovery (%)
1	None	107.04	
	1:2	55.32	103.4
	1:4	28.53	106.6
	1:8	15.23	113.8
	1:16	7.02	105.0
2	None	51.73	
	1:2	28.56	110.4
	1:4	13.14	101.6
	1:8	6.16	95.2
	1:16	3.67	113.5
3	None	33.37	
	1:2	16.24	97.4
	1:4	9.16	109.8
	1:8	4.66	111.8
	1:16	2.31	110.9

## 9 LIMITATIONS OF USE

### 9.1 Interfering Substances

Any improper handling of samples or modification of this test might influence the results.

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

The assay contains reagents to minimize interference of HAMA and heterophilic antibodies. However, extremely high titers of HAMA or heterophilic antibodies may interfere with the test results.

### 9.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of CA 15-3 in a sample.

### 9.3 High-Dose-Hook Effect

No hook effect was observed in this test up to 3000 U/mL of CA 15-3.

**Revised 6 Sept. 2006 (Vers. 1.0)****NOT FOR SALE IN THE USA****10 LEGAL ASPECTS****10.1 Reliability of Results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG<sup>®</sup>.

**10.2 Therapeutical Consequences**

Therapeutical consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 10.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutical consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutical consequences.

**10.3 Liability**

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 10.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

**Revised 6 Sept. 2006 (Vers. 1.0)****NOT FOR SALE IN THE USA****11 REFERENCES**

1. Cancer Facts & Figures - American Cancer Society 1992.
2. Parker S, et al.: Cancer Statistics. CA 1996; 65: 5-27.
3. Chittoor S, Swain S: Adjuvant Therapy in Early Breast Cancer. Am Fam Physician 1991;44:453-462.
4. Patterson AGH, Lees AW, Hanson J, et al.: Impact of chemotherapy on survival in breast cancer. Lancet 1980;2:312.
5. Fey MF, Brunner KW, Sonntag RW: Prognostic factors in metastatic breast cancer. Cancer Clin Trials 1981;4:237-247.
6. Hilkens, J., et al.: Monoclonal antibodies against human milkfat globule membranes. Prot. Biol. Fluids 1981;29:813-816.
7. Hilkens, J., et al.: Monoclonal antibodies against human milkfat globule membranes useful in carcinoma research. Prot. Biol. Fluids 1984;31:1013-1016.
8. Hilkens, J., et al.: MAM-6 antigen, a new serum marker for breast cancer monitoring. Cancer Research 1986;46:2582-2587.
9. Hilkens, J., et al.: Monoclonal antibodies against human milkfat globule membranes detecting differentiation antigens of the mammary gland and its tumors. Int. J. Cancer 1984;34:197-206.
10. Kufe, D., et al.: Differential reactivity of a novel monoclonal antibody (DF3) with human malignant versus benign breast tumor. Hybridoma 1984;3(3):223-232.
11. Sekine, H., et al.: Purification and characterization of a high molecular weight glycoprotein detectable in human milk and breast carcinomas. J. Immunol. 1985;135(5):3610-3615.
12. Fujino, N., et al.: Clinical evaluation of an immunoradiometric assay for CA 15-3 antigen associated with human mammary carcinomas: Comparison with carcinoembryonic antigen. Jpn. J. Clin. Oncol. 1986;46:335-345.
13. Colomer, R., et al.: CA 15-3 Early results of a new breast cancer marker. Anticancer Research 1986;6:683-684.
14. Hayes, D.F., et al.: Comparison of circulating CA 15-3 and carcinoembryonic antigen levels in patients with breast cancer. J. Clin. Oncol. 1986;4(10):1542-1550.
15. Pons-Anticet, D.F.M., et al.: Value of CA 15-3 in the follow-up of breast cancer patients. Brit. J. Cancer 1987;55(5):567-569.