1,25-(OH)₂-Vitamin D ImmuTube[®] ELISA Kit

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Hersteller


1,25-(OH)₂-Vitamin D ImmuTube® ELISA Kit

*For the in vitro determination of 1,25-(OH)₂ vitamin D in
plasma and serum*

Valid from 21.07.2011

K 2113

Immundiagnostik AG

Content

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	18
6. SAMPLE PREPARATION	19
7. ASSAY PROCEDURE	19
<i>Principle of the test</i>	19
<i>Extraction of 1,25-(OH)₂-Vitamin D₃</i>	20
<i>Pre-incubation</i>	21
<i>Test procedure</i>	21
8. RESULTS	22
9. QUALITY CONTROL	23
<i>Expected values</i>	23
10. PRECAUTIONS	23
11. TECHNICAL HINTS	24
12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24
13. REFERENCES	24

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of 1,25-dihydroxy vitamin D in serum and plasma. Before the ELISA test, column chromatography by "ImmuTube®" is used for separation and concentration of 1,25-(OH)₂-Vitamin D in the samples. The test kit is for *in vitro* diagnostic use only.

2. INTRODUCTION

Vitamin D is either produced in the skin (under the influence of UV light) or taken up from nourishment. The storage type of vitamin D, namely 25-hydroxy vitamin D, is formed in the liver. The hormone 1,25-dihydroxy vitamin D (D hormone) is formed in a second hydroxylation step in the kidney. The responsible enzyme, the kidney 1 α -hydroxylase, is subjected to a rigid control through hormones (especially parathyroid hormone) and its activity is influenced by the serum concentrations of calcium and phosphate.

The serum concentration of 1,25-dihydroxy vitamin D normally re-adjusts itself to the demands of metabolism. Deviations from the normal range of 1,25-dihydroxy vitamin D must therefore always be interpreted in the context of the remaining parameters of the calcium metabolism. The serum concentration of 1,25-dihydroxy vitamin D decreases only in seldom cases of vitamin D deficiency. For the diagnosis of vitamin D deficiency the precursor metabolite, 25-hydroxyvitamin D, should be measured.

The reason for a non-physiological deficiency of 1,25-dihydroxy vitamin D can be found in metabolic disturbances, caused either by genetic defects of the enzyme 1 α -hydroxylase (rare) or kidney malfunctions (more common). Even a slightly impaired kidney function can lead to a decrease of the 1,25-dihydroxy vitamin D concentration.

Since 1,25-dihydroxy vitamin D has important functions in calcium metabolism as well as supplementing secretion of parathyroid hormone from the parathyroid glands, increasing kidney malfunctioning leads to development of renal osteopathy, which is characterized by osteomalacia and osteitis fibrosa.

Treatment of renal osteopathy consists of the administration of 1,25-dihydroxy vitamin D (calcitriol) or the prohormone 1 α -hydroxy vitamin D. In renal tubules malfunctions decreased or relatively low levels of 1,25-dihydroxy vitamin D (e.g. diabetes insipidus, Fanconi-syndrom) are found. A non-physiological over-production of 1,25-dihydroxy vitamin D arises in granulomatosis (e.g. sarcoidosis), where extra-renal synthesis of 1,25-dihydroxy vitamin D occurs. This can lead to hypercalcaemia. Also in idiopathic hypercalciuria a relatively high level of 1,25-dihydroxy vitamin D is found. Increased concentrations of 1,25-dihydroxy vitamin D can be seen in case of non-functional vitamin D receptors (rare), during calcium deficient nutrition, as well as a result from overproduction of parathyroid hormone (primary hyperthyroidism).

Most of the 1,25-(OH)₂-Vitamin D does not circulate in the free form. Like 25-(OH)-Vitamin D, 1,25-(OH)₂-Vitamin D is bound and transported by the Vitamin D Binding Protein and by Albumin. It is nearly insoluble in water.

The following binding percentages of 1,25-(OH)₂-Vitamin D and of 25-(OH)-Vitamin D by Vitamin D Binding Protein (VDBP) and by Albumin are mentioned in literature:

	1,25-(OH)₂-Vitamin D	25-(OH)-Vitamin D
VDBP	60 %	85 %
Albumin	38 %	14.6 %
free	2 %	0.4 %

The serum concentrations of 25-(OH)-Vitamin D and of 1,25-(OH)₂-Vitamin D are quite different. In the literature a reference value range between 30-100 ng/ml for 25-(OH)-Vitamin D and 17-65 pg/ml for 1,25-(OH)₂-Vitamin D are considered.

Indications:

- Defect of kidney functions
Chronic kidney failure
Haemodialysis following kidney transplantation
- Renal osteopathy
- Osteomalacia from various types of vitamin D metabolism disturbances
- Kidney tubules function disturbances (diabetes insipidus, Fanconi-Syndrom)
- Monitoring of therapy with active vitamin D metabolites
- Ideopathic hypercalciuria
- Hypercalcaemia

3. MATERIAL SUPPLIED

Catalogue No.	Content	Kit Components	Quantity
K2113MTP	PLATE	One holder with precoated strips	12x8 Wells
K2113WP	WASHBUF	ELISA wash buffer concentrate 10x	100 ml
K2113E	ETHANOL	Ethanol, ready-to-use	1.5 ml
K2113A1	AB	Detection antibody, anti 1,25-(OH) ₂ vitamin D, stock solution	100 µl
K2113ST	STD	Standard, ready-to-use (for range see specification)	6 x 1.5 ml
K2113KO	CTRL	Controls, ready-to-use (for range see specification)	2 x 1.5 ml
K2113K	CONJ	Conjugate, polyclonal peroxidase-labeled antibody, ready-to-use	25 ml
K2113TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	2 x 15 ml
K2113AC	STOP	ELISA stop solution, ready-to-use	15 ml
K2113CWB	COLBUF	Wash buffer for ImmuTube® columns	2 x 65 ml
K2113ER	ELUREAG	Elution reagent for ImmuTube® columns	30 ml
K2113ABBUF	ABBUF	Antibody dilution buffer, ready-to-use	2 x 15 ml
K2113FOL	FOL	Foil to cover the microtiter plate	2 x 1

ADDITIONALLY REQUIRED AND ONLY IN COMBINATION WITH K2113

K2113.SI	COLUMNS	ImmuTube® columns for separation of 1,25-(OH) ₂ vitamin D in the samples	96 columns
----------	---------	---	------------

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest., **sterile and pyrogen-free**)
- 75 x 12 mm glass tubes (no plastic)
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge
- Vortex-Mixer
- Vacuum centrifuge or nitrogen distributor
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)
- "Mixing-rotator" (end-over-end rotation) or bottle roller

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used for analysis of 88 samples, 6 standards and 2 controls, when single values are obtained. However, we recommend to run the test in duplicate.
- The ELISA **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for two weeks**.
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.
- Dilute the **AB** (detection antibody) **1:1000** in **ABBUF** (antibody dilution buffer). In order to use the kit in two separate runs, 2 **ABBUF**-vials are provided with the kit package. Transfer **15 µl** from the **AB**-vial into one **15 ml ABBUF**-vial and label it correspondingly. The so **diluted antibody** is stable for **4 weeks** when stored at 2-8°C. The **AB** (detection antibody) is stable at **2-8°C** until the expiry date stated on the label.

6. SAMPLE PREPARATION

Serum/plasma samples

Fresh collected blood should be centrifuged within one hour. Store samples at -20 °C if not assayed on the same day. Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.

Sample volume of **300 µl** is needed. Sample volume < 300 µl leads to wrong results.

Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive Enzyme-Immuno-Assay (EIA) technique with a selected monoclonal antibody recognizing 1,25-dihydroxy vitamin D.

For a reliable determination of 1,25-(OH)₂-Vitamin D, it is necessary to separate the 1,25-(OH)₂-Vitamin D from the 1000 fold excess of 25-(OH)-Vitamin D and to concentrate it. This is performed by the use of a column with immobilized monoclonal anti 1,25 (OH)₂-Vitamin D mouse antibodies.

Standards, controls and patient samples which are assayed for 1,25-dihydroxy vitamin D are incubated after the extraction step with the detection antibody. The pre-incubated solution is then transferred to the microplate coated with 1,25-dihydroxy vitamin D. During this incubation step, 1,25-dihydroxy vitamin D in the sample and a fixed amount of 1,25-dihydroxy vitamin D bound to the microtiter well compete for the binding of the detection antibodies. Then a peroxidase-conjugated anti-mouse antibody is added into each microplate well and a complex of 1,25-dihydroxy vitamin D - detection antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 1,25-dihydroxy vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 1,25-dihydroxy vitamin D in the samples is determined from this curve.

Extraction of 1,25-(OH)₂-Vitamin D₃

1. Prior to use in the assay allow all samples and reagents to come to room temperature (18-26 °C). Before use, mix well samples and reagents.
2. Place Imm Tubes® in a suitable rack and make sure that no suspension remained on the Imm Tubes® cover
3. Label the covers of Imm Tubes®, open Imm Tubes®, add quickly 300 µl of STD/ SAMPLE/CTRL (Standard/Sample/Control), close and mix gently .
4. "Mix-rotate" (end-over-end rotation) intensively for 90 min at room temperature (18-26 °C). Let the remaining separation material on the inner side of the cover flow down
5. Insert Imm Tubes® in plastic reagent vials (min. volume of 5 ml), centrifuge for 1 min at 1000g .
6. Open the Imm Tubes® by breaking the outlets and unscrewing the covers. Centrifuge for 2 min at 1000g to dryness .
7. Add 600 µl of COLBUF (Imm Tube® Washing Buffer) and allow to interact for 10 min . Centrifuge for 2 min at 1000g to dryness.
7a. Add 600 µl of COLBUF (Imm Tube® Washing Buffer) and immediately centrifuge for 2 min at 1000g .
8. Add 600 µl of sterile aqua bidest. and centrifuge for 2 min at 1000g to dryness.
9. Repeat step 8: add 600 µl of sterile aqua bidest. and centrifuge for 2 min at 1000g to dryness. Discard the through-flow.
10. Label fresh glass tubes, transfer Imm Tubes® in the labeled glass tubes.
11. Add 250 µl of ELUREAG (Elution Reagent), wait 1 min , centrifuge for 1 min at 1000g and collect the eluate with the 1,25-(OH) ₂ -Vitamin D in the glass tubes.
12. Evaporate the 250 µl eluate under a nitrogen stream at 37 °C or in a vacuum centrifuge.

Pre-incubation

Note, always dilute the **AB** (detection antibody) approximately 2 hours before use in the assay (see P. 5., Preparation and storage of reagents).

1. Add 10 µl of ETHANOL (ethanol) into each glass tube. Immediately after adding ethanol gently vortex each tube to avoid any possible evaporation.
2. Add 225 µl diluted detection antibody into each glass tube. The antibody solution is viscous. Mix thoroughly.
3. Cover glass tubes and incubate for 1 hour at room temperature (18-26°C).

Test procedure

1. Mark the positions of STD/SAMPLE/CTRL (Standards/Sample/Control) on a protocol sheet.
2. Take microtiter strips out of the microtiter module. Unused strips must be covered with the enclosed foil, stored at 2-8° C and used within 12 weeks.
3. Add 200 µl of STD/SAMPLE/CTRL into respective well. All these solutions are viscous. Please pipette slowly and carefully. We recommend to wet the pipette tip before using it to transfer the pre-incubate.
4. Cover the plate with the enclosed foil tightly and incubate for 16-22 hours at 4-8°C* .
5. Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
6. Add 200 µl CONJ (conjugate) into each well.
7. Cover the plate tightly with the enclosed foil and incubate for exactly 1 hour at room temperature (18-26°C) on a horizontal mixer.

- | |
|--|
| 8. Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper. |
| 9. Add 200 µl of SUB (substrate) into each well. |
| 10. Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark**. |
| 11. Add 50 µl of STOP (stop solution) into each well, mix thoroughly. |
| 12. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference. |

* As with any competitive immunoassay, consistent incubation times and temperature are essential for accurate plate-to-plate comparisons. Fluctuations in overnight incubation can lead to increased inter-assay CV's. It is therefore recommended to use always the same incubation time, i.e. 18 hours.

** The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the „4-Parameter-algorithm“.

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. Point-to-Point-algorithm

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

9. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Reference range (plasma or serum):

Suckling babes and children up to 12:	30 - 180 pg/ml
Teenager and young adults (age 13-21):	30 - 100 pg/ml
Healthy adults (age 21-60) :	17 - 65 pg/ml
Persons older than 61:	10 - 50 pg/ml
Pregnant women (8-42 week):	ca. 60% higher values than adults

The normal range is independent of the season.

We recommend each laboratory to establish its own norm concentration range.

10. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

11. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

13. REFERENCES

1. Durham B.W. et al. (1995) Ann. Clin. Biochem., 32, 77
2. Hollis B.W. (1996) Calcif. Tissue Int., 58, 4
3. Hollis B.W. (1995) Clinical Chemistry, 41, 1313
4. Iqbal S.J. et al. (1996) Clinical Chemistry, 42, 112
5. Schmidt-Gayk H, Armbruster FP, Bouillon R (Eds.) (1990) Calcium Regulating Hormones, Vitamin D Metabolites, and Cyclic AMP. Assays and Their Clinical Application. Springer-Verlag Berlin, Heidelberg, New York, London, Paris, Tokyo, Hong Kong, Barcelona
6. Armbruster FP et al. (2000) Clin Lab 46(3-4):165-66
7. Erben RG et al. (2002) Mol Endocrinol 16(7):1524-37
8. Peterlik M. and Cross H. S. (2005) Review: Vitamin D and calcium deficits predispose for multiple chronic diseases. European Journal of Clinical Investigation, 35, 290-304

Used Symbols:

Store at

Catalog Number

In Vitro Diagnostic Device

No. of tests

Manufacturer

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49 (0) 62 51/70 19 00

Fax: +49 (0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com