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Urotensin II ELISA Kit

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
Urotensin II ELISA Kit

For the in vitro determination of Urotensin II in serum/plasma and cell culture supernatant

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Table of contents	Page
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	15
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
TEST PROCEDURE	17
10. RESULTS	20
TYPICAL CALIBRATION CURVE	20
11. LIMITATIONS	21
12. QUALITY CONTROL	21
13. LITERATURE	21
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	21

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Manual


Urotensin II ELISA Kit

*For the in vitro determination of Urotensin II in serum/plasma
and cell culture supernatant*

Valid from 22.07.2008

REF

K 1013

**RUO**

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Urotensin II** in serum/plasma and cell culture supernatant. For research use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Urotensin II can be found in the spinal cord tissue from many different species (fish, vertebrae, primaten) and probably plays a role there as a hormone in the neuro-secretory-system. It acts like a peptide, consisting of 11 amino acids and the cyclic region (CFWKYC), which is responsible for biological activity.

In 1999, Urotensin II was first defined as a ligand for the GPR14 receptor, which belongs to the Somatostatin/Opioid receptor family and among other things, is expressed in the endothelial cells of the blood vessels.

Just as In-Vitro diagnostics can show isolated blood vessels, Urotensin II acts here as a vaso-constrictive "Somatostatin-similar" peptide, nevertheless, the vaso-constrictive capacity far exceeds that of the aforementioned vaso-active substances eg. Endotheline-1, Serotonin, Noradrenalin or Angiotensin II. Urotensin II shows overall functional and regional differences in its effects depending on the species. At the moment, there is still no generally agreed statement about the physiological role of Urotensin II on the heart-circulatory system.

3. PRINCIPLE OF THE TEST

This test kit works on the basis of a competitive enzyme immuno-assay. The Urotensin II in the sample competes with the Urotensin-II peptide coated to the plate.

In the first step the sample or calibrator will be incubated with the specific antibody (rabbit anti-Urotensin II). After adding the mixture to the plate, the Urotensin-II in calibrator or sample competes with the Urotensin-II coated to the plate. The more analytic concentration in the sample, or calibrator, the less antibody will be bound to the immobilized peptide.

The detection and quantifying follows a peroxidase-defining secondary antibody and the appropriate reactant conversion. The formed chromogen complex quantity is inversely proportional to the analytic content. TMB is used as substrate. The resulting chromogen compound can be photographically measured at 450nm.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 1013MTP	PLATE	One holder with precoated strips	96
K 1013WP	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K 1013AP	ASYBUF	Assay buffer, ready-to-use	1 x 100 ml
K 1013A2	AB	Antibody, ready to use	1x 8 ml
K 1013K	CONJ	Conjugate (goat anti rabbit, HRP-labeled), ready to use	1 x 15 ml
K 1013ST	STD	Standards, lyophilized (0; 32; 160; 800; 4000; 20000 pgl/ml)	6 vials
K 1013TMB	SUB	TMB Substrate (Tetramethylbenzidin), ready-to-use	1 x 15 ml
K 1013AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

On request we will send you 3 calibrator sets free of charge. Any further calibrator set will be charged.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Distilled water
- Bidistilled (aqua bidest.) or deionized water
- Deep freezer -20 °C
- Precision pipettors calibrated to deliver 10-1000 µl
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Water bath or heating block
- Standard laboratory glass or plastic vials, cups, etc. (one time products)
- Microtiter plate reader 450 nm
(reference wave length 620 or 690 nm)
- Seal cover for microtiter plates

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.**
- Reagents with a **volume less than 100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. The **WASHBUF** (wash buffer concentrate) is stable **at 2-8°C** until the expiry date stated on the label. **Diluted buffer solution can be stored in a closed flask at 2-8°C for one month.**
- The lyophilized **STD** (standards) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards and control are not stable.
- The **AB** (antibody) is to remove from the kit and is to be stored at -20 °C. The **AB** is stable until the expiry date stated on the label of kit.
- All other test reagents are ready to use. Test reagents are stable at 2-8°C until the expiry date stated on the label of kit.

7. PRECAUTIONS

- For research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date stated on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum/Plasma

Samples have to be stored at -20°C . Dilute the samples **1:10** in **ASYBUF** (assay buffer) before use.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- The assay should always be performed according the enclosed manual.

Pre-incubation

1. Add **150 µl AB** (Antibody) and **150 µl STD** (Standard) and prediluted sample into a 1.5 ml test tube (Eppendorf).
2. Incubate for **1 hour** at 37 °C while shaking.

Test procedure

1.	Bring all reagents and samples to room temperature (18-26 °C) and mix well
2.	Mark the positions of STD (Standards) and samples in duplicate on a protocol sheet
3.	Take as many microtiter strips as needed from kit. Store unused strips covered at 2-8° C. Strips are stable until expiry date stated on the label
4.	Wash each well 5x with 250 µl diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
5.	Add 100 µl pre-incubated STD (Standard) and sample in duplicates into the respective wells
6.	Cover and incubate for 2 hours shaking on a horizontal mixer at 37 °C **
7.	Aspirate the contents of each well. Wash each well 5x with 250 µl diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper

8.	Add 100 µl CONJ (Peroxidase-labeled antibody)
9.	Cover and incubate for 45 min shaking on a horizontal mixer at room temperature (18-26°C)
10.	Aspirate the contents of each well. Wash each well 5x with 250 µl diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
11.	Add 100 µl SUB (TMB substrate)
12.	Incubate for 5 - 15 min shaking on a horizontal mixer in the dark at room temperature (18-26°C) *
13.	Add 50 µl STOP (stop solution) and mix shortly
14.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

*The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.


**The above incubation steps at 37 °C on a horizontal mixer are recommended by the producer. If there is no possibility to incubate at 37 °C, while shaking, we recommend to incubate at 37 °C without any shaking.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

Typical calibration curve

Calibrator	OD mean	Concentration [pg/ml]
1	1.186	0
2	0.939	32
3	0.625	160
4	0.406	800
5	0.214	4.000
6	0.107	20.000

The data is for demonstration only and cannot be used in place of data generations at the time of the assay.

Serum/Plasma

For the calculation of the sample values, the results must be multiplied by the dilution factor.

11. LIMITATIONS

Samples with Urotensin II levels greater than the highest calibrator should be diluted and re-assayed.

12. QUALITY CONTROL

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

13. LITERATURE

1. Ames et al. (1999). Nature 40: 282-286
2. Coulouarn et al. (1998). PNAS 22;95 (26):15803-8

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For reserch use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate

protective clothing. Any spills should be wiped out immediately with copious quantities of water.

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number