Osteonectin EIA Kit

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1. INTENDED USE	12
2. SUMMARY AND EXPLANATION OF THE TEST	12
3. PRINCIPLE OF THE TEST	12
4. MATERIAL SUPPLIED	13
5. MATERIAL REQUIRED BUT NOT SUPPLIED	13
6. PREPARATION AND STORAGE OF REAGENTS	14
7. PRECAUTIONS	14
8. SPECIMEN COLLECTION AND PREPARATION	15
9. ASSAY PROCEDURE	15
PROCEDURAL NOTES	15
TEST PROCEDURE	16
10. RESULTS	17
11. LIMITATIONS	17
12. QUALITY CONTROL	17
EXPECTED VALUES	18
13. REFERENCES	18
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	18

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Manual

Osteonectin EIA Kit


For the in vitro determination of Osteonectin in serum

For research use only

Valid from 20.04.2010

K 4231

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim
Tel.: ++49 6251 70190-0
Fax: ++ 49 6251 849430
e.mail: Info@immundiagnostik.com
www.immundiagnostik.com

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Osteonectin** in serum. It is for research use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Osteonectin (SPARC, BM-40) is a 32,700 molecular weight, Ca^{2+} -binding glyco-protein which is a synthetic product of osteoblasts, endothelial cells, and megakaryocytes. Osteonectin is involved in tissue remodeling both from the perspective of bone metabolism and endothelial cell proliferation and repair. Therefore, the quantitative analysis of this protein could serve as a useful marker for the researcher studying the biochemical processes governing vascular wound repair, platelet activation, and skeletal metabolism.

3. PRINCIPLE OF THE TEST

This enzyme immuno assay (EIA) can be used for the determination of **Osteonectin** in serum.

The test principle is based on a competition between antigen in the sample or standards and the antigen coated on the wells of microplate. A peroxidase-conjugated antibody is used for detection and quantification, and tetramethylbenzidine (TMB) as a peroxidase substrate. The enzymatic reaction is terminated by acidic stop solution. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from standard. Osteonectin present in the patient samples is determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 4231MTP	MTP	Microtiter plate, precoated	12 x 8 wells
K 4231WB	WASHBUF	ELISA wash concentrate 10 x	100 ml
K 4231 PV	SAMPLEBUF	Sample dilution buffer	100 ml
K 4231 A	AB	Antibody (Anti-Osteonectin antibody)	12 ml
K 4231K	CONJ	Conjugate (Peroxidase-labeled), ready to use	22 ml
K 4231ST	STD	Calibrators, ready to use (0; 0.04; 0.12; 0.37; 1.11; 3.33 µg/ml)	1 x 6 vials
K 4231KO1	CTRL	Control, ready to use	1 vial
K 4231KO2	CTRL	Control, ready to use	1 vial
K 4231TMB	SUB	TMB substrate (Tetramethylbenzidine)	2 x 15 ml
K 4231AC	STOP	ELISA stop solution, ready to use	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 5-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **STD** (standards), **CTRL** (controls) and the **AB** (antibody) can be stored with a temperature of **2-8° C** for two weeks. For a storage for more than two weeks the components must be stored with a temperature of **-20° C**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum

Serum samples must be diluted **1:10** before performing the assay, e.g.

25 µl serum + **225 µl** SAMPLEBUF (sample dilution buffer), mix well.

100 µl of the diluted sample per well are used in the test.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component within the same assay.
- The quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

Wash the precoated microtiter plate **5 x with 250 µl** ELISA wash buffer.
Carry out the tests in duplicate.

1. Add **100 µl** of **STD** (standards), **CTRL** (controls) and diluted **samples** into each well in duplicate
2. Add **100 µl** of **AB** (anti-osteonectin antibody solution) into each well
3. Incubate over night at 4°C shaking on a horizontal mixer
4. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer
5. Add **200 µl** of **CONJ** (conjugate) into each well
6. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature
7. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer
8. Add **200 µl** of **SUB** (substrate solution) into each well
9. Incubate for **10-20 minutes** at room temperature
10. Add **50 µl** of **STOP** (stop solution) into each well and mix shortly
11. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Serum

Multiply the result by **10** to get the real concentration.

11. LIMITATIONS

Samples with Osteonectin concentration greater than the calibrator with the highest concentration (with the lowest OD), should be diluted and re-assayed.

12. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Reference range

We recommend each laboratory to establish its own baseline values.

13. REFERENCES

1. Whitson, S.W., et al., J. Cell Biol. 99, 607 (1984).
2. Otsuka, K., et al., J. Biol. Chem. 259, 9805 (1984).
3. Kuwata, F., et al., J. Biol. Chem. 260, 6993 (1985).
4. Sage, H., et al., J. Biol. Chem. 259, 3993 (1984).
5. Sage, H., et al., J. Cell. Physiol. 127, 373 (1986).
6. Kelm, Jr., R.J., et al., Blood 80, 3112 (1992).
7. Sage, E.H. and Bornstein, P. J. Biol. Chem. 266, 14831 (1992).

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The test components which are made of human serum are tested for HVB and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are for research use only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any alterations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number