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Retinol-bindendes Protein (RBP)/ RBP4 ELISA Kit

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
For the in vitro determination des RBP/RBP4 in plasma, serum
and urine

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1. INTENDED USE	15
2. SUMMARY AND EXPLANATION OF THE TEST	15
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	17
7. PRECAUTIONS	18
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
10. RESULTS	21
11. LIMITATIONS	21
12. QUALITY CONTROL	22
EXPECTED VALUES	22
13. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
SENSITIVITY	23
SAMPLE DILUTION	23
14. REFERENCES	24
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24

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In-Vitro-Diagnostikum

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Hersteller

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
Retinol-binding protein (RBP)/ RBP4 ELISA Kit

For the in vitro determination of RBP/RBP4 in plasma, serum and
urine

Gültig ab / Valid from 06.04.2010

K 6120

1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of **free Retinol-binding protein (RBP)/RBP4** as well as **RBP4 complexed with transthyretin** in plasma, serum and urine. For in vitro diagnostic use only. of free as well as complexed

2. SUMMARY AND EXPLANATION OF THE TEST

Retinol-binding protein (RBP)/RBP4 is a small (21kD) transport protein for vitamin A which forms a complex with prealbumin in blood but loses its affinity for prealbumin once the vitamin has been delivered to the target cells. The free RBP/RBP4 molecule is rapidly filtered at the glomerulus and catabolized in the renal tubules after resorption by the proximal tubular cells (like other small molecules e.g. β -2 microglobulin). In kidney disease with prevailing tubular changes these proteins are not reabsorbed and appear in the urine.

As published by Yang et al. (2005) the retinol-binding protein (**RBP)/RBP4** seems to play a key role in the development of insulin resistance. The fat cell derived peptide RBP/RBP4 also modulates the glucose homeostasis and impairs the insulin sensitivity as well as insulin resistance. The elevation of serum RBP/RBP4 causes systemic insulin resistance, whereas its reduction improves the insulin action.

As a conclusion from the results, the authors suggest that **RBP/RBP4** alters insulin sensitivity in part by affecting insulin signalling in muscle through alterations in the amount of tyrosine-phosphorylated IRS-1 and PI(3)K activation. Thus, RBP/RBP4 may contribute to the pathogenesis of type 2 diabetes, and lowering RBP/RBP4 could be a new strategy for treating type 2 diabetes.

Indications

- Early detection of tubular proteinuria
- Chronic liver diseases
- Cadmium poisoning
- Studies of insulin resistance

3. PRINCIPLE OF THE TEST

This Enzyme-Linked Immunosorbent Assay (ELISA) can be used for quantitative determination of **Retinol-binding protein (RBP)/RBP4** in plasma, serum and urine. In a first incubation step, RBP/RBP4 in the samples is bound to polyclonal rabbit anti RBP/RBP4 antibodies, immobilized on the microtitre plate. A peroxidase-conjugated anti RBP/RBP4 antibody is used for detection and quantification, and tetramethylbenzidine (TMB) as a peroxidase substrate. Samples are quantified by referring their optical density to a Lot-dependant master calibration curve and the use of a Calibrator that is run with each test.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6120MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6120WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6120K	CONJ	Conjugate, (rabbit anti RBP/RBP4, peroxidase-labeled)	200 µl
K 6120PV	SAMPLEBUF	Sample dilution buffer, ready-to-use	100 ml
K 6120CAL	CAL	Calibrator, lyophilized	2 x 1 vial
K 6120KO1	CTRL	Control, lyophilized	2 vials
K 6120KO2	CTRL	Control, lyophilized	2 vials
K 6120TMB	SUB	TMB Substrate (Tetramethylbenzidin), ready-to-use	15 ml
K 6120AC	STOP	ELISA Stop solution, ready-to-use	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 5-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The ELISA **WASHBUF** (wash buffer concentrate) must be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- Use **100 µl of diluted wash buffer as a BLANK**. Pipette into the respective well.
- Lyophilized **CAL** (calibrator) and **CTRL** (controls) must be reconstituted with **500 µl** aqua bidest. Allow the vial to stand for minimum 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted Calibrator and Controls can be stored for two weeks at 2 - 8 °C.

- **CONJ** (conjugate) (POD-Antibody) must be diluted **1:100** in wash buffer (100 µl CONJ + 10 ml wash buffer). The antibody is stable at 2 -8 °C until expiry date stated on the label. **Diluted antibody solution is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date stated on the label of test package when stored at **2-8°C**.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop Solution consists of diluted Sulfuric Acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Plasma and serum

Samples can be stored for two weeks at 4°C. For longer storage, freeze at or below -20°C.

Dilute samples **1:5000** in **SAMPLEBUF** (sample dilution buffer) before use.

Dilution I:

20 µl sample + 980 µl SAMPLEBUF = 1:50

Dilution II:

10 µl Dilution I + 990 µl SAMPLEBUF = 1:100

To avoid some potential for error in RBP/RBP4 estimation due to small volumes, we recommend the following dilution steps:

Dilution I:

20 µl Sample + 980 µl SAMPLEBUF = 1:50

Dilution II:

50 µl Dilution I + 450 µl SAMPLEBUF = 1:10

Dilution III:

50 µl Dilution II + 450 µl SAMPLEBUF = 1:10

Urine

Adjust the urine to pH between 6 and 8 with 1 N NaOH. Samples are stable at 2-8°C for 2 weeks. For longer storage, freeze at or below -20°C.

Before use, dilute urine **1:10 in SAMPLEBUF** (sample dilution buffer), e.g.
100 µl urine + 900 µl SAMPLEBUF

Urine with a **RBP4 concentration > 330 µg/l** must be diluted **1:100**, e.g.
10 µl urine + 990 µl SAMPLEBUF

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

Wash the precoated **PLATE** (microtiter plate) **5 x with 250 µl** ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl of BLANK/CAL/CTRL/SAMPLE** (Blank / Calibrator / Controls / Sample) in duplicate into respective well.
2. Incubate for **1 hour at room temperature** shaking on a horizontal mixer.
3. Decant the content of the plate and wash the wells **5 x with 250 µl** of washing buffer.
4. Add **100 µl** of diluted **CONJ** (conjugate) into each well.
5. Incubate for **1 hour** at room temperature, shaking on a horizontal mixer.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** of washing buffer.
7. Add **100 µl** of **SUB** (TMB substrate solution).
8. Incubate for **10-20 minutes** at room temperature, shaking slightly until color differences are sufficient.
9. Add **50 µl** of **STOP** (stop solution) and mix shortly.
10. Determine absorption immediately with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

10. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the Calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A, B, C and D. In both cases the optical density of the calibrator (CAL) is essential.

Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Caution: Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of replicate values should always be examined after automatic evaluation of the results. If this option is not available with the used program, a control of the replicate values should be done manually by the operator.

Serum or Plasma

Multiply the result by **5000** to get the real concentration.

Urine

Multiply the result by the **dilution factor** to get the real concentration.

11. LIMITATIONS

Samples with RBP/RBP4 levels greater than the highest calibrator value, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercially control samples as internal quality control.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal range:

Plasma or Serum:

Adults 20 – 75 mg/l

Newborn 11 – 34 mg/l

Age 6 months 18 – 50 mg/l

Urine: 0.01 – 0.54 mg/l

It is recommended for each laboratory to establish its own normal range.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik RBP/RBP4 ELISA test was calculated from 16 determinations on each of two samples.

Intra-Assay CV n= 16

Sample	RBP/RBP4 Mean value [µg/l]	Intra-Assay CV [%]
1	24.1	5
2	11.1	5

The total precision (inter-assay variation) of the Immundiagnostik RBP/RBP4 ELISA test was calculated from data on 2 samples obtained by different technicians on different days.

Inter-Assay CV n= 25

Sample	RBP/RBP4 Mean value [µg/l]	Inter-Assay CV [%]
1	4.4	9.8
2	6.9	9.7

Sensitivity

The detection limit was defined as $B_0 + 2SD$ and determined to be 0.9 µg/l.

Sample dilution

Linearity n= 1

One patient sample was diluted. The results are shown below:

Sample	Dilution	Expected [µg/l]	Measured [µg/l]
A	1:7000	4.8	4.8
	1:14000	2.8	2.4
	1:28000	1.2	1.2
	1:56000	0.6	0.8

14. REFERENCES

1. Graham TE, Wason CJ, Blüher M, Kahn BB (2007) Shortcomings in methodology complicate measurements of serum retinol binding protein (RBP4) in insulin-resistant human subjects. Diabetologia DOI 10.1007/s00125-006-0557-0
2. Graham TE, Yang Q, Bluher M, Hammarstedt A, Ciaraldi TP, Henry RR, Wason CJ, Oberbach A, Jansson PA, Smith U, Kahn BB (2006) Retinol-binding protein 4 and insulin resistance in lean, obese, and diabetic subjects. N Engl J Med 354(24):2552-63
3. Yang et al. (2005), Nature 436:356-62
4. Blumsohn A et al. (1990) Clinica Chimica Acta 195 :133-38
5. Bernard AM et al. (1982) Clinica Chimica Acta 126 :1-7

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number