Adiponectin (total) ELISA Kit

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15. LITERATUR

1. Blüher (Blueher), M., et al. (2007) Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care* 30(2): p. 280-5
2. Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelty in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
3. Fernandez-Real, J.M., et al. (2004) Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care* 27(3): p. 739-45
4. Halperin F, Beckman JA, Patti ME, Trijillo ME, Garvin M, Craeger MA, Scherer PE and Goldfine AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154
5. Higashiura, K., et al. (2004) Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)* 61(6): p. 753-9
6. Koenig W, Khuseynova N, Baumert J, Meisinger C and Löwel H (2006) Serum Concentrations of Adiponectin and Risk of Type 2 Diabetes Mellitus and Coronary Heart Disease in Apparently Healthy Middle-Aged Men: Results From the 18-Year Follow-Up of a Large Cohort From Southern Germany. *JACC* 48:1369-1377

7. Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. Diabetes 55:249-259
8. Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. Clin Chem 50:1511-1525
9. Nakamura, Y., et al. (2004) Implications of plasma concentrations of adiponectin in patients with coronary artery disease. Heart 90(5): p. 528-33
10. Pajvani, U.B., et al. (2003) Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. J Biol Chem 278(11): p. 9073-85
11. Shibata, R., et al. (2004) Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. J Biol Chem 279(27): p. 28670-4

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
Adiponectin (total) ELISA Kit

For the in vitro determination of human *adiponectin* in
serum and plasma

Valid from 26.05.2011

K 6250

Content

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	16
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. PREPARATION AND STORAGE OF REAGENTS	17
6. SAMPLE PREPARATION	17
Serum and Plasma	17
7. ASSAY PROCEDURE	18
Principle of the test	18
Test procedure	18
8. RESULTS	19
9. LIMITATIONS	20
10. QUALITY CONTROL	20
Expected values	20
11. PERFORMANCE CHARACTERISTICS	21
Precision and reproducibility	21
Recovery	21
Analytical Sensitivity	22
Linearity	22
12. PRECAUTIONS	22
13. TECHNICAL HINTS	23
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
15. REFERENCES	24

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) Kit is intended for the quantitative determination of human adiponectin in serum and plasma. It is for in vitro diagnostic use only.

2. INTRODUCTION

Adiponectin consists of 244 amino acids and has an approximate molecular weight of 30 kDa. The adiponectin concentration in blood is 5-30 µg/ml that accounts for about 0.01% of the total plasma proteins. Adiponectin has four domains: a signal peptide, a variable domain, a collagen-like N-terminal domain and a globular C-terminal domain. Adiponectin exists in different oligomer forms in vivo. Beside the trimer and dimer also high molecular multimers are present. Adiponectin is mainly expressed by adipocytes, but also muscle cells and hepatocytes can synthesize it. The only, until now known, natural inducer of the adiponectin synthesis is IGF-I. Two different adiponectin receptors are known. Both of them are ubiquitarily expressed, though their distribution in the tissues varies: Adiponectin Receptor 1 (AdipoR1) is synthesized in muscle and AdipoR2 in liver tissue. The adiponectin significance for the human organism is not completely clear until now. Studies show, that adiponectin correlates negatively with BMI and could participate at the energy metabolism through the regulation of fatty acid oxidation. Adiponectin influences further physiological processes like angiogenesis. It is also associated with glucose and lipid metabolism. Adiponectin levels are associated with insulin resistance and accordingly linked with Type II Diabetes. Furthermore, it is involved in inflammatory processes and therewith of importance for development of arteriosclerosis and coronary artery diseases. Blüher et al. (2007) showed a correlation of total serum adiponectin to insulin sensitivity and the ability of total serum adiponectin levels to predict insulin resistance and impaired glucose tolerance. Low adiponectin concentrations result in inhibition of fatty acid oxidation and are associated with insulin resistance and metabolic syndrome as well as arteriosclerosis.

Indications

- Energy metabolism and body weight regulation
- Metabolic syndrome
- Type 2 diabetes
- Coronary artery disease
- Atherosclerosis

3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6250MTP	PLATE	Precoated microtiter plate	12 x 8 wells
K 6250WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6250VP	SAMLEBUF	Sample dilution buffer	2 x 100 ml
K 6250ST	STD	Adiponectin standards, lyophilized (0; 1.4; 5.5; 22; 88 ng/ml)	2 x 5 vials
K 6250KO1	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K 6250KO2	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K 6250K	CONJ	Conjugate, ready to use	15 ml
K 6250TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	15 ml
K 6250AC	STOP	ELISA stop solution, ready to use	15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker with 37 °C incubator
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA WASHBUF** (wash buffer concentrate) must be diluted with aqua bidist. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidist.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2–8°C for one month**.
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at **2–8°C** until the expiry date stated on the label. Reconstitution details are given in the data sheet.
- All other test reagents are ready to use. The test reagents are stable until the expiry date given on the label when stored at **2–8°C**.

6. SAMPLE PREPARATION

Serum and Plasma

Samples are diluted **1:1000** in **SAMPLEBUF** (sample dilution buffer) before use. E. g.:

Dilution I:

25 µl sample + 975 µl **SAMPLEBUF** = **1:40**

Dilution II:

40 µl of dilution I + 960 µl **SAMPLEBUF** = **1:25**

For analysis, pipette **100 µl** of **dilution II** per well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the two-site "sandwich" technique with one monoclonal and one polyclonal antibody that bind to human Adiponectin.

Standards, controls and diluted patient samples which are assayed for human Adiponectin are added to wells of microplate coated with a high affine monoclonal anti-human Adiponectin antibody. During the first incubation step, Adiponectin in the samples is bound by the immobilized antibody. Then a peroxidase labeled conjugate is added to each well and the following complex is formed: capture antibody - human Adiponectin – Peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the Adiponectin concentration of sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standard. Adiponectin present in the patient samples, is determined directly from this curve.

Test procedure

1. Bring all reagents and samples to room temperature (18-26 °C) and mix well
2. Mark the positions of STD /SAMPLE/CTRL (Standards/Sample/Controls) in duplicate on a protocol sheet
3. Take as many microtiter strips as needed from kit. Store unused strips covered at 2-8° C. Strips are stable until expiry date stated on the label
4. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
5. Add 100 µl of STD/SAMPLE/CTRL (Standard/Sample/Controls) in duplicate into respective well
6. Cover plate tightly and incubate for 1 hour at room temperature (18-26 °C) on a horizontal mixer
7. Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
8. Add 100 µl CONJ (conjugate) into each well

9.	Cover plate tightly and incubate for 1 hour at room temperature (18-26 °C) on a horizontal mixer
10.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
11.	Add 100 µl of SUB (substrate) into each well
12.	Incubate for 10–20 minutes at room temperature (18-26°C) in the dark*
13.	Add 50 µl of STOP (stop solution) into each well, mix thoroughly
14.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

* The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the „4-Parameter-algorithm“.

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Serum and plasma samples

To obtain the final Adiponectin concentration of the samples, multiply the estimated value by the dilution factor of **1000**.

9. LIMITATIONS

Stool samples with Adiponectin levels greater than the highest standard value, should be diluted with SAMPLEBUF (sample dilution buffer), and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal ranges

A significant correlation between adiponectin serum values and age as well as gender of the probands is reported, in turn the correlation between the respective BMI seems to be less significant.

Company	Adiponectin mean value* [µg/mL]		
	Adiponectin total - Kit	Normal glucose tolerance	Impaired glucose tolerance
LINCO		8.95 ± 0.55	3.38 ± 0.26
Mediagnost		8.81 ± 3.43	3.51 ± 1.47
			Type 2 diabetes
			3.48 ± 0.42
			3.82 ± 2.16

*Values from Reference 1.

We recommend each laboratory to establish its own norm concentration range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) was calculated from 30 replicate determinations on each one of two samples.

Intra-Assay (n = 30)		
Sample	Adiponectin [µg/mL]	CV [%]
1	5.25	3.2
2	5.92	3.4

Inter-Assay (n = 12)		
Sample	Adiponectin [µg/ml]	CV [%]
1	6.9	6.3
2	8.4	6.1

Recovery

Two samples were spiked with 3 different Adiponectin standards and measured using this assay.

n=2

Sample [ng/mL]	Spike [ng/mL]	Adiponectin expected [ng/mL]	Adiponectin measured [ng/mL]
3.8	2	5.8	6.1
3.8	4	7.8	7.6
3.8	6	9.8	9.2
2.1	5	7.1	8.1
2.1	7	9.1	9.6
2.1	9	11.1	10.8

Analytical Sensitivity

The sensitivity was set as B0 + 2 SD. The zero-standard was measured 20 times.

Sample	Adiponectin mean value [OD]	Standard variation [SD]	Detection limit [ng/mL]
1	0.006	0.001	0.338

Linearity

Two patient samples were diluted with wash buffer and analyzed. The results are shown below:

n= 2

Sample	Dilution	Expected [µg/ml]	Measured [µg/ml]
A	1 : 600	7.65	7.65
	1 : 1200	3.83	4.14
	1 : 2400	1.91	2.33
	1 : 4800	0.96	1.23
B	1 : 300	7.60	7.60
	1 : 600	3.80	3.87
	1 : 1200	1.90	2.11
	1 : 2400	0.95	1.17

12. PRECAUTIONS

- For in vitro diagnostic use only.
- Quality control guidelines should be observed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thi-

merosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- Stop solution contains sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

15. REFERENCES

1. Blüher (Blueher), M., et al. (2007) Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care* 30(2): p. 280-5
2. Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelties in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
3. Fernandez-Real, J.M., et al. (2004) Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care* 27(3): p. 739-45
4. Halperin F, Beckman JA, Patti ME, Trijillo ME, Garvin M, Craeger MA, Scherer PE and Goldfine AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154
5. Higashiura, K., et al. (2004) Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)* 61(6): p. 753-9
6. Koenig W, Khuseynova N, Baumert J, Meisinger C and Löwel H (2006) Serum Concentrations of Adiponectin and Risk of Type 2 Diabetes Mellitus and Coronary Heart Disease in Apparently Healthy Middle-Aged Men: Results From the 18-Year Follow-Up of a Large Cohort From Southern Germany. *JACC* 48:1369-1377
7. Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. *Diabetes* 55:249-259
8. Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. *Clin Chem* 50:1511-1525
9. Nakamura, Y., et al. (2004) Implications of plasma concentrations of adiponectin in patients with coronary artery disease. *Heart* 90(5): p. 528-33
10. Pajvani, U.B., et al. (2003) Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. *J Biol Chem* 278(11): p. 9073-85
11. Shibata, R., et al. (2004) Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. *J Biol Chem* 279(27): p. 28670-4

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49 (0) 62 51/70 19 00

Fax: +49 (0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com