α_2 -Makroglobulin ELISA Kit

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TABLE OF CONTENTS	PAGE
1. INTENDED USE	16
2. SUMMARY AND EXPLANATION OF THE TEST	16
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	17
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	18
7. PRECAUTIONS	18
8. SPECIMEN COLLECTION AND PREPARATION	19
9. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
10. RESULTS	21
11. LIMITATIONS	21
12. QUALITY CONTROL	22
EXPECTED VALUES	22
13. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
RECOVERY	23
SENSITIVITY	24
SAMPLE DILUTION	25
14. REFERENCES	26
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26

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Manual


α_2 -Macroglobulin ELISA Kit

For the in vitro determination of α_2 -Macroglobulin in urine,
serum and plasma

Valid from 31.03.2010

K 6610A

1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of **α 2-Macroglobulin** in urine, serum and plasma. For in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Alpha 2-Macroglobulin is one of the biggest plasma proteins with a molecular weight of 650-900kDa, dependent on the degree of glycosylation. It consists of 4 subunits. **Alpha 2-Macroglobulin** functions as protease inhibitor of all known classes of endopeptidase (serine-, cysteine-, aspartat- and metalloproteases).

The measurement of urinary proteins in the urine allows the diagnosis of proteinuria. In case of a microalbuminuria (selective proteinuria) increased levels (> 150 mg protein/24h) of albumin and transferrin (marker proteins of the tubulus) could be measured. In case of severe diseases of the kidney with pathological changes of the basal membran in the glomerulus (Chronic Glomerulonephritis, Tubulopathy, Nephrotic Syndrom) the protein level will rise up to > 3 g/24h and also the proteintyp will change to high molecular proteins like β -Lipoproteins and **α 2-Macroglobulin** (non-selective proteinuria).

Additionally previous investigations have shown that the determination of the two acute-phase proteins in the urins, C-reactive protein (CRP) and **α 2-Macroglobulin**, allows a noninvasive diagnosis of acute renal transplant dysfunction and rejection.

Indication

- Marker for glomerulonephropathy
- Monitoring of renal transplant dysfunction and rejection (together with CRP)

3. PRINCIPLE OF THE TEST

In a first incubation step, the **α 2-Macroglobulin** in the samples is bound to polyclonal rabbit antibodies (in excess), which are immobilized to the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labeled anti **α 2-Macroglobulin** (POD-Antibody) antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of **α 2-Macroglobulin** in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the results obtained from the calibrators. **α 2-Macroglobulin**, present in the patient samples, is determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K6610AMTP	PLATE	One holder with pre-coated strips	12 x 8 wells
K6610AWP	WASHBUF	ELISA wash buffer concentrate (10x)	1 x 100 ml
K6610AK	CONJ	Conjugate, (rabbit anti α -2-Macroglobulin peroxidase-labeled)	1 x 200 μ l
K6610APV	SAMPLEBUF	Sample dilution buffer, ready-to-use	2 x 100 ml
K6610AST	STD	Calibrators, lyophilized	2 x 6 vials
K6610AKO1	CTRL	Control, lyophilized	2 x 1 vial
K6610AKO2	CTRL	Control, lyophilized	2 x 1 vial
K6610ATMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K6610AAC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with Aqua dest. **1:10** before use (add 450 ml A. dest. to 50 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The reagent tubes with the lower amounts should be centrifuged before use, to guarantee that the required amount is available for use.
- The **Conjugate** (CONJ; POD antibody) must be diluted **1:101** in ELISA wash buffer (100 µl CONJ + 10 ml ELISA wash buffer). The undiluted antibody is stable at 2 -8 °C until expiry date given on the label. **Diluted antibody solution is not stable and cannot be stored.**
- The lyophilized **calibrators** (STD) and **controls** (CTRL) are stable at 2-8°C until the expiry date. The calibrators (STD) and controls (CTRL) must be reconstituted with **500 µl** of aqua bidest. Allow the vial to stand for 10 minutes and then mix thoroughly by gentle inversion to ensure complete reconstitution. Reconstituted calibrators and controls can be stored at **-20°C for 4 weeks and used once after thawing.**

7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop Solution consists of diluted Sulfuric Acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Plasma and serum

Plasma and sera are stable at 2-8°C for about 14 days. For long time storage we recommend -20°C:

Plasma and sera must be diluted **1:50 000** with SAMPLEBUF (sample dilution buffer). Dilution in three steps is recommended.

For example:

50 µl serum/plasma+ **950 µl** SAMPLEBUF, mix well (**1:20**)

20 µl of the 1:20 dilution + **980 µl** SAMPLEBUF, mix well (**1:1000**)

20 µl of the 1:1000 dilution + **980 µl** SAMPLEBUF, mix well (**1:50 000**)

For analysis, pipette **100 µl** of the **final dilution** per well.

Urine

Urine samples must be diluted before the assay **1:5** in SAMPLEBUF (sample dilution buffer).

For example:

100 µl urine+ **400 µl** SAMPLEBUF, mix well (**1:5**).

For analysis, pipette **100 µl** of the **diluted urine** per well.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

Wash the precoated microtiter plate **5 x with 250 μ l** ELISA wash buffer (WASHBUF). After the final washing step, remove residual buffer by tapping the plate on absorbent paper. Carry out the tests in duplicate.

1.	Add 100 μl of each STD (standard), CTRL (controls) and SAMPLE (patient samples) in the corresponding well.
2.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) shaking on a horizontal mixer or for 2 hours at room temperature (18-26°C) without any shaking .
3.	Decant the content of the plate and wash the wells 5 x with 250μl ELISA WASHBUF (wash buffer). After the final washing step, remove residual buffer by tapping the plate on absorbent paper.
4.	Add 100 μl of diluted CONJ (Conjugate).
5.	Incubate for 1 hour , shaking on a horizontal mixer, at room temperature.
6.	Decant the content of the plate and wash the wells 5 x with 250μl ELISA WASHBUF. After the final washing step, remove residual buffer by tapping the plate on absorbent paper.
7.	Add 100 μl of SUB (TMB substrate) solution.
8.	Incubate for 10-20 minutes at room temperature.
9.	Add 50 μl of STOP (stop solution) and mix shortly.
10.	Determine absorption with an ELISA reader at 450 nm against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-Parameter-algorithm".

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Plasma and serum samples

The estimated plasma and serum concentration must be multiplied by the dilution factor of **50 000**.

Urine samples

The estimated concentration must be multiplied by the dilution factor of **5**.

11. LIMITATIONS

Samples with α 2-Macroglobulin levels greater than the highest calibrator should be diluted in SAMPLEBUF (sample dilution buffer) and re-assayed.

12. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Plasma and serum: 1.3 – 3.0 g/l

Urine: < 0.18 mg/l; corresponds to 180 ng/ml

We recommend each laboratory to establish its own normal range.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay-Variation

The intra-assay variation was calculated from 20 determinations on each one of two urine or serum samples measured by one person.

Intra-Assay VK n= 20

Urine

Sample	α_2 -Macroglobulin mean value [ng/mL]	Intra-Assay CV [%]
1	163.98	2.53
2	84.17	3.27

Serum

Sample	α_2 -Macroglobulin mean value [mg/mL]	Intra-Assay CV [%]
1	3.92	2.61
2	2.18	6.37

Inter-Assay-Variation

The inter-assay variation was calculated from data on 2 urine or two serum samples obtained by different technicians in 10 different assays on different days.

Intra-Assay VK n= 10

Urine

Sample	α_2 -Macroglobulin mean value [ng/mL]	Intra-Assay CV [%]
1	153.46	3.67
2	83.65	3.43

Serum

Sample	α_2 -Macroglobulin mean value [mg/mL]	Intra-Assay CV [%]
1	4.01	3.31
2	2.00	3.99

Recovery

Two α_2 -Macroglobulin-containing urine or serum samples were spiked with 3 different α_2 -Macroglobulin-concentrations and measured.

Recovery n=3

Urine

Sample [ng/ml]	Spike [ng/ml]	α_2 -Macroglobulin expected [ng/ml]	α_2 -Macroglobulin measured [ng/ml]
19.1	200	219.1	199.3
19.1	100	119.1	118.4
19.1	50	69.1	64.02
34.7	200	234.7	223.8
34.7	100	134.7	118.6
34.7	50	84.7	80.46

Serum

Sample [ng/ml]	Spike [ng/ml]	α_2 -Macroglobulin expected [ng/ml]	α_2 -Macroglobulin measured [ng/ml]
17.3	200	217.3	192.0
17.3	100	117.3	107.0
17.3	50	67.3	64.0
31.8	200	231.8	227.5
31.8	100	131.8	123.3
31.8	50	81.8	80.4

Sensitivity

The sensitivity was set as $B_0 + 2 \text{ SD}$. The zero-standard was measured 20 times.

Sample	α_2 -Macroglobulin mean value [OD]	Standard variation	Detection limit [ng/mL]
1	0.005	0.0013	1.68

Sample dilution

n= 2

Two patient urine or serum samples were diluted with sample dilution buffer and measured with the assay. The results are shown below:

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15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent; these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number