α_1 -Antitrypsin ELISA Kit

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Hersteller

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
α_1 -Antitrypsin ELISA KIT

For the in vitro determination of α_1 -Antitrypsin in stool

Valid from 31.03.2011

K 6760

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of α_1 -Antitrypsin in stool. It is for in vitro diagnostic use only.

2. INTRODUCTION

The α_1 -antitrypsin acts as a primary inhibitor of elastase from polymorph nuclear neutrophilic granulocytes (PMN) and is released during inflammatory processes in order to reduce the proteolytic activity of the PMN-elastase in the inflammation region. In addition, it inhibits via complex formation a series of serine proteinases, like blood-clotting proteinases, trypsin, chymotrypsin, etc. Thus, the α_1 -antitrypsin plays an important regulatory as well as anti-inflammatory role.

The α_1 -antitrypsin is a linear glycoprotein with a molecular weight of ca. 52 kDa (394 amino acid residues), a free cysteine residue and three carbohydrate side chains. It is predominantly synthesized in the liver but also by intestinal macrophages, monocytes and epithelial cells. Although α_1 -antitrypsin is the main serine proteinase inhibitor in human plasma, the proof of fecal α_1 -antitrypsin has become an important marker for intestinal protein loss and permeability, as it is able to resist degradation in the gut due to its anti-proteolytic activity. It will therefore stay intact and it is possible to detect it in the faeces using an immunoassay. Moreover, the measurement of fecal α_1 -antitrypsin-concentration is used to evaluate and monitor chronic inflammatory intestinal diseases. In the clinical routine, the α_1 -antitrypsin-clearance (ratio of the α_1 -antitrypsin-ELISA-values of stool and serum samples) has been established along with the sole determination of the 24h- α_1 -antitrypsin-secretion in stool. Thus, the group of J. S. Fordtran (Strygler et al. 1990) reports, that the sole determination of the α_1 -antitrypsin-concentration in stool yielded false positive or false negative results in 21% of the patients compared to the α_1 -antitrypsin-clearance-measurements.

In a comparative study with the radial immunodiffusion (RID), routinely used in the clinical diagnostics, the α_1 -antitrypsin-ELISA developed by Immundiagnostik demonstrated significant advantages in the analysis of serum, stool and Caco-2-cell culture supernatants (Faust et al. 2001):

- 1) The α_1 -antitrypsin-concentrations obtained by the ELISA were on average about 30 % higher than the corresponding values from the radial immunodiffusion measurements.

- 2) Only the ELISA-system detected α_1 -antitrypsin in the cell culture supernatants.

The results clearly demonstrate that our α_1 -antitrypsin-ELISA-test is more sensitive than other routinely used methods and that it **recognizes the hepatic as well as the enteral α_1 -antitrypsin form**. This newly developed test represents a promising alternative to the use of current clinical routine methods. It is superior to the radial immunodiffusion especially in cases with extremely high protein loss through the gut. The combination of two specific antibodies eliminates to a large extent the possibility of false negative results guaranteeing reliable diagnoses. This newly developed test is a non-invasive, simple test for the detection of protein loss through the gut.

Indications

- Enteric protein loss and intestinal permeability
- Morbus Crohn
- Necrotic enterocolitis
- Viral, bacterial or allergic inflammation

3. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 6760MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6760WB	WASHBUF	ELISA wash concentrate 10x	2 x 100 ml
K 6760K	CONJ	Conjugate, (sheep anti α_1 -antitrypsin, peroxidase-labelled)	200 μ l
K 6760CAL	CAL	Calibrator, ready-to-use	1 vial
K 6760KO1	CTRL	Control, ready-to-use	1 vial
K 6760KO2	CTRL	Control, ready-to-use	1 vial
K 6760TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 6760AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

The used calibrator was calibrated on the WHO reference material CRM 470.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μ l** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C using a water bath before dilution. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- Use **100 μ l of diluted wash buffer as a BLANK**. Pipette into the respective well.
- The **CAL** (calibrator) and the **CTRL** (controls) are stable at **2-8°C** until the expiry date stated on the label.
- The **CONJ** (conjugate) must be diluted **1: 100** in wash buffer (100 μ l CONJ + 9900 μ l wash buffer). The undiluted conjugate is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

6. SAMPLE PREPARATION

Extraction of the stool sample

The **wash buffer** is used as a sample extraction buffer. We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample has to be suspended in 5 ml buffer.

Constant buffer volume: 5 ml

Constant dilution factor: 1:50

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample not depending on the sample amount.

Constant buffer volume: 5 ml

The dilution factor varies depending on the sample amount which has to be considered in the subsequent calculations as shown below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor remains constant.

Variable buffer volume

Constant dilution factor: 1:50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results.

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix stool sample and buffer; vortex for at least 30 sec. depending on the stool consistency.

Transfer approx. 1 ml stool suspension to an Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm (\approx 13000 g).

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Dilution of samples

Stool samples

The supernatant is centrifuged at 13.000 rpm for 2 min. and diluted **1:500 in wash buffer**.

For example:

40 μ l supernatant + **960 μ l** wash buffer, mix well (**dilution I**) (1:25)

50 μ l of this dilution I + **950 μ l** wash buffer, mix well (**dilution II**) (1:20)

For analysis, pipette **100 μ l** of the supernatant of **dilution step II** per well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the "sandwich" technique with two selected polyclonal antibodies that bind to human α -1-antitrypsin.

Calibrator, controls and prediluted samples which are assayed for human α -1-antitrypsin are added into the wells of a micro plate coated with a high affine polyclonal anti-human α -1-antitrypsin antibody. During the first incubation step, α -1-antitrypsin is bound by the immobilized antibody. Then a peroxidase-conjugated polyclonal anti-human α -1-antitrypsin antibody is added into each microtiter well and a "sandwich" of capture antibody - human α -1-antitrypsin – peroxidase-conjugate is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of α -1-antitrypsin. Samples are quantified by referring their optical density to a Lot-dependant master calibration curve and the use of a Calibrator that is run with each test.

Test procedure

1.	Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well
2.	Mark the positions of BLANK/CAL/CTRL/SAMPLE (Blank / Calibrator / Controls/ Sample) in duplicate on a protocol sheet
3.	Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label
4.	Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
5.	Add 100 µl of BLANK/CAL/CTRL/SAMPLE (Blank / Calibrator / Controls / Sample) in duplicate into respective well
6.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer
7.	Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
8.	Add 100 µl CONJ (conjugate) into each well
9.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer

- | |
|---|
| 10. Discard the contents of each well. Wash 5 times by dispensing 250 μl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper |
| 11. Add 100 μl of SUB (substrate) into each well |
| 12. Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark* |
| 13. Add 50 μl of STOP (stop solution) into each well, mix thoroughly |
| 14. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference |

*The intensity of the colour change is temperature sensitive. We recommend to observe the colour change and to stop the reaction upon good differentiation.

8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the Calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A, B, C and D. In both cases the optical density of the calibrator (CAL) is essential.

Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Caution: Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of replicate values should always be examined after automatic evaluation of the results. If this option is not available with the used program, a control of the replicate values should be done manually by the operator.

Stool samples

To obtain the α_1 -antitrypsin concentration in stool samples, use the corresponding dilution factor according to the sample preparation:

For the α_1 -Antitrypsin concentration of fecal samples, calculate as described in the following example:

weight: 80 mg (1ml stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 500

dilution factor: 62,5 x 500 = 31250

Multiply the result with **31250** to get the real concentration. **The dilution factor depends on the weight of the faeces.**

9. LIMITATIONS

Stool samples with α -1-antitrypsin concentration greater than the highest standard value should be further diluted with wash buffer and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Cut off-value for healthy people (Stool): < 26,8 mg/dl

The cut off-value should be used as a guideline only. It is recommended that each laboratory establishes its own cut off-value for its patient population.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay-Variation

The precision (intra-assay variation) of the Immundiagnostik α -1-Antitrypsin ELISA test was calculated from 20 replicate determinations on each one of two samples.

Intra-Assay CV n= 20

Sample	α -1-Antitrypsin mean value [mg/dl]	Intra-Assay CV [%]
1	15,2	4,5
2	42,2	13,1

Inter-Assay-Variation

The total precision (inter-assay variation) of the Immundiagnostik α_1 -Antitrypsin ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

Sample	α_1 -Antitrypsin mean value [mg/dl]	Inter-Assay CV [%]
1	16,15	9,8
2	54,46	14,8

Recovery

Two samples were spiked with different α_1 -antitrypsin calibrator amounts and measured with the assay.

Recovery n=4

Sample [μ g/l]	Spike [μ g/l]	α_1 -Antitrypsin expected [μ g/l]	α_1 -Antitrypsin measured [μ g/l]
6.3	1.65	7.95	7.2
6.3	5	11.3	11.4
6.3	15	21.3	20.9
6.3	45	51.3	46.7
5.6	1.65	7.3	7.0
5.6	5	10.6	10.9
5.6	15	20.6	17.5
5.6	45	50.6	46.7

Sensitivity

The sensitivity limit was set as $B_0 + 2 \text{ SD}$. The Zero-standard was measured 20 times.

n=20

Sample	α -1-Antitrypsin mean value [OD]	Standard variation	Detection limit [mg/dl]
1	0,117	0,015	1,8

Linearity

Two patient serum samples were diluted with wash buffer. The results are shown below:

n= 2

Sample	Dilution	Expected [mg/dl]	Measured [mg/dl]
A	1:12500	48	48.88
	1:25000	24.5	23.25
	1:50000	12.3	11.9
	1:100000	6.1	6
B	1:12500	158.4	158.4
	1:25000	79.3	99
	1:50000	39.6	33
	1:100000	19.8	22.1

Cross reactivity

No cross reactivity with other plasma proteins in stool was observed.

No cross reactivity with alpha-1-antitrypsin in mouse serum was observed.

12. PRECAUTIONS

- For in vitro diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

15. REFERENCES

1. Arndt et al., 1993, Crohn; Klin. Lab. **11**: 867-876
2. Stein, J., 1996, 3. Post-graduertenkurs der DGVS
3. Arndt et al., 1992, The Lancet, Vol **340**: OCT, 24
4. Karbach et al., 1989, Gastroenterol. **27**: 362
5. Strygler B et al., 1990, Gastroenterol. **99**(5):1380-7
6. Faust D et al., 2001, Gastroenterol. **39**(9):769-74
7. Faust D et al., 2002, Clin Exp Immunol. **128**(2):279-84

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number

For research use only