# S100A12 ELISA Kit

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**Manual**


# S100A12 ELISA Kit

*For the in vitro determination of S100A12 in stool*

Valid from 08.11.2010

**K 6938**

## 1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) Kit is intended for the quantitative determination of S100A12 in stool. It is for *in vitro* diagnostic use only.

## 2. INTRODUCTION

Alternative names: Calgranulin C, EN-RAGE

Members of the S100 protein family are low molecular mass acidic proteins characterized by cell-type-specific expression and the presence of 2 EF-hand calcium-binding domains. The S100 proteins or calgranulins are expressed in activated granulocytes under conditions of chronic inflammation. Hofmann et al. (1999) reported that RAGE is a central cell surface receptor for S100A12, which they referred to as EN-RAGE (Extracellular Newly identified RAGE-binding protein), and related members of the S100/calgranulin super family. Interaction of S100A12 with cellular RAGE on endothelium, mononuclear phagocytes, and lymphocytes triggered cellular activation, with generation of key pro-inflammatory mediators. In murine models, blockade of S100A12 quenched delayed-type hypersensitivity and inflammatory colitis by arresting activation of central signaling pathways and expression of inflammatory mediators.

Kaiser et al. (2007) determined the S100A12 concentration in stool samples of 171 patients suffering from infectious gastroenteritis, Crohn's disease, ulcerative colitis or irritable bowel syndrome and 24 healthy controls. Faecal S100A12 correlated with the degree of intestinal inflammation. Patients with inflammatory bowel disease had significantly higher S100A12 levels compared to these with irritable bowel syndrome or healthy controls. The S100A12 content in stool correlated with the degree of intestinal inflammation.

Faecal S100A12 is a novel non-invasive marker distinguishing inflammatory bowel disease (IBD) from irritable bowel syndrome (IBS) or healthy individuals with a high sensitivity and specificity. S100A12 stool concentrations can be used as an indicator of inflammation activity in Morbus Crohn's disease and Colitis ulcerosa. Since S100A12 is secreted by activated neutrophils and acts as a pro-inflammatory mediator, it can be used as marker for diagnosis and therapeutic control of inflammatory bowel diseases.

### Indication

- Marker for acute inflammation
- Estimation of gastrointestinal inflammation degree
- Parameter for monitoring Morbus Crohn's disease and Colitis ulcerosa

### 3. MATERIAL SUPPLIED

| Cat. No   | Content   | Kit Components                                                              | Quantity     |
|-----------|-----------|-----------------------------------------------------------------------------|--------------|
| K 6938MTP | PLATE     | One holder with precoated strips                                            | 12 x 8 wells |
| K 6938WP  | WASHBUF   | ELISA wash buffer concentrate 10x                                           | 1 x 100 ml   |
| K 6938EP  | EXBUF     | Extraction buffer concentrate 2,5x                                          | 2 x 100 ml   |
| K 6938PV  | SAMPLEBUF | Sample dilution buffer, ready-to-use                                        | 1 x 100 ml   |
| K 6938ST  | STD       | S100A12 standards, lyophilized<br>(0; 0.22; 0.66; 2.0; 6.0; 18.0; 54 ng/ml) | 2 x 7 vials  |
| K 6938KO1 | CTRL      | Control, lyophilized<br>(see specification for range)                       | 2 x 1 vial   |
| K 6938KO2 | CTRL      | Control, lyophilized<br>(see specification for range)                       | 2 x 1 vial   |
| K 6938K   | CONJ      | Conjugate                                                                   | 200 µl       |
| K 6938TMB | SUB       | TMB substrate (Tetramethylbenzidine),<br>ready to use                       | 15 ml        |
| K 6938AC  | STOP      | ELISA stop solution, ready to use                                           | 15 ml        |

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 13000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm  
(reference wave length 620 or 690 nm)

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) must be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml a. bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. **The buffer concentrate is stable at 2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask **at 2-8°C for one month.**
- The **EXBUF** (extraction buffer concentrate) must be diluted with aqua bidist. **1:2.5** before use (100 ml EXBUF + 150 ml aqua bidist.), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37°C in a water bath. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for three months.**
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable **at 2-8°C** until the expiry date stated on the label. The **STD** (standards) and **CTRL** (controls) must be reconstituted with **aqua bidest. (reconstitution volume, see product specification)**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards are not stable and can not be stored.**
- The **CONJ** (conjugate) must be diluted **1:100** in diluted wash buffer (100 µl CONJ + 9.9 ml diluted wash buffer). **Undiluted conjugate** is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date given on the label when stored at **2-8°C**.

## 6. SAMPLE PREPARATION

### *Extraction of the stool sample*

We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing **100 mg of stool sample** (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 10 745804 332). The stool sample must be suspended in **5 ml EXBUF** (extraction buffer).

**Constant buffer volume: 5 ml**

**Constant dilution factor: 1:50**

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!
  - a. Add **5 ml EXBUF** to the stool sample not depending on the sample amount.

**Constant buffer volume: 5 ml**

**The dilution factor varies** depending on the sample amount which has to be considered in the subsequent calculations. Please have a look at the table for the correction factors below:

| Weight [mg] | Dilution factor |
|-------------|-----------------|
| 80          | 62.5            |
| 82          | 60.9            |
| 84          | 59.5            |
| 86          | 58.1            |
| 88          | 56.8            |
| 90          | 55.6            |
| 92          | 54.3            |
| 94          | 53.2            |
| 96          | 52.1            |
| 98          | 51.0            |
| <b>100</b>  | <b>50</b>       |

| Weight [mg] | Dilution factor |
|-------------|-----------------|
| 102         | 49.0            |
| 104         | 48.1            |
| 106         | 47.2            |
| 108         | 46.3            |
| 110         | 45.5            |
| 112         | 44.6            |
| 114         | 43.9            |
| 116         | 43.1            |
| 118         | 42.4            |
| 120         | 41.6            |

b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor remains constant.

### Variable buffer volume

#### Constant dilution factor: 1:50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results .

| Weight [mg] | Buffer Volume [ml] |
|-------------|--------------------|
| 80          | 4.0                |
| 82          | 4.1                |
| 84          | 4.2                |
| 86          | 4.3                |
| 88          | 4.4                |
| 90          | 4.5                |
| 92          | 4.6                |
| 94          | 4.7                |
| 96          | 4.8                |
| 98          | 4.9                |
| <b>100</b>  | <b>5.0</b>         |

| Weight [mg] | Buffer Volume [ml] |
|-------------|--------------------|
| 102         | 5.1                |
| 104         | 5.2                |
| 106         | 5.3                |
| 108         | 5.4                |
| 110         | 5.5                |
| 112         | 5.6                |
| 114         | 5.7                |
| 116         | 5.8                |
| 118         | 5.9                |
| 120         | 6.0                |

After  
buffer  
deposition  
consistency.

Transfer ca. 1 ml stool suspension (**dilution step I**) to an Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm (= 13000 g).

## Dilution of samples

### Stool samples

The supernatant of the extraction (dilution step I) is diluted **1:50 with SAMPLEBUF** (Sample dilution buffer). For example:

20 µl supernatant (dilution I) + 980 µl SAMPLEBUF = 1:50 (**dilution step II**)

### Final dilution: 1:2500

For analysis, pipette **100 µl** of the supernatant of **dilution step II** per well.

We recommend to store the extracted samples (dilution step I) for no longer than 24 h at 2-8 °C. The extracted samples are not stable and cannot be frozen and stored for later use.

## **7. ASSAY PROCEDURE**

### *Principle of the test*

The S100A12 ELISA Kit employs the quantitative sandwich enzyme immunoassay technique.

A monoclonal antibody specific for S100A12 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and the immobilized antibody binds any S100A12 present. After washing away any unbound substances, an HRP-conjugated polyclonal antibody specific for S100A12 is added to the wells. Following a wash to remove any unbound antibody HRP-conjugate, the remaining conjugate is allowed to react with the HRP substrate tetramethylbenzidine. The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured at 450 nm. The absorbance is proportional to the concentration of S100A12. A standard curve is constructed by plotting absorbance values versus S100A12 concentrations of calibrators, and concentrations of unknown samples are determined using this standard curve.

### Test procedure

|    |                                                                                                                                                                                                                      |
|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. | Bring all <b>reagents and samples to room temperature</b> (18-26 °C) and mix well                                                                                                                                    |
| 2. | Mark the <b>positions of STD /SAMPLE/CTRL</b> (Standards/Sample/Controls) in duplicate on a protocol sheet                                                                                                           |
| 3. | <b>Take</b> as many <b>microtiter strips as needed</b> from kit. Store unused strips covered at 2-8° C. Strips are stable until expiry date stated on the label                                                      |
| 4. | Wash each well <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper                                     |
| 5. | Add <b>100 µl of STD/SAMPLE/CTRL</b> (Standards/Sample/Controls) in duplicate into respective well                                                                                                                   |
| 6. | Cover plate tightly and <b>incubate for 1 hour at room temperature</b> (18-26 °C) on a horizontal mixer                                                                                                              |
| 7. | Aspirate the contents of each well. Wash each well <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper |
| 8. | Add <b>100 µl CONJ</b> (conjugate) into each well                                                                                                                                                                    |
| 9. | Cover plate tightly and <b>incubate for 1 hour at room temperature</b> (18-26 °C) on a horizontal mixer                                                                                                              |

- |                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 10. Aspirate the contents of each well. Wash each well <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper                                                                                                                                  |
| 11. Add <b>100 µl of SUB</b> (substrate) into each well                                                                                                                                                                                                                                                                                                   |
| 12. Incubate for <b>10 - 20 minutes at room temperature</b> (18-26°C) in the dark*                                                                                                                                                                                                                                                                        |
| 13. Add <b>50 µl of STOP</b> (stop solution) into each well, mix thoroughly                                                                                                                                                                                                                                                                               |
| 14. Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference |

\*The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

### 1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

### 2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

## Stool samples

To obtain the S100A12 concentration in stool samples, multiply the estimated value by the dilution factor according to the sample preparation:

**Sample preparation 1 and 2b:      dilution factor constant :      1 : 50**

Multiply the obtained result by **2500** (dilution step I x dilution step II) to get the final concentration.

**Sample preparation 2a:                      dilution factor is variable**

The corresponding dilution factor is taken from the table and should be multiplied by 50 (dilution step II) to get the final concentration.

## 9. LIMITATIONS

Stool samples with S100A12 levels greater than the highest standard value, should be diluted with SAMPLEBUF, and be re-assayed.

## 10. QUALITY CONTROL

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

We recommend each laboratory to establish its own norm concentration range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

| Intra-Assay (n=16) |                 |        |
|--------------------|-----------------|--------|
| Sample             | S100A12 [ng/ml] | VK [%] |
| 1                  | 3718.4          | 7.1    |
| 2                  | 918.3           | 3.5    |

| Inter-Assay (n=22) |                 |        |
|--------------------|-----------------|--------|
| Sample             | S100A12 [ng/ml] | VK [%] |
| 1                  | 3003.64         | 15.42  |
| 2                  | 1363.49         | 20.18  |

### *Recovery*

Two samples were spiked with 2 different S100A12 standards and measured using this assay.

n=2

| Sample | Spike [ng/ml] | S100A12 expected [ng/ml] | S100A12 measured [ng/ml] |
|--------|---------------|--------------------------|--------------------------|
| 1      | 0.56          | 0.892                    | 0.803                    |
| 1      | 0.28          | 0.602                    | 0.56                     |
| 2      | 1.125         | 1.642                    | 1.701                    |
| 2      | 0.56          | 1.077                    | 0.98                     |

## Sensitivity

The sensitivity was set as  $B_0 + 1SD$ . The zero-standard was measured 10 times.

| Sample | S100A12 mean value [OD] | Standard variation (SD) | Detection limit [ng/ml] |
|--------|-------------------------|-------------------------|-------------------------|
| 1      | 0.011                   | 0.0008                  | 0.016*                  |

\* Corresponds to a concentration of 39.6 ng/l at a stool sample dilution of 1 : 2500.

## Linearity

Two patient samples were diluted with wash buffer and analyzed. The results are shown below:

n= 2

| Sample | Dilution | Expected [ng/ml] | Measured [ng/ml] |
|--------|----------|------------------|------------------|
| P1     | 1:2500   | 8797.3           | 8797.3           |
|        | 1:5000   | 4399             | 4239.8           |
|        | 1:10000  | 2199.5           | 1986             |
| P2     | 1:2500   | 2016.1           | 2016.1           |
|        | 1:5000   | 1008             | 932.2            |
|        | 1:10000  | 504              | 362              |

## 12. REFERENCES

1. Hofmann MA, Drury S; Fu C, Qu W, Taguchi A, Lu Y, Avila C, Kambham N, Bierhaus A, Nawroth P, Neurath MF, Slattery T, Beach, D, McClary J, Nagashima M, Morser J, Stern D, Schmidt AM. RAGE mediates a novel proinflammatory axis: a central cell surface receptor for S100/calgranulin polypeptides. *Cell*. **97**: 889-901, 1999
2. Karl J, Wild N, Tacke M, Andres H, Garczarek U, Rollinger W, Zolg W. Improved diagnosis of colorectal cancer using a combination of fecal occult blood and novel fecal protein markers. *Clin Gastroenterol Hepatol*. 2008 Oct;6(10):1122-8.
3. Foell D, Wittkowski H, Ren Z, Turton J, Pang G, Daebritz J, Ehrchen J, Heidemann J, Borody T, Roth J, Clancy R. Phagocyte-specific S100 proteins are released from affected mucosa and promote immune responses during inflammatory bowel disease. *J Pathol*. 2008 Oct;216(2):183-92.
4. Sidler MA, Leach ST, Day AS. Fecal S100A12 and fecal calprotectin as noninvasive markers for inflammatory bowel disease in children. *Inflamm Bowel Dis*. 2008 Mar;14(3):359-66.
5. Kaiser T, Langhorst J, Wittkowski H, Becker K, Friedrich AW, Rueffer A, Dobos GJ, Roth J, Foell D. Faecal S100A12 as a non-invasive marker distinguishing inflammatory bowel disease from irritable bowel syndrome. *Gut* 2007;56:1706-1713
6. de Jong Naomi SH, Leach ST, Day AS. Fecal S100A12: A Novel Noninvasive Marker in Children with Crohn`s Disease. *Inflamm Bowel Dis*. July 2006;12(7):566–572

## 13. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Quality control guidelines should be observed.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

## 14. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number