PhiCal® Calprotectin ELISA Kit

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Hersteller

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
PhiCal[®] Calprotectin ELISA Kit

***For the in vitro determination of PhiCal[®] Calprotectin
(MRP 8/14) in stool***

Valid from 02.08.2011

K 6967

[PhiCal[®]: registered German trademark of Immundiagnostik AG]

[Not sold in the USA and Norway]

Content

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	19
6. SAMPLE PREPARATION	20
<i>Extraction of the stool sample</i>	20
<i>Dilution of samples</i>	21
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Test procedure</i>	22
8. RESULTS	23
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Expected values</i>	23
<i>Precision and reproducibility</i>	24
<i>Recovery</i>	24
<i>Analytical Sensitivity</i>	25
<i>Linearity</i>	25
<i>Cross-reactivity</i>	25
12. PRECAUTIONS	26
13. TECHNICAL HINTS	26
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	27

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) Kit is intended for the quantitative determination of **PhiCal®** Calprotectin (MRP 8/14) in stool. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Alternative names:

Calgranulin A: MRP8, S100A8, CP-10 (in mouse)

Calgranulin B: MRP14, S100A9,

MRP8/14: L1, (p8,14), p34

Calprotectin is a calcium-binding protein secreted predominantly by neutrophils and monocytes. Fecal Calprotectin is a marker for neoplastic and inflammatory gastrointestinal diseases. It is often difficult to distinguish between irritable bowel syndrome and chronic inflammatory bowel disease. This leads in many cases to extensive and unnecessary colonoscopic examinations. The Calprotectin test allows clear differentiation between the two patient groups. Fecal Calprotectin levels correlate significantly with histologic and endoscopic assessment of disease activity in *Morbus Crohn's* disease and ulcerative colitis as well as with the fecal excretion of indium-111-labelled neutrophilic granulocytes that has been suggested as the "gold standard" of disease activity in inflammatory bowel disease. However, measuring 111-indium-labeled granulocytes is very costly (patient's hospitalization, analysis and disposal of isotopic material) and is connected with radioactive exposition of the patients. For this reason, a repeated application to children and pregnant women is not recommended.

Elevated levels of Calprotectin are a much better predictor of relapse than standard inflammatory markers (CRP, ESR HB). Comparing this marker with standard fecal occult blood screening in colorectal cancer demonstrates clearly the diagnostic advantages of the fecal Calprotectin test. The parameter is of a high diagnostic value: if the Calprotectin level in stool is low, the probability is high that no organic intestinal disease exists.

Indications

- Marker for acute inflammation
- Estimation of gastrointestinal inflammation degree
- Parameter for monitoring *Morbus Crohn's* disease, *Colitis ulcerosa* or the patient's status after removal of polyps.
- Discrimination between patients with inflammatory bowel disease (acute *Morbus Crohn's* disease and ulcerative colitis) and irritable bowel syndrome when using a fecal test system.

3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6967MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6967WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6967EP	EXBUF	Extraction buffer concentrate 2.5x	2 x 90 ml
K 6937PV	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 100 ml
K 6967CAL	CAL	Calprotectin calibrator, lyophilized	2 vials
K 6967KO1	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K 6967KO2	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K 6967K	CONJ	Conjugate, ready to use	15 ml
K 6967TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	15 ml
K 6967AC	STOP	ELISA stop solution, ready to use	15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25°C (≤18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA WASHBUF** (wash buffer concentrate) must be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be re-dissolved at 37°C in a water bath before dilution. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. Diluted **buffer solution (wash buffer)** can be stored in a closed flask at **2–8°C for one month**.
- Use **100 µl of SAMPLEBUF** (sample dilution buffer) **as a BLANK**. Pipette into the respective well.
- The **EXBUF** (extraction buffer concentrate) must be diluted with ultra pure water **1:2.5** before use (90 ml EXBUF + 135 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37°C in a water bath. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2–8°C for three months**.
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at **2–8°C** until the expiry date stated on the label. The **STD** (standards) and **CTRL** (controls) must be re-constituted with **500 µl** of ultra pure water. After the reconstitution, vortex the vial for 4 seconds and allow the pellet to dissolve by a subsequent incubation for 10 minutes to insure complete reconstitution. **Reconstituted standards and controls can be stored at 2-8°C for four weeks**.
- All other test reagents are ready to use. The test reagents are stable until the expiry date given on the label when stored at **2–8°C**.

6. SAMPLE PREPARATION

Extraction of the stool sample

Diluted extraction buffer is used as a sample extraction buffer. We recommend the following sample preparation:

1a. Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

SAS with 0.75 ml Buffer:

Applied amount of stool:	15 mg
Buffer Volume:	0.75 ml
Dilution Factor:	1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw Stool Sample has to be thawed. For remarkably inhomogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube with 0.75 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off and leave 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.
- e) Allow sample to stand for app. 10 minutes until sediment has settled down. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure, the sediment will not be dispersed again.

1b. Sample preparation kit from Roche Diagnostics, Mannheim, Germany (Cat. No. 10 745 804 322)

Alternatively, other stool sample preparation kits (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany) can be used. In the Roche sample preparation kit, 100 mg of stool sample are suspended in 5 ml of extraction buffer using a vibrator mixer (e.g. Vortex mixer). Centrifugation of the suspension is recommended.

Dilution I (1a. or 1b.)**1:50***Dilution of samples***Stool samples**

The supernatant of the extraction (dilution step I) is diluted **1:50 with SAMPLEBUF** (sample dilution buffer). For example:

20 µl supernatant (dilution I) + 980 µl **SAMPLEBUF** = 1:50 (**dilution step II**)

For analysis, pipette **100 µl** of the supernatant of **dilution step II** per well.

Calprotectin in stool is described to be stable for approximately 6 days. Nevertheless, we recommend to store the samples for not more than 48 h at 2-8 °C. Long term storage is recommended at -20 °C. Allow frozen samples to thaw slowly, preferably at 2-8° C over night and warm the samples to room temperature before analysis.

(Poullis A et al. (2002) Aliment Pharmacol Thr 16:675-681)

7. ASSAY PROCEDURE*Principle of the test*

The assay utilizes the two-site "sandwich" technique with two selected monoclonal antibodies that bind to human Calprotectin.

Calibrator, controls and diluted patient samples which are assayed for human Calprotectin are added to wells of microplate coated with a high affine monoclonal anti-human Calprotectin antibody. During the first incubation step, Calprotectin in the samples is bound by the immobilized antibody. Then a peroxidase labeled conjugate is added to each well and the following complex is formed: capture antibody - human Calprotectin - Peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the Calprotectin concentration of sample. Samples are quantified by referring their optical density to a Lot-dependant master calibration curve and the use of a Calibrator that is run with each test.

Test procedure

1.	Bring all reagents and samples to room temperature (18-26 °C) and mix well
2.	Mark the positions of BLANK/CAL/CTRL/SAMPLE (Blank/Calibrator/Controls/ Sample) in duplicate on a protocol sheet t
3.	Take as many microtiter strips as needed from kit. Store unused strips covered at 2-8° C. Strips are stable until expiry date stated on the label
4.	Add 100 µl of BLANK/CAL/CTRL/SAMPLE (Blank/Calibrator/Controls/ Sample) in duplicate into respective well
5.	Cover plate tightly and incubate for 30 minutes at room temperature (18-26°C)
6.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
7.	Add 100 µl CONJ (conjugate) into each well
8.	Cover plate tightly and incubate for 30 minutes at room temperature (18-26°C)
9.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper
10.	Add 100 µl of SUB (substrate) into each well
11.	Incubate for 10–20 minutes at room temperature (18-26°C) in the dark*
12.	Add 50 µl of STOP (stop solution) into each well, mix thoroughly
13.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used

* The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.

8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the Calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A,B,C and D. In both cases the optical density of the calibrator (CAL) is essential.

Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Caution: Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of replicate values should always be examined after automatic evaluation of the results. If this option is not available with the used program, a control of the replicate values should be done manually by the operator.

Stool samples

To obtain the calprotectin concentration in the stool samples, multiply the obtained result by **2500** (dilution step I x dilution step II).

9. LIMITATIONS

Stool samples with Calprotectin levels greater than the highest standard value, should be diluted with SAMPLEBUF (sample dilution buffer), and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal ranges

1 g stool is equivalent to 1ml

- **The median value in healthy adults is about 25 mg/kg.**
- **Samples giving values above 50 mg/kg are regarded as positive.**

We recommend each laboratory to establish its own norm concentration range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)		
Sample	Calprotectin [ng/mL]	CV [%]
1	89.2	5.6
2	229.1	3.2

Inter-Assay (n = 12)		
Sample	Calprotectin [ng/mL]	CV [%]
1	107.8	4.4
2	476.1	8.9

Recovery

Two samples were spiked with different Calprotectin standards and measured using this assay. n=2

Sample [ng/mL]	Spike [ng/mL]	Calprotectin expected [ng/mL]	Calprotectin measured [ng/mL]
1	10.5	28.5	29.6
1	17.5	35.3	35.6
1	40.5	57.9	59.3
1	63.3	79.9	83.2
1	173.2	188.1	188.7
2	10.5	31.2	33.3
2	17.5	38.1	41.0
2	40.5	60.6	62.7
2	63.3	82.4	90.8

Analytical Sensitivity

The sensitivity was set as $B0 + 6 \cdot SD$. The zero-standard was measured 22 times.

Sample	Calprotectin mean value [OD]	Standard variation [SD]	Detection limit [ng/mL]
1	0.021	0.00188	2.099

* Corresponds to 5.2475 mg/L at 1:2500 dilution of the stool sample.

Linearity

Two patient samples were diluted with wash buffer and analyzed. The results are shown below:

n= 2

Sample	Dilution	Expected [ng/mL]	Measured [ng/mL]
A	2500	120	120
	5000	60	61,8
	10000	30	32,4
B	2500	116	116
	5000	58	59,4
	10000	29	28,9

Cross-reactivity

No cross-reactivity was observed to the following plasma proteins:

Lysozyme - 0%

PMN-Elastase - 0%

Myeloperoxidase - 0%

Lactoferrin - 0%

12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Quality control guidelines should be observed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution contains sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

15. REFERENCES

1. Poullis A et al. (2002) *Aliment Pharmacol Ther* 16:675-681
2. Fagerhol et al. (2000) *The Lancet* 356:1783-1784
3. Tibble et al. (2000) *Gut* 47:506-513
4. Tibble et al. (2000) *Gastro* 119:15-22

Publications based on Immundiagnostik's calprotectin ELISA

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Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49 (0) 62 51/70 19 00

Fax: +49 (0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com

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[PhiCal®: registered German trademark of Immundiagnostik AG]